


Validation Data Sheet

| dPCR CNV Probe Assay | FRYL |
|------------------------|---------------------------------|
| Assay Catalog No. | DCC0000368 |
| dPCR Instrument | QIAcuity digital PCR instrument |
| dPCR Master Mix | QIAcuity Probe PCR Master Mix |
| Experimental Sample(s) | genomic DNA, XbaI-digested |

1D Scatter Plot: Single-well data of target assay FRYL (Green Channel) and reference assay TERT (Yellow Channel) in different gDNA input amounts (A2,B2,C2: 6ng/rxn, D2,E2,F2: 4ng/rxn, G2,H2: NTCs)



2D Scatter Plot: Single-well data of target assay FRYL (Green Channel, X-axis) and reference assay TERT (Yellow Channel, Y-axis) in 4ng gDNA input amount.

