



Validation Data Sheet

dPCR CNV Probe Assay	PLEKHF1
Assay Catalog No.	DCC0000430
dPCR Instrument	QIAcuity digital PCR instrument
dPCR Master Mix	QIAcuity Probe PCR Master Mix
Experimental Sample(s)	genomic DNA, XbaI-digested

1D Scatter Plot: Single-well data of target assay PLEKHF1 (Green Channel) and reference assay TERT (Yellow Channel) in different gDNA input amounts (A8,B8,C8: 6ng/rxn, D8,E8,F8: 4ng/rxn, G8,H8: NTCs)



2D Scatter Plot: Single-well data of target assay PLEKHF1 (Green Channel, X-axis) and reference assay TERT (Yellow Channel, Y-axis) in 4ng gDNA input amount.

