



# Validation Data Sheet

<b>dPCR CNV Probe Assay</b>	<b>ANKRD30B</b>
Assay Catalog No.	DCC0000510
dPCR Instrument	QIAcuity digital PCR instrument
dPCR Master Mix	QIAcuity Probe PCR Master Mix
Experimental Sample(s)	genomic DNA, XbaI-digested

**1D Scatter Plot:** Single-well data of target assay SMN2 (Green Channel) and reference assay TERT (Yellow Channel) in different gDNA input amounts (A8,B8,C8: 6ng/rxn, D8,E8,F8: 4ng/rxn, G8,H8: NTCs)



**2D Scatter Plot:** Single-well data of target assay SMN2 (Green Channel, X-axis) and reference assay TERT (Yellow Channel, Y-axis) in 4ng gDNA input amount.

