

# Validation Data Sheet

<b>dPCR CNV Probe Assay</b>	<b>NFATC2</b>
Assay Catalog No.	DCG0000011
dPCR Instrument	QIAcuity digital PCR instrument
dPCR Master Mix	QIAcuity Probe PCR Master Mix
Experimental Sample(s)	genomic DNA, XbaI-digested

**1D Scatter Plot:** Single-well data of target assay NFATC2 (Green Channel) and reference assay AP3B1 (Yellow Channel) in different gDNA input amounts (A3,B3,C3: 6ng/rxn, D3,E3,F3: 4ng/rxn, G3,H3: NTCs)



**2D Scatter Plot:** Single-well data of target assay NFATC2 (Green Channel, X-axis) and reference assay AP3B1 (Yellow Channel, Y-axis) in 4ng gDNA input amount.

