



# Validation Data Sheet

dPCR CNV Probe Assay	CD44
Assay Catalog No.	DCG0000231
dPCR Instrument	QIAcuity digital PCR instrument
dPCR Master Mix	QIAcuity Probe PCR Master Mix
Experimental Sample(s)	genomic DNA, XbaI-digested

**1D Scatter Plot:** Single-well data of target assay CD44 (Green Channel) and reference assay AP3B1 (Yellow Channel) in different gDNA input amounts (A9,B9,C9: 6ng/rxn, D9,E9,F9: 4ng/rxn, G9,H9: NTCs)



**2D Scatter Plot:** Single-well data of target assay CD44 (Green Channel, X-axis) and reference assay AP3B1 (Yellow Channel, Y-axis) in 4ng gDNA input amount.

