


# Validation Data Sheet

dPCR CNV Probe Assays	WNT1
Assay Catalog No.	DCG0000396
dPCR Instrument	QIAcuity digital PCR instrument
dPCR Master Mix	QIAcuity Probe PCR Master Mix
Experimental Sample(s)	genomic DNA, XbaI-digested

**1D Scatter Plot:** Single-well data of target assay WNT1 (Green Channel) and reference assay TERT (Yellow Channel) in different gDNA input amounts (A3,B3,C3: 6ng/rxn, D3,E3,F3: 4ng/rxn, G3,H3: NTCs)



**2D Scatter Plot:** Single-well data of target assay WNT1 (Green Channel, X-axis) and reference assay TERT(Yellow Channel, Y-axis) in 4ng gDNA input amount.

