


# Validation Data Sheet

| dPCR CNV Probe Assay   | URI1                            |
|------------------------|---------------------------------|
| Assay Catalog No.      | DCG0000579                      |
| dPCR Instrument        | QIAcuity digital PCR instrument |
| dPCR Master Mix        | QIAcuity Probe PCR Master Mix   |
| Experimental Sample(s) | genomic DNA, XbaI-digested      |

**1D Scatter Plot:** Single-well data of target assay URI1 (Green Channel) and reference assay TERT (Yellow Channel) in different gDNA input amounts (A8,B8,C8: 4ng/rxn, D8,E8,F8: 6ng/rxn, G8,H8: NTCs)



**2D Scatter Plot:** Single-well data of target assay URI1 (Green Channel, X-axis) and reference assay TERT (Yellow Channel, Y-axis) in 4ng gDNA input amount.

