

dPCR CNV Probe Assay

Validation Data Sheet

dPCR CNV Probe Assay

Hk2

Assay Catalog No.	DCG0000607
dPCR Instrument	QIAcuity® digital PCR instrument
dPCR Master Mix	QIAcuity Probe PCR Master Mix
Experimental Sample(s)	genomic DNA, <i>EcoRI</i> -digested




Figure 1. 1D Scatter Plot. Single-well data of target assay Hk2 (Green Channel) and reference assay TERT (Yellow Channel) in different gDNA input amounts (A12, B12, C12: 6 ng/rxn; D12, E12, F12: 4 ng/rxn; G12, H12: NTCs).




Figure 2. 2D Scatter Plot. Single-well data of target assay Hk2 (Green Channel, X-axis) and reference assay TERT (Yellow Channel, Y-axis) in 4 ng gDNA input amount.