



# Validation Data Sheet

dPCR CNV Probe Assays	AP3B1
Assay Catalog No.	DCR0000238
dPCR Instrument	QIAcuity digital PCR instrument
dPCR Master Mix	QIAcuity Probe PCR Master Mix
Experimental Sample(s)	genomic DNA, XbaI-digested

**1D Scatter Plot:** Single-well data of target assay AP3B1 (Green Channel) and reference assay TERT (Yellow Channel) in different gDNA input amounts (A2,B2,C2: 6ng/rxn, D2,E2,F2: 4ng/rxn, G2,H2: NTCs)



**2D Scatter Plot:** Single-well data of target assay AP3B1 (Green Channel, X-axis) and reference assay TERT(Yellow Channel, Y-axis) in 4ng gDNA input amount.

