

RT² IncRNA PCR Array (96-Well Format and 384-Well [4 x 96] Format)

Human Cell Development & Differentiation

Cat. no. 330721 LAHS-003Z

For IncRNA expression analysis by pathway and disease using laboratory-verified SYBR[®] Green qPCR assays

Format	For use with the following real-time cyclers
RT ² IncRNA PCR Array, Format A	Applied Biosystems [®] models 5700, 7000, 7300, 7500, 7700, 7900HT, QuantStudio [®] 6K Flex and QuantStudio 12K Flex (96-well standard block), ViiA [™] 7 (96-well block); Bio-Rad [®] models iCycler [®] , iQ [™] 5, MyiQ [™] , MyiQ2; Bio-Rad/MJ Research Chromo4 [™] ; Eppendorf [®] Mastercycler [®] ep realplex models 2, 2s, 4, 4s; Stratagene [®] models Mx3005P [®] , Mx3000P [®] ; Takara TP-800, TP-900 (Dice [®])
RT ² IncRNA PCR Array, Format C	Applied Biosystems models 7500 (Fast block), 7900HT (fast block), QuantStudio 6K/12K Flex (96-well fast block), StepOnePlus [™] , ViiA7 (fast block)
RT ² IncRNA PCR Array, Format D	Bio-Rad CFX96 [™] , CFX Connect [™] , CFX96 Touch [™] Deep Well; Bio-Rad/MJ Research models DNA Engine Opticon [®] , DNA Engine Opticon 2; Stratagene
RT ² IncRNA PCR Array, Format E	Applied Biosystems models 7900HT (384-well block), ViiA 7 (384-well block); Bio-Rad CFX384 [™]
RT ² IncRNA PCR Array, Format F	Roche [®] LightCycler [®] 480 (96-well block)
RT ² IncRNA PCR Array, Format G	Roche LightCycler 480 (384-well block)
RT ² IncRNA PCR Array, Format H	Fluidigm [®] BioMark [™]
RT ² IncRNA PCR Array, Format R	Rotor-Gene [®] Q

Description

The Human Cell Development & Differentiation RT² lncRNA PCR Array profiles the expression of 84 long noncoding RNAs (lncRNAs) differentially expressed during cellular differentiation and organism development. This array provides researchers with a convenient and accurate way to analyze the lncRNAs most relevant to cell fate and cell lineage decisions. These lncRNAs have been carefully selected based on results published in peer-reviewed journals that suggest a correlation with different stages of development from stem cells to terminal differentiation. The profiling results from this array can potentially serve as a useful molecular marker for specific stem cells, may yield insights into the molecular mechanisms behind specific stem cell differentiation processes, and/or possibly help identify lncRNAs associated with the function of specific mature cells or tissue types. A set of controls present on this array enables data analysis using the $\Delta\Delta\text{CT}$ method of relative quantification, assessment of reverse transcription performance, and assessment of PCR performance. Using SYBR Green-based real-time PCR, the expression of a focused panel of lncRNAs related to development and differentiation can be easily and reliably analyzed with this RT² lncRNA PCR Array.

For further details, consult the *RT² lncRNA PCR Array Handbook*.

Shipping and storage

RT² lncRNA PCR Arrays in formats A, C, D, E, F, G, and R are shipped at ambient temperature, on dry ice, or on blue ice packs, depending on destination and accompanying products. RT² lncRNA PCR Arrays in format H are shipped on dry ice or blue ice packs.

For long term storage, keep plates at -20°C .

Note: Ensure that you have the correct RT² lncRNA PCR Array format for your real-time cycler (see table above).

Note: Open the package and store the products appropriately immediately on receipt.

Array layout (96-well)

For 384-well 4 x 96 PCR arrays, genes are present in a staggered format. Refer to the *RT² lncRNA PCR Array Handbook* for layout.

	1	2	3	4	5	6	7	8	9	10	11	12
A	ADAMTS9-AS2	AIRN	ANKRD20A5P	BACE1-AS	BDNF-AS	CEBPA-AS1	CRNDE	DACT3-AS1	DANCR	DBH-AS1	DIO3OS	DLX6-AS1
B	DNAJC27-AS1	DPP10-AS1	EGOT	EMX2OS	FENDRR	FOXD2-AS1	FRMD6-AS1	GAS5	GATA3-AS1	H19	HAR1A	HAR1B
C	HOTAIR	HOTAIRM1	HOTTIP	HOXA11-AS	HOXA-AS2	HOXA-AS3	IGF2-AS	IPW	JPX	KCND3-AS1	KCNIP4-IT1	KCNQ1OT1
D	LAMP5-AS1	LINC00458	LINC00568	LINC00599	LINC01057	LINC01108	LINC01109	LINC01159	LINC-ROR	LOC101060553	LOC101926895	LOC143666
E	LOC150051	LOC91450	MALAT1	MEG3	MEG9	MIAT	NCBP2-AS2	NEAT1	NR2F1-AS1	OIP5-AS1	PRKCCQ-AS1	RMST
F	RP11-104J23.1	RP11-127I20.7	RP11-292E2.1	RP11-5P4.1	RP11-5P4.2	RP11-669E14.4	RP11-932A10.1	RP11-991C1.1	SNAP25-AS1	SNHG7	SNHG8	SOC32-AS1
G	SOX2-OT	SSSCA1-AS1	TINCR	TRIM52-AS1	TRPC7-AS1	TSIX	TUG1	TUNAR	UCA1	WT1-AS	XIST	ZFAS1
H	ACTB	B2M	RPLP0	RN7SK	SNORA73A	HGDC	RTC	RTC	RTC	PPC	PPC	PPC

Gene table: RT² lncRNA PCR Array

Position	UniGene	GenBank	Symbol	Description
A01	N/A	ENST00000460833	ADAMTS9-AS2	ADAMTS9 antisense RNA 2
A02	N/A	ENST00000601203	AIRN	Antisense of IGF2R non-protein coding RNA
A03	Hs.632352	NR_040113	ANKRD20A5P	Ankyrin repeat domain 20 family, member A5, pseudogene
A04	N/A	NR_037803	BACE1-AS	BACE1 antisense RNA
A05	Hs.675323	NR_002832	BDNF-AS	BDNF antisense RNA
A06	Hs.287521	NR_026887	CEBPA-AS1	Hypothetical LOC80054
A07	N/A	ENST00000501177	CRNDE	Colorectal neoplasia differentially expressed (non-protein coding)
A08	N/A	ENST00000500689	DACT3-AS1	DACT3 antisense RNA 1
A09	Hs.744077	NR_024031	DANCR	KIAA0114
A10	N/A	ENST00000425189	DBH-AS1	DBH antisense RNA 1
A11	Hs.745484	NR_002770	DIO3OS	DIO3 opposite strand (non-protein coding)
A12	Hs.34969	NR_015448	DLX6-AS1	DLX6 antisense RNA 1
B01	Hs.436366	NR_034113	DNAJC27-AS1	Hypothetical LOC729723
B02	Hs.97540	NR_036580	DPP10-AS1	Hypothetical LOC389023
B03	N/A	ENST00000414938	EGOT	Eosinophil granule ontogeny transcript (non-protein coding)
B04	Hs.312592	NR_002791	EMX2OS	EMX2 opposite strand (non-protein coding)
B05	Hs.448825	NR_033925	FENDRR	Hypothetical LOC400550
B06	N/A	ENST00000445551	FOXD2-AS1	FOXD2 antisense RNA 1 (head to head)
B07	Hs.645410	NR_037676	FRMD6-AS1	Chromosome 14 open reading frame 82
B08	Hs.736055	NR_002578	GAS5	Growth arrest-specific 5 (non-protein coding)
B09	Hs.568835	NR_024255	GATA3-AS1	Hypothetical LOC399717
B10	Hs.533566	NR_002196	H19	H19, imprinted maternally expressed transcript (non-protein coding)
B11	Hs.549606	NR_003244	HAR1A	Highly accelerated region 1A (non-protein coding)
B12	N/A	ENST00000447910	HAR1B	Highly accelerated region 1B (non-protein coding)
C01	Hs.197076	NR_003716	HOTAIR	Hox transcript antisense RNA (non-protein coding)
C02	N/A	ENST00000425358	HOTAIRM1	HOXA transcript antisense RNA, myeloid-specific 1
C03	N/A	ENST00000421733	HOTTIP	HOXA distal transcript antisense RNA
C04	Hs.587427	NR_002795	HOXA11-AS	HOXA11 antisense RNA 1 (non-protein coding)
C05	N/A	ENST00000517	HOXA-AS2	HOXA cluster antisense RNA 2

Position	UniGene	GenBank	Symbol	Description
		550		
C06	N/A	ENST00000518848	HOXA-AS3	HOXA cluster antisense RNA 3
C07	Hs.716962	NR_028043	IGF2-AS	Insulin-like growth factor 2 antisense
C08	N/A	NR_023915	IPW	Imprinted in Prader-Willi syndrome (non-protein coding)
C09	N/A	ENST00000414209	JPX	JPX transcript, XIST activator (non-protein coding)
C10	N/A	ENST00000419258	KCND3-AS1	KCND3 antisense RNA 1 [Source:HGNC Symbol;Acc:40317]
C11	N/A	NR_002813	KCNIP4-IT1	KCNIP4 intronic transcript 1 (non-protein coding)
C12	Hs.741312	NR_002728	KCNQ1OT1	KCNQ1 overlapping transcript 1 (non-protein coding)
D01	N/A	ENST00000443469	LAMP5-AS1	LAMP5 antisense RNA 1
D02	N/A	ENST00000423442	LINC00458	Long intergenic non-protein coding RNA 458
D03	N/A	ENST00000416894	LINC00568	Long intergenic non-protein coding RNA 568
D04	Hs.12513	NR_024281	LINC00599	Hypothetical LOC157627
D05	N/A	ENST00000414374	LINC01057	Long intergenic non-protein coding RNA 1057
D06	N/A	ENST00000427276	LINC01108	Long intergenic non-protein coding RNA 1108 [Source:HGNC Symbol;Acc:49234]
D07	N/A	ENST00000603837	LINC01109	Long intergenic non-protein coding RNA 1109 [Source:HGNC Symbol;Acc:49235]
D08	N/A	ENST00000433433	LINC01159	Long intergenic non-protein coding RNA 1159 [Source:HGNC Symbol;Acc:49514]
D09	N/A	ENST00000553704	LINC-ROR	Long intergenic non-protein coding RNA, regulator of reprogramming
D10	N/A	XR_243076	LOC101060553	Uncharacterized LOC101060553
D11	N/A	ENST00000558105	LOC101926895	Uncharacterized LOC101926895
D12	N/A	NR_026967	LOC143666	Uncharacterized LOC143666
E01	Hs.742325	XR_109661	LOC150051	Hypothetical LOC150051
E02	N/A	NR_026998	LOC91450	Uncharacterized LOC91450
E03	Hs.642877	NR_002819	MALAT1	Metastasis associated lung adenocarcinoma transcript 1 (non-protein coding)
E04	Hs.654863	NR_002766	MEG3	Maternally expressed 3 (non-protein coding)
E05	N/A	ENST00000429368	MEG9	Maternally expressed 9 (non-protein coding)
E06	Hs.517502	NR_003491	MIAT	Myocardial infarction associated transcript (non-protein coding)
E07	Hs.118820	NR_024388	NCBP2-AS2	Hypothetical LOC152217
E08	N/A	NR_028272	NEAT1	Nuclear paraspeckle assembly transcript 1 (non-protein coding)
E09	Hs.457407	NR_015369	NR2F1-AS1	Hypothetical LOC441094
E10	N/A	ENST00000500949	OIP5-AS1	OIP5 antisense RNA 1
E11	N/A	ENST00000414894	PRKCC-AS1	PRKCC antisense RNA 1
E12	N/A	NR_024037	RMST	Rhabdomyosarcoma 2 associated transcript (non-protein coding)
F01	N/A	ENST00000588294	RP11-104J23.1	RP11-104J23.1
F02	N/A	ENST00000588099	RP11-127I20.7	RP11-127I20.7
F03	N/A	ENST00000498457	RP11-292E2.1	RP11-292E2.1
F04	N/A	ENST00000420285	RP11-5P4.1	RP11-5P4.1
F05	N/A	ENST00000418091	RP11-5P4.2	RP11-5P4.2
F06	N/A	ENST00000574252	RP11-669E14.4	RP11-669E14.4
F07	N/A	ENST00000556602	RP11-932A10.1	RP11-932A10.1
F08	N/A	ENST00000553435	RP11-991C1.1	RP11-991C1.1
F09	N/A	ENST00000421143	SNAP25-AS1	SNAP25 antisense RNA 1
F10	Hs.741305	NR_003672	SNHG7	Small nucleolar RNA host gene 7 (non-protein coding)
F11	Hs.535762	NR_003584	SNHG8	Small nucleolar RNA host gene 8 (non-protein coding)
F12	N/A	ENST00000499137	SOCS2-AS1	SOCS2 antisense RNA 1
G01	N/A	ENST00000410	SOX2-OT	SOX2 overlapping transcript (non-protein coding)

Position	UniGene	GenBank	Symbol	Description
		534		
G02	N/A	ENST00000567594	SSSCA1-AS1	SSSCA1 antisense RNA 1 (head to head)
G03	Hs.515575	NR_027064	TINCR	Placenta-specific 2 (non-protein coding)
G04	N/A	ENST00000433265	TRIM52-AS1	TRIM52 antisense RNA 1 (head to head)
G05	N/A	ENST00000514459	TRPC7-AS1	TRPC7 antisense RNA 1 [Source:HGNC Symbol;Acc:40936]
G06	Hs.529901	NR_003255	TSIX	TSIX transcript, XIST antisense RNA (non-protein coding)
G07	Hs.554829	NR_002323	TUG1	Taurine upregulated 1 (non-protein coding)
G08	N/A	ENST00000503525	TUNAR	TCL1 upstream neural differentiation-associated RNA
G09	Hs.644234	NR_015379	UCA1	Urothelial cancer associated 1 (non-protein coding)
G10	Hs.567499	NR_023920	WT1-AS	WT1 antisense RNA (non-protein coding)
G11	Hs.529901	NR_001564	XIST	X (inactive)-specific transcript (non-protein coding)
G12	Hs.356766	NR_003604	ZFAS1	ZNF1 antisense RNA 1
H01	Hs.520640	NM_001101	ACTB	Actin, beta
H02	Hs.534255	NM_004048	B2M	Beta-2-microglobulin
H03	Hs.546285	NM_001002	RPLP0	Ribosomal protein, large, P0
H04	N/A	NR_001445	RN7SK	RNA, 7SK small nuclear
H05	N/A	NR_002907	SNORA73A	Small nucleolar RNA, H/ACA box 73A
H06	N/A	SA_00105	HGDC	Human Genomic DNA Contamination
H07	N/A	SA_00104	RTC	Reverse Transcription Control
H08	N/A	SA_00104	RTC	Reverse Transcription Control
H09	N/A	SA_00104	RTC	Reverse Transcription Control
H10	N/A	SA_00103	PPC	Positive PCR Control
H11	N/A	SA_00103	PPC	Positive PCR Control
H12	N/A	SA_00103	PPC	Positive PCR Control

Functional Gene Grouping:

Adult Neural Stem Cells: AIRN, FOXD2-AS1, KCND3-AS1, SIX3-AS1.

Cell-Fate Programming & Reprogramming: DLX6-AS1, FENDRR, CEBPA-AS1, HOTAIR, HOTAIRM1, HOTTIP, JPX, LINC01159, LINC-ROR, MEG3, SIX3-AS1, TINCR, TSIX, TUNAR, XIST.

Dendritic Cell Differentiation: ANKRD20A5P, RP11-292E2.1.

Embryonic & Induced Pluripotent Stem Cells: LINC01109, LINC-ROR, SOX2-OT, TSIX.

Embryonic Development: GAS5, LOC101926895, NCBP2-AS2, NR2F1-AS1, RP11-104J23.1, RP11-127I20.7, RP11-5P4.1, RP11-5P4.2, RP11-669E14.4, RP11-932A10.1, RP11-991C1.1, TRPC7-AS1, UCA1.

Embryonic Stem Cell Pluripotency & Differentiation: MALAT1, NEAT1, SNHG7, TSIX, TUG1.

Mesenchymal Stem Cells & Osteoblast Differentiation: MEG3.

Neurogenesis: ADAMTS9-AS2, BACE1-AS, BDNF-AS, DACT3-AS1, DBH-AS1, DIO3OS, DLX6-AS1, DNAJC27-AS1, DPP10-AS1, EMX2OS, FRMD6-AS1, H19, HAR1A, HAR1B, KCNIP4-IT1, KCNQ1OT1, LAMP5-AS1, LINC00458, LINC00599, LINC01057, LINC01108, LINC01109, LOC101060553, LOC143666, LOC150051, LOC91450, MALAT1, MIAT, NEAT1, OIP5-AS1, PRKCQ-AS1, SIX3-AS1, SNAP25-AS1, SNHG8, SOCS2-AS1, SSSCA1-AS1, TRIM52-AS1, TUG1, RMST.

Skeletal Muscle Development: H19, MALAT1, MEG3, NEAT1.

Known or Predicted Associations with Cell Differentiation & Development -Related miRNA: COX10-AS1, CROCCP2, EPB41L4A-AS1, LINC00116, LINC00493, LINC00657, MIR17HG, RP11-1134I14.8, SENP3-EIF4A1, SNHG15, TTC28-AS1.

Other Differentiation & Development Related lncRNAs: CRNDE, DANCR, EGOT, GATA3-AS1, HOXA11-AS, HOXA-AS2, HOXA-AS3, IGF2-AS, IPW, LINC00568, MEG9, WT1-AS, ZFAS1.

Related products

For optimal performance, RT² lncRNA PCR Arrays should be used together with the RT² First Strand Kit for cDNA synthesis and RT² SYBR[®] Green qPCR Mastermixes for PCR.

Product	Contents	Cat. no.
RT ² First Strand Kit (12)	Enzymes and reagents for cDNA synthesis	330401
RT ² First Strand Kit (50)	Enzymes and reagents for cDNA synthesis	330404
RT ² SYBR Green qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with real-time cyclers that do not require a reference dye, including: Bio-Rad models CFX96, CFX384, DNA Engine Opticon 2; Bio-Rad/MJ Research Chromo4; Roche LightCycler 480 (96-well and 384-well); all other cyclers	330500
RT ² lncRNA qPCR Assays	Laboratory-verified qPCR assays for lncRNA expression	Varies
RT ² SYBR Green ROX™ qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the following real-time cyclers: Applied Biosystems models 5700, 7000, 7300, 7500 [Standard and FAST], 7700, 7900HT 96-well block [Standard and FAST] and 384-well block, StepOnePlus; Eppendorf Mastercycler ep realplex models 2, 2S, 4, 4S; Stratagene models Mx3000P, Mx3005P, Mx4000; Takara TP-800	330520
RT ² SYBR Green Fluor qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the following real-time cyclers: Bio-Rad models iCycler, iQ5, MyiQ, MyiQ2	330510

* Larger kit sizes available; please inquire.

RT² lncRNA PCR Array products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Find out more at www.qiagen.com/search/rt2-lncRNA-pcr-arrays!

Trademarks: QIAGEN[®] (QIAGEN Group); Applied Biosystems[®], ViiA[™], StepOnePlus[™], ROX[™] (Applied Biosystems Corporation or its subsidiaries); Bio-Rad[®], iCycler[®], iQ[™], MyiQ[™], Chromo4[™], CFX96[™], DNA Engine Opticon[®], CFX384[™] (Bio-Rad Laboratories, Inc.); Stratagene[®], Mx3005P[®], Mx3000P[®], Mx4000[®] (Stratagene); Eppendorf[®], Mastercycler[®] (Eppendorf AG); Roche[®], LightCycler[®] (Roche Group); Fluidigm[®] BioMark[™] (Fluidigm Corporation); SYBR[®] (Molecular Probes, Inc.).

1086804 08/2014 © 2014 QIAGEN, all rights reserved.

www.qiagen.com

Australia ■ 1-800-243-800
Austria ■ 0800-281011
Belgium ■ 0800-79612
Brazil ■ 0800-557779
Canada ■ 800-572-9613
China ■ 800-988-0325
Denmark ■ 80-885945

Finland ■ 0800-914416
France ■ 01-60-920-930
Germany ■ 02103-29-12000
Hong Kong ■ 800 933 965
India ■ 1-800-102-4114
Ireland ■ 1800 555 049
Italy ■ 800-787980

Japan ■ 03-6890-7300
Korea (South) ■ 080-000-7145
Luxembourg ■ 8002 2076
Malaysia ■ 603-7981-5510
Mexico ■ 01-800-7742-436
The Netherlands ■ 0800-0229592
Norway ■ 800-18859

Singapore ■ 1800-742-4368
Spain ■ 91-630-7050
Sweden ■ 020-790282
Switzerland ■ 055-254-22-11
Taiwan ■ 0800-665-957
UK ■ 0808-234-3665
USA ■ 800-426-8157

