

# RT<sup>2</sup> lncRNA PreAMP Primer Mix

## Human Cell Differentiation & Development

Cat. no. 330741 LBH-003Z and LBH-003ZH

### For cDNA synthesis prior to real-time PCR

**Note:** Each RT<sup>2</sup> lncRNA PreAMP Primer Mix is specific for a particular pathway or disease RT<sup>2</sup> lncRNA PCR Array. Check the label to verify that the correct pathway-specific RT<sup>2</sup> lncRNA PreAMP Primer Mix is used for the RT<sup>2</sup> lncRNA PCR Array. Verify that the lot number of the RT<sup>2</sup> lncRNA PreAMP Primer Mix is compatible with that of the RT<sup>2</sup> lncRNA PCR Array to be used. If the RT<sup>2</sup> lncRNA PreAMP Primer Mix and the RT<sup>2</sup> lncRNA PCR Array have been purchased at different times, check with QIAGEN Technical Service to ensure their compatibility. For purchase of an RT<sup>2</sup> lncRNA PCR Array to use with this RT<sup>2</sup> lncRNA PreAMP Primer Mix, check with QIAGEN Technical Services before ordering to ensure the available RT<sup>2</sup> lncRNA PCR Arrays are compatible.

For further details, consult the *RT<sup>2</sup> lncRNA PreAMP cDNA Synthesis Handbook*.

### Shipping and storage

RT<sup>2</sup> lncRNA PreAMP Primer Mix is shipped frozen or at ambient temperature. Store at –20°C. Avoid repeated freezing and thawing. Under these conditions, RT<sup>2</sup> lncRNA PreAMP Primer Mix can be kept for at least 6 months from the date of receipt without any reduction in performance.

**Note:** Open the package and store the products appropriately immediately on receipt.

### Kit Contents

<b>Catalog number</b>	LBH-003Z	LBH-003ZH
<b>Product format</b>	RT <sup>2</sup> lncRNA PreAMP Primer Mix	RT <sup>2</sup> lncRNA PreAMP Primer Mix, Format H
<b>Number of reactions</b>	For 12 × 25 µl preamplification reactions	For 480 µl preamplification reactions
<b>Contents</b>	90 µl RT <sup>2</sup> lncRNA PreAMP Primer Mix	1.66 ml RT <sup>2</sup> lncRNA PreAMP Primer Mix



## Related products

Product	Contents	Cat. no.
RT <sup>2</sup> PreAMP cDNA Synthesis Kit (12)	For preamplification of cDNA from small or FFPE samples of RNA for analysis on RT <sup>2</sup> lncRNA PCR Arrays	330451
RT <sup>2</sup> lncRNA PCR Arrays	For lncRNA expression analysis by pathway and disease using laboratory-verified SYBR® Green qPCR assays	330721
RT <sup>2</sup> SYBR® Green qPCR Mastermix (12)	Contains 12 x 1.35 ml tubes: for 12 x 96-well RT <sup>2</sup> PCR Arrays, 8 x 384-well RT <sup>2</sup> PCR Arrays or 1200 x 25 µl reactions or 3000 x 10 µl reactions	330500
RT <sup>2</sup> SYBR® Green ROX qPCR Mastermix (12)	Contains 12 x 1.35 ml tubes: for 12 x 96-well RT <sup>2</sup> PCR Arrays, 8 x 384-well RT <sup>2</sup> PCR Arrays or 1200 x 25 µl reactions or 3000 x 10 µl reactions	330520
RT <sup>2</sup> SYBR® Green Fluor qPCR Mastermix (12)	Contains 12 x 1.35 ml tubes: for 12 x 96-well RT <sup>2</sup> PCR Arrays, 8 x 384-well RT <sup>2</sup> PCR Arrays or 1200 x 25 µl reactions or 3000 x 10 µl reactions	330510
RT <sup>2</sup> SYBR® Green FAST Mastermix (12)	Contains 12 x 1.35 ml tubes: for 12 x 100-well ring discs on the QIAGEN Rotor-Gene Q or 1500 x 20 µl reactions	330620

RT<sup>2</sup> lncRNA PreAMP Primer Mix products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at [www.qiagen.com](http://www.qiagen.com) or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, Rotor-Gene® (QIAGEN Group); SYBR® (Molecular Probes, Inc.); Roche®, LightCycler® (Roche Group); Applied Biosystems®, ROX™, StepOnePlus™, ViiA™ (Applied Biosystems Corporation or its subsidiaries); Eppendorf®, Mastercycler® (Eppendorf AG); Stratagene®, Mx3005P®, Mx3000P®, Mx4000® (Stratagene); Bio-Rad®, iCycler®, Chromo4™, CFX96™, DNA Engine Opticon®, CFX384™, iQ™, MyiQ™ (Bio-Rad Laboratories, Inc.).

1086828 08/2014 © 2014 QIAGEN, all rights reserved.

[www.qiagen.com](http://www.qiagen.com)

Canada ■ 800-572-9613

Ireland ■ 1800 555 049

Norway ■ 800-18859

China ■ 8621-3865-3865

Italy ■ 800-787980

Singapore ■ 1800-742-4368

Denmark ■ 80-885945

Japan ■ 03-6890-7300

Spain ■ 91-630-7050

Australia ■ 1-800-243-800

Finland ■ 0800-914416

Korea (South) ■ 080-000-7145

Sweden ■ 020-790282

Austria ■ 0800/281010

France ■ 01-60-920-930

Luxembourg ■ 8002 2076

Switzerland ■ 055-254-22-11

Belgium ■ 0800-79612

Germany ■ 02103-29-12000

Mexico ■ 01-800-7742-436

UK ■ 01293-422-911

Brazil ■ 0800-557779

Hong Kong ■ 800 933 965

The Netherlands ■ 0800 0229592

USA ■ 800-426-8157



Sample & Assay Technologies