

Zebrafish RT² qPCR Primer gDNA Control

For detection of zebrafish genomic DNA contamination

The RT² qPCR Primer Assay for zebrafish genomic DNA contamination specifically detects non-transcribed, zebrafish genomic DNA contamination with a high level of sensitivity. For protocol details, refer to the *RT² qPCR Primer Assay Handbook*.

Contents

Product	Contents	Cat. no.
Zebrafish RT ² qPCR Primer gDNA Control (200)	RT ² qPCR Primer Assay for zebrafish genomic DNA contamination; One zebrafish genomic DNA RT ² qPCR Primer Assay consisting of 2 primers for 200 PCR reactions; 10 μM each primer	330011 PA-038

Shipping and storage

The Zebrafish RT² qPCR Primer gDNA Control is shipped at ambient temperature. For long-term storage, keep tubes at -20°C. If stored under these conditions, RT² qPCR Primer gDNA Controls are stable for 6 months after receipt.

Note: Open the package and store the products appropriately immediately on receipt.



RT² qPCR Primer Assay products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN® (QIAGEN Group).

1074388 11/2012 © 2012 QIAGEN, all rights reserved.

www.qiagen.com

Canada = 800-572-9613

Ireland = 1800 555 049

Norway = 800-18859

China = 8621-3865-3865

Italy = 800-787980

Singapore = 1800-742-4368

Denmark = 80-885945

Japan = 03-6890-7300

Spain = 91-630-7050

Australia = 1-800-243-800

Finland = 0800-914416

Korea (South) = 080-000-7145

Sweden = 020-790282

Austria = 0800/281010

France = 01-60-920-930

Luxembourg = 8002 2076

Switzerland = 055-254-22-11

Belgium = 0800-79612

Germany = 02103-29-12000

Mexico = 01-800-7742-436

UK = 01293-422-911

Brazil = 0800-557779

Hong Kong = 800 933 965

The Netherlands = 0800 0229592

USA = 800-426-8157



Sample & Assay Technologies