

RT² Profiler PCR Array (96-Well Format and 384-Well [4 x 96] Format)

Human Cellular Stress Responses

Cat. no. 330231 PAHS-019ZA

For pathway expression analysis

Format	For use with the following real-time cyclers
RT ² Profiler PCR Array, Format A	Applied Biosystems [®] models 5700, 7000, 7300, 7500, 7700, 7900HT, ViiA™ 7 (96-well block); Bio-Rad [®] models iCycler [®] , iQ™ 5, MyiQ™, MyiQ2; Bio-Rad/MJ Research Chromo4™; Eppendorf [®] Mastercycler [®] ep realplex models 2, 2s, 4, 4s; Stratagene [®] models Mx3005P [®] , Mx3000P [®] ; Takara TP-800
RT ² Profiler PCR Array, Format C	Applied Biosystems models 7500 (Fast block), 7900HT (Fast block), StepOnePlus™, ViiA 7 (Fast block)
RT ² Profiler PCR Array, Format D	Bio-Rad CFX96™; Bio-Rad/MJ Research models DNA Engine Opticon [®] , DNA Engine Opticon 2; Stratagene Mx4000 [®]
RT ² Profiler PCR Array, Format E	Applied Biosystems models 7900HT (384-well block), ViiA 7 (384-well block); Bio-Rad CFX384™
RT ² Profiler PCR Array, Format F	Roche [®] LightCycler [®] 480 (96-well block)
RT ² Profiler PCR Array, Format G	Roche LightCycler 480 (384-well block)
RT ² Profiler PCR Array, Format H	Fluidigm [®] BioMark™



Description

The Human Cellular Stress Responses RT² Profiler PCR Array profiles the expression of 84 genes involved in cellular stress response. This array includes genes encoding antioxidant and pro-oxidant enzymes, molecular chaperones and proteins involved in xenobiotic metabolism. Representative heat shock proteins and cytochrome P450 genes are also included. Using real-time PCR, you can easily and reliably analyze expression of a focused panel of genes related to cellular stress response with this array.

For further details, consult the *RT² Profiler PCR Array Handbook*.

Shipping and storage

RT² Profiler PCR Arrays in formats A, C, D, E, F, and G are shipped at ambient temperature, on dry ice, or blue ice packs depending on destination and accompanying products. RT² Profiler PCR Arrays in format H are shipped on dry ice or blue ice packs.

For long term storage, keep plates at –20°C.

Note: Ensure that you have the correct RT² Profiler PCR Array format for your real-time cycler (see table above).

Note: Open the package and store the products appropriately immediately on receipt.

Array layout (96-well)

For 384-well 4 x 96 PCR arrays, genes are present in a staggered format. Refer to the *RT² Profiler PCR Array Handbook* for layout.

	1	2	3	4	5	6	7	8	9	10	11	12
A	ACADSB	BAG1	CALR	CANX	CAT	CCT2	CCT3	CCT4	CCT5	CCT7	CCT8	CES1
B	CLU	CRYAA	CRYAB	CYP11A1	CYP11B2	CYP17A1	CYP1A1	CYP1A2	CYP1B1	CYP2C19	CYP2C9	CYP2D6
C	CYP2E1	CYP2F1	CYP3A4	CYP4A11	CYP4B1	CYP7A1	CYP7B1	DNAJA1	DNAJA2	DNAJA3	DNAJA4	DNAJB1
D	DNAJB11	DNAJB12	DNAJB2	DNAJB4	DNAJB5	DNAJB6	DNAJB9	DNAJC1	DNAJC4	DNAJC5	DNAJC6	DNAJC7
E	DNAJC8	DNAJC9	EPHX2	FMO1	FMO4	FMO5	GPX1	GPX2	GSR	GSTA1	GSTA5	GSTM5
F	HMOX1	HMOX2	HOPX	HSPA12A	HSPA1A	HSPA1L	HSPA2	HSPA4	HSPA5	HSPA8	HSPA9	HSPB1
G	HSPB2	HSPB3	HSPD1	HSPF1	PPID	SERPINH1	SOD1	SOD2	SOD3	ST13	TCP1	XDH
H	ACTB	B2M	GAPDH	HPRT1	RPLP0	HGDC	RTC	RTC	RTC	PPC	PPC	PPC

Gene table: RT² Profiler PCR Array

Position	UniGene	GenBank	Symbol	Description
A01	Hs.81934	NM_001609	ACADSB	Acyl-CoA dehydrogenase, short/branched chain
A02	Hs.377484	NM_004323	BAG1	BCL2-associated athanogene
A03	Hs.515162	NM_004343	CALR	Calreticulin
A04	Hs.699155	NM_001746	CANX	Calnexin
A05	Hs.502302	NM_001752	CAT	Catalase
A06	Hs.189772	NM_006431	CCT2	Chaperonin containing TCP1, subunit 2 (beta)
A07	Hs.491494	NM_005998	CCT3	Chaperonin containing TCP1, subunit 3 (gamma)
A08	Hs.421509	NM_006430	CCT4	Chaperonin containing TCP1, subunit 4 (delta)
A09	Hs.1600	NM_012073	CCT5	Chaperonin containing TCP1, subunit 5 (epsilon)
A10	Hs.368149	NM_006429	CCT7	Chaperonin containing TCP1, subunit 7 (eta)
A11	Hs.125113	NM_006585	CCT8	Chaperonin containing TCP1, subunit 8 (theta)
A12	Hs.558865	NM_001266	CES1	Carboxylesterase 1
B01	Hs.436657	NM_001831	CLU	Clusterin
B02	Hs.184085	NM_000394	CRYAA	Crystallin, alpha A
B03	Hs.408767	NM_001885	CRYAB	Crystallin, alpha B
B04	Hs.303980	NM_000781	CYP11A1	Cytochrome P450, family 11, subfamily A, polypeptide 1
B05	Hs.632054	NM_000498	CYP11B2	Cytochrome P450, family 11, subfamily B, polypeptide 2
B06	Hs.438016	NM_000102	CYP17A1	Cytochrome P450, family 17, subfamily A, polypeptide 1
B07	Hs.72912	NM_000499	CYP1A1	Cytochrome P450, family 1, subfamily A, polypeptide 1
B08	Hs.1361	NM_000761	CYP1A2	Cytochrome P450, family 1, subfamily A, polypeptide 2
B09	Hs.154654	NM_000104	CYP1B1	Cytochrome P450, family 1, subfamily B, polypeptide 1
B10	Hs.282409	NM_000769	CYP2C19	Cytochrome P450, family 2, subfamily C, polypeptide 19
B11	Hs.282624	NM_000771	CYP2C9	Cytochrome P450, family 2, subfamily C, polypeptide 9
B12	Hs.648256	NM_000106	CYP2D6	Cytochrome P450, family 2, subfamily D, polypeptide 6
C01	Hs.12907	NM_000773	CYP2E1	Cytochrome P450, family 2, subfamily E, polypeptide 1
C02	Hs.558318	NM_000774	CYP2F1	Cytochrome P450, family 2, subfamily F, polypeptide 1
C03	Hs.654391	NM_017460	CYP3A4	Cytochrome P450, family 3, subfamily A, polypeptide 4
C04	Hs.1645	NM_000778	CYP4A11	Cytochrome P450, family 4, subfamily A, polypeptide 11
C05	Hs.436317	NM_000779	CYP4B1	Cytochrome P450, family 4, subfamily B, polypeptide 1
C06	Hs.1644	NM_000780	CYP7A1	Cytochrome P450, family 7, subfamily A, polypeptide 1
C07	Hs.667720	NM_004820	CYP7B1	Cytochrome P450, family 7, subfamily B, polypeptide 1
C08	Hs.445203	NM_001539	DNAJA1	DnaJ (Hsp40) homolog, subfamily A, member 1
C09	Hs.368078	NM_005880	DNAJA2	DnaJ (Hsp40) homolog, subfamily A, member 2
C10	Hs.459779	NM_005147	DNAJA3	DnaJ (Hsp40) homolog, subfamily A, member 3
C11	Hs.513053	NM_018602	DNAJA4	DnaJ (Hsp40) homolog, subfamily A, member 4
C12	Hs.515210	NM_006145	DNAJB1	DnaJ (Hsp40) homolog, subfamily B, member 1
D01	Hs.317192	NM_016306	DNAJB11	DnaJ (Hsp40) homolog, subfamily B, member 11
D02	Hs.696014	NM_017626	DNAJB12	DnaJ (Hsp40) homolog, subfamily B, member 12
D03	Hs.77768	NM_006736	DNAJB2	DnaJ (Hsp40) homolog, subfamily B, member 2
D04	Hs.13852	NM_007034	DNAJB4	DnaJ (Hsp40) homolog, subfamily B, member 4
D05	Hs.237506	NM_012266	DNAJB5	DnaJ (Hsp40) homolog, subfamily B, member 5
D06	Hs.490745	NM_005494	DNAJB6	DnaJ (Hsp40) homolog, subfamily B, member 6
D07	Hs.6790	NM_012328	DNAJB9	DnaJ (Hsp40) homolog, subfamily B, member 9
D08	Hs.499000	NM_022365	DNAJC1	DnaJ (Hsp40) homolog, subfamily C, member 1
D09	Hs.172847	NM_005528	DNAJC4	DnaJ (Hsp40) homolog, subfamily C, member 4

Position	UniGene	GenBank	Symbol	Description
D10	Hs.164419	NM_025219	DNAJC5	DnaJ (Hsp40) homolog, subfamily C, member 5
D11	Hs.647643	NM_014787	DNAJC6	DnaJ (Hsp40) homolog, subfamily C, member 6
D12	Hs.500156	NM_003315	DNAJC7	DnaJ (Hsp40) homolog, subfamily C, member 7
E01	Hs.433540	NM_014280	DNAJC8	DnaJ (Hsp40) homolog, subfamily C, member 8
E02	Hs.654694	NM_015190	DNAJC9	DnaJ (Hsp40) homolog, subfamily C, member 9
E03	Hs.212088	NM_001979	EPHX2	Epoxide hydrolase 2, cytoplasmic
E04	Hs.1424	NM_002021	FMO1	Flavin containing monooxygenase 1
E05	Hs.386502	NM_002022	FMO4	Flavin containing monooxygenase 4
E06	Hs.642706	NM_001461	FMO5	Flavin containing monooxygenase 5
E07	Hs.76686	NM_000581	GPX1	Glutathione peroxidase 1
E08	Hs.2704	NM_002083	GPX2	Glutathione peroxidase 2 (gastrointestinal)
E09	Hs.271510	NM_000637	GSR	Glutathione reductase
E10	Hs.446309	NM_145740	GSTA1	Glutathione S-transferase alpha 1
E11	Hs.646984	NM_153699	GSTA5	Glutathione S-transferase alpha 5
E12	Hs.75652	NM_000851	GSTM5	Glutathione S-transferase mu 5
F01	Hs.517581	NM_002133	HMOX1	Heme oxygenase (decycling) 1
F02	Hs.284279	NM_002134	HMOX2	Heme oxygenase (decycling) 2
F03	Hs.654864	NM_139211	HOPX	HOP homeobox
F04	Hs.654682	NM_025015	HSPA12A	Heat shock 70kDa protein 12A
F05	Hs.728810	NM_005345	HSPA1A	Heat shock 70kDa protein 1A
F06	Hs.690634	NM_005527	HSPA1L	Heat shock 70kDa protein 1-like
F07	Hs.728938	NM_021979	HSPA2	Heat shock 70kDa protein 2
F08	Hs.90093	NM_002154	HSPA4	Heat shock 70kDa protein 4
F09	Hs.716396	NM_005347	HSPA5	Heat shock 70kDa protein 5 (glucose-regulated protein, 78kDa)
F10	Hs.702021	NM_006597	HSPA8	Heat shock 70kDa protein 8
F11	Hs.184233	NM_004134	HSPA9	Heat shock 70kDa protein 9 (mortalin)
F12	Hs.520973	NM_001540	HSPB1	Heat shock 27kDa protein 1
G01	Hs.709660	NM_001541	HSPB2	Heat shock 27kDa protein 2
G02	Hs.41707	NM_006308	HSPB3	Heat shock 27kDa protein 3
G03	Hs.595053	NM_002156	HSPD1	Heat shock 60kDa protein 1 (chaperonin)
G04	Hs.1197	NM_002157	HSPE1	Heat shock 10kDa protein 1 (chaperonin 10)
G05	Hs.581725	NM_005038	PPID	Peptidylprolyl isomerase D
G06	Hs.596449	NM_001235	SERPINH1	Serpin peptidase inhibitor, clade H (heat shock protein 47), member 1, (collagen binding protein 1)
G07	Hs.443914	NM_000454	SOD1	Superoxide dismutase 1, soluble
G08	Hs.487046	NM_000636	SOD2	Superoxide dismutase 2, mitochondrial
G09	Hs.2420	NM_003102	SOD3	Superoxide dismutase 3, extracellular
G10	Hs.707999	NM_003932	ST13	Suppression of tumorigenicity 13 (colon carcinoma) (Hsp70 interacting protein)
G11	Hs.363137	NM_030752	TCP1	T-complex 1
G12	Hs.250	NM_000379	XDH	Xanthine dehydrogenase
H01	Hs.520640	NM_001101	ACTB	Actin, beta
H02	Hs.534255	NM_004048	B2M	Beta-2-microglobulin
H03	Hs.592355	NM_002046	GAPDH	Glyceraldehyde-3-phosphate dehydrogenase
H04	Hs.412707	NM_000194	HPRT1	Hypoxanthine phosphoribosyltransferase 1
H05	Hs.546285	NM_001002	RPLP0	Ribosomal protein, large, P0
H06	N/A	SA_00105	HGDC	Human Genomic DNA Contamination
H07	N/A	SA_00104	RTC	Reverse Transcription Control
H08	N/A	SA_00104	RTC	Reverse Transcription Control
H09	N/A	SA_00104	RTC	Reverse Transcription Control
H10	N/A	SA_00103	PPC	Positive PCR Control
H11	N/A	SA_00103	PPC	Positive PCR Control
H12	N/A	SA_00103	PPC	Positive PCR Control

Related products

For optimal performance, RT² Profiler PCR Arrays should be used together with the RT² First Strand Kit for cDNA synthesis and RT² SYBR[®] Green qPCR Mastermixes for PCR.

Product	Contents	Cat. no.
RT ² First Strand Kit (12)	Enzymes and reagents for cDNA synthesis	330401
RT ² SYBR Green qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with real-time cyclers that do not require a reference dye, including: Bio-Rad models CFX96, CFX384, DNA Engine Opticon 2; Bio-Rad/MJ Research Chromo4; Roche LightCycler 480 (96-well and 384-well); all other cyclers	330500
RT ² SYBR Green ROX [™] qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the following real-time cyclers: Applied Biosystems models 5700, 7000, 7300, 7500 [Standard and FAST], 7700, 7900HT 96-well block [Standard and FAST] and 384-well block, StepOnePlus; Eppendorf Mastercycler ep realplex models 2, 2S, 4, 4S; Stratagene models Mx3000P, Mx3005P, Mx4000; Takara TP-800	330520
RT ² SYBR Green Fluor qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the following real-time cyclers: Bio-Rad models iCycler, iQ5, MyiQ, MyiQ2	330510

* Larger kit sizes available; please inquire.

RT² Profiler PCR Array products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN® (QIAGEN Group); Applied Biosystems®, ViiA™, StepOnePlus™, ROX™ (Applied Biosystems Corporation or its subsidiaries); Bio-Rad®, iCycler®, iQ™, MyiQ™, Chromo4™, CFX96™, DNA Engine Opticon®, CFX384™ (Bio-Rad Laboratories, Inc.); Stratagene®, Mx3005P®, Mx3000P®, Mx4000® (Stratagene); Eppendorf®, Mastercycler® (Eppendorf AG); Roche®, LightCycler® (Roche Group); Fluidigm® BioMark™ (Fluidigm Corporation); SYBR® (Molecular Probes, Inc.).

1066029 03/2011 © 2011 QIAGEN, all rights reserved.

www.qiagen.com

Canada ■ 800-572-9613

Ireland ■ 1800 555 049

Norway ■ 800-18859

China ■ 8621-3865-3865

Italy ■ 800-787980

Singapore ■ 1800-742-4368

Denmark ■ 80-885945

Japan ■ 03-6890-7300

Spain ■ 91-630-7050

Australia ■ 1-800-243-800

Finland ■ 0800-914416

Korea (South) ■ 080-000-7145

Sweden ■ 020-790282

Austria ■ 0800/281010

France ■ 01-60-920-930

Luxembourg ■ 8002 2076

Switzerland ■ 055-254-22-11

Belgium ■ 0800-79612

Germany ■ 02103-29-12000

Mexico ■ 01-800-7742-436

UK ■ 01293-422-911

Brazil ■ 0800-557779

Hong Kong ■ 800 933 965

The Netherlands ■ 0800 0229592

USA ■ 800-426-8157



Sample & Assay Technologies