

RT² Profiler PCR Array (96-Well Format and 384-Well [4 x 96] Format)

Human MAP Kinase Signaling Pathway

Cat. no. 330231 PAHS-061ZA

For pathway expression analysis

Format	For use with the following real-time cyclers
RT ² Profiler PCR Array, Format A	Applied Biosystems® models 5700, 7000, 7300, 7500, 7700, 7900HT, ViiA™ 7 (96-well block); Bio-Rad® models iCycler®, iQ™ 5, MyiQ™, MyiQ2; Bio-Rad/MJ Research Chromo4™; Eppendorf® Mastercycler® ep realplex models 2, 2s, 4, 4s; Stratagene® models Mx3005P®, Mx3000P®; Takara TP-800
RT ² Profiler PCR Array, Format C	Applied Biosystems models 7500 (Fast block), 7900HT (Fast block), StepOnePlus™, ViiA 7 (Fast block)
RT ² Profiler PCR Array, Format D	Bio-Rad CFX96™; Bio-Rad/MJ Research models DNA Engine Opticon®, DNA Engine Opticon 2; Stratagene Mx4000®
RT ² Profiler PCR Array, Format E	Applied Biosystems models 7900HT (384-well block), ViiA 7 (384-well block); Bio-Rad CFX384™
RT ² Profiler PCR Array, Format F	Roche® LightCycler® 480 (96-well block)
RT ² Profiler PCR Array, Format G	Roche LightCycler 480 (384-well block)
RT ² Profiler PCR Array, Format H	Fluidigm® BioMark™



Sample & Assay Technologies

Description

The Human MAP Kinase Signaling Pathway RT² Profiler PCR Array profiles the expression of 84 genes related to the MAP kinase (MAPK) signaling pathway. Members of the MKKK, MKK, and MAPK families are represented on this array. Transcription factors and genes whose expression is induced by MAP kinase signaling are included as well. Raf regulating proteins and MEKK1 interacting proteins are found on this array along with cell cycle proteins regulated by the Erk1/2 pathway. Genes related to scaffolding and anchoring are also included. Using real-time PCR, you can easily and reliably analyze expression of a focused panel of genes related to the MAPK signaling pathway with this array.

For further details, consult the *RT² Profiler PCR Array Handbook*.

Shipping and storage

RT² Profiler PCR Arrays in formats A, C, D, E, F, and G are shipped at ambient temperature, on dry ice, or blue ice packs depending on destination and accompanying products. RT² Profiler PCR Arrays in format H are shipped on dry ice or blue ice packs.

For long term storage, keep plates at -20°C .

Note: Ensure that you have the correct RT² Profiler PCR Array format for your real-time cycler (see table above).

Note: Open the package and store the products appropriately immediately on receipt.

Array layout (96-well)

For 384-well 4 x 96 PCR arrays, genes are present in a staggered format. Refer to the *RT² Profiler PCR Array Handbook* for layout.

	1	2	3	4	5	6	7	8	9	10	11	12
A	ARAF	ATF2	BRAF	CCNA1	CCNA2	CCNB1	CCNB2	CCND1	CCND2	CCND3	CCNE1	CDC42
B	CDK2	CDK4	CDK6	CDKN1A	CDKN1B	CDKN1C	CDKN2A	CDKN2B	CDKN2C	CDKN2D	CHUK	COL1A1
C	CREB1	CREBBP	DLK1	E2F1	EGFR	EGR1	ELK1	ETS1	ETS2	FOS	GRB2	HRAS
D	HSPA5	HSPB1	JUN	KRAS	KSR1	LAMTOR3	MAP2K1	MAP2K2	MAP2K3	MAP2K4	MAP2K5	MAP2K6
E	MAP2K7	MAP3K1	MAP3K2	MAP3K3	MAP3K4	MAP4K1	MAPK1	MAPK10	MAPK11	MAPK12	MAPK13	MAPK14
F	MAPK3	MAPK6	MAPK7	MAPK8	MAPK8IP2	MAPK9	MAPKAPK2	MAPKAPK3	MAX	MEF2C	MKNK1	MOS
G	MST1	MYC	NFATC4	NRAS	PAK1	PRDX6	RAC1	RAF1	RB1	SFN	SMAD4	TP53
H	ACTB	B2M	GAPDH	HPRT1	RPLP0	HGDC	RTC	RTC	RTC	PPC	PPC	PPC

Gene table: RT² Profiler PCR Array

Position	UniGene	GenBank	Symbol	Description
A01	Hs.446641	NM_001654	ARAF	V-raf murine sarcoma 3611 viral oncogene homolog
A02	Hs.592510	NM_001880	ATF2	Activating transcription factor 2
A03	Hs.550061	NM_004333	BRAF	V-raf murine sarcoma viral oncogene homolog B1
A04	Hs.417050	NM_003914	CCNA1	Cyclin A1
A05	Hs.58974	NM_001237	CCNA2	Cyclin A2
A06	Hs.23960	NM_031966	CCNB1	Cyclin B1
A07	Hs.194698	NM_004701	CCNB2	Cyclin B2
A08	Hs.523852	NM_053056	CCND1	Cyclin D1
A09	Hs.376071	NM_001759	CCND2	Cyclin D2
A10	Hs.534307	NM_001760	CCND3	Cyclin D3
A11	Hs.244723	NM_001238	CCNE1	Cyclin E1
A12	Hs.690198	NM_001791	CDC42	Cell division cycle 42 (GTP binding protein, 25kDa)
B01	Hs.19192	NM_001798	CDK2	Cyclin-dependent kinase 2
B02	Hs.95577	NM_000075	CDK4	Cyclin-dependent kinase 4
B03	Hs.119882	NM_001259	CDK6	Cyclin-dependent kinase 6
B04	Hs.370771	NM_000389	CDKN1A	Cyclin-dependent kinase inhibitor 1A (p21, Cip1)
B05	Hs.238990	NM_004064	CDKN1B	Cyclin-dependent kinase inhibitor 1B (p27, Kip1)
B06	Hs.106070	NM_000076	CDKN1C	Cyclin-dependent kinase inhibitor 1C (p57, Kip2)
B07	Hs.512599	NM_000077	CDKN2A	Cyclin-dependent kinase inhibitor 2A (melanoma, p16, inhibits CDK4)
B08	Hs.72901	NM_004936	CDKN2B	Cyclin-dependent kinase inhibitor 2B (p15, inhibits CDK4)
B09	Hs.728783	NM_078626	CDKN2C	Cyclin-dependent kinase inhibitor 2C (p18, inhibits CDK4)
B10	Hs.435051	NM_001800	CDKN2D	Cyclin-dependent kinase inhibitor 2D (p19, inhibits CDK4)
B11	Hs.198998	NM_001278	CHUK	Conserved helix-loop-helix ubiquitous kinase
B12	Hs.172928	NM_000088	COL1A1	Collagen, type I, alpha 1
C01	Hs.516646	NM_004379	CREB1	CAMP responsive element binding protein 1
C02	Hs.459759	NM_004380	CREBBP	CREB binding protein
C03	Hs.533717	NM_003836	DLK1	Delta-like 1 homolog (Drosophila)
C04	Hs.654393	NM_005225	E2F1	E2F transcription factor 1
C05	Hs.488293	NM_005228	EGFR	Epidermal growth factor receptor
C06	Hs.326035	NM_001964	EGR1	Early growth response 1
C07	Hs.181128	NM_005229	ELK1	ELK1, member of ETS oncogene family
C08	Hs.369438	NM_005238	ETS1	V-ets erythroblastosis virus E26 oncogene homolog 1 (avian)
C09	Hs.644231	NM_005239	ETS2	V-Ets erythroblastosis virus E26 oncogene homolog 2 (avian)
C10	Hs.728789	NM_005252	FOS	FBJ murine osteosarcoma viral oncogene homolog
C11	Hs.444356	NM_002086	GRB2	Growth factor receptor-bound protein 2
C12	Hs.37003	NM_005343	HRAS	V-Ha-ras Harvey rat sarcoma viral oncogene homolog
D01	Hs.716396	NM_005347	HSPA5	Heat shock 70kDa protein 5 (glucose-regulated protein, 78kDa)
D02	Hs.520973	NM_001540	HSPB1	Heat shock 27kDa protein 1
D03	Hs.714791	NM_002228	JUN	Jun proto-oncogene
D04	Hs.505033	NM_004985	KRAS	V-Ki-ras2 Kirsten rat sarcoma viral oncogene homolog
D05	Hs.133534	NM_014238	KSR1	Kinase suppressor of ras 1
D06	Hs.716375	NM_021970	LAMTOR3	Late endosomal/lysosomal adaptor, MAPK and MTOR activator 3
D07	Hs.145442	NM_002755	MAP2K1	Mitogen-activated protein kinase kinase 1
D08	Hs.465627	NM_030662	MAP2K2	Mitogen-activated protein kinase kinase 2
D09	Hs.514012	NM_002756	MAP2K3	Mitogen-activated protein kinase kinase 3

Position	UniGene	GenBank	Symbol	Description
D10	Hs.514681	NM_003010	MAP2K4	Mitogen-activated protein kinase kinase 4
D11	Hs.114198	NM_002757	MAP2K5	Mitogen-activated protein kinase kinase 5
D12	Hs.463978	NM_002758	MAP2K6	Mitogen-activated protein kinase kinase 6
E01	Hs.531754	NM_145185	MAP2K7	Mitogen-activated protein kinase kinase 7
E02	Hs.657756	NM_005921	MAP3K1	Mitogen-activated protein kinase kinase kinase 1
E03	Hs.145605	NM_006609	MAP3K2	Mitogen-activated protein kinase kinase kinase 2
E04	Hs.29282	NM_002401	MAP3K3	Mitogen-activated protein kinase kinase kinase 3
E05	Hs.390428	NM_005922	MAP3K4	Mitogen-activated protein kinase kinase kinase 4
E06	Hs.95424	NM_007181	MAP4K1	Mitogen-activated protein kinase kinase kinase kinase 1
E07	Hs.431850	NM_002745	MAPK1	Mitogen-activated protein kinase 1
E08	Hs.125503	NM_002753	MAPK10	Mitogen-activated protein kinase 10
E09	Hs.57732	NM_002751	MAPK11	Mitogen-activated protein kinase 11
E10	Hs.432642	NM_002969	MAPK12	Mitogen-activated protein kinase 12
E11	Hs.178695	NM_002754	MAPK13	Mitogen-activated protein kinase 13
E12	Hs.485233	NM_001315	MAPK14	Mitogen-activated protein kinase 14
F01	Hs.861	NM_002746	MAPK3	Mitogen-activated protein kinase 3
F02	Hs.411847	NM_002748	MAPK6	Mitogen-activated protein kinase 6
F03	Hs.150136	NM_002749	MAPK7	Mitogen-activated protein kinase 7
F04	Hs.138211	NM_002750	MAPK8	Mitogen-activated protein kinase 8
F05	Hs.558180	NM_012324	MAPK8IP2	Mitogen-activated protein kinase 8 interacting protein 2
F06	Hs.484371	NM_002752	MAPK9	Mitogen-activated protein kinase 9
F07	Hs.643566	NM_004759	MAPKAPK2	Mitogen-activated protein kinase-activated protein kinase 2
F08	Hs.234521	NM_004635	MAPKAPK3	Mitogen-activated protein kinase-activated protein kinase 3
F09	Hs.285354	NM_002382	MAX	MYC associated factor X
F10	Hs.653394	NM_002397	MEF2C	Myocyte enhancer factor 2C
F11	Hs.371594	NM_003684	MKNK1	MAP kinase interacting serine/threonine kinase 1
F12	Hs.533432	NM_005372	MOS	V-mos Moloney murine sarcoma viral oncogene homolog
G01	Hs.512587	NM_020998	MST1	Macrophage stimulating 1 (hepatocyte growth factor-like)
G02	Hs.202453	NM_002467	MYC	V-myc myelocytomatosis viral oncogene homolog (avian)
G03	Hs.77810	NM_004554	NFATC4	Nuclear factor of activated T-cells, cytoplasmic, calcineurin-dependent 4
G04	Hs.486502	NM_002524	NRAS	Neuroblastoma RAS viral (v-ras) oncogene homolog
G05	Hs.435714	NM_002576	PAK1	P21 protein (Cdc42/Rac)-activated kinase 1
G06	Hs.120	NM_004905	PRDX6	Peroxiredoxin 6
G07	Hs.413812	NM_006908	RAC1	Ras-related C3 botulinum toxin substrate 1 (rho family, small GTP binding protein Rac1)
G08	Hs.159130	NM_002880	RAF1	V-raf-1 murine leukemia viral oncogene homolog 1
G09	Hs.408528	NM_000321	RB1	Retinoblastoma 1
G10	Hs.523718	NM_006142	SFN	Stratifin
G11	Hs.75862	NM_005359	SMAD4	SMAD family member 4
G12	Hs.654481	NM_000546	TP53	Tumor protein p53
H01	Hs.520640	NM_001101	ACTB	Actin, beta
H02	Hs.534255	NM_004048	B2M	Beta-2-microglobulin
H03	Hs.592355	NM_002046	GAPDH	Glyceraldehyde-3-phosphate dehydrogenase
H04	Hs.412707	NM_000194	HPRT1	Hypoxanthine phosphoribosyltransferase 1
H05	Hs.546285	NM_001002	RPLP0	Ribosomal protein, large, P0
H06	N/A	SA_00105	HGDC	Human Genomic DNA Contamination
H07	N/A	SA_00104	RTC	Reverse Transcription Control
H08	N/A	SA_00104	RTC	Reverse Transcription Control
H09	N/A	SA_00104	RTC	Reverse Transcription Control
H10	N/A	SA_00103	PPC	Positive PCR Control
H11	N/A	SA_00103	PPC	Positive PCR Control
H12	N/A	SA_00103	PPC	Positive PCR Control

Related products

For optimal performance, RT² Profiler PCR Arrays should be used together with the RT² First Strand Kit for cDNA synthesis and RT² SYBR[®] Green qPCR Mastermixes for PCR.

Product	Contents	Cat. no.
RT ² First Strand Kit (12)	Enzymes and reagents for cDNA synthesis	330401
RT ² SYBR Green qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with real-time cyclers that do not require a reference dye, including: Bio-Rad models CFX96, CFX384, DNA Engine Opticon 2; Bio-Rad/MJ Research Chromo4; Roche LightCycler 480 (96-well and 384-well); all other cyclers	330500
RT ² SYBR Green ROX [™] qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the following real-time cyclers: Applied Biosystems models 5700, 7000, 7300, 7500 [Standard and FAST], 7700, 7900HT 96-well block [Standard and FAST] and 384-well block, StepOnePlus; Eppendorf Mastercycler ep realplex models 2, 2S, 4, 4S; Stratagene models Mx3000P, Mx3005P, Mx4000; Takara TP-800	330520
RT ² SYBR Green Fluor qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the following real-time cyclers: Bio-Rad models iCycler, iQ5, MyiQ, MyiQ2	330510

* Larger kit sizes available; please inquire.

RT² Profiler PCR Array products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN® (QIAGEN Group); Applied Biosystems®, ViiA™, StepOnePlus™, ROX™ (Applied Biosystems Corporation or its subsidiaries); Bio-Rad®, iCycler®, iQ™, MyiQ™, Chromo4™, CFX96™, DNA Engine Opticon®, CFX384™ (Bio-Rad Laboratories, Inc.); Stratagene®, Mx3005P®, Mx3000P®, Mx4000® (Stratagene); Eppendorf®, Mastercycler® (Eppendorf AG); Roche®, LightCycler® (Roche Group); Fluidigm® BioMark™ (Fluidigm Corporation); SYBR® (Molecular Probes, Inc.).

1066029 03/2011 © 2011 QIAGEN, all rights reserved.

www.qiagen.com

Canada ■ 800-572-9613

Ireland ■ 1800 555 049

Norway ■ 800-18859

China ■ 8621-3865-3865

Italy ■ 800-787980

Singapore ■ 1800-742-4368

Denmark ■ 80-885945

Japan ■ 03-6890-7300

Spain ■ 91-630-7050

Australia ■ 1-800-243-800

Finland ■ 0800-914416

Korea (South) ■ 080-000-7145

Sweden ■ 020-790282

Austria ■ 0800/281010

France ■ 01-60-920-930

Luxembourg ■ 8002 2076

Switzerland ■ 055-254-22-11

Belgium ■ 0800-79612

Germany ■ 02103-29-12000

Mexico ■ 01-800-7742-436

UK ■ 01293-422-911

Brazil ■ 0800-557779

Hong Kong ■ 800 933 965

The Netherlands ■ 0800 0229592

USA ■ 800-426-8157



Sample & Assay Technologies