# RT<sup>2</sup> Profiler PCR Array (Rotor-Gene® Format) Human Dopamine & Serotonin Pathway

Cat. no. 330231 PAHS-158ZR

#### For pathway expression analysis

Format	For use with the following real-time cyclers
RT <sup>2</sup> Profiler PCR Array,	Rotor-Gene Q, other Rotor-Gene cyclers
Format R	

#### **Description**

The Human Dopamine & Serotonin Pathway RT<sup>2</sup> Profiler PCR array profiles the expression of 84 genes associated with the dopamine and serotonin systems. Dopamine and serotonin are 2 of the major neurotransmitter systems in the mammalian nervous system. Dopamine affects brain processes that control both motor and emotional behavior and plays a role in the brain's reward mechanism. Serotonin is critical in temperature regulation, sensory perception, locomotion, sleep, and psychosis. Pharmacological agents targeting dopaminergic/serotoninergic neurotransmission have been clinically used to manage several neurological and psychiatric disorders including Parkinson's disease, schizophrenia, bipolar disorder, depression, attention deficit and hyperactivity disorder (ADHD), and addiction. Besides significant progress in understanding their structural, genetic and pharmacological properties, recent studies have uncovered the complexity, intricacy, and plasticity of intracellular signaling mechanisms involved in dopamine and serotonin receptor function. These receptors act through diverse G-protein coupled and G-protein independent mechanisms that trigger downstream intracellular signal transduction events involving the cAMP/PKA, PI-3Kinase/AKT, phospholipase A2 (PLA2), and phospholipase C (PLC) pathways. These pathways in turn regulate various functions including synthesis, transport and degradation of dopamine and serotonin as well as the transcriptional regulation key genes linked to multiple neuropathological conditions. Using real time PCR, research studies can easily and reliably analyze the expression of a focused panel of genes associated with and regulated by the dopamine and serotonin systems with this array.

For further details, consult the RT<sup>2</sup> Profiler PCR Array Handbook.

#### Shipping and storage

RT<sup>2</sup> Profiler PCR Arrays in the Rotor-Gene format are shipped at ambient temperature, on dry ice, or blue ice packs depending on destination and accompanying products.

For long term storage, keep plates at -20°C.

**Note**: Ensure that you have the correct RT<sup>2</sup> Profiler PCR Array format for your real-time cycler (see table above).

Note: Open the package and store the products appropriately immediately on receipt.



## **Array layout**

The 96 real-time assays in the Rotor-Gene format are located in wells 1–96 of the Rotor-Disc<sup>™</sup> (plate A1–A12=Rotor-Disc 1–12, plate B1–B12=Rotor-Disc 13–24, etc.). To maintain data analysis compatibility, wells 97–100 do not contain real-time assays but will contain master mix to account for weight balance.

## Gene table: RT<sup>2</sup> Profiler PCR Array

Position	UniGene	GenBank	Symbol	Description	
A01	Hs.192215	NM_021116	ADCY1	Adenylate cyclase 1 (brain)	
A02	Hs.481545	NM 020546	ADCY2	Adenylate cyclase 2 (brain)	
A03	Hs.708074	NM_004036	ADCY3	Adenylate cyclase 3	
A04	Hs.655144	NM 183357	ADCY5	Adenylate cyclase 5	
A05	Hs.99913	NM 000684	ADRB1	Adrenergic, beta-1-, receptor	
A06	Hs.591251	NM 000024	ADRB2	Adrenergic, beta-2-, receptor, surface	
A07	Hs.83636	NM_001619	ADRBK1	Adrenergic, beta, receptor kinase 1	
A08	Hs.657494	NM 005160	ADRBK2	Adrenergic, beta, receptor kinase 2	
A09	Hs.525622	NM 005163	AKT1	V-akt murine thymoma viral oncogene homolog 1	
A10	Hs.631535	NM 001626	AKT2	V-akt murine thymoma viral oncogene homolog 2	
A11	Hs.498292	NM 005465	AKT3	V-akt murine thymoma viral oncogene homolog 3 (protein kinase B, gamma	
A12	Hs.654431	NM 000697	ALOX12	Arachidonate 12-lipoxygenase	
B01	Hs.434980	NM 000484	APP	Amyloid beta (A4) precursor protein	
B02	Hs.503284	NM 004041	ARRB1	Arrestin, beta 1	
B03	Hs.435811	NM 004313	ARRB2	Arrestin, beta 2	
B04	Hs.502182	NM 001709	BDNF	Brain-derived neurotrophic factor	
B05	Hs.141125	NM_004346	CASP3	Caspase 3, apoptosis-related cysteine peptidase	
B06		NM 004935	CASP3 CDK5		
B06 B07	Hs.647078		COMT	Cyclin-dependent kinase 5	
	Hs.370408	NM_000754		Catechol-O-methyltransferase	
B08	Hs.516646	NM_004379	CREB1	CAMP responsive element binding protein 1	
B09	Hs.648256	NM_000106	CYP2D6	Cytochrome P450, family 2, subfamily D, polypeptide 6	
B10	Hs.591890	NM_000787	DBH	Dopamine beta-hydroxylase (dopamine beta-monooxygenase)	
B11	Hs.359698	NM_000790	DDC	Dopa decarboxylase (aromatic L-amino acid decarboxylase)	
B12	Hs.2624	NM_000794	DRD1	Dopamine receptor D1	
C01	Hs.73893	NM_000795	DRD2	Dopamine receptor D2	
C02	Hs.121478	NM_000796	DRD3	Dopamine receptor D3	
C03	Hs.99922	NM_000797	DRD4	Dopamine receptor D4	
C04	Hs.380681	NM_000798	DRD5	Dopamine receptor D5	
C05	Hs.171695	NM_004417	DUSP1	Dual specificity phosphatase 1	
C06	Hs.116092	NM_004441	EPHB1	EPH receptor B1	
C07	Hs.728789	NM_005252	FOS	FBJ murine osteosarcoma viral oncogene homolog	
C08	Hs.248114	NM_000514	GDNF	Glial cell derived neurotrophic factor	
C09	Hs.514227	NM_002055	GFAP	Glial fibrillary acidic protein	
C10	Hs.32959	NM 001004056	GRK4	G protein-coupled receptor kinase 4	
C11	Hs.524625	NM 005308	GRK5	G protein-coupled receptor kinase 5	
C12	Hs.235116	NM 002082	GRK6	G protein-coupled receptor kinase 6	
D01	Hs.466828	NM 019884	GSK3A	Glycogen synthase kinase 3 alpha	
D02	Hs.445733	NM 002093	GSK3B	Glycogen synthase kinase 3 beta	
D03	Hs.247940	NM 000524	HTR1A	5-hydroxytryptamine (serotonin) receptor 1A	
D04	Hs.123016	NM 000863	HTR1B	5-hydroxytryptamine (serotonin) receptor 1B	
D05	Hs.121482	NM 000864	HTR1D	5-hydroxytryptamine (serotonin) receptor 1D	
D06	Hs.1611	NM 000865	HTR1E	5-hydroxytryptamine (serotonin) receptor 1E	
D07	Hs.248136	NM 000866	HTR1F	5-hydroxytryptamine (serotonin) receptor 1F	
D07	Hs.654586	NM 000621	HTR2A	5-hydroxytryptamine (serotonin) receptor 2A	
D09	Hs.421649	NM 000867	HTR2B	5-hydroxytryptamine (serotonin) receptor 2B	
D10	Hs.149037	NM 000868	HTR2C	5-hydroxytryptamine (serotonin) receptor 2C	
D10	Hs.413899	NM 000869	HTR3A	5-hydroxytryptamine (serotonin) receptor 3A	
D11	Hs.241377	NM 006028	HTR3B	5-hydroxytryptamine (serotonin) receptor 3A 5-hydroxytryptamine (serotonin) receptor 3B	
E01	Hs.483773	NM 000870	HTR4	5-hydroxytryptamine (serotonin) receptor 4	
E02		NM 024012	HTR5A	, , , , , , , , ,	
	Hs.65791			5-hydroxytryptamine (serotonin) receptor 5A	
E03	Hs.22180	NM_000871	HTR6	5-hydroxytryptamine (serotonin) receptor 6	
E04	Hs.73739	NM_000872	HTR7	5-hydroxytryptamine (serotonin) receptor 7 (adenylate cyclase-coupled)	
E05	Hs.567295	NM_002222	ITPR1	Inositol 1,4,5-trisphosphate receptor, type 1	
E06	Hs.183109	NM_000240	MAOA	Monoamine oxidase A	
E07	Hs.654473	NM_000898	MAOB	Monoamine oxidase B	
E08	Hs.431850	NM_002745	MAPK1	Mitogen-activated protein kinase 1	
E09	Hs.524430	NM_002135	NR4A1	Nuclear receptor subfamily 4, group A, member 1	

Position	UniGene	GenBank	Symbol	Description	
E10	Hs.279522	NM_006981	NR4A3	Nuclear receptor subfamily 4, group A, member 3	
E11	Hs.348762	NM_006661	PDE10A	Phosphodiesterase 10A	
E12	Hs.89901	NM_006202	PDE4A	Phosphodiesterase 4A, cAMP-specific	
F01	Hs.198072	NM_002600	PDE4B	Phosphodiesterase 4B, cAMP-specific	
F02	Hs.132584	NM_000923	PDE4C	Phosphodiesterase 4C, cAMP-specific	
F03	Hs.117545	NM_006203	PDE4D	Phosphodiesterase 4D, cAMP-specific	
F04	Hs.22584	NM_024411	PDYN	Prodynorphin	
F05	Hs.553498	NM_006218	PIK3CA	Phosphoinositide-3-kinase, catalytic, alpha polypeptide	
F06	Hs.32942	NM_002649	PIK3CG	Phosphoinositide-3-kinase, catalytic, gamma polypeptide	
F07	Hs.319438	NM_000929	PLA2G5	Phospholipase A2, group V	
F08	Hs.431173	NM_015192	PLCB1	Phospholipase C, beta 1 (phosphoinositide-specific)	
F09	Hs.355888	NM_004573	PLCB2	Phospholipase C, beta 2	
F10	Hs.591953	NM_000932	PLCB3	Phospholipase C, beta 3 (phosphatidylinositol-specific)	
F11	Hs.286192	NM_181505	PPP1R1B	Protein phosphatase 1, regulatory (inhibitor) subunit 1B	
F12	Hs.631630	NM_002730	PRKACA	Protein kinase, cAMP-dependent, catalytic, alpha	
G01	Hs.196384	NM_000963	PTGS2	Prostaglandin-endoperoxide synthase 2 (prostaglandin G/H synthase and	
	113.170004		11032	cyclooxygenase)	
G02	Hs.158322	NM_003053	SLC18A1	Solute carrier family 18 (vesicular monoamine), member 1	
G03	Hs.654476	NM_003054	SLC18A2	Solute carrier family 18 (vesicular monoamine), member 2	
G04	Hs.406	NM_001044	SLC6A3	Solute carrier family 6 (neurotransmitter transporter, dopamine), member 3	
G05	Hs.591192	NM_001045	SLC6A4	Solute carrier family 6 (neurotransmitter transporter, serotonin), member 4	
G06	Hs.271771	NM_000345	SNCA	Synuclein, alpha (non A4 component of amyloid precursor)	
G07	Hs.426463	NM_005460	SNCAIP	Synuclein, alpha interacting protein	
G08	Hs.445503	NM_003178	SYN2	Synapsin II	
G09	Hs.183671	NM_005651	TDO2	Tryptophan 2,3-dioxygenase	
G10	Hs.435609	NW_000360	TH	Tyrosine hydroxylase	
G11	Hs.591999	NM_004179	TPH1	Tryptophan hydroxylase 1	
G12	Hs.376337	NM_173353	TPH2	Tryptophan hydroxylase 2	
H01	Hs.520640	NM_001101	ACTB	Actin, beta	
H02	Hs.534255	NM_004048	B2M	Beta-2-microglobulin	
H03	Hs.592355	NM_002046	GAPDH	Glyceraldehyde-3-phosphate dehydrogenase	
H04	Hs.412707	NM_000194	HPRT1	Hypoxanthine phosphoribosyltransferase 1	
H05	Hs.546285	NM_001002	RPLP0	Ribosomal protein, large, PO	
H06	N/A	SA_00105	HGDC	Human Genomic DNA Contamination	
H07	N/A	SA_00104	RTC	Reverse Transcription Control	
H08	N/A	SA_00104	RTC	Reverse Transcription Control	
H09	N/A	SA_00104	RTC	Reverse Transcription Control	
H10	N/A	SA_00103	PPC	Positive PCR Control	
H11	N/A	SA_00103	PPC	Positive PCR Control	
H12	N/A	SA 00103	PPC	Positive PCR Control	

# **Related products**

For optimal performance, RT<sup>2</sup> Profiler PCR Arrays should be used together with the RT<sup>2</sup> First Strand Kit for cDNA synthesis and RT<sup>2</sup> SYBR<sup>®</sup> Green qPCR Mastermixes for PCR.

Product	Contents	Cat. no.
RT <sup>2</sup> First Strand Kit (12)	Enzymes and reagents for cDNA synthesis	330401
RT <sup>2</sup> SYBR Green ROX <sup>™</sup> FAST Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the Rotor-Gene Q and other Rotor-Gene cyclers	330620

<sup>\*</sup> Larger kit sizes available; please inquire.

RT<sup>2</sup> Profiler PCR Array products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen. com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, Rotor-Gene®, Rotor-Disc™ (QIAGEN Group); ROX™ (Applera Corporation or its subsidiaries); SYBR® (Molecular Probes, Inc.).

1067688 03/2011 © 2011 QIAGEN, all rights reserved.

www.aiaaen.com Australia • 1-800-243-800 Austria • 0800/281010 Belgium • 0800-79612 Brazil • 0800-557779

Canada • 800-572-9613 China • 8621-3865-3865 Denmark ■ 80-885945 Finland • 0800-914416 France • 01-60-920-930 Germany • 02103-29-12000 Hong Kong • 800 933 965

Ireland • 1800 555 049 Italy • 800-787980 Japan ■ 03-6890-7300 Korea (South) • 080-000-7145 Luxembourg ■ 8002 2076 Mexico • 01-800-7742-436 The Netherlands • 0800 0229592 USA • 800-426-8157

Norway • 800-18859 Singapore ■ 1800-742-4368 Spain ■ 91-630-7050 Sweden • 020-790282 Switzerland ■ 055-254-22-11 UK • 01293-422-911

