

RT² Profiler PCR Array (Rotor-Gene[®] Format)

Mouse Insulin Signaling Pathway

Cat. no. 330231 PAMM-030ZR

For pathway expression analysis

Format	For use with the following real-time cyclers
RT ² Profiler PCR Array, Format R	Rotor-Gene Q, other Rotor-Gene cyclers

Description

The Mouse Insulin Signaling Pathway RT² Profiler PCR Array profiles the expression of 84 genes related to insulin-responsive genes. This array represents genes associated with the insulin receptor and target genes for insulin signaling. Genes involved in the metabolism of carbohydrates, lipids, and proteins and other related biological responses are also included. This array contains members of the PI3K and MAPK Pathways as well. Using real-time PCR, you can easily and reliably analyze expression of a focused panel of genes related to insulin-responsive with this array.

For further details, consult the *RT² Profiler PCR Array Handbook*.

Shipping and storage

RT² Profiler PCR Arrays in the Rotor-Gene format are shipped at ambient temperature, on dry ice, or blue ice packs depending on destination and accompanying products.

For long term storage, keep plates at -20°C.

Note: Ensure that you have the correct RT² Profiler PCR Array format for your real-time cycler (see table above).

Note: Open the package and store the products appropriately immediately on receipt.



Array layout

The 96 real-time assays in the Rotor-Gene format are located in wells 1–96 of the Rotor-Disc™ (plate A1–A12=Rotor-Disc 1–12, plate B1–B12=Rotor-Disc 13–24, etc.). To maintain data analysis compatibility, wells 97–100 do not contain real-time assays but will contain master mix to account for weight balance.

Gene table: RT² Profiler PCR Array

Position	UniGene	GenBank	Symbol	Description
A01	Mm.31374	NM_133360	Acaca	Acetyl-Coenzyme A carboxylase alpha
A02	Mm.356689	NM_015729	Acox1	Acyl-Coenzyme A oxidase 1, palmitoyl
A03	Mm.389380	NM_013460	Adra1d	Adrenergic receptor, alpha 1d
A04	Mm.4665	NM_009636	Aebp1	AE binding protein 1
A05	Mm.6645	NM_009652	Akt1	Thymoma viral proto-oncogene 1
A06	Mm.177194	NM_007434	Akt2	Thymoma viral proto-oncogene 2
A07	Mm.235194	NM_011785	Akt3	Thymoma viral proto-oncogene 3
A08	Mm.220946	NM_009703	Araf	V-raf murine sarcoma 3611 viral oncogene homolog
A09	Mm.238213	NM_009743	Bcl2l1	Bcl2-like 1
A10	Mm.245513	NM_139294	Braf	Braf transforming gene
A11	Mm.8687	NM_007598	Cap1	CAP, adenylate cyclase-associated protein 1 (yeast)
A12	Mm.266871	NM_007619	Cbl	Casitas B-lineage lymphoma
B01	Mm.349667	NM_007678	Cebpa	CCAAT/enhancer binding protein (C/EBP), alpha
B02	Mm.439656	NM_009883	Cebpb	CCAAT/enhancer binding protein (C/EBP), beta
B03	Mm.4407	NM_013459	Cfd	Complement factor D (adipsin)
B04	Mm.156	NM_010070	Dok1	Docking protein 1
B05	Mm.243323	NM_010071	Dok2	Docking protein 2
B06	Mm.33910	NM_013739	Dok3	Docking protein 3
B07	Mm.240885	NM_019819	Dusp14	Dual specificity phosphatase 14
B08	Mm.28839	NM_145371	Eif2b1	Eukaryotic translation initiation factor 2B, subunit 1 (alpha)
B09	Mm.6700	NM_007918	Eif4ebp1	Eukaryotic translation initiation factor 4E binding protein 1
B10	Mm.280913	NM_007948	Ercc1	Excision repair cross-complementing rodent repair deficiency, complementation group 1
B11	Mm.423078	NM_019395	Fbp1	Fructose bisphosphatase 1
B12	Mm.246513	NM_010234	Fos	FBJ osteosarcoma oncogene
C01	Mm.135965	NM_177798	Frs2	Fibroblast growth factor receptor substrate 2
C02	Mm.392206	NM_144939	Frs3	Fibroblast growth factor receptor substrate 3
C03	Mm.18064	NM_008061	G6pc	Glucose-6-phosphatase, catalytic
C04	Mm.140768	NM_021331	G6pc2	Glucose-6-phosphatase, catalytic, 2
C05	Mm.277409	NM_021356	Gab1	Growth factor receptor bound protein 2-associated protein 1
C06	Mm.220358	NM_010292	Gck	Glucokinase
C07	Mm.252391	NM_010271	Gpd1	Glycerol-3-phosphate dehydrogenase 1 (soluble)
C08	Mm.273117	NM_010345	Grb10	Growth factor receptor bound protein 10
C09	Mm.439649	NM_008163	Grb2	Growth factor receptor bound protein 2
C10	Mm.394930	NM_019827	Gsk3b	Glycogen synthase kinase 3 beta
C11	Mm.255848	NM_013820	Hk2	Hexokinase 2
C12	Mm.334313	NM_008284	Hras1	Harvey rat sarcoma virus oncogene 1
D01	Mm.275742	NM_010513	Igf1r	Insulin-like growth factor I receptor
D02	Mm.3862	NM_010514	Igf2	Insulin-like growth factor 2
D03	Mm.21300	NM_008341	Igfbp1	Insulin-like growth factor binding protein 1
D04	Mm.46269	NM_008386	Ins1	Insulin I
D05	Mm.476857	NM_013564	InsI3	Insulin-like 3
D06	Mm.4952	NM_010570	Irs1	Insulin receptor substrate 1
D07	Mm.407207	NM_001081212	Irs2	Insulin receptor substrate 2
D08	Mm.275071	NM_010591	Jun	Jun oncogene
D09	Mm.4292	NM_013692	Klf10	Kruppel-like factor 10
D10	Mm.383182	NM_021284	Kras	V-Ki-ras2 Kirsten rat sarcoma viral oncogene homolog
D11	Mm.3213	NM_010700	Ldlr	Low density lipoprotein receptor
D12	Mm.277072	NM_008493	Lep	Lepin
E01	Mm.248907	NM_008927	Map2k1	Mitogen-activated protein kinase kinase 1
E02	Mm.196581	NM_011949	Mapk1	Mitogen-activated protein kinase 1
E03	Mm.21158	NM_020009	Mtor	Mechanistic target of rapamycin (serine/threonine kinase)
E04	Mm.181485	NM_010878	Nck1	Non-catalytic region of tyrosine kinase adaptor protein 1
E05	Mm.2893	NM_010927	Nos2	Nitric oxide synthase 2, inducible
E06	Mm.154796	NM_023456	Npy	Neuropeptide Y
E07	Mm.29856	NM_028994	Pck2	Phosphoenolpyruvate carboxykinase 2 (mitochondrial)
E08	Mm.10504	NM_011062	Pdpk1	3-phosphoinositide dependent protein kinase 1

Position	UniGene	GenBank	Symbol	Description
E09	Mm.260521	NM_008839	Pik3ca	Phosphatidylinositol 3-kinase, catalytic, alpha polypeptide
E10	Mm.213128	NM_029094	Pik3cb	Phosphatidylinositol 3-kinase, catalytic, beta polypeptide
E11	Mm.259333	NM_001024955	Pik3r1	Phosphatidylinositol 3-kinase, regulatory subunit, polypeptide 1 (p85 alpha)
E12	Mm.12945	NM_008841	Pik3r2	Phosphatidylinositol 3-kinase, regulatory subunit, polypeptide 2 (p85 beta)
F01	Mm.383180	NM_013631	Pklr	Pyruvate kinase liver and red blood cell
F02	Mm.3020	NM_011146	Pparg	Peroxisome proliferator activated receptor gamma
F03	Mm.1970	NM_031868	Ppp1ca	Protein phosphatase 1, catalytic subunit, alpha isoform
F04	Mm.7980	NM_011102	Prkcc	Protein kinase C, gamma
F05	Mm.291554	NM_008857	Prki	Protein kinase C, iota
F06	Mm.28561	NM_008860	Prkcz	Protein kinase C, zeta
F07	Mm.1270	NM_011164	Prl	Prolactin
F08	Mm.277916	NM_011201	Ptfn1	Protein tyrosine phosphatase, non-receptor type 1
F09	Mm.29855	NM_011213	Ptprf	Protein tyrosine phosphatase, receptor type, F
F10	Mm.184163	NM_029780	Raf1	V-raf-leukemia viral oncogene 1
F11	Mm.1181	NM_022984	Retn	Resistin
F12	Mm.301827	NM_009097	Rps6ka1	Ribosomal protein S6 kinase polypeptide 1
G01	Mm.389894	NM_009101	Rras	Harvey rat sarcoma oncogene, subgroup R
G02	Mm.276572	NM_025846	Rras2	Related RAS viral (r-ras) oncogene homolog 2
G03	Mm.250422	NM_008871	Serpine1	Serine (or cysteine) peptidase inhibitor, clade E, member 1
G04	Mm.86595	NM_011368	Shc1	Src homology 2 domain-containing transforming protein C1
G05	Mm.330113	NM_011989	Slc27a4	Solute carrier family 27 (fatty acid transporter), member 4
G06	Mm.21002	NM_011400	Slc2a1	Solute carrier family 2 (facilitated glucose transporter), member 1
G07	Mm.210815	NM_009166	Sorbs1	Sorbin and SH3 domain containing 1
G08	Mm.434583	NM_009231	Sos1	Son of sevenless homolog 1 (Drosophila)
G09	Mm.278701	NM_011480	Srebf1	Sterol regulatory element binding transcription factor 1
G10	Mm.441333	NM_009375	Tg	Thyroglobulin
G11	Mm.4177	NM_009463	Ucp1	Uncoupling protein 1 (mitochondrial, proton carrier)
G12	Mm.282184	NM_009505	Vegfa	Vascular endothelial growth factor A
H01	Mm.328431	NM_007393	Actb	Actin, beta
H02	Mm.163	NM_009735	B2m	Beta-2 microglobulin
H03	Mm.343110	NM_008084	Gapdh	Glyceraldehyde-3-phosphate dehydrogenase
H04	Mm.3317	NM_010368	Gusb	Glucuronidase, beta
H05	Mm.2180	NM_008302	Hsp90ab1	Heat shock protein 90 alpha (cytosolic), class B member 1
H06	N/A	SA_00106	MGDC	Mouse Genomic DNA Contamination
H07	N/A	SA_00104	RTC	Reverse Transcription Control
H08	N/A	SA_00104	RTC	Reverse Transcription Control
H09	N/A	SA_00104	RTC	Reverse Transcription Control
H10	N/A	SA_00103	PPC	Positive PCR Control
H11	N/A	SA_00103	PPC	Positive PCR Control
H12	N/A	SA_00103	PPC	Positive PCR Control

Related products

For optimal performance, RT² Profiler PCR Arrays should be used together with the RT² First Strand Kit for cDNA synthesis and RT² SYBR[®] Green qPCR Mastermixes for PCR.

Product	Contents	Cat. no.
RT ² First Strand Kit (12)	Enzymes and reagents for cDNA synthesis	330401
RT ² SYBR Green ROX [™] FAST Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the Rotor-Gene Q and other Rotor-Gene cyclers	330620

* Larger kit sizes available; please inquire.

RT² Profiler PCR Array products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN[®], Rotor-Gene[®], Rotor-Disc[™] (QIAGEN Group); ROX[™] (Applied Biosystems or its subsidiaries); SYBR[®] (Molecular Probes, Inc.).

1067688 03/2011 © 2011 QIAGEN, all rights reserved.

www.qiagen.com

Canada ■ 800-572-9613

Ireland ■ 1800 555 049

Norway ■ 800-18859

China ■ 8621-3865-3865

Italy ■ 800-787980

Singapore ■ 1800-742-4368

Denmark ■ 80-885945

Japan ■ 03-6890-7300

Spain ■ 91-630-7050

Australia ■ 1-800-243-800

Finland ■ 0800-914416

Korea (South) ■ 080-000-7145

Sweden ■ 020-790282

Austria ■ 0800/281010

France ■ 01-60-920-930

Luxembourg ■ 8002 2076

Switzerland ■ 055-254-22-11

Belgium ■ 0800-79612

Germany ■ 02103-29-12000

Mexico ■ 01-800-7742-436

UK ■ 01293-422-911

Brazil ■ 0800-557779

Hong Kong ■ 800 933 965

The Netherlands ■ 0800 0229592

USA ■ 800-426-8157



Sample & Assay Technologies