

RT² Profiler PCR Array (96-Well Format and 384-Well [4 x 96] Format)

Mouse Gap Junctions

Cat. no. 330231 PAMM-144ZA

For pathway expression analysis

Format	For use with the following real-time cyclers
RT ² Profiler PCR Array, Format A	Applied Biosystems® models 5700, 7000, 7300, 7500, 7700, 7900HT, ViiA™ 7 (96-well block); Bio-Rad® models iCycler®, iQ™ 5, MyiQ™, MyiQ2; Bio-Rad/MJ Research Chromo4™; Eppendorf® Mastercycler® ep realplex models 2, 2s, 4, 4s; Stratagene® models Mx3005P®, Mx3000P®; Takara TP-800
RT ² Profiler PCR Array, Format C	Applied Biosystems models 7500 (Fast block), 7900HT (Fast block), StepOnePlus™, ViiA 7 (Fast block)
RT ² Profiler PCR Array, Format D	Bio-Rad CFX96™; Bio-Rad/MJ Research models DNA Engine Opticon®, DNA Engine Opticon 2; Stratagene Mx4000®
RT ² Profiler PCR Array, Format E	Applied Biosystems models 7900HT (384-well block), ViiA 7 (384-well block); Bio-Rad CFX384™
RT ² Profiler PCR Array, Format F	Roche® LightCycler® 480 (96-well block)
RT ² Profiler PCR Array, Format G	Roche LightCycler 480 (384-well block)
RT ² Profiler PCR Array, Format H	Fluidigm® BioMark™



Sample & Assay Technologies

Description

The Mouse Gap Junctions RT² Profiler PCR Array profiles the expression of 84 key genes encoding components, interactors, and regulators of gap junctions. Gap junctions are comprised of complexes of innexins and connexins whose extracellular domains dimerize with similar complexes in an adjacent cell. These intercellular complexes form channels, each with its own conductance and molecular permeability. The channels connect cytoplasm, allowing ions and small molecules to pass through and mediating communication between the adjacent cells. Cell surface receptors for neurotransmitters, cytokines, growth factors, and lysophosphatidic acid activate protein kinase, G-protein, and secondary messenger signaling pathways to regulate gap junctions via changes in connexin phosphorylation and membrane potential. The connexins directly bind tubulins, and their downstream signaling pathways regulate and recruit microtubules to help define cell shape and mediate intracellular transport. A wide variety of cell types express gap junctions including cardiomyocytes, keratinocytes, astrocytes, endothelial cells, and smooth muscle cells. Gap junctions regulate many biological processes such as cellular growth and differentiation, embryonic development, excitable cell contraction, immune responses, neural activity, tissue homeostasis, and metabolic transport. Mutations in connexin genes or other disruptions of gap junction function contribute to the pathophysiology of cardiovascular disease, neurological disorders, and developmental abnormalities. Profiling the expression of gap junction components and regulators may lead to a better understanding of molecular mechanisms behind gap-junction-mediated cell biology. Using real-time PCR, research studies can easily and reliably analyze the expression of a focused panel of genes involved in gap junctions with this array.

For further details, consult the *RT² Profiler PCR Array Handbook*.

Shipping and storage

RT² Profiler PCR Arrays in formats A, C, D, E, F, and G are shipped at ambient temperature, on dry ice, or blue ice packs depending on destination and accompanying products. RT² Profiler PCR Arrays in format H are shipped on dry ice or blue ice packs.

For long term storage, keep plates at -20°C .

Note: Ensure that you have the correct RT² Profiler PCR Array format for your real-time cycler (see table above).

Note: Open the package and store the products appropriately immediately on receipt.

Array layout (96-well)

For 384-well 4 x 96 PCR arrays, genes are present in a staggered format. Refer to the *RT² Profiler PCR Array Handbook* for layout.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Adcy1	Adcy2	Adcy3	Adrb2	Cav1	Cdk1	Csnk1d	Ctnnb1	Dbn1	Egfr	Gja1	Gja3
B	Gja4	Gja5	Gja6	Gja8	Gjb1	Gjb2	Gjb3	Gjb4	Gjb5	Gjb6	Gjc2	Gjc3
C	Gjd2	Gje1	Gnai1	Grb2	Grm1	Gucy1a3	Gucy1b3	Hras1	Htr2a	Itpr1	Itpr2	Kras
D	Lpar1	Map2k1	Map2k2	Map2k5	Map3k2	Mapk1	Mapk3	Mapk7	Nov	Nras	Panx1	Panx2
E	Panx3	Pdgfra	Pdgfrb	Plcb1	Plcb2	Prkaca	Prkacb	Prkca	Prkcb	Prkcc	Prkg1	Prkg2
F	Raf1	Sos1	Sos2	Src	Tjap1	Tjp1	Tjp2	Tuba1a	Tuba1b	Tuba1c	Tuba3a	Tuba4a
G	Tuba8	Tubb1	Tubb2a	Tubb2b	Tubb2c	Tubb3	Tubb4	Tubb5	Tubb6	Tubd1	Tube1	Tubg1
H	Actb	B2m	Gapdh	Gusb	Hsp90ab1	MGDC	RTC	RTC	RTC	PPC	PPC	PPC

Gene table: RT² Profiler PCR Array

Position	UniGene	GenBank	Symbol	Description
A01	Mm.259733	NM_009622	Adcy1	Adenylate cyclase 1
A02	Mm.390617	NM_153534	Adcy2	Adenylate cyclase 2
A03	Mm.70546	NM_138305	Adcy3	Adenylate cyclase 3
A04	Mm.5598	NM_007420	Adrb2	Adrenergic receptor, beta 2
A05	Mm.28278	NM_007616	Cav1	Caveolin 1, caveolae protein
A06	Mm.281367	NM_007659	Cdk1	Cyclin-dependent kinase 1
A07	Mm.216227	NM_139059	Csnk1d	Casein kinase 1, delta
A08	Mm.291928	NM_007614	Ctnnb1	Catenin (cadherin associated protein), beta 1
A09	Mm.19016	NM_019813	Dbn1	Drebrin 1
A10	Mm.8534	NM_007912	Egfr	Epidermal growth factor receptor
A11	Mm.378921	NM_010288	Gja1	Gap junction protein, alpha 1
A12	Mm.57207	NM_016975	Gja3	Gap junction protein, alpha 3
B01	Mm.24615	NM_008120	Gja4	Gap junction protein, alpha 4
B02	Mm.281816	NM_008121	Gja5	Gap junction protein, alpha 5
B03	Mm.123113	NM_001001496	Gja6	Gap junction protein, alpha 6
B04	Mm.468160	NM_008123	Gja8	Gap junction protein, alpha 8
B05	Mm.21198	NM_008124	Gjb1	Gap junction protein, beta 1
B06	Mm.390683	NM_008125	Gjb2	Gap junction protein, beta 2
B07	Mm.90003	NM_008126	Gjb3	Gap junction protein, beta 3
B08	Mm.56906	NM_008127	Gjb4	Gap junction protein, beta 4
B09	Mm.26859	NM_010291	Gjb5	Gap junction protein, beta 5
B10	Mm.25652	NM_008128	Gjb6	Gap junction protein, beta 6
B11	Mm.40016	NM_175452	Gjc2	Gap junction protein, gamma 2
B12	Mm.377888	NM_080450	Gjc3	Gap junction protein, gamma 3
C01	Mm.389394	NM_010290	Gjd2	Gap junction protein, delta 2
C02	Mm.332771	NM_029722	Gje1	Gap junction protein, epsilon 1
C03	Mm.254629	NM_010305	Gnai1	Guanine nucleotide binding protein (G protein), alpha inhibiting 1
C04	Mm.439649	NM_008163	Grb2	Growth factor receptor bound protein 2
C05	Mm.391904	NM_016976	Grm1	Glutamate receptor, metabotropic 1
C06	Mm.143831	NM_021896	Gucy1a3	Guanylate cyclase 1, soluble, alpha 3
C07	Mm.9445	NM_017469	Gucy1b3	Guanylate cyclase 1, soluble, beta 3
C08	Mm.334313	NM_008284	Hras1	Harvey rat sarcoma virus oncogene 1
C09	Mm.214351	NM_172812	Htr2a	5-hydroxytryptamine (serotonin) receptor 2A
C10	Mm.227912	NM_010585	Itpr1	Inositol 1,4,5-trisphosphate receptor 1
C11	Mm.393003	NM_019923	Itpr2	Inositol 1,4,5-trisphosphate receptor 2
C12	Mm.383182	NM_021284	Kras	V-Ki-ras2 Kirsten rat sarcoma viral oncogene homolog
D01	Mm.4772	NM_010336	Lpar1	Lysophosphatidic acid receptor 1
D02	Mm.248907	NM_008927	Map2k1	Mitogen-activated protein kinase kinase 1
D03	Mm.275436	NM_023138	Map2k2	Mitogen-activated protein kinase kinase 2
D04	Mm.325746	NM_011840	Map2k5	Mitogen-activated protein kinase kinase 5
D05	Mm.211762	NM_011946	Map3k2	Mitogen-activated protein kinase kinase kinase 2
D06	Mm.196581	NM_011949	Mapk1	Mitogen-activated protein kinase 1
D07	Mm.8385	NM_011952	Mapk3	Mitogen-activated protein kinase 3
D08	Mm.38172	NM_011841	Mapk7	Mitogen-activated protein kinase 7
D09	Mm.5167	NM_010930	Nov	Nephroblastoma overexpressed gene

Position	UniGene	GenBank	Symbol	Description
D10	Mm.400954	NM_010937	Nras	Neuroblastoma ras oncogene
D11	Mm.142253	NM_019482	Panx1	Pannexin 1
D12	Mm.483059	NM_001002005	Panx2	Pannexin 2
E01	Mm.217159	NM_172454	Panx3	Pannexin 3
E02	Mm.221403	NM_011058	Pdgfra	Platelet derived growth factor receptor, alpha polypeptide
E03	Mm.4146	NM_008809	Pdgfrb	Platelet derived growth factor receptor, beta polypeptide
E04	Mm.330607	NM_019677	Plcb1	Phospholipase C, beta 1
E05	Mm.215156	NM_177568	Plcb2	Phospholipase C, beta 2
E06	Mm.19111	NM_008854	Prkaca	Protein kinase, cAMP dependent, catalytic, alpha
E07	Mm.16766	NM_011100	Prkacb	Protein kinase, cAMP dependent, catalytic, beta
E08	Mm.222178	NM_011101	Prkca	Protein kinase C, alpha
E09	Mm.207496	NM_008855	Prkcb	Protein kinase C, beta
E10	Mm.7980	NM_011102	Prkcc	Protein kinase C, gamma
E11	Mm.381172	NM_011160	Prkg1	Protein kinase, cGMP-dependent, type I
E12	Mm.263002	NM_008926	Prkg2	Protein kinase, cGMP-dependent, type II
F01	Mm.184163	NM_029780	Raf1	V-raf-leukemia viral oncogene 1
F02	Mm.434583	NM_009231	Sos1	Son of sevenless homolog 1 (Drosophila)
F03	Mm.3770	NM_001135559	Sos2	Son of sevenless homolog 2 (Drosophila)
F04	Mm.22845	NM_009271	Src	Rous sarcoma oncogene
F05	Mm.284594	NM_028751	Tjap1	Tight junction associated protein 1
F06	Mm.4342	NM_009386	Tjp1	Tight junction protein 1
F07	Mm.104744	NM_011597	Tjp2	Tight junction protein 2
F08	Mm.439690	NM_011653	Tuba1a	Tubulin, alpha 1A
F09	Mm.392113	NM_011654	Tuba1b	Tubulin, alpha 1B
F10	Mm.88212	NM_009448	Tuba1c	Tubulin, alpha 1C
F11	Mm.287784	NM_009446	Tuba3a	Tubulin, alpha 3A
F12	Mm.1155	NM_009447	Tuba4a	Tubulin, alpha 4A
G01	Mm.32884	NM_017379	Tuba8	Tubulin, alpha 8
G02	Mm.45285	NM_001080971	Tubb1	Tubulin, beta 1
G03	Mm.422827	NM_009450	Tubb2a	Tubulin, beta 2a
G04	Mm.379227	NM_023716	Tubb2b	Tubulin, beta 2b
G05	Mm.227260	NM_146116	Tubb2c	Tubulin, beta 2c
G06	Mm.40068	NM_023279	Tubb3	Tubulin, beta 3
G07	Mm.7420	NM_009451	Tubb4	Tubulin, beta 4
G08	Mm.273538	NM_011655	Tubb5	Tubulin, beta 5
G09	Mm.181860	NM_026473	Tubb6	Tubulin, beta 6
G10	Mm.383207	NM_019756	Tubd1	Tubulin, delta 1
G11	Mm.297956	NM_028006	Tube1	Epsilon-tubulin 1
G12	Mm.479145	NM_134024	Tubg1	Tubulin, gamma 1
H01	Mm.328431	NM_007393	Actb	Actin, beta
H02	Mm.163	NM_009735	B2m	Beta-2 microglobulin
H03	Mm.343110	NM_008084	Gapdh	Glyceraldehyde-3-phosphate dehydrogenase
H04	Mm.3317	NM_010368	Gusb	Glucuronidase, beta
H05	Mm.2180	NM_008302	Hsp90ab1	Heat shock protein 90 alpha (cytosolic), class B member 1
H06	N/A	SA_00106	MGDC	Mouse Genomic DNA Contamination
H07	N/A	SA_00104	RTC	Reverse Transcription Control
H08	N/A	SA_00104	RTC	Reverse Transcription Control
H09	N/A	SA_00104	RTC	Reverse Transcription Control
H10	N/A	SA_00103	PPC	Positive PCR Control
H11	N/A	SA_00103	PPC	Positive PCR Control
H12	N/A	SA_00103	PPC	Positive PCR Control

Related products

For optimal performance, RT² Profiler PCR Arrays should be used together with the RT² First Strand Kit for cDNA synthesis and RT² SYBR[®] Green qPCR Mastermixes for PCR.

Product	Contents	Cat. no.
RT ² First Strand Kit (12)	Enzymes and reagents for cDNA synthesis	330401
RT ² SYBR Green qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with real-time cyclers that do not require a reference dye, including: Bio-Rad models CFX96, CFX384, DNA Engine Opticon 2; Bio-Rad/MJ Research Chromo4; Roche LightCycler 480 (96-well and 384-well); all other cyclers	330500
RT ² SYBR Green ROX™ qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the following real-time cyclers: Applied Biosystems models 5700, 7000, 7300, 7500 [Standard and FAST], 7700, 7900HT 96-well block [Standard and FAST] and 384-well block, StepOnePlus; Eppendorf Mastercycler ep realplex models 2, 2S, 4, 4S; Stratagene models Mx3000P, Mx3005P, Mx4000; Takara TP-800	330520
RT ² SYBR Green Fluor qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the following real-time cyclers: Bio-Rad models iCycler, iQ5, MyiQ, MyiQ2	330510

* Larger kit sizes available; please inquire.

RT² Profiler PCR Array products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN® (QIAGEN Group); Applied Biosystems®, ViiA™, StepOnePlus™, ROX™ (Applied Biosystems Corporation or its subsidiaries); Bio-Rad®, iCycler®, iQ™, MyiQ™, Chromo4™, CFX96™, DNA Engine Opticon®, CFX384™ (Bio-Rad Laboratories, Inc.); Stratagene®, Mx3005P®, Mx3000P®, Mx4000® (Stratagene); Eppendorf®, Mastercycler® (Eppendorf AG); Roche®, LightCycler® (Roche Group); Fluidigm® BioMark™ (Fluidigm Corporation); SYBR® (Molecular Probes, Inc.).

1066029 03/2011 © 2011 QIAGEN, all rights reserved.

www.qiagen.com

Canada ■ 800-572-9613

Ireland ■ 1800 555 049

Norway ■ 800-18859

China ■ 8621-3865-3865

Italy ■ 800-787980

Singapore ■ 1800-742-4368

Denmark ■ 80-885945

Japan ■ 03-6890-7300

Spain ■ 91-630-7050

Australia ■ 1-800-243-800

Finland ■ 0800-914416

Korea (South) ■ 080-000-7145

Sweden ■ 020-790282

Austria ■ 0800/281010

France ■ 01-60-920-930

Luxembourg ■ 8002 2076

Switzerland ■ 055-254-22-11

Belgium ■ 0800-79612

Germany ■ 02103-29-12000

Mexico ■ 01-800-7742-436

UK ■ 01293-422-911

Brazil ■ 0800-557779

Hong Kong ■ 800 933 965

The Netherlands ■ 0800 0229592

USA ■ 800-426-8157



Sample & Assay Technologies