

RT² Profiler PCR Array (Rotor-Gene® Format)

Mouse Osteoporosis

Cat. no. 330231 PAMM-170ZR

For pathway expression analysis

Format	For use with the following real-time cyclers
RT ² Profiler PCR Array, Format R	Rotor-Gene Q, other Rotor-Gene cyclers

Description

The Mouse Osteoporosis RT² Profiler PCR Array profiles the expression of 84 genes involved in pathogenesis of osteoporosis (OP). Advanced age, gender, and immobilization are major risk factors for developing OP, and additional contributing factors include diminished sex steroid production in post-menopausal women. OP is a metabolic disorder of the bones characterized by low bone mineral density (BMD) and increased incidence of fractures due to disruption of bone remodeling — the balance between bone resorption and bone formation. Bone remodeling is conducted by osteoclasts (cells responsible for bone resorption) and by osteoblasts (cells responsible for bone formation). Osteoblasts have a central role in bone metabolism and are responsible for bone matrix synthesis and mineralization, synthesis of growth factors and hormones, and regulation of osteoclastogenesis for bone resorption. In OP, a pathological imbalance in the bone remodeling process is typically linked to a disrupted RANKL/OPG signaling equilibrium wherein elevated RANKL levels favor resorption through osteoclast formation, function, and survival with lowered BMD. Recent evidence also suggests that inflammation plays a significant role in disrupting osteoclast–osteoblast equilibrium, which affects BMD. Enormous research efforts are underway to determine the molecular mechanisms of pathogenesis of OP with the aim of obtaining novel targets for its treatment and prevention as well as the identification of early diagnostic markers. The genes profiled with this array are associated with osteoblast and osteoclast activity including WNT and BMP signaling pathways, ECM and bone matrix remodeling, and cytokines and growth factors currently associated with OP molecular pathogenesis. A set of controls present on each array enables data analysis using the $\Delta\Delta\text{CT}$ method of relative quantification and assessment of reverse transcription performance, genomic DNA contamination, and PCR performance. Using real-time PCR, research studies can easily and reliably analyze the expression of a

focused panel of genes involved in osteoporosis with this array.

For further details, consult the *RT² Profiler PCR Array Handbook*.

Shipping and storage

RT² Profiler PCR Arrays in the Rotor-Gene format are shipped at ambient temperature, on dry ice, or blue ice packs depending on destination and accompanying products.

For long term storage, keep plates at –20°C.

Note: Ensure that you have the correct RT² Profiler PCR Array format for your real-time cycler (see table above).

Note: Open the package and store the products appropriately immediately on receipt.



Array layout

The 96 real-time assays in the Rotor-Gene format are located in wells 1–96 of the Rotor-Disc™ (plate A1–A12=Rotor-Disc 1–12, plate B1–B12=Rotor-Disc 13–24, etc.). To maintain data analysis compatibility, wells 97–100 do not contain real-time assays but will contain master mix to account for weight balance.

Gene table: RT² Profiler PCR Array

Position	UniGene	GenBank	Symbol	Description
A01	Mm.46354	NM_007388	Acp5	Acid phosphatase 5, tartrate resistant
A02	Mm.66952	NM_173029	Adcy10	Adenylate cyclase 10
A03	Mm.12286	NM_007440	Alox12	Arachidonate 12-lipoxygenase
A04	Mm.4584	NM_009660	Alox15	Arachidonate 15-lipoxygenase
A05	Mm.41072	NM_009662	Alox5	Arachidonate 5-lipoxygenase
A06	Mm.288186	NM_007431	Alpl	Alkaline phosphatase, liver/bone/kidney
A07	Mm.39005	NM_013476	Ar	Androgen receptor
A08	Mm.389459	NM_007541	Bglap	Bone gamma carboxyglutamate protein
A09	Mm.103205	NM_007553	Bmp2	Bone morphogenetic protein 2
A10	Mm.595	NM_007557	Bmp7	Bone morphogenetic protein 7
A11	Mm.4361	NM_007587	Calca	Calcitonin/calcitonin-related polypeptide, alpha
A12	Mm.4642	NM_007588	Calcr	Calcitonin receptor
B01	Mm.1186	NM_009801	Car2	Carbonic anhydrase 2
B02	Mm.103619	NM_013803	Casr	Calcium-sensing receptor
B03	Mm.271833	NM_011611	Cd40	CD40 antigen
B04	Mm.270587	NM_011930	Clcn7	Chloride channel 7
B05	Mm.297251	NM_009924	Cnr2	Cannabinoid receptor 2 (macrophage)
B06	Mm.277735	NM_007742	Col1a1	Collagen, type I, alpha 1
B07	Mm.277792	NM_007743	Col1a2	Collagen, type I, alpha 2
B08	Mm.100940	NM_007744	Comt	Catechol-O-methyltransferase
B09	Mm.20904	NM_019922	Crtap	Cartilage associated protein
B10	Mm.272085	NM_007802	Ctsk	Cathepsin K
B11	Mm.1262	NM_007809	Cyp17a1	Cytochrome P450, family 17, subfamily a, polypeptide 1
B12	Mm.5199	NM_007810	Cyp19a1	Cytochrome P450, family 19, subfamily a, polypeptide 1
C01	Mm.378235	NM_016974	Dbp	D site albumin promoter binding protein
C02	Mm.214717	NM_010051	Dkk1	Dickkopf homolog 1 (Xenopus laevis)
C03	Mm.27254	NM_008813	Enpp1	Ectonucleotide pyrophosphatase/phosphodiesterase 1
C04	Mm.9213	NM_007956	Esr1	Estrogen receptor 1 (alpha)
C05	Mm.2561	NM_010157	Esr2	Estrogen receptor 2 (beta)
C06	Mm.386776	NM_007953	Esrra	Estrogen related receptor, alpha
C07	Mm.265716	NM_010206	Fgfr1	Fibroblast growth factor receptor 1
C08	Mm.16340	NM_010207	Fgfr2	Fibroblast growth factor receptor 2
C09	Mm.389327	NM_010285	Ghrh	Growth hormone releasing hormone
C10	Mm.28328	NM_008288	Hsd11b1	Hydroxysteroid 11-beta dehydrogenase 1
C11	Mm.268521	NM_010512	Igf1	Insulin-like growth factor 1
C12	Mm.141936	NM_008342	Igfbp2	Insulin-like growth factor binding protein 2
D01	Mm.4392	NM_008357	Il15	Interleukin 15
D02	Mm.1019	NM_031168	Il6	Interleukin 6
D03	Mm.2856	NM_010559	Il6ra	Interleukin 6 receptor, alpha
D04	Mm.482186	NM_001033228	Itga1	Integrin alpha 1
D05	Mm.87150	NM_016780	Itgb3	Integrin beta 3
D06	Mm.277072	NM_008493	Lep	Leptin
D07	Mm.27961	NM_019782	Lepre1	Leprean 1
D08	Mm.271854	NM_008512	Lrp1	Low density lipoprotein receptor-related protein 1
D09	Mm.274581	NM_008513	Lrp5	Low density lipoprotein receptor-related protein 5
D10	Mm.321990	NM_008514	Lrp6	Low density lipoprotein receptor-related protein 6
D11	Mm.87787	NM_010735	Lta	Lymphotoxin A
D12	Mm.3900	NM_013589	Ltbp2	Latent transforming growth factor beta binding protein 2
E01	Mm.389466	NM_011839	Mab2112	Mab-21-like 2 (C. elegans)
E02	Mm.29564	NM_008610	Mmp2	Matrix metalloproteinase 2
E03	Mm.3514	NM_010834	Mstn	Myostatin
E04	Mm.89959	NM_010840	Mthfr	5,10-methylenetetrahydrofolate reductase
E05	Mm.329560	NM_016791	Nfatc1	Nuclear factor of activated T-cells, cytoplasmic, calcineurin-dependent 1
E06	Mm.135266	NM_008711	Nog	Noggin
E07	Mm.258415	NM_008713	Nos3	Nitric oxide synthase 3, endothelial cell
E08	Mm.154796	NM_023456	Npy	Neuropeptide Y
E09	Mm.129481	NM_008173	Nr3c1	Nuclear receptor subfamily 3, group C, member 1

Position	UniGene	GenBank	Symbol	Description
E10	Mm.42026	NM_011027	P2rx7	Purinergic receptor P2X, ligand-gated ion channel, 7
E11	Mm.79983	NM_011961	Plod2	Procollagen lysine, 2-oxoglutarate 5-dioxygenase 2
E12	Mm.1270	NM_011164	Prl	Prolactin
F01	Mm.284168	NM_020623	Pth	Parathyroid hormone
F02	Mm.3542	NM_011199	Pth1r	Parathyroid hormone 1 receptor
F03	Mm.28440	NM_008970	Pthlh	Parathyroid hormone-like peptide
F04	Mm.391013	NM_009820	Runx2	Runt related transcription factor 2
F05	Mm.281691	NM_013834	Sfrp1	Secreted frizzled-related protein 1
F06	Mm.42095	NM_016687	Sfrp4	Secreted frizzled-related protein 4
F07	Mm.1431	NM_011367	Shbg	Sex hormone binding globulin
F08	Mm.265602	NM_024449	Sost	Sclerostin
F09	Mm.291442	NM_009242	Sparc	Secreted acidic cysteine rich glycoprotein
F10	Mm.288474	NM_009263	Spp1	Secreted phosphoprotein 1
F11	Mm.277406	NM_009283	Stat1	Signal transducer and activator of transcription 1
F12	Mm.248380	NM_011577	Tgfb1	Transforming growth factor, beta 1
G01	Mm.206505	NM_011594	Timp2	Tissue inhibitor of metalloproteinase 2
G02	Mm.116683	NM_009397	Tnfaip3	Tumor necrosis factor, alpha-induced protein 3
G03	Mm.6251	NM_009399	Tnfrsf11a	Tumor necrosis factor receptor superfamily, member 11a
G04	Mm.15383	NM_008764	Tnfrsf11b	Tumor necrosis factor receptor superfamily, member 11b (osteoprotegerin)
G05	Mm.235328	NM_011610	Tnfrsf1b	Tumor necrosis factor receptor superfamily, member 1b
G06	Mm.249221	NM_011613	Tnfsf11	Tumor necrosis factor (ligand) superfamily, member 11
G07	Mm.173847	NM_011648	Tshr	Thyroid stimulating hormone receptor
G08	Mm.3280	NM_011658	Twist1	Twist homolog 1 (Drosophila)
G09	Mm.245084	NM_009504	Vdr	Vitamin D receptor
G10	Mm.282184	NM_009505	Vegfa	Vascular endothelial growth factor A
G11	Mm.4709	NM_011718	Wnt10b	Wingless related MMTV integration site 10b
G12	Mm.1367	NM_009522	Wnt3a	Wingless-related MMTV integration site 3A
H01	Mm.328431	NM_007393	Actb	Actin, beta
H02	Mm.163	NM_009735	B2m	Beta-2 microglobulin
H03	Mm.343110	NM_008084	Gapdh	Glyceraldehyde-3-phosphate dehydrogenase
H04	Mm.3317	NM_010368	Gusb	Glucuronidase, beta
H05	Mm.2180	NM_008302	Hsp90ab1	Heat shock protein 90 alpha (cytosolic), class B member 1
H06	N/A	SA_00106	MGDC	Mouse Genomic DNA Contamination
H07	N/A	SA_00104	RTC	Reverse Transcription Control
H08	N/A	SA_00104	RTC	Reverse Transcription Control
H09	N/A	SA_00104	RTC	Reverse Transcription Control
H10	N/A	SA_00103	PPC	Positive PCR Control
H11	N/A	SA_00103	PPC	Positive PCR Control
H12	N/A	SA_00103	PPC	Positive PCR Control

Related products

For optimal performance, RT² Profiler PCR Arrays should be used together with the RT² First Strand Kit for cDNA synthesis and RT² SYBR[®] Green qPCR Mastermixes for PCR.

Product	Contents	Cat. no.
RT ² First Strand Kit (12)	Enzymes and reagents for cDNA synthesis	330401
RT ² SYBR Green ROX [™] FAST Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the Rotor-Gene Q and other Rotor-Gene cyclers	330620

* Larger kit sizes available; please inquire.

RT² Profiler PCR Array products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, Rotor-Gene®, Rotor-Disc™ (QIAGEN Group); ROX™ (Applied Biosystems Corporation or its subsidiaries); SYBR® (Molecular Probes, Inc.).

1067688 03/2011 © 2011 QIAGEN, all rights reserved.

www.qiagen.com

Australia ■ 1-800-243-800

Austria ■ 0800/281010

Belgium ■ 0800-79612

Brazil ■ 0800-557779

Canada ■ 800-572-9613

China ■ 8621-3865-3865

Denmark ■ 80-885945

Finland ■ 0800-914416

France ■ 01-60-920-930

Germany ■ 02103-29-12000

Hong Kong ■ 800 933 965

Ireland ■ 1800 555 049

Italy ■ 800-787980

Japan ■ 03-6890-7300

Korea (South) ■ 080-000-7145

Luxembourg ■ 8002 2076

Mexico ■ 01-800-7742-436

The Netherlands ■ 0800 0229592

Norway ■ 800-18859

Singapore ■ 1800-742-4368

Spain ■ 91-630-7050

Sweden ■ 020-790282

Switzerland ■ 055-254-22-11

UK ■ 01293-422-911

USA ■ 800-426-8157



Sample & Assay Technologies