

RT² Profiler PCR Array (96-Well Format and 384-Well [4 x 96] Format)

Mouse miR-210 Targets

Cat. no. 330231 PAMM-6009ZA

For pathway expression analysis

Format	For use with the following real-time cyclers
RT ² Profiler PCR Array, Format A	Applied Biosystems® models 5700, 7000, 7300, 7500, 7700, 7900HT, ViiA™ 7 (96-well block); Bio-Rad® models iCycler®, iQ™ 5, MyiQ™, MyiQ2; Bio-Rad/MJ Research Chromo4™; Eppendorf® Mastercycler® ep realplex models 2, 2s, 4, 4s; Stratagene® models Mx3005P®, Mx3000P®; Takara TP-800
RT ² Profiler PCR Array, Format C	Applied Biosystems models 7500 (Fast block), 7900HT (Fast block), StepOnePlus™, ViiA 7 (Fast block)
RT ² Profiler PCR Array, Format D	Bio-Rad CFX96™; Bio-Rad/MJ Research models DNA Engine Opticon®, DNA Engine Opticon 2; Stratagene Mx4000®
RT ² Profiler PCR Array, Format E	Applied Biosystems models 7900HT (384-well block), ViiA 7 (384-well block); Bio-Rad CFX384™
RT ² Profiler PCR Array, Format F	Roche® LightCycler® 480 (96-well block)
RT ² Profiler PCR Array, Format G	Roche LightCycler 480 (384-well block)
RT ² Profiler PCR Array, Format H	Fluidigm® BioMark™



Sample & Assay Technologies

Description

The Mouse miR-210 Targets RT² Profiler PCR Array profiles the expression of 84 mmu-miR-210-3p target genes. This panel of 84 genes includes currently known experimentally verified plus bioinformatically predicted target genes regulated by mmu-miR-210-3p . miRNA target gene expression analysis provides further insight into the function of these specific miRNAs. Using real-time PCR, research studies can easily and reliably analyze the expression of a focused panel of genes likely to be regulated by miR-210 with this array.

For further details, consult the *RT² Profiler PCR Array Handbook*.

Shipping and storage

RT² Profiler PCR Arrays in formats A, C, D, E, F, and G are shipped at ambient temperature, on dry ice, or blue ice packs depending on destination and accompanying products. RT² Profiler PCR Arrays in format H are shipped on dry ice or blue ice packs.

For long term storage, keep plates at –20°C.

Note: Ensure that you have the correct RT² Profiler PCR Array format for your real-time cycler (see table above).

Note: Open the package and store the products appropriately immediately on receipt.

Array layout (96-well)

For 384-well 4 x 96 PCR arrays, genes are present in a staggered format. Refer to the *RT² Profiler PCR Array Handbook* for layout.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Acvr1b	Adams5	Adcy1	Adcy9	Aifm3	Aldh5a1	Ankrd34a	Armc1	Atg7	B4gal5	BC057022	Bdnf
B	Casp8ap2	Ccbp2	Ccdc97	Chn1	Chrn1	Coro2b	Cox10	Cpeb2	Ddah1	Dimt1	Dnajc15	Drd5
C	Dsc3	Dtx1	E2f3	Efn1	Efn3	Efn2	Enpp5	Fam116a	Fam73b	Foxa1	Gpd1l	Gsr
D	Hhip	Hif3a	Hoxa1	Hoxa9	Htr7	Iscu	Kcmf1	Kctd11	Klf12	Lair1	Lamp1	Mef2a
E	Mid1ip1	Mll2	Mnt	Ncam1	Ndufa4	Necab1	Neurod2	Nptx1	Nr1d2	Pptc7	Ptpn1	Rad52
F	Rcan2	Rgma	Rtn1	Scrt1	Sema3g	Sept8	Setd8	Slc25a26	Smg5	Syngap1	Tcf7l2	Thsd7a
G	Tmem151b	Tppp	Txndc9	Txnip	Xist	Xpnp3	Zadhf2	Zfand3	Zfp148	Zhye28	Zmat3	Zmiz1
H	Actb	B2m	Gapdh	Gusb	Hsp90ab1	MGDC	RTC	RTC	RTC	PPC	PPC	PPC

Gene table: RT² Profiler PCR Array

Position	UniGene	GenBank	Symbol	Description
A01	Mm.308467	NM_007395	Acvr1b	Activin A receptor, type 1B
A02	Mm.112933	NM_011782	Adams5	A disintegrin-like and metallopeptidase (reprolysin type) with thrombospondin type 1 motif, 5 (aggrecanase-2)
A03	Mm.259733	NM_009622	Adcy1	Adenylate cyclase 1
A04	Mm.439750	NM_009624	Adcy9	Adenylate cyclase 9
A05	Mm.40038	NM_175178	Aifm3	Apoptosis-inducing factor, mitochondrion-associated 3
A06	Mm.393311	NM_172532	Aldh5a1	Aldehyde dehydrogenase family 5, subfamily A1
A07	Mm.267664	NM_001024851	Ankrd34a	Ankyrin repeat domain 34A
A08	Mm.316007	NM_028840	Armc1	Armadillo repeat containing 1
A09	Mm.275332	NM_028835	Atg7	Autophagy-related 7 (yeast)
A10	Mm.200886	NM_019835	B4gal5	UDP-Gal:betaGlcNAc beta 1,4-galactosyltransferase, polypeptide 5
A11	Mm.194212	NM_001004180	BC057022	CDNA sequence BC057022
A12	Mm.1442	NM_007540	Bdnf	Brain derived neurotrophic factor
B01	Mm.22279	NM_011997	Casp8ap2	Caspase 8 associated protein 2
B02	Mm.258105	NM_021609	Ccbp2	Chemokine binding protein 2
B03	Mm.278798	NM_028771	Ccdc97	Coiled-coil domain containing 97
B04	Mm.475464	NM_175752	Chn1	Chimerin (chimaerin) 1
B05	Mm.86425	NM_009601	Chrn1	Cholinergic receptor, nicotinic, beta polypeptide 1 (muscle)
B06	Mm.335229	NM_175484	Coro2b	Coronin, actin binding protein, 2B
B07	Mm.340211	NM_178379	Cox10	COX10 homolog, cytochrome c oxidase assembly protein, heme A: farnesyltransferase (yeast)
B08	Mm.7233	NM_175937	Cpeb2	Cytoplasmic polyadenylation element binding protein 2
B09	Mm.234247	NM_026993	Ddah1	Dimethylarginine dimethylaminohydrolase 1
B10	Mm.9563	NM_025447	Dimt1	DIM1 dimethyladenosine transferase 1-like (S. cerevisiae)
B11	Mm.248046	NM_025384	Dnajc15	DnaJ (Hsp40) homolog, subfamily C, member 15
B12	Mm.167154	NM_013503	Drd5	Dopamine receptor D5
C01	Mm.89935	NM_007882	Dsc3	Desmocollin 3
C02	Mm.1645	NM_008052	Dtx1	Deltex 1 homolog (Drosophila)
C03	Mm.268356	NM_010093	E2f3	E2F transcription factor 3
C04	Mm.15675	NM_010107	Efn1	Ephrin A1
C05	Mm.331159	NM_010108	Efn3	Ephrin A3
C06	Mm.323188	NM_183141	Efn2	Leucine rich repeat and fibronectin type III, extracellular 2
C07	Mm.30145	NM_032003	Enpp5	Ectonucleotide pyrophosphatase/phosphodiesterase 5
C08	Mm.260512	NM_145969	Fam116a	Family with sequence similarity 116, member A
C09	Mm.254974	NM_175392	Fam73b	Family with sequence similarity 73, member B
C10	Mm.4578	NM_008259	Foxa1	Forkhead box A1
C11	Mm.38198	NM_175380	Gpd1l	Glycerol-3-phosphate dehydrogenase 1-like
C12	Mm.283573	NM_010344	Gsr	Glutathione reductase
D01	Mm.254493	NM_020259	Hhip	Hedgehog-interacting protein
D02	Mm.135110	NM_016868	Hif3a	Hypoxia inducible factor 3, alpha subunit
D03	Mm.197	NM_010449	Hoxa1	Homeobox A1
D04	Mm.4694	NM_010456	Hoxa9	Homeobox A9
D05	Mm.254266	NM_008315	Htr7	5-hydroxytryptamine (serotonin) receptor 7
D06	Mm.29497	NM_025526	Iscu	IscU iron-sulfur cluster scaffold homolog (E. coli)
D07	Mm.29194	NM_019715	Kcmf1	Potassium channel modulatory factor 1

Position	UniGene	GenBank	Symbol	Description
D08	Mm.239498	NM_153143	Kctd11	Potassium channel tetramerisation domain containing 11
D09	Mm.396499	NM_010636	Klf12	Kruppel-like factor 12
D10	Mm.290880	NM_178611	Lair1	Leukocyte-associated Ig-like receptor 1
D11	Mm.16716	NM_010684	Lamp1	Lysosomal-associated membrane protein 1
D12	Mm.132788	NM_001033713	Mef2a	Myocyte enhancer factor 2A
E01	Mm.29429	NM_026524	Mid1ip1	Mid1 interacting protein 1 (gastrulation specific G12-like (zebrafish))
E02	Mm.264889	NM_001033276	Mil2	Myeloid/lymphoid or mixed-lineage leukemia 2
E03	Mm.3759	NM_010813	Mnt	Max binding protein
E04	Mm.4974	NM_010875	Ncam1	Neural cell adhesion molecule 1
E05	Mm.415865	NM_010886	Ndufa4	NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 4
E06	Mm.389031	NM_178617	Necab1	N-terminal EF-hand calcium binding protein 1
E07	Mm.4814	NM_010895	Neurod2	Neurogenic differentiation 2
E08	Mm.5142	NM_008730	Nptx1	Neuronal pentraxin 1
E09	Mm.26587	NM_011584	Nr1d2	Nuclear receptor subfamily 1, group D, member 2
E10	Mm.485447	NM_177242	Ptpc7	PTC7 protein phosphatase homolog (S. cerevisiae)
E11	Mm.277916	NM_011201	Ptpn1	Protein tyrosine phosphatase, non-receptor type 1
E12	Mm.149	NM_011236	Rad52	RAD52 homolog (S. cerevisiae)
F01	Mm.251242	NM_030598	Rcan2	Regulator of calcineurin 2
F02	Mm.333943	NM_177740	Rgma	RGM domain family, member A
F03	Mm.221275	NM_153457	Rtn1	Reticulon 1
F04	Mm.218747	NM_130893	Scr1	Scratch homolog 1, zinc finger protein (Drosophila)
F05	Mm.27812	NM_001025379	Sema3g	Sema domain, immunoglobulin domain (Ig), short basic domain, secreted, (semaphorin) 3G
F06	Mm.274399	NM_033144	Sept8	Septin 8
F07	Mm.137966	NM_030241	Setd8	SET domain containing (lysine methyltransferase) 8
F08	Mm.280725	NM_026255	Slc25a26	Solute carrier family 25 (mitochondrial carrier, phosphate carrier), member 26
F09	Mm.439830	NM_178246	Smg5	Smg-5 homolog, nonsense mediated mRNA decay factor (C. elegans)
F10	Mm.328888	XM_915205	Syngap1	Synaptic Ras GTPase activating protein 1 homolog (rat)
F11	Mm.139815	NM_009333	Tcf7l2	Transcription factor 7-like 2, T-cell specific, HMG-box
F12	Mm.11989	NM_001164805	Thsd7a	Thrombospondin, type 1, domain containing 7A
G01	Mm.329064	NM_001013749	Tmem151b	Transmembrane protein 151B
G02	Mm.39752	NM_182839	Tppp	Tubulin polymerization promoting protein
G03	Mm.28438	NM_172054	Txndc9	Thioredoxin domain containing 9
G04	Mm.410189	NM_023719	Txnip	Thioredoxin interacting protein
G05	Mm.274770	NR_001570	Xist	Inactive X specific transcripts
G06	Mm.390833	NM_177310	Xpnpep3	X-prolyl aminopeptidase (aminopeptidase P) 3, putative
G07	Mm.147052	NM_146090	Zadh2	Zinc binding alcohol dehydrogenase, domain containing 2
G08	Mm.220972	NM_148926	Zfand3	Zinc finger, AN1-type domain 3
G09	Mm.392667	NM_011749	Zfp148	Zinc finger protein 148
G10	Mm.113388	NM_001015039	Zfyve28	Zinc finger, FYVE domain containing 28
G11	Mm.35705	NM_009517	Zmat3	Zinc finger matrin type 3
G12	Mm.227484	NM_183208	Zmiz1	Zinc finger, MIZ-type containing 1
H01	Mm.328431	NM_007393	Actb	Actin, beta
H02	Mm.163	NM_009735	B2m	Beta-2 microglobulin
H03	Mm.343110	NM_008084	Gapdh	Glyceraldehyde-3-phosphate dehydrogenase
H04	Mm.3317	NM_010368	Gusb	Glucuronidase, beta
H05	Mm.2180	NM_008302	Hsp90ab1	Heat shock protein 90 alpha (cytosolic), class B member 1
H06	N/A	SA_00106	MGDC	Mouse Genomic DNA Contamination
H07	N/A	SA_00104	RTC	Reverse Transcription Control
H08	N/A	SA_00104	RTC	Reverse Transcription Control
H09	N/A	SA_00104	RTC	Reverse Transcription Control
H10	N/A	SA_00103	PPC	Positive PCR Control
H11	N/A	SA_00103	PPC	Positive PCR Control
H12	N/A	SA_00103	PPC	Positive PCR Control

Related products

For optimal performance, RT² Profiler PCR Arrays should be used together with the RT² First Strand Kit for cDNA synthesis and RT² SYBR[®] Green qPCR Mastermixes for PCR.

Product	Contents	Cat. no.
RT ² First Strand Kit (12)	Enzymes and reagents for cDNA synthesis	330401
RT ² SYBR Green qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with real-time cyclers that do not require a reference dye, including: Bio-Rad models CFX96, CFX384, DNA Engine Opticon 2; Bio-Rad/MJ Research Chromo4; Roche LightCycler 480 (96-well and 384-well); all other cyclers	330500
RT ² SYBR Green ROX [™] qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the following real-time cyclers: Applied Biosystems models 5700, 7000, 7300, 7500 [Standard and FAST], 7700, 7900HT 96-well block [Standard and FAST] and 384-well block, StepOnePlus; Eppendorf Mastercycler ep realplex models 2, 2S, 4, 4S; Stratagene models Mx3000P, Mx3005P, Mx4000; Takara TP-800	330520
RT ² SYBR Green Fluor qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the following real-time cyclers: Bio-Rad models iCycler, iQ5, MyiQ, MyiQ2	330510

* Larger kit sizes available; please inquire.

RT² Profiler PCR Array products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN® (QIAGEN Group); Applied Biosystems®, ViiA™, StepOnePlus™, ROX™ (Applied Biosystems Corporation or its subsidiaries); Bio-Rad®, iCycler®, iQ™, MyiQ™, Chromo4™, CFX96™, DNA Engine Opticon®, CFX384™ (Bio-Rad Laboratories, Inc.); Stratagene®, Mx3005P®, Mx3000P®, Mx4000® (Stratagene); Eppendorf®, Mastercycler® (Eppendorf AG); Roche®, LightCycler® (Roche Group); Fluidigm® BioMark™ (Fluidigm Corporation); SYBR® (Molecular Probes, Inc.).

1066029 03/2011 © 2011 QIAGEN, all rights reserved.

www.qiagen.com

Canada ■ 800-572-9613

China ■ 8621-3865-3865

Denmark ■ 80-885945

Finland ■ 0800-914416

France ■ 01-60-920-930

Germany ■ 02103-29-12000

Hong Kong ■ 800 933 965

Ireland ■ 1800 555 049

Italy ■ 800-787980

Japan ■ 03-6890-7300

Korea (South) ■ 080-000-7145

Luxembourg ■ 8002 2076

Mexico ■ 01-800-7742-436

The Netherlands ■ 0800 0229592

Norway ■ 800-18859

Singapore ■ 1800-742-4368

Spain ■ 91-630-7050

Sweden ■ 020-790282

Switzerland ■ 055-254-22-11

UK ■ 01293-422-911

USA ■ 800-426-8157

Australia ■ 1-800-243-800

Austria ■ 0800/281010

Belgium ■ 0800-79612

Brazil ■ 0800-557779



Sample & Assay Technologies