

# RT<sup>2</sup> Profiler PCR Array (Rotor-Gene® Format)

## Rat Housekeeping Genes

Cat. no. 330231 PARN-000ZR

For pathway expression analysis

Format	For use with the following real-time cyclers
RT <sup>2</sup> Profiler PCR Array, Format R	Rotor-Gene Q, other Rotor-Gene cyclers

### Description

The Rat Housekeeping Genes RT<sup>2</sup> Profiler PCR Array profiles the expression of eight (8) identical sets of twelve (12) commonly used housekeeping genes. This array can be easily used to identify genes with a constant level of expression among your different experimental conditions for use in normalizing your relative gene expression profiling experiment. Housekeeping genes encode proteins that are usually essential for the maintenance of cellular function. Their expression often, but not always, remains constant under most experimental conditions. Some of these genes may be regulated under certain circumstances and may vary with cell type. Using real-time PCR, you can easily and reliably evaluate the expression of these housekeeping genes among your different experiment conditions with this array.

For further details, consult the *RT<sup>2</sup> Profiler PCR Array Handbook*.

### Shipping and storage

RT<sup>2</sup> Profiler PCR Arrays in the Rotor-Gene format are shipped at ambient temperature, on dry ice, or blue ice packs depending on destination and accompanying products.

For long term storage, keep plates at –20°C.

**Note:** Ensure that you have the correct RT<sup>2</sup> Profiler PCR Array format for your real-time cycler (see table above).

**Note:** Open the package and store the products appropriately immediately on receipt.



## Array layout

The 96 real-time assays in the Rotor-Gene format are located in wells 1–96 of the Rotor-Disc™ (plate A1–A12=Rotor-Disc 1–12, plate B1–B12=Rotor-Disc 13–24, etc.). To maintain data analysis compatibility, wells 97–100 do not contain real-time assays but will contain master mix to account for weight balance.

## Gene table: RT<sup>2</sup> Profiler PCR Array

Position	UniGene	GenBank	Symbol	Description
A01	Rn.94978	NM_031144	Actb	Actin, beta
A02	Rn.1868	NM_012512	B2m	Beta-2 microglobulin
A03	Rn.47	NM_012583	Hprt1	Hypoxanthine phosphoribosyltransferase 1
A04	Rn.98667	NM_001004082	Hsp90ab1	Heat shock protein 90 alpha (cytosolic), class B member 1
A05	Rn.107896	NM_017025	Ldha	Lactate dehydrogenase A
A06	Rn.8381	NM_001012356	Nono	Non-POU domain containing, octamer-binding
A07	Rn.108127	NM_053291	Pgk1	Phosphoglycerate kinase 1
A08	Rn.129133	XM_345576	Ppih	Peptidylprolyl isomerase H (cyclophilin H)
A09	Rn.973	NM_001007604	Rplp1	Ribosomal protein, large, P1
A10	Rn.101725	NM_130428	Sdha	Succinate dehydrogenase complex, subunit A, flavoprotein (Fp)
A11	Rn.22712	NM_001004198	Tbp	TATA box binding protein
A12	Rn.98672	NM_022712	Tfrc	Transferrin receptor
B01	Rn.94978	NM_031144	Actb	Actin, beta
B02	Rn.1868	NM_012512	B2m	Beta-2 microglobulin
B03	Rn.47	NM_012583	Hprt1	Hypoxanthine phosphoribosyltransferase 1
B04	Rn.98667	NM_001004082	Hsp90ab1	Heat shock protein 90 alpha (cytosolic), class B member 1
B05	Rn.107896	NM_017025	Ldha	Lactate dehydrogenase A
B06	Rn.8381	NM_001012356	Nono	Non-POU domain containing, octamer-binding
B07	Rn.108127	NM_053291	Pgk1	Phosphoglycerate kinase 1
B08	Rn.129133	XM_345576	Ppih	Peptidylprolyl isomerase H (cyclophilin H)
B09	Rn.973	NM_001007604	Rplp1	Ribosomal protein, large, P1
B10	Rn.101725	NM_130428	Sdha	Succinate dehydrogenase complex, subunit A, flavoprotein (Fp)
B11	Rn.22712	NM_001004198	Tbp	TATA box binding protein
B12	Rn.98672	NM_022712	Tfrc	Transferrin receptor
C01	Rn.94978	NM_031144	Actb	Actin, beta
C02	Rn.1868	NM_012512	B2m	Beta-2 microglobulin
C03	Rn.47	NM_012583	Hprt1	Hypoxanthine phosphoribosyltransferase 1
C04	Rn.98667	NM_001004082	Hsp90ab1	Heat shock protein 90 alpha (cytosolic), class B member 1
C05	Rn.107896	NM_017025	Ldha	Lactate dehydrogenase A
C06	Rn.8381	NM_001012356	Nono	Non-POU domain containing, octamer-binding
C07	Rn.108127	NM_053291	Pgk1	Phosphoglycerate kinase 1
C08	Rn.129133	XM_345576	Ppih	Peptidylprolyl isomerase H (cyclophilin H)
C09	Rn.973	NM_001007604	Rplp1	Ribosomal protein, large, P1
C10	Rn.101725	NM_130428	Sdha	Succinate dehydrogenase complex, subunit A, flavoprotein (Fp)
C11	Rn.22712	NM_001004198	Tbp	TATA box binding protein
C12	Rn.98672	NM_022712	Tfrc	Transferrin receptor
D01	Rn.94978	NM_031144	Actb	Actin, beta
D02	Rn.1868	NM_012512	B2m	Beta-2 microglobulin
D03	Rn.47	NM_012583	Hprt1	Hypoxanthine phosphoribosyltransferase 1
D04	Rn.98667	NM_001004082	Hsp90ab1	Heat shock protein 90 alpha (cytosolic), class B member 1
D05	Rn.107896	NM_017025	Ldha	Lactate dehydrogenase A
D06	Rn.8381	NM_001012356	Nono	Non-POU domain containing, octamer-binding
D07	Rn.108127	NM_053291	Pgk1	Phosphoglycerate kinase 1
D08	Rn.129133	XM_345576	Ppih	Peptidylprolyl isomerase H (cyclophilin H)
D09	Rn.973	NM_001007604	Rplp1	Ribosomal protein, large, P1
D10	Rn.101725	NM_130428	Sdha	Succinate dehydrogenase complex, subunit A, flavoprotein (Fp)
D11	Rn.22712	NM_001004198	Tbp	TATA box binding protein
D12	Rn.98672	NM_022712	Tfrc	Transferrin receptor
E01	Rn.94978	NM_031144	Actb	Actin, beta
E02	Rn.1868	NM_012512	B2m	Beta-2 microglobulin
E03	Rn.47	NM_012583	Hprt1	Hypoxanthine phosphoribosyltransferase 1
E04	Rn.98667	NM_001004082	Hsp90ab1	Heat shock protein 90 alpha (cytosolic), class B member 1
E05	Rn.107896	NM_017025	Ldha	Lactate dehydrogenase A
E06	Rn.8381	NM_001012356	Nono	Non-POU domain containing, octamer-binding
E07	Rn.108127	NM_053291	Pgk1	Phosphoglycerate kinase 1
E08	Rn.129133	XM_345576	Ppih	Peptidylprolyl isomerase H (cyclophilin H)
E09	Rn.973	NM_001007604	Rplp1	Ribosomal protein, large, P1

Position	UniGene	GenBank	Symbol	Description
E10	Rn.101725	NM_130428	Sdha	Succinate dehydrogenase complex, subunit A, flavoprotein (Fp)
E11	Rn.22712	NM_001004198	Tbp	TATA box binding protein
E12	Rn.98672	NM_022712	Tfrc	Transferrin receptor
F01	Rn.94978	NM_031144	Actb	Actin, beta
F02	Rn.1868	NM_012512	B2m	Beta-2 microglobulin
F03	Rn.47	NM_012583	Hprt1	Hypoxanthine phosphoribosyltransferase 1
F04	Rn.98667	NM_001004082	Hsp90ab1	Heat shock protein 90 alpha (cytosolic), class B member 1
F05	Rn.107896	NM_017025	Ldha	Lactate dehydrogenase A
F06	Rn.8381	NM_001012356	Nono	Non-POU domain containing, octamer-binding
F07	Rn.108127	NM_053291	Pgk1	Phosphoglycerate kinase 1
F08	Rn.129133	XM_345576	Ppih	Peptidylprolyl isomerase H (cyclophilin H)
F09	Rn.973	NM_001007604	Rplp1	Ribosomal protein, large, P1
F10	Rn.101725	NM_130428	Sdha	Succinate dehydrogenase complex, subunit A, flavoprotein (Fp)
F11	Rn.22712	NM_001004198	Tbp	TATA box binding protein
F12	Rn.98672	NM_022712	Tfrc	Transferrin receptor
G01	Rn.94978	NM_031144	Actb	Actin, beta
G02	Rn.1868	NM_012512	B2m	Beta-2 microglobulin
G03	Rn.47	NM_012583	Hprt1	Hypoxanthine phosphoribosyltransferase 1
G04	Rn.98667	NM_001004082	Hsp90ab1	Heat shock protein 90 alpha (cytosolic), class B member 1
G05	Rn.107896	NM_017025	Ldha	Lactate dehydrogenase A
G06	Rn.8381	NM_001012356	Nono	Non-POU domain containing, octamer-binding
G07	Rn.108127	NM_053291	Pgk1	Phosphoglycerate kinase 1
G08	Rn.129133	XM_345576	Ppih	Peptidylprolyl isomerase H (cyclophilin H)
G09	Rn.973	NM_001007604	Rplp1	Ribosomal protein, large, P1
G10	Rn.101725	NM_130428	Sdha	Succinate dehydrogenase complex, subunit A, flavoprotein (Fp)
G11	Rn.22712	NM_001004198	Tbp	TATA box binding protein
G12	Rn.98672	NM_022712	Tfrc	Transferrin receptor
H01	Rn.94978	NM_031144	Actb	Actin, beta
H02	Rn.1868	NM_012512	B2m	Beta-2 microglobulin
H03	Rn.47	NM_012583	Hprt1	Hypoxanthine phosphoribosyltransferase 1
H04	Rn.98667	NM_001004082	Hsp90ab1	Heat shock protein 90 alpha (cytosolic), class B member 1
H05	Rn.107896	NM_017025	Ldha	Lactate dehydrogenase A
H06	Rn.8381	NM_001012356	Nono	Non-POU domain containing, octamer-binding
H07	Rn.108127	NM_053291	Pgk1	Phosphoglycerate kinase 1
H08	Rn.129133	XM_345576	Ppih	Peptidylprolyl isomerase H (cyclophilin H)
H09	Rn.973	NM_001007604	Rplp1	Ribosomal protein, large, P1
H10	Rn.101725	NM_130428	Sdha	Succinate dehydrogenase complex, subunit A, flavoprotein (Fp)
H11	Rn.22712	NM_001004198	Tbp	TATA box binding protein
H12	Rn.98672	NM_022712	Tfrc	Transferrin receptor

## Related products

For optimal performance, RT<sup>2</sup> Profiler PCR Arrays should be used together with the RT<sup>2</sup> First Strand Kit for cDNA synthesis and RT<sup>2</sup> SYBR<sup>®</sup> Green qPCR Mastermixes for PCR.

Product	Contents	Cat. no.
RT <sup>2</sup> First Strand Kit (12)	Enzymes and reagents for cDNA synthesis	330401
RT <sup>2</sup> SYBR Green ROX <sup>™</sup> FAST Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the Rotor-Gene Q and other Rotor-Gene cyclers	330620

\* Larger kit sizes available; please inquire.

RT<sup>2</sup> Profiler PCR Array products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at [www.qiagen.com](http://www.qiagen.com) or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN<sup>®</sup>, Rotor-Gene<sup>®</sup>, Rotor-Disc<sup>™</sup> (QIAGEN Group); ROX<sup>™</sup> (Applied Biosystems or its subsidiaries); SYBR<sup>®</sup> (Molecular Probes, Inc.).

1067688 03/2011 © 2011 QIAGEN, all rights reserved.

[www.qiagen.com](http://www.qiagen.com)

Australia ■ 1-800-243-800

Austria ■ 0800/281010

Belgium ■ 0800-79612

Brazil ■ 0800-557779

Canada ■ 800-572-9613

China ■ 8621-3865-3865

Denmark ■ 80-885945

Finland ■ 0800-914416

France ■ 01-60-920-930

Germany ■ 02103-29-12000

Hong Kong ■ 800 933 965

Ireland ■ 1800 555 049

Italy ■ 800-787980

Japan ■ 03-6890-7300

Korea (South) ■ 080-000-7145

Luxembourg ■ 8002 2076

Mexico ■ 01-800-7742-436

The Netherlands ■ 0800 0229592

Norway ■ 800-18859

Singapore ■ 1800-742-4368

Spain ■ 91-630-7050

Sweden ■ 020-790282

Switzerland ■ 055-254-22-11

UK ■ 01293-422-911

USA ■ 800-426-8157



Sample & Assay Technologies