

# RT<sup>2</sup> Profiler PCR Array (96-Well Format and 384-Well [4 x 96] Format)

## Rat Inflammatory Cytokines & Receptors

Cat. no. 330231 PARN-011ZA

For pathway expression analysis

| Format                                       | For use with the following real-time cyclers                                                                                                                                                                                                                                          |
|----------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| RT <sup>2</sup> Profiler PCR Array, Format A | Applied Biosystems® models 5700, 7000, 7300, 7500, 7700, 7900HT, ViiA™ 7 (96-well block); Bio-Rad® models iCycler®, iQ™ 5, MyiQ™, MyiQ2; Bio-Rad/MJ Research Chromo4™; Eppendorf® Mastercycler® ep realplex models 2, 2s, 4, 4s; Stratagene® models Mx3005P®, Mx3000P®; Takara TP-800 |
| RT <sup>2</sup> Profiler PCR Array, Format C | Applied Biosystems models 7500 (Fast block), 7900HT (Fast block), StepOnePlus™, ViiA 7 (Fast block)                                                                                                                                                                                   |
| RT <sup>2</sup> Profiler PCR Array, Format D | Bio-Rad CFX96™; Bio-Rad/MJ Research models DNA Engine Opticon®, DNA Engine Opticon 2; Stratagene Mx4000®                                                                                                                                                                              |
| RT <sup>2</sup> Profiler PCR Array, Format E | Applied Biosystems models 7900HT (384-well block), ViiA 7 (384-well block); Bio-Rad CFX384™                                                                                                                                                                                           |
| RT <sup>2</sup> Profiler PCR Array, Format F | Roche® LightCycler® 480 (96-well block)                                                                                                                                                                                                                                               |
| RT <sup>2</sup> Profiler PCR Array, Format G | Roche LightCycler 480 (384-well block)                                                                                                                                                                                                                                                |
| RT <sup>2</sup> Profiler PCR Array, Format H | Fluidigm® BioMark™                                                                                                                                                                                                                                                                    |



Sample & Assay Technologies

## Description

The Rat Inflammatory Cytokines & Receptors RT<sup>2</sup> Profiler PCR Array profiles the expression of 84 key genes mediating the inflammatory response. Acute inflammation occurs in response to cell damage due to infection or injury. During this process, cellular and plasma derived factors encourage extravasation, the recruitment of circulating immune cells into the affected tissue. These immune cells in turn up-regulate the expression of inflammatory cytokines that recruit additional immune cells to mount an immune response to any invading organisms and to further promote and eventually resolve the inflammatory response. Chronic inflammation, or expression of these cytokines and receptors at low levels over long periods of time, promotes various pathological conditions including allergies and asthma, cardiovascular system disorders (atherosclerosis), central nervous system disorders (Alzheimer's disease), fibrosis, and rheumatoid arthritis. This array contains genes involved in mediating immune cascade reactions during inflammation. The chemokines, cytokines, and interleukins involved in the inflammatory response are represented as well as their receptors. Profiling the expression of these inflammatory cytokine and receptor genes helps determine the state of and the mechanisms behind inflammation in your experimental model system. Using real-time PCR, you can easily and reliably analyze expression of a focused panel of genes related to inflammation with this array.

For further details, consult the *RT<sup>2</sup> Profiler PCR Array Handbook*.

## Shipping and storage

RT<sup>2</sup> Profiler PCR Arrays in formats A, C, D, E, F, and G are shipped at ambient temperature, on dry ice, or blue ice packs depending on destination and accompanying products. RT<sup>2</sup> Profiler PCR Arrays in format H are shipped on dry ice or blue ice packs.

For long term storage, keep plates at –20°C.

**Note:** Ensure that you have the correct RT<sup>2</sup> Profiler PCR Array format for your real-time cycler (see table above).

**Note:** Open the package and store the products appropriately immediately on receipt.

---

## Array layout (96-well)

For 384-well 4 x 96 PCR arrays, genes are present in a staggered format. Refer to the *RT<sup>2</sup> Profiler PCR Array Handbook* for layout.

|   | 1      | 2                        | 3     | 4     | 5         | 6        | 7        | 8        | 9         | 10       | 11      | 12    |
|---|--------|--------------------------|-------|-------|-----------|----------|----------|----------|-----------|----------|---------|-------|
| A | Aimp1  | Bmp2                     | Ccl11 | Ccl12 | Ccl17     | Ccl19    | Ccl2     | Ccl20    | Ccl22     | Ccl24    | Ccl3    | Ccl4  |
| B | Ccl5   | Ccl6                     | Ccl7  | Ccl9  | Ccr1      | Ccr10    | Ccr2     | Ccr3     | Ccr4      | Ccr5     | Ccr6    | Ccr8  |
| C | Cd40lg | Csf1                     | Csf2  | Csf3  | Cx3cl1    | Cx3cr1   | Cxcl1    | Cxcl10   | Cxcl11    | Cxcl12   | Cxcl2   | Cxcl5 |
| D | Cxcl9  | Cxcr2                    | Cxcr3 | Cxcr5 | Faslg     | Ifng     | Il10ra   | Il11     | Il13      | Il15     | Il16    | Il17a |
| E | Il17b  | Il17f                    | Il1a  | Il1b  | Il1r1     | Il1rn    | Il21     | Il27     | Il2rb     | Il2rg    | Il3     | Il33  |
| F | Il4    | Il5                      | Il5ra | Il6r  | Il6st     | Il7      | Il8ra    | Lta      | Ltb       | Mif      | Nampt   | Osm   |
| G | Pf4    | RGD1561905<br>_predicted | Spp1  | Tnf   | Tnfrsf11b | Tnfrsf10 | Tnfrsf11 | Tnfrsf13 | Tnfrsf13b | Tnfrsf14 | Tnfrsf4 | Vegfa |
| H | Actb   | B2m                      | Hprt1 | LdhA  | Rplp1     | RGDC     | RTC      | RTC      | RTC       | PPC      | PPC     | PPC   |

## Gene table: RT<sup>2</sup> Profiler PCR Array

| Position | UniGene   | GenBank      | Symbol | Description                                                                    |
|----------|-----------|--------------|--------|--------------------------------------------------------------------------------|
| A01      | Rn.3192   | NM_053757    | Aimp1  | Aminoacyl tRNA synthetase complex-interacting multifunctional protein 1        |
| A02      | Rn.90931  | NM_017178    | Bmp2   | Bone morphogenetic protein 2                                                   |
| A03      | Rn.10632  | NM_019205    | Ccl11  | Chemokine (C-C motif) ligand 11                                                |
| A04      | Rn.137780 | NM_001105822 | Ccl12  | Chemokine (C-C motif) ligand 12                                                |
| A05      | Rn.46679  | NM_057151    | Ccl17  | Chemokine (C-C motif) ligand 17                                                |
| A06      | Rn.12445  | NM_001108661 | Ccl19  | Chemokine (C-C motif) ligand 19                                                |
| A07      | Rn.4772   | NM_031530    | Ccl2   | Chemokine (C-C motif) ligand 2                                                 |
| A08      | Rn.10722  | NM_019233    | Ccl20  | Chemokine (C-C motif) ligand 20                                                |
| A09      | Rn.48727  | NM_057203    | Ccl22  | Chemokine (C-C motif) ligand 22                                                |
| A10      | Rn.138822 | NM_001013045 | Ccl24  | Chemokine (C-C motif) ligand 24                                                |
| A11      | Rn.10139  | NM_013025    | Ccl3   | Chemokine (C-C motif) ligand 3                                                 |
| A12      | Rn.37880  | NM_053858    | Ccl4   | Chemokine (C-C motif) ligand 4                                                 |
| B01      | Rn.8019   | NM_031116    | Ccl5   | Chemokine (C-C motif) ligand 5                                                 |
| B02      | Rn.7857   | NM_001004202 | Ccl6   | Chemokine (C-C motif) ligand 6                                                 |
| B03      | Rn.26815  | NM_001007612 | Ccl7   | Chemokine (C-C motif) ligand 7                                                 |
| B04      | Rn.145427 | NM_001012357 | Ccl9   | Chemokine (C-C motif) ligand 9                                                 |
| B05      | Rn.34673  | NM_020542    | Ccr1   | Chemokine (C-C motif) receptor 1                                               |
| B06      | Rn.49607  | NM_001108836 | Ccr10  | Chemokine (C-C motif) receptor 10                                              |
| B07      | Rn.211983 | NM_021866    | Ccr2   | Chemokine (C-C motif) receptor 2                                               |
| B08      | Rn.9890   | NM_053958    | Ccr3   | Chemokine (C-C motif) receptor 3                                               |
| B09      | Rn.81076  | NM_133532    | Ccr4   | Chemokine (C-C motif) receptor 4                                               |
| B10      | Rn.10736  | NM_053960    | Ccr5   | Chemokine (C-C motif) receptor 5                                               |
| B11      | Rn.161767 | NM_001013145 | Ccr6   | Chemokine (C-C motif) receptor 6                                               |
| B12      | N/A       | XM_236704    | Ccr8   | Chemokine (C-C motif) receptor 8                                               |
| C01      | Rn.44218  | NM_053353    | Cd40lg | CD40 ligand                                                                    |
| C02      | Rn.83632  | NM_023981    | Csf1   | Colony stimulating factor 1 (macrophage)                                       |
| C03      | Rn.44285  | XM_340799    | Csf2   | Colony stimulating factor 2 (granulocyte-macrophage)                           |
| C04      | Rn.53973  | NM_017104    | Csf3   | Colony stimulating factor 3 (granulocyte)                                      |
| C05      | Rn.107266 | NM_134455    | Cx3cl1 | Chemokine (C-X3-C motif) ligand 1                                              |
| C06      | Rn.10482  | NM_133534    | Cx3cr1 | Chemokine (C-X3-C motif) receptor 1                                            |
| C07      | Rn.10907  | NM_030845    | Cxcl1  | Chemokine (C-X-C motif) ligand 1 (melanoma growth stimulating activity, alpha) |
| C08      | Rn.10584  | NM_139089    | Cxcl10 | Chemokine (C-X-C motif) ligand 10                                              |
| C09      | Rn.13664  | NM_182952    | Cxcl11 | Chemokine (C-X-C motif) ligand 11                                              |
| C10      | Rn.54439  | NM_022177    | Cxcl12 | Chemokine (C-X-C motif) ligand 12 (stromal cell-derived factor 1)              |
| C11      | Rn.10230  | NM_053647    | Cxcl2  | Chemokine (C-X-C motif) ligand 2                                               |
| C12      | Rn.44449  | NM_022214    | Cxcl5  | Chemokine (C-X-C motif) ligand 5                                               |
| D01      | Rn.7391   | NM_145672    | Cxcl9  | Chemokine (C-X-C motif) ligand 9                                               |
| D02      | Rn.90347  | NM_017183    | Cxcr2  | Chemokine (C-X-C motif) receptor 2                                             |
| D03      | Rn.24787  | NM_053415    | Cxcr3  | Chemokine (C-X-C motif) receptor 3                                             |
| D04      | Rn.127121 | NM_053303    | Cxcr5  | Chemokine (C-X-C motif) receptor 5                                             |
| D05      | Rn.9725   | NM_012908    | Faslg  | Fas ligand (TNF superfamily, member 6)                                         |
| D06      | Rn.10795  | NM_138880    | Ifng   | Interferon gamma                                                               |
| D07      | Rn.64632  | NM_057193    | Il10ra | Interleukin 10 receptor, alpha                                                 |

| Position | UniGene   | GenBank      | Symbol                   | Description                                            |
|----------|-----------|--------------|--------------------------|--------------------------------------------------------|
| D08      | Rn.198483 | NM_133519    | Il11                     | Interleukin 11                                         |
| D09      | Rn.9921   | NM_053828    | Il13                     | Interleukin 13                                         |
| D10      | Rn.2490   | NM_013129    | Il15                     | Interleukin 15                                         |
| D11      | Rn.104665 | NM_001105749 | Il16                     | Interleukin 16                                         |
| D12      | Rn.218513 | NM_001106897 | Il17a                    | Interleukin 17A                                        |
| E01      | Rn.50003  | NM_053789    | Il17b                    | Interleukin 17B                                        |
| E02      | Rn.34022  | NM_001015011 | Il17f                    | Interleukin 17F                                        |
| E03      | Rn.12300  | NM_017019    | Il1a                     | Interleukin 1 alpha                                    |
| E04      | Rn.9869   | NM_031512    | Il1b                     | Interleukin 1 beta                                     |
| E05      | Rn.9758   | NM_013123    | Il1r1                    | Interleukin 1 receptor, type I                         |
| E06      | Rn.162640 | NM_022194    | Il1rn                    | Interleukin 1 receptor antagonist                      |
| E07      | Rn.214949 | NM_001108943 | Il21                     | Interleukin 21                                         |
| E08      | N/A       | XM_344962    | Il27                     | Interleukin 27                                         |
| E09      | Rn.5832   | NM_013195    | Il2rb                    | Interleukin 2 receptor, beta                           |
| E10      | Rn.14508  | NM_080889    | Il2rg                    | Interleukin 2 receptor, gamma                          |
| E11      | Rn.10652  | NM_031513    | Il3                      | Interleukin 3                                          |
| E12      | Rn.106849 | NM_001014166 | Il33                     | Interleukin 33                                         |
| F01      | Rn.108255 | NM_201270    | Il4                      | Interleukin 4                                          |
| F02      | Rn.44227  | NM_021834    | Il5                      | Interleukin 5                                          |
| F03      | Rn.67068  | NM_053645    | Il5ra                    | Interleukin 5 receptor, alpha                          |
| F04      | Rn.1716   | NM_017020    | Il6r                     | Interleukin 6 receptor                                 |
| F05      | Rn.12138  | NM_001008725 | Il6st                    | Interleukin 6 signal transducer                        |
| F06      | Rn.10793  | NM_013110    | Il7                      | Interleukin 7                                          |
| F07      | Rn.138115 | NM_019310    | Il8ra                    | Interleukin 8 receptor, alpha                          |
| F08      | Rn.160577 | NM_080769    | Lta                      | Lymphotoxin alpha (TNF superfamily, member 1)          |
| F09      | Rn.203016 | NM_212507    | Ltb                      | Lymphotoxin beta (TNF superfamily, member 3)           |
| F10      | Rn.2661   | NM_031051    | Mif                      | Macrophage migration inhibitory factor                 |
| F11      | Rn.203508 | NM_177928    | Nampt                    | Nicotinamide phosphoribosyltransferase                 |
| F12      | Rn.127158 | NM_001006961 | Osm                      | Oncostatin M                                           |
| G01      | Rn.6800   | NM_001007729 | Pf4                      | Platelet factor 4                                      |
| G02      | Rn.145200 | XM_345342    | RGD156190<br>5_predicted | Complement component 5                                 |
| G03      | Rn.8871   | NM_012881    | Spp1                     | Secreted phosphoprotein 1                              |
| G04      | Rn.2275   | NM_012675    | Tnf                      | Tumor necrosis factor (TNF superfamily, member 2)      |
| G05      | Rn.202973 | NM_012870    | Tnfrsf11b                | Tumor necrosis factor receptor superfamily, member 11b |
| G06      | Rn.83627  | NM_145681    | Tnfsf10                  | Tumor necrosis factor (ligand) superfamily, member 10  |
| G07      | Rn.64517  | NM_057149    | Tnfsf11                  | Tumor necrosis factor (ligand) superfamily, member 11  |
| G08      | Rn.19955  | NM_001009623 | Tnfsf13                  | Tumor necrosis factor (ligand) superfamily, member 13  |
| G09      | Rn.219918 | NM_001109112 | Tnfsf13b                 | Tumor necrosis factor (ligand) superfamily, member 13b |
| G10      | N/A       | XM_236794    | Tnfsf14                  | Tumor necrosis factor (ligand) superfamily, member 14  |
| G11      | Rn.30043  | NM_053552    | Tnfsf4                   | Tumor necrosis factor (ligand) superfamily, member 4   |
| G12      | Rn.1923   | NM_031836    | Vegfa                    | Vascular endothelial growth factor A                   |
| H01      | Rn.94978  | NM_031144    | Actb                     | Actin, beta                                            |
| H02      | Rn.1868   | NM_012512    | B2m                      | Beta-2 microglobulin                                   |
| H03      | Rn.47     | NM_012583    | Hprt1                    | Hypoxanthine phosphoribosyltransferase 1               |
| H04      | Rn.107896 | NM_017025    | Ldha                     | Lactate dehydrogenase A                                |
| H05      | Rn.973    | NM_001007604 | Rplp1                    | Ribosomal protein, large, P1                           |
| H06      | N/A       | U26919       | RGDC                     | Rat Genomic DNA Contamination                          |
| H07      | N/A       | SA_00104     | RTC                      | Reverse Transcription Control                          |
| H08      | N/A       | SA_00104     | RTC                      | Reverse Transcription Control                          |
| H09      | N/A       | SA_00104     | RTC                      | Reverse Transcription Control                          |
| H10      | N/A       | SA_00103     | PPC                      | Positive PCR Control                                   |
| H11      | N/A       | SA_00103     | PPC                      | Positive PCR Control                                   |
| H12      | N/A       | SA_00103     | PPC                      | Positive PCR Control                                   |

## Related products

For optimal performance, RT<sup>2</sup> Profiler PCR Arrays should be used together with the RT<sup>2</sup> First Strand Kit for cDNA synthesis and RT<sup>2</sup> SYBR<sup>®</sup> Green qPCR Mastermixes for PCR.

| Product                                                         | Contents                                                                                                                                                                                                                                                                                                                                                            | Cat. no. |
|-----------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|
| RT <sup>2</sup> First Strand Kit (12)                           | Enzymes and reagents for cDNA synthesis                                                                                                                                                                                                                                                                                                                             | 330401   |
| RT <sup>2</sup> SYBR Green qPCR Mastermix (2)*                  | For 2 x 96 assays in 96-well plates; suitable for use with real-time cyclers that do not require a reference dye, including: Bio-Rad models CFX96, CFX384, DNA Engine Opticon 2; Bio-Rad/MJ Research Chromo4; Roche LightCycler 480 (96-well and 384-well); all other cyclers                                                                                       | 330500   |
| RT <sup>2</sup> SYBR Green ROX <sup>™</sup> qPCR Mastermix (2)* | For 2 x 96 assays in 96-well plates; suitable for use with the following real-time cyclers: Applied Biosystems models 5700, 7000, 7300, 7500 [Standard and FAST], 7700, 7900HT 96-well block [Standard and FAST] and 384-well block, StepOnePlus; Eppendorf Mastercycler ep realplex models 2, 2S, 4, 4S; Stratagene models Mx3000P, Mx3005P, Mx4000; Takara TP-800 | 330520   |
| RT <sup>2</sup> SYBR Green Fluor qPCR Mastermix (2)*            | For 2 x 96 assays in 96-well plates; suitable for use with the following real-time cyclers: Bio-Rad models iCycler, iQ5, MyiQ, MyiQ2                                                                                                                                                                                                                                | 330510   |

\* Larger kit sizes available; please inquire.

RT<sup>2</sup> Profiler PCR Array products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at [www.qiagen.com](http://www.qiagen.com) or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN® (QIAGEN Group); Applied Biosystems®, ViiA™, StepOnePlus™, ROX™ (Applied Biosystems Corporation or its subsidiaries); Bio-Rad®, iCycler®, iQ™, MyiQ™, Chromo4™, CFX96™, DNA Engine Opticon®, CFX384™ (Bio-Rad Laboratories, Inc.); Stratagene®, Mx3005P®, Mx3000P®, Mx4000® (Stratagene); Eppendorf®, Mastercycler® (Eppendorf AG); Roche®, LightCycler® (Roche Group); Fluidigm® BioMark™ (Fluidigm Corporation); SYBR® (Molecular Probes, Inc.).

1066029 03/2011 © 2011 QIAGEN, all rights reserved.

[www.qiagen.com](http://www.qiagen.com)

Canada ■ 800-572-9613

China ■ 8621-3865-3865

Denmark ■ 80-885945

Finland ■ 0800-914416

France ■ 01-60-920-930

Germany ■ 02103-29-12000

Hong Kong ■ 800 933 965

Ireland ■ 1800 555 049

Italy ■ 800-787980

Japan ■ 03-6890-7300

Korea (South) ■ 080-000-7145

Luxembourg ■ 8002 2076

Mexico ■ 01-800-7742-436

The Netherlands ■ 0800 0229592

Norway ■ 800-18859

Singapore ■ 1800-742-4368

Spain ■ 91-630-7050

Sweden ■ 020-790282

Switzerland ■ 055-254-22-11

UK ■ 01293-422-911

USA ■ 800-426-8157

Australia ■ 1-800-243-800

Austria ■ 0800/281010

Belgium ■ 0800-79612

Brazil ■ 0800-557779



Sample & Assay Technologies