

RT² Profiler PCR Array (Rotor-Gene® Format)

Rat Signal Transduction PathwayFinder

Cat. no. 330231 PARN-014ZR

For pathway expression analysis

Format	For use with the following real-time cyclers
RT ² Profiler PCR Array, Format R	Rotor-Gene Q, other Rotor-Gene cyclers

Description

The Rat Signal Transduction PathwayFinder RT² Profiler PCR Array profiles the expression of 84 key genes responsive to signal transduction pathway activation or inhibition. Cellular signaling forms a complex network of gene interactions involving multiple signal transduction pathways. Each pathway ultimately increases or decreases the expression of its target genes resulting in alteration of cellular processes. Changes in target gene expression suggest signaling pathway activation or inhibition. However, gene expression results in the same pathway vary widely, depending on model systems and experimental conditions. Therefore, multiple target genes from each pathway should be examined to ensure accurate signaling pathway identification in a variety of model systems. In addition, analyzing multiple pathways simultaneously tests for signaling pathway crosstalk. This array includes target genes for 10 commonly studied signal transduction pathways, including pathways important for developmental, immunological, metabolic, and stress-activated processes. Results obtained with this array can suggest pathways that are potentially activated or inhibited by an experimental stimulus for further follow-up studies. Using real-time PCR, research studies can easily and reliably analyze the expression of a focused panel of genes involved in cellular signaling with this array.

For further details, consult the *RT² Profiler PCR Array Handbook*.

Shipping and storage

RT² Profiler PCR Arrays in the Rotor-Gene format are shipped at ambient temperature, on dry ice, or blue ice packs depending on destination and accompanying products.

For long term storage, keep plates at –20°C.

Note: Ensure that you have the correct RT² Profiler PCR Array format for your real-time cycler (see table above).

Note: Open the package and store the products appropriately immediately on receipt.



Sample & Assay Technologies

Array layout

The 96 real-time assays in the Rotor-Gene format are located in wells 1–96 of the Rotor-Disc™ (plate A1–A12=Rotor-Disc 1–12, plate B1–B12=Rotor-Disc 13–24, etc.). To maintain data analysis compatibility, wells 97–100 do not contain real-time assays but will contain master mix to account for weight balance.

Gene table: RT² Profiler PCR Array

Position	UniGene	GenBank	Symbol	Description
A01	Rn.54820	NM_057107	Acs13	Acyl-CoA synthetase long-chain family member 3
A02	Rn.87821	NM_053623	Acs14	Acyl-CoA synthetase long-chain family member 4
A03	Rn.105862	NM_053607	Acs15	Acyl-CoA synthetase long-chain family member 5
A04	Rn.10232	NM_012715	Adm	Adrenomedullin
A05	Rn.10520	NM_012780	Arnt	Aryl hydrocarbon receptor nuclear translocator
A06	Rn.2423	NM_024403	Atf4	Activating transcription factor 4 (tax-responsive enhancer element B67)
A07	Rn.162212	NM_024355	Axin2	Axin 2
A08	Rn.10668	NM_017059	Bax	Bcl2-associated X protein
A09	Rn.25176	NM_173837	Bbc3	Bcl-2 binding component 3
A10	Rn.9996	NM_016993	Bcl2	B-cell CLL/lymphoma 2
A11	Rn.19770	NM_133416	Bcl2a1d	B-cell leukemia/lymphoma 2 related protein A1d
A12	Rn.10323	NM_031535	Bcl2l1	Bcl2-like 1
B01	Rn.64578	NM_023987	Birc3	Baculoviral IAP repeat-containing 3
B02	Rn.90931	NM_017178	Bmp2	Bone morphogenetic protein 2
B03	Rn.10318	NM_012827	Bmp4	Bone morphogenetic protein 4
B04	Rn.27923	NM_017259	Btg2	BTG family, member 2
B05	Rn.162391	NM_001107956	Car9	Carbonic anhydrase 9
B06	Rn.8019	NM_031116	Ccl5	Chemokine (C-C motif) ligand 5
B07	Rn.22279	NM_171992	Ccnd1	Cyclin D1
B08	Rn.96083	NM_022267	Ccnd2	Cyclin D2
B09	Rn.10089	NM_080782	Cdkn1a	Cyclin-dependent kinase inhibitor 1A
B10	Rn.29897	NM_031762	Cdkn1b	Cyclin-dependent kinase inhibitor 1B
B11	Rn.202620	NM_013154	Cebpd	CCAAT/enhancer binding protein (C/EBP), delta
B12	Rn.11389	NM_012930	Cpt2	Carnitine palmitoyltransferase 2
C01	Rn.83632	NM_023981	Csf1	Colony stimulating factor 1 (macrophage)
C02	Rn.161784	NM_024159	Dab2	Disabled homolog 2 (Drosophila)
C03	Rn.37227	NM_031507	Egfr	Epidermal growth factor receptor
C04	Rn.19723	NM_012843	Emp1	Epithelial membrane protein 1
C05	Rn.11365	NM_017001	Epo	Erythropoietin
C06	Rn.36412	NM_012556	Fabp1	Fatty acid binding protein 1, liver
C07	Rn.162521	NM_139194	Fas	Fas (TNF receptor superfamily, member 6)
C08	Rn.10326	NM_133550	Fcer2	Fc fragment of IgE, low affinity II, receptor for (CD23)
C09	Rn.11306	NM_012953	Fosl1	Fos-like antigen 1
C10	Rn.54447	NM_012848	Fth1	Ferritin, heavy polypeptide 1
C11	Rn.10250	NM_024127	Gadd45a	Growth arrest and DNA-damage-inducible, alpha
C12	Rn.35886	NM_001008321	Gadd45b	Growth arrest and DNA-damage-inducible, beta
D01	Rn.92350	NM_133293	Gata3	GATA binding protein 3
D02	Rn.8365	NM_012815	Gclc	Glutamate-cysteine ligase, catalytic subunit
D03	Rn.2460	NM_017305	Gclm	Glutamate cysteine ligase, modifier subunit
D04	Rn.19721	NM_053906	Gsr	Glutathione reductase
D05	Rn.4028	NM_053523	Herpud1	Homocysteine-inducible, endoplasmic reticulum stress-inducible, ubiquitin-like domain member 1
D06	Rn.19727	NM_024360	Hes1	Hairy and enhancer of split 1 (Drosophila)
D07	Rn.22422	NM_024383	Hes5	Hairy and enhancer of split 5 (Drosophila)
D08	Rn.98202	XM_342216	Hey1	Hairy/enhancer-of-split related with YRPW motif 1
D09	Rn.58672	NM_130417	Hey2	Hairy/enhancer-of-split related with YRPW motif 2
D10	Rn.22252	NM_001107977	Heyl	Hairy/enhancer-of-split related with YRPW motif-like
D11	Rn.3160	NM_012580	Hmox1	Heme oxygenase (decycling) 1
D12	Rn.12	NM_012967	Icam1	Intercellular adhesion molecule 1
E01	Rn.2113	NM_012797	Id1	Inhibitor of DNA binding 1
E02	Rn.10795	NM_138880	Ifng	Interferon gamma
E03	Rn.3723	NM_019242	Irf1	Interferon-related developmental regulator 1
E04	Rn.6396	NM_012591	Irf1	Interferon regulatory factor 1
E05	Rn.88804	NM_019147	Jag1	Jagged 1
E06	Rn.107896	NM_017025	Ldha	Lactate dehydrogenase A
E07	Rn.127809	NM_133393	Lfng	LFNG O-fucosylpeptide 3-beta-N-acetylglucosaminyltransferase
E08	Rn.17847	NM_001009717	Lrg1	Leucine-rich alpha-2-glycoprotein 1

Position	UniGene	GenBank	Symbol	Description
E09	Rn.129914	NM_021846	Mcl1	Myeloid cell leukemia sequence 1
E10	Rn.10282	NM_012864	Mmp7	Matrix metalloproteinase 7
E11	Rn.12072	NM_012603	Myc	Myelocytomatosis oncogene
E12	Rn.25046	NM_001105721	Notch1	Notch homolog 1, translocation-associated (Drosophila)
F01	Rn.11234	NM_017000	Nqo1	NAD(P)H dehydrogenase, quinone 1
F02	Rn.87449	NM_133306	Olr1	Oxidized low density lipoprotein (lectin-like) receptor 1
F03	Rn.223	NM_022381	Pcna	Proliferating cell nuclear antigen
F04	Rn.96181	NM_013141	Ppard	Peroxisome proliferator-activated receptor delta
F05	Rn.105585	NM_053566	Ptch1	Patched homolog 1 (Drosophila)
F06	Rn.55115	NM_017045	Rb1	Retinoblastoma 1
F07	Rn.29367	NM_012620	Serpine1	Serpin peptidase inhibitor, clade E (nexin, plasminogen activator inhibitor type 1), member 1
F08	Rn.145068	XM_231115	Slc27a4	Solute carrier family 27 (fatty acid transporter), member 4
F09	Rn.3205	NM_138827	Slc2a1	Solute carrier family 2 (facilitated glucose transporter), member 1
F10	Rn.127801	NM_053565	Socs3	Suppressor of cytokine signaling 3
F11	Rn.110441	XM_001066536	Sorbs1	Sorbin and SH3 domain containing 1
F12	Rn.107103	NM_181550	Sqstm1	Sequestosome 1
G01	Rn.33229	NM_032612	Stat1	Signal transducer and activator of transcription 1
G02	Rn.2275	NM_012675	Tnf	Tumor necrosis factor (TNF superfamily, member 2)
G03	Rn.83627	NM_145681	Tnfsf10	Tumor necrosis factor (ligand) superfamily, member 10
G04	Rn.29777	NM_053800	Txn1	Thioredoxin 1
G05	Rn.67581	NM_031614	Txnrd1	Thioredoxin reductase 1
G06	Rn.1923	NM_031836	Vegfa	Vascular endothelial growth factor A
G07	Rn.63486	NM_031716	Wisp1	WNT1 inducible signaling pathway protein 1
G08	Rn.138108	NM_001105714	Wnt1	Wingless-type MMTV integration site family, member 1
G09	Rn.17219	XM_342308	Wnt2b	Wingless-type MMTV integration site family, member 2B
G10	Rn.218621	XM_220546	Wnt3a	Wingless-type MMTV integration site family, member 3A
G11	Rn.48749	NM_022631	Wnt5a	Wingless-type MMTV integration site family, member 5A
G12	Rn.22149	NM_001108226	Wnt6	Wingless-type MMTV integration site family, member 6
H01	Rn.94978	NM_031144	Actb	Actin, beta
H02	Rn.1868	NM_012512	B2m	Beta-2 microglobulin
H03	Rn.47	NM_012583	Hprt1	Hypoxanthine phosphoribosyltransferase 1
H04	Rn.107896	NM_017025	Ldha	Lactate dehydrogenase A
H05	Rn.973	NM_001007604	Rplp1	Ribosomal protein, large, P1
H06	N/A	U26919	RGDC	Rat Genomic DNA Contamination
H07	N/A	SA_00104	RTC	Reverse Transcription Control
H08	N/A	SA_00104	RTC	Reverse Transcription Control
H09	N/A	SA_00104	RTC	Reverse Transcription Control
H10	N/A	SA_00103	PPC	Positive PCR Control
H11	N/A	SA_00103	PPC	Positive PCR Control
H12	N/A	SA_00103	PPC	Positive PCR Control

Related products

For optimal performance, RT² Profiler PCR Arrays should be used together with the RT² First Strand Kit for cDNA synthesis and RT² SYBR[®] Green qPCR Mastermixes for PCR.

Product	Contents	Cat. no.
RT ² First Strand Kit (12)	Enzymes and reagents for cDNA synthesis	330401
RT ² SYBR Green ROX™ FAST Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the Rotor-Gene Q and other Rotor-Gene cyclers	330620

* Larger kit sizes available; please inquire.

RT² Profiler PCR Array products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN[®], Rotor-Gene[®], Rotor-Disc[™] (QIAGEN Group); ROX[™] (Applied Biosystems Corporation or its subsidiaries); SYBR[®] (Molecular Probes, Inc.).

1067688 03/2011 © 2011 QIAGEN, all rights reserved.

www.qiagen.com

Canada ■ 800-572-9613

China ■ 8621-3865-3865

Denmark ■ 80-885945

Finland ■ 0800-914416

France ■ 01-60-920-930

Germany ■ 02103-29-12000

Hong Kong ■ 800 933 965

Ireland ■ 1800 555 049

Italy ■ 800-787980

Japan ■ 03-6890-7300

Korea (South) ■ 080-000-7145

Luxembourg ■ 8002 2076

Mexico ■ 01-800-7742-436

The Netherlands ■ 0800 0229592

Norway ■ 800-18859

Singapore ■ 1800-742-4368

Spain ■ 91-630-7050

Sweden ■ 020-790282

Switzerland ■ 055-254-22-11

UK ■ 01293-422-911

USA ■ 800-426-8157

Australia ■ 1-800-243-800

Austria ■ 0800/281010

Belgium ■ 0800-79612

Brazil ■ 0800-557779



Sample & Assay Technologies