

RT² Profiler PCR Array (96-Well Format and 384-Well [4 x 96] Format)

Rat Obesity

Cat. no. 330231 PARN-017ZA

For pathway expression analysis

Format	For use with the following real-time cyclers
RT ² Profiler PCR Array, Format A	Applied Biosystems® models 5700, 7000, 7300, 7500, 7700, 7900HT, ViiA™ 7 (96-well block); Bio-Rad® models iCycler®, iQ™ 5, MyiQ™, MyiQ2; Bio-Rad/MJ Research Chromo4™; Eppendorf® Mastercycler® ep realplex models 2, 2s, 4, 4s; Stratagene® models Mx3005P®, Mx3000P®; Takara TP-800
RT ² Profiler PCR Array, Format C	Applied Biosystems models 7500 (Fast block), 7900HT (Fast block), StepOnePlus™, ViiA 7 (Fast block)
RT ² Profiler PCR Array, Format D	Bio-Rad CFX96™; Bio-Rad/MJ Research models DNA Engine Opticon®, DNA Engine Opticon 2; Stratagene Mx4000®
RT ² Profiler PCR Array, Format E	Applied Biosystems models 7900HT (384-well block), ViiA 7 (384-well block); Bio-Rad CFX384™
RT ² Profiler PCR Array, Format F	Roche® LightCycler® 480 (96-well block)
RT ² Profiler PCR Array, Format G	Roche LightCycler 480 (384-well block)
RT ² Profiler PCR Array, Format H	Fluidigm® BioMark™



Sample & Assay Technologies

Description

The Rat Obesity RT² Profiler PCR Array profiles the expression of 84 genes related to obesity. The control of a relative stable body weight depends on a complex interaction among the hormone axes in the peripheral nervous system and neurotransmitter signals in the central nervous system that play a crucial role in maintaining a balanced energy metabolism. This PCR Array includes obesity-related genes that are directly involved in the regulation of energy intake and expenditure. The genes encode for orexigenic peptides, hormones, and receptors; anorectic peptides, hormones, receptors; and central and peripheral signaling molecules involved in energy expenditure. Using real-time PCR, you can easily and reliably analyze expression of a focused panel of genes related to Obesity with this array.

For further details, consult the *RT² Profiler PCR Array Handbook*.

Shipping and storage

RT² Profiler PCR Arrays in formats A, C, D, E, F, and G are shipped at ambient temperature, on dry ice, or blue ice packs depending on destination and accompanying products. RT² Profiler PCR Arrays in format H are shipped on dry ice or blue ice packs.

For long term storage, keep plates at –20°C.

Note: Ensure that you have the correct RT² Profiler PCR Array format for your real-time cycler (see table above).

Note: Open the package and store the products appropriately immediately on receipt.

Array layout (96-well)

For 384-well 4 x 96 PCR arrays, genes are present in a staggered format. Refer to the *RT² Profiler PCR Array Handbook* for layout.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Adcyap1	Adcyap1r1	Adipoq	Adipor1	Adipor2	Adra2b	Adrb1	Agrp	Apoa4	Atrn	Bdnf	Brs3
B	C3	Calca	Calcr	Cartpt	Cck	Cckar	Clps	Cnr1	Cntf	Cntfr	Crh	Crhr1
C	Drd1a	Drd2	Gal	Galr1	Gcg	Gcgr	Gh1	Ghr	Ghrl	Ghsr	Glp1r	Grp
D	Grpr	Hcrt	Hcrt1	Hrh1	Htr2c	Iapp	Il1a	Il1b	Il1r1	Il6	Il6r	Ins1
E	Ins2	Insr	Lep	Lep1	Mc3r	Mchr1	Nmb	Nmbr	Nmu	Nmur1	Npy	Npy1r
F	Nr3c1	Ntrk1	Nts	Ntsr1	Oprk1	Oprm1	Pomc	Ppara	Pparg	Ppargc1a	Prlhr	Ptpn1
G	Pyy	Ramp3	Sigmar1	Sort1	Sst	Sstr1	Thrb	Tnf	Trh	Trhr	Ucn	Ucp1
H	Actb	B2m	Hprt1	Ldha	Rplp1	RGDC	RTC	RTC	RTC	PPC	PPC	PPC

Gene table: RT² Profiler PCR Array

Position	UniGene	GenBank	Symbol	Description
A01	Rn.202559	NM_016989	Adcyap1	Adenylate cyclase activating polypeptide 1
A02	Rn.88408	NM_133511	Adcyap1r1	Adenylate cyclase activating polypeptide 1 receptor 1
A03	Rn.24299	NM_144744	Adipoq	Adiponectin, C1Q and collagen domain containing
A04	Rn.104556	NM_207587	Adipor1	Adiponectin receptor 1
A05	Rn.101984	NM_001037979	Adipor2	Adiponectin receptor 2
A06	Rn.10296	NM_138505	Adra2b	Adrenergic, alpha-2B-, receptor
A07	Rn.87064	NM_012701	Adrb1	Adrenergic, beta-1-, receptor
A08	Rn.137597	NM_033650	Agrp	Agouti related protein homolog (mouse)
A09	Rn.15739	NM_012737	Apoa4	Apolipoprotein A-IV
A10	Rn.53846	NM_031351	Atrn	Attractin
A11	Rn.11266	NM_012513	Bdnf	Brain-derived neurotrophic factor
A12	Rn.86415	NM_152845	Brs3	Bombesin-like receptor 3
B01	Rn.11378	NM_016994	C3	Complement component 3
B02	Rn.90085	NM_017338	Calca	Calcitonin-related polypeptide alpha
B03	Rn.10062	NM_053816	Calcr	Calcitonin receptor
B04	Rn.89164	NM_017110	Cartpt	CART prepropeptide
B05	Rn.9781	NM_012829	Cck	Cholecystokinin
B06	Rn.10184	NM_012688	Cckar	Cholecystokinin A receptor
B07	Rn.6714	NM_013139	Clps	Colipase, pancreatic
B08	Rn.89774	NM_012784	Cnr1	Cannabinoid receptor 1 (brain)
B09	Rn.6067	NM_013166	Cntf	Ciliary neurotrophic factor
B10	Rn.55036	NM_001003929	Cntfr	Ciliary neurotrophic factor receptor
B11	Rn.10349	NM_031019	Crh	Corticotropin releasing hormone
B12	Rn.10499	NM_030999	Crhr1	Corticotropin releasing hormone receptor 1
C01	Rn.24039	NM_012546	Drd1a	Dopamine receptor D1A
C02	Rn.87299	NM_012547	Drd2	Dopamine receptor D2
C03	Rn.8929	NM_033237	Gal	Galanin prepropeptide
C04	Rn.10213	NM_012958	Galr1	Galanin receptor 1
C05	Rn.54383	NM_012707	Gcg	Glucagon
C06	Rn.11225	NM_172092	Gcgr	Glucagon receptor
C07	Rn.146351	NM_001034848	Gh1	Growth hormone 1
C08	Rn.2178	NM_017094	Ghr	Growth hormone receptor
C09	Rn.42103	NM_021669	Ghrl	Ghrelin/obestatin prepropeptide
C10	Rn.74241	NM_032075	Ghsr	Growth hormone secretagogue receptor
C11	Rn.11408	NM_012728	Glp1r	Glucagon-like peptide 1 receptor
C12	Rn.10930	NM_133570	Grp	Gastrin releasing peptide
D01	Rn.10316	NM_012706	Grpr	Gastrin releasing peptide receptor
D02	Rn.7628	NM_013179	Hcrt	Hypocretin
D03	Rn.88262	NM_013064	Hcrt1	Hypocretin (orexin) receptor 1
D04	Rn.81032	NM_017018	Hrh1	Histamine receptor H 1
D05	Rn.9935	NM_012765	Htr2c	5-hydroxytryptamine (serotonin) receptor 2C
D06	Rn.11394	NM_012586	Iapp	Islet amyloid polypeptide
D07	Rn.12300	NM_017019	Il1a	Interleukin 1 alpha
D08	Rn.9869	NM_031512	Il1b	Interleukin 1 beta
D09	Rn.9758	NM_013123	Il1r1	Interleukin 1 receptor, type 1

Position	UniGene	GenBank	Symbol	Description
D10	Rn.9873	NM_012589	Il6	Interleukin 6
D11	Rn.1716	NM_017020	Il6r	Interleukin 6 receptor
D12	Rn.962	NM_019129	Ins1	Insulin 1
E01	Rn.989	NM_019130	Ins2	Insulin 2
E02	Rn.9876	NM_017071	Insr	Insulin receptor
E03	Rn.44444	NM_013076	Lep	Leptin
E04	Rn.9891	NM_012596	Lepr	Leptin receptor
E05	Rn.215838	NM_001025270	Mc3r	Melanocortin 3 receptor
E06	Rn.10822	NM_031758	Mchr1	Melanin-concentrating hormone receptor 1
E07	Rn.130371	NM_001109149	Nmb	Neuromedin B
E08	Rn.89709	NM_012799	Nmbr	Neuromedin B receptor
E09	Rn.47720	NM_022239	Nmu	Neuromedin U
E10	Rn.54457	NM_023100	Nmur1	Neuromedin U receptor 1
E11	Rn.9714	NM_012614	Npy	Neuropeptide Y
E12	Rn.11642	NM_001113357	Npy1r	Neuropeptide Y receptor Y1
F01	Rn.90070	NM_012576	Nr3c1	Nuclear receptor subfamily 3, group C, member 1
F02	Rn.39098	NM_021589	Ntrk1	Neurotrophic tyrosine kinase, receptor, type 1
F03	Rn.60814	NM_001102381	Nts	Neurotensin
F04	Rn.200149	NM_001108967	Ntsr1	Neurotensin receptor 1
F05	Rn.89571	NM_017167	Oprk1	Opioid receptor, kappa 1
F06	Rn.10118	NM_013071	Oprm1	Opioid receptor, mu 1
F07	Rn.108195	NM_139326	Pomc	Proopiomelanocortin
F08	Rn.9753	NM_013196	Ppara	Peroxisome proliferator activated receptor alpha
F09	Rn.23443	NM_013124	Pparg	Peroxisome proliferator-activated receptor gamma
F10	Rn.19172	NM_031347	Ppargc1a	Peroxisome proliferator-activated receptor gamma, coactivator 1 alpha
F11	Rn.138127	NM_139193	Prlhr	Prolactin releasing hormone receptor
F12	Rn.11317	NM_012637	Plpn1	Protein tyrosine phosphatase, non-receptor type 1
G01	Rn.13173	NM_001034080	Pyy	Peptide YY (mapped)
G02	Rn.48672	NM_020100	Ramp3	Receptor (G protein-coupled) activity modifying protein 3
G03	Rn.1129	NM_030996	Sigmar1	Sigma non-opioid intracellular receptor 1
G04	Rn.11286	NM_031767	Sort1	Sortilin 1
G05	Rn.34418	NM_012659	Sst	Somatostatin
G06	Rn.42915	NM_012719	Sstr1	Somatostatin receptor 1
G07	Rn.88692	NM_012672	Thrb	Thyroid hormone receptor beta
G08	Rn.2275	NM_012675	Tnf	Tumor necrosis factor (TNF superfamily, member 2)
G09	Rn.22	NM_013046	Trh	Thyrotropin releasing hormone
G10	Rn.9962	NM_013047	Trhr	Thyrotropin releasing hormone receptor
G11	Rn.10190	NM_019150	Ucn	Urocortin
G12	Rn.10281	NM_012682	Ucp1	Uncoupling protein 1 (mitochondrial, proton carrier)
H01	Rn.94978	NM_031144	Actb	Actin, beta
H02	Rn.1868	NM_012512	B2m	Beta-2 microglobulin
H03	Rn.47	NM_012583	Hprt1	Hypoxanthine phosphoribosyltransferase 1
H04	Rn.107896	NM_017025	Ldha	Lactate dehydrogenase A
H05	Rn.973	NM_001007604	Rplp1	Ribosomal protein, large, P1
H06	N/A	U26919	RGDC	Rat Genomic DNA Contamination
H07	N/A	SA_00104	RTC	Reverse Transcription Control
H08	N/A	SA_00104	RTC	Reverse Transcription Control
H09	N/A	SA_00104	RTC	Reverse Transcription Control
H10	N/A	SA_00103	PPC	Positive PCR Control
H11	N/A	SA_00103	PPC	Positive PCR Control
H12	N/A	SA_00103	PPC	Positive PCR Control

Related products

For optimal performance, RT² Profiler PCR Arrays should be used together with the RT² First Strand Kit for cDNA synthesis and RT² SYBR[®] Green qPCR Mastermixes for PCR.

Product	Contents	Cat. no.
RT ² First Strand Kit (12)	Enzymes and reagents for cDNA synthesis	330401
RT ² SYBR Green qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with real-time cyclers that do not require a reference dye, including: Bio-Rad models CFX96, CFX384, DNA Engine Opticon 2; Bio-Rad/MJ Research Chromo4; Roche LightCycler 480 (96-well and 384-well); all other cyclers	330500
RT ² SYBR Green ROX [™] qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the following real-time cyclers: Applied Biosystems models 5700, 7000, 7300, 7500 [Standard and FAST], 7700, 7900HT 96-well block [Standard and FAST] and 384-well block, StepOnePlus; Eppendorf Mastercycler ep realplex models 2, 2S, 4, 4S; Stratagene models Mx3000P, Mx3005P, Mx4000; Takara TP-800	330520
RT ² SYBR Green Fluor qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the following real-time cyclers: Bio-Rad models iCycler, iQ5, MyiQ, MyiQ2	330510

* Larger kit sizes available; please inquire.

RT² Profiler PCR Array products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN® (QIAGEN Group); Applied Biosystems®, ViiA™, StepOnePlus™, ROX™ (Applied Biosystems Corporation or its subsidiaries); Bio-Rad®, iCycler®, iQ™, MyiQ™, Chromo4™, CFX96™, DNA Engine Opticon®, CFX384™ (Bio-Rad Laboratories, Inc.); Stratagene®, Mx3005P®, Mx3000P®, Mx4000® (Stratagene); Eppendorf®, Mastercycler® (Eppendorf AG); Roche®, LightCycler® (Roche Group); Fluidigm® BioMark™ (Fluidigm Corporation); SYBR® (Molecular Probes, Inc.).

1066029 03/2011 © 2011 QIAGEN, all rights reserved.

www.qiagen.com

Canada ■ 800-572-9613

China ■ 8621-3865-3865

Denmark ■ 80-885945

Finland ■ 0800-914416

France ■ 01-60-920-930

Germany ■ 02103-29-12000

Hong Kong ■ 800 933 965

Ireland ■ 1800 555 049

Italy ■ 800-787980

Japan ■ 03-6890-7300

Korea (South) ■ 080-000-7145

Luxembourg ■ 8002 2076

Mexico ■ 01-800-7742-436

The Netherlands ■ 0800 0229592

Norway ■ 800-18859

Singapore ■ 1800-742-4368

Spain ■ 91-630-7050

Sweden ■ 020-790282

Switzerland ■ 055-254-22-11

UK ■ 01293-422-911

USA ■ 800-426-8157

Australia ■ 1-800-243-800

Austria ■ 0800/281010

Belgium ■ 0800-79612

Brazil ■ 0800-557779



Sample & Assay Technologies