

RT² Profiler PCR Array (96-Well Format and 384-Well [4 x 96] Format)

Rat DNA Damage Signaling Pathway

Cat. no. 330231 PARN-029ZA

For pathway expression analysis

Format	For use with the following real-time cyclers
RT ² Profiler PCR Array, Format A	Applied Biosystems® models 5700, 7000, 7300, 7500, 7700, 7900HT, ViiA™ 7 (96-well block); Bio-Rad® models iCycler®, iQ™ 5, MyiQ™, MyiQ2; Bio-Rad/MJ Research Chromo4™; Eppendorf® Mastercycler® ep realplex models 2, 2s, 4, 4s; Stratagene® models Mx3005P®, Mx3000P®; Takara TP-800
RT ² Profiler PCR Array, Format C	Applied Biosystems models 7500 (Fast block), 7900HT (Fast block), StepOnePlus™, ViiA 7 (Fast block)
RT ² Profiler PCR Array, Format D	Bio-Rad CFX96™; Bio-Rad/MJ Research models DNA Engine Opticon®, DNA Engine Opticon 2; Stratagene Mx4000®
RT ² Profiler PCR Array, Format E	Applied Biosystems models 7900HT (384-well block), ViiA 7 (384-well block); Bio-Rad CFX384™
RT ² Profiler PCR Array, Format F	Roche® LightCycler® 480 (96-well block)
RT ² Profiler PCR Array, Format G	Roche LightCycler 480 (384-well block)
RT ² Profiler PCR Array, Format H	Fluidigm® BioMark™



Sample & Assay Technologies

Description

The Rat DNA Damage Signaling RT² Profiler PCR Array profiles the expression of 84 genes involved in DNA damage signaling pathways. The genes featured are those associated with the ATR/ATM signaling network and transcriptional targets of DNA damage response. DNA damage can result in cell cycle arrest, apoptosis, and the stabilization and repair of the cellular genome. Using real-time PCR, you can easily and reliably analyze expression of a focused panel of genes related to DNA Damage Signaling with this array.

For further details, consult the *RT² Profiler PCR Array Handbook*.

Shipping and storage

RT² Profiler PCR Arrays in formats A, C, D, E, F, and G are shipped at ambient temperature, on dry ice, or blue ice packs depending on destination and accompanying products. RT² Profiler PCR Arrays in format H are shipped on dry ice or blue ice packs.

For long term storage, keep plates at -20°C .

Note: Ensure that you have the correct RT² Profiler PCR Array format for your real-time cycler (see table above).

Note: Open the package and store the products appropriately immediately on receipt.

Array layout (96-well)

For 384-well 4 x 96 PCR arrays, genes are present in a staggered format. Refer to the *RT² Profiler PCR Array Handbook* for layout.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Abl1	Apex1	Atm	Atrx	Bard1	Bax	Bbc3	Blm	Brca1	Brca2	Cdc25a	Cdc25c
B	Cdkn1a	Chek1	Chek2	Csnk2a2	Dclre1a	Ddb2	Ddit3	Ercc1	Ercc2	Exo1	Fanca	Fancc
C	Fancd2	Fancg	Fen1	Gadd45a	Gadd45g	Hus1	Lig1	Mbd4	Mgmt	Mif	Mlh1	Mlh3
D	Mpg	Mre11a	Msh2	Msh3	Nbn	Nthl1	Ogg1	Parp1	Parp2	Pcna	Pms1	Pms2
E	Pold3	Pole	Polh	Poli	Ppm1d	Ppp1r15a	Prkdc	Ptfg1	Rad1	Rad17	Rad18	Rad21
F	Rad50	Rad51	Rad51c	Rad51ll	Rad52	Rad9	Rev1	Rnf8	Rpa1	Sir1	Smc1a	Smc3
G	Sumo1	Terf1	Topbp1	Tp53	Tp53bp1	Ung	Wrn	Wrnip1	Xpc	Xrcc1	Xrcc2	Xrcc6
H	Actb	B2m	Hprt1	Ldha	Rplp1	RGDC	RTC	RTC	RTC	PPC	PPC	PPC

Gene table: RT² Profiler PCR Array

Position	UniGene	GenBank	Symbol	Description
A01	Rn.3105	NM_001100850	Abl1	C-abl oncogene 1, receptor tyrosine kinase
A02	Rn.5949	NM_024148	Apex1	APEX nuclease (multifunctional DNA repair enzyme) 1
A03	Rn.214048	NM_001106821	Atm	Ataxia telangiectasia mutated homolog (human)
A04	Rn.107838	NM_001105757	Atrx	Alpha thalassemia/mental retardation syndrome X-linked (RAD54 homolog, <i>S. cerevisiae</i>)
A05	Rn.48735	NM_022622	Bard1	BRCA1 associated RING domain 1
A06	Rn.10668	NM_017059	Bax	Bcl2-associated X protein
A07	Rn.25176	NM_173837	Bbc3	Bcl-2 binding component 3
A08	Rn.211198	NM_001107526	Blm	Bloom syndrome, RecQ helicase-like
A09	Rn.48840	NM_012514	Brca1	Breast cancer 1
A10	Rn.103225	NM_031542	Brca2	Breast cancer 2
A11	Rn.11390	NM_133571	Cdc25a	Cell division cycle 25 homolog A (<i>S. pombe</i>)
A12	Rn.162298	NM_001107396	Cdc25c	Cell division cycle 25 homolog C (<i>S. pombe</i>)
B01	Rn.10089	NM_080782	Cdkn1a	Cyclin-dependent kinase inhibitor 1A
B02	Rn.33267	NM_080400	Chek1	CHK1 checkpoint homolog (<i>S. pombe</i>)
B03	Rn.163213	NM_053677	Chek2	CHK2 checkpoint homolog (<i>S. pombe</i>)
B04	Rn.24013	NM_001107409	Csnk2a2	Casein kinase 2, alpha prime polypeptide
B05	Rn.41244	NM_001106201	Dclre1a	DNA cross-link repair 1A, PSO2 homolog (<i>S. cerevisiae</i>)
B06	Rn.154614	XM_242065	Ddb2	Damage specific DNA binding protein 2
B07	Rn.11183	NM_024134	Ddit3	DNA-damage inducible transcript 3
B08	Rn.7320	NM_001106228	Ercc1	Excision repair cross-complementing rodent repair deficiency, complementation group 1
B09	Rn.74906	NM_001172809	Ercc2	Excision repair cross-complementing rodent repair deficiency, complementation group 2
B10	Rn.52638	NM_001107198	Exo1	Exonuclease 1
B11	Rn.163415	NM_001108455	Fanca	Fanconi anemia, complementation group A
B12	Rn.10798	NM_012557	Fancc	Fanconi anemia, complementation group C
C01	Rn.203979	NM_001001719	Fancd2	Fanconi anemia, complementation group D2
C02	Rn.198713	XM_001076851	Fancg	Fanconi anemia, complementation group G
C03	Rn.16664	NM_053430	Fen1	Flap structure-specific endonuclease 1
C04	Rn.10250	NM_024127	Gadd45a	Growth arrest and DNA-damage-inducible, alpha
C05	Rn.16950	NM_001077640	Gadd45g	Growth arrest and DNA-damage-inducible, gamma
C06	Rn.12812	NM_001109092	Hus1	HUS1 checkpoint homolog (<i>S. pombe</i>)
C07	Rn.20467	NM_030855	Lig1	Ligase I, DNA, ATP-dependent
C08	Rn.179073	XM_001059437	Mbd4	Methyl-CpG binding domain protein 4
C09	Rn.9836	NM_012861	Mgmt	O-6-methylguanine-DNA methyltransferase
C10	Rn.2661	NM_031051	Mif	Macrophage migration inhibitory factor
C11	Rn.20391	NM_031053	Mlh1	MutL homolog 1 (<i>E. coli</i>)
C12	Rn.119254	NM_001108043	Mlh3	MutL homolog 3 (<i>E. coli</i>)
D01	Rn.11241	NM_012601	Mpg	N-methylpurine-DNA glycosylase
D02	Rn.209040	NM_022279	Mre11a	MRE11 meiotic recombination 11 homolog A (<i>S. cerevisiae</i>)
D03	Rn.3174	NM_031058	Msh2	MutS homolog 2 (<i>E. coli</i>)
D04	Rn.162551	XM_001065837	Msh3	MutS homolog 3 (<i>E. coli</i>)
D05	Rn.25214	NM_138873	Nbn	Nibrin
D06	Rn.14632	NM_001105728	Nthl1	Nth (endonuclease III)-like 1 (<i>E.coli</i>)

Position	UniGene	GenBank	Symbol	Description
D07	Rn.22623	NM_030870	Ogg1	8-oxoguanine DNA glycosylase
D08	Rn.11327	NM_013063	Parp1	Poly (ADP-ribose) polymerase 1
D09	Rn.22730	NM_001106030	Parp2	Poly (ADP-ribose) polymerase 2
D10	Rn.223	NM_022381	Pcna	Proliferating cell nuclear antigen
D11	Rn.47945	NM_001009535	Pms1	Postmeiotic segregation increased 1 (S. cerevisiae)
D12	Rn.102072	NM_001105908	Pms2	PMS2 postmeiotic segregation increased 2 (S. cerevisiae)
E01	Rn.162469	NM_001024750	Pold3	Polymerase (DNA-directed), delta 3, accessory subunit
E02	Rn.40977	NM_001107152	Pole	Polymerase (DNA directed), epsilon
E03	Rn.53711	NM_001108204	Polh	Polymerase (DNA directed), eta
E04	Rn.136855	NM_001106137	Poli	Polymerase (DNA directed), iota
E05	Rn.15540	NM_001105825	Ppm1d	Protein phosphatase 1D magnesium-dependent, delta isoform
E06	Rn.2232	NM_133546	Ppp1r15a	Protein phosphatase 1, regulatory (inhibitor) subunit 15A
E07	Rn.24110	NM_001108327	Prkdc	Protein kinase, DNA activated, catalytic polypeptide
E08	Rn.271	NM_022391	Pttg1	Pituitary tumor-transforming 1
E09	Rn.27018	NM_001106419	Rad1	RAD1 homolog (S. pombe)
E10	Rn.154275	NM_001024778	Rad17	RAD17 homolog (S. pombe)
E11	Rn.22793	NM_001077673	Rad18	RAD18 homolog (S. cerevisiae)
E12	Rn.3991	NM_001025701	Rad21	RAD21 homolog (S. pombe)
F01	Rn.51136	NM_022246	Rad50	RAD50 homolog (S. cerevisiae)
F02	Rn.214052	NM_001109204	Rad51	RAD51 homolog (RecA homolog, E. coli) (S. cerevisiae)
F03	Rn.160309	NM_001129777	Rad51c	Rad51 homolog c (S. cerevisiae)
F04	Rn.124221	XM_576058	Rad51l1	RAD51-like 1 (S. cerevisiae)
F05	Rn.8154	NM_001106617	Rad52	RAD52 homolog (S. cerevisiae)
F06	Rn.145016	XM_219684	Rad9	RAD9 homolog (S. pombe)
F07	Rn.94641	NM_001108213	Rev1	REV1 homolog (S. cerevisiae)
F08	Rn.12463	NM_001025727	Rnf8	Ring finger protein 8
F09	Rn.14841	NM_001047843	Rpa1	Replication protein A1
F10	Rn.219976	NM_001107627	Sirt1	Sirtuin (silent mating type information regulation 2 homolog) 1 (S. cerevisiae)
F11	Rn.11763	NM_031683	Smc1a	Structural maintenance of chromosomes 1A
F12	Rn.11074	NM_031583	Smc3	Structural maintenance of chromosomes 3
G01	Rn.1221	NM_001009672	Sumo1	SMT3 suppressor of mif two 3 homolog 1 (S. cerevisiae)
G02	Rn.33853	NM_001012464	Terf1	Telomeric repeat binding factor (NIMA-interacting) 1
G03	Rn.204161	XM_236578	Topbp1	Topoisomerase (DNA) II binding protein 1
G04	Rn.54443	NM_030989	Tp53	Tumor protein p53
G05	Rn.221896	NM_001106501	Tp53bp1	Tumor protein p53 binding protein 1
G06	Rn.104379	NM_001013124	Ung	Uracil-DNA glycosylase
G07	Rn.144144	XM_214361	Wrrn	Werner syndrome
G08	Rn.8783	NM_172332	Wrrnip1	Werner helicase interacting protein 1
G09	Rn.22820	NM_001107874	Xpc	Xeroderma pigmentosum, complementation group C
G10	Rn.13754	NM_053435	Xrcc1	X-ray repair complementing defective repair in Chinese hamster cells 1
G11	Rn.8224	NM_001109215	Xrcc2	X-ray repair complementing defective repair in Chinese hamster cells 2
G12	Rn.161996	NM_139080	Xrcc6	X-ray repair complementing defective repair in Chinese hamster cells 6
H01	Rn.94978	NM_031144	Actb	Actin, beta
H02	Rn.1868	NM_012512	B2m	Beta-2 microglobulin
H03	Rn.47	NM_012583	Hprt1	Hypoxanthine phosphoribosyltransferase 1
H04	Rn.107896	NM_017025	Ldha	Lactate dehydrogenase A
H05	Rn.973	NM_001007604	Rplp1	Ribosomal protein, large, P1
H06	N/A	U26919	RGDC	Rat Genomic DNA Contamination
H07	N/A	SA_00104	RTC	Reverse Transcription Control
H08	N/A	SA_00104	RTC	Reverse Transcription Control
H09	N/A	SA_00104	RTC	Reverse Transcription Control
H10	N/A	SA_00103	PPC	Positive PCR Control
H11	N/A	SA_00103	PPC	Positive PCR Control
H12	N/A	SA_00103	PPC	Positive PCR Control

Related products

For optimal performance, RT² Profiler PCR Arrays should be used together with the RT² First Strand Kit for cDNA synthesis and RT² SYBR[®] Green qPCR Mastermixes for PCR.

Product	Contents	Cat. no.
RT ² First Strand Kit (12)	Enzymes and reagents for cDNA synthesis	330401
RT ² SYBR Green qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with real-time cyclers that do not require a reference dye, including: Bio-Rad models CFX96, CFX384, DNA Engine Opticon 2; Bio-Rad/MJ Research Chromo4; Roche LightCycler 480 (96-well and 384-well); all other cyclers	330500
RT ² SYBR Green ROX™ qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the following real-time cyclers: Applied Biosystems models 5700, 7000, 7300, 7500 [Standard and FAST], 7700, 7900HT 96-well block [Standard and FAST] and 384-well block, StepOnePlus; Eppendorf Mastercycler ep realplex models 2, 2S, 4, 4S; Stratagene models Mx3000P, Mx3005P, Mx4000; Takara TP-800	330520
RT ² SYBR Green Fluor qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the following real-time cyclers: Bio-Rad models iCycler, iQ5, MyiQ, MyiQ2	330510

* Larger kit sizes available; please inquire.

RT² Profiler PCR Array products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN® (QIAGEN Group); Applied Biosystems®, ViiA™, StepOnePlus™, ROX™ (Applied Biosystems Corporation or its subsidiaries); Bio-Rad®, iCycler®, iQ™, MyiQ™, Chromo4™, CFX96™, DNA Engine Opticon®, CFX384™ (Bio-Rad Laboratories, Inc.); Stratagene®, Mx3005P®, Mx3000P®, Mx4000® (Stratagene); Eppendorf®, Mastercycler® (Eppendorf AG); Roche®, LightCycler® (Roche Group); Fluidigm® BioMark™ (Fluidigm Corporation); SYBR® (Molecular Probes, Inc.).

1066029 03/2011 © 2011 QIAGEN, all rights reserved.

www.qiagen.com

Canada ■ 800-572-9613

Ireland ■ 1800 555 049

Norway ■ 800-18859

China ■ 8621-3865-3865

Italy ■ 800-787980

Singapore ■ 1800-742-4368

Denmark ■ 80-885945

Japan ■ 03-6890-7300

Spain ■ 91-630-7050

Australia ■ 1-800-243-800

Finland ■ 0800-914416

Korea (South) ■ 080-000-7145

Sweden ■ 020-790282

Austria ■ 0800/281010

France ■ 01-60-920-930

Luxembourg ■ 8002 2076

Switzerland ■ 055-254-22-11

Belgium ■ 0800-79612

Germany ■ 02103-29-12000

Mexico ■ 01-800-7742-436

UK ■ 01293-422-911

Brazil ■ 0800-557779

Hong Kong ■ 800 933 965

The Netherlands ■ 0800 0229592

USA ■ 800-426-8157



Sample & Assay Technologies