

# RT<sup>2</sup> Profiler PCR Array (Rotor-Gene® Format)

## Rat Growth Factors

Cat. no. 330231 PARN-041ZR

For pathway expression analysis

Format	For use with the following real-time cyclers
RT <sup>2</sup> Profiler PCR Array, Format R	Rotor-Gene Q, other Rotor-Gene cyclers

### Description

The Rat Growth Factors RT<sup>2</sup> Profiler PCR Array profiles the expression of 84 genes related to growth factors. Growth factors play a vital role in various normal biological processes such as embryogenesis, wound healing and inflammation. This array contains angiogenic growth factors and regulators of apoptosis. Genes involved in cell differentiation are included as well. Also represented are genes related to embryonic development as well as genes involved in tissue-specific development. Using real-time PCR, you can easily and reliably analyze expression of a focused panel of genes related to the growth factors with this array.

For further details, consult the *RT<sup>2</sup> Profiler PCR Array Handbook*.

### Shipping and storage

RT<sup>2</sup> Profiler PCR Arrays in the Rotor-Gene format are shipped at ambient temperature, on dry ice, or blue ice packs depending on destination and accompanying products.

For long term storage, keep plates at –20°C.

**Note:** Ensure that you have the correct RT<sup>2</sup> Profiler PCR Array format for your real-time cycler (see table above).

**Note:** Open the package and store the products appropriately immediately on receipt.



---

Sample & Assay Technologies

## Array layout

The 96 real-time assays in the Rotor-Gene format are located in wells 1–96 of the Rotor-Disc™ (plate A1–A12=Rotor-Disc 1–12, plate B1–B12=Rotor-Disc 13–24, etc.). To maintain data analysis compatibility, wells 97–100 do not contain real-time assays but will contain master mix to account for weight balance.

## Gene table: RT<sup>2</sup> Profiler PCR Array

Position	UniGene	GenBank	Symbol	Description
A01	Rn.19121	NM_012902	Amh	Anti-Mullerian hormone
A02	Rn.154711	NM_053397	Artn	Artemin
A03	Rn.11266	NM_012513	Bdnf	Brain-derived neurotrophic factor
A04	Rn.9305	NM_031323	Bmp1	Bone morphogenetic protein 1
A05	Rn.118010	NM_001031824	Bmp10	Bone morphogenetic protein 10
A06	Rn.90931	NM_017178	Bmp2	Bone morphogenetic protein 2
A07	Rn.208979	NM_017105	Bmp3	Bone morphogenetic protein 3
A08	Rn.10318	NM_012827	Bmp4	Bone morphogenetic protein 4
A09	Rn.135799	NM_001108168	Bmp5	Bone morphogenetic protein 5
A10	Rn.40476	NM_013107	Bmp6	Bone morphogenetic protein 6
A11	Rn.18030	XM_342591	Bmp7	Bone morphogenetic protein 7
A12	Rn.126011	NM_001109432	Bmp8a	Bone morphogenetic protein 8a
B01	Rn.130193	NM_207615	Clcf1	Cardiotrophin-like cytokine factor 1
B02	Rn.83632	NM_023981	Csf1	Colony stimulating factor 1 (macrophage)
B03	Rn.44285	XM_340799	Csf2	Colony stimulating factor 2 (granulocyte-macrophage)
B04	Rn.53973	NM_017104	Csf3	Colony stimulating factor 3 (granulocyte)
B05	Rn.10907	NM_030845	Cxcl1	Chemokine (C-X-C motif) ligand 1 (melanoma growth stimulating activity, alpha)
B06	Rn.54439	NM_022177	Cxcl12	Chemokine (C-X-C motif) ligand 12 (stromal cell-derived factor 1)
B07	Rn.6075	NM_012842	Egf	Epidermal growth factor
B08	Rn.42897	NM_021689	Ereg	Epiregulin
B09	Rn.88013	NM_012846	Fgf1	Fibroblast growth factor 1
B10	Rn.44439	NM_012951	Fgf10	Fibroblast growth factor 10
B11	Rn.73083	NM_130816	Fgf11	Fibroblast growth factor 11
B12	Rn.15709	NM_053428	Fgf13	Fibroblast growth factor 13
C01	Rn.44476	NM_022223	Fgf14	Fibroblast growth factor 14
C02	Rn.81230	NM_130753	Fgf15	Fibroblast growth factor 15
C03	Rn.44353	NM_019198	Fgf17	Fibroblast growth factor 17
C04	Rn.44478	NM_019199	Fgf18	Fibroblast growth factor 18
C05	Rn.31808	NM_019305	Fgf2	Fibroblast growth factor 2
C06	Rn.198777	NM_130751	Fgf22	Fibroblast growth factor 22
C07	Rn.81236	NM_130817	Fgf3	Fibroblast growth factor 3
C08	Rn.81223	NM_053809	Fgf4	Fibroblast growth factor 4
C09	Rn.44445	NM_022211	Fgf5	Fibroblast growth factor 5
C10	Rn.81222	NM_131908	Fgf6	Fibroblast growth factor 6
C11	Rn.98842	NM_022182	Fgf7	Fibroblast growth factor 7
C12	Rn.73565	NM_133286	Fgf8	Fibroblast growth factor 8
D01	Rn.25174	NM_012952	Fgf9	Fibroblast growth factor 9
D02	Rn.10796	NM_031761	Figf	C-fos induced growth factor
D03	Rn.8707	NM_024375	Gdf10	Growth differentiation factor 10
D04	Rn.51292	XM_343148	Gdf11	Growth differentiation factor 11
D05	Rn.127822	XM_001066344	Gdf5	Growth differentiation factor 5
D06	Rn.53970	NM_019139	Gdnf	Glial cell derived neurotrophic factor
D07	Rn.10468	NM_017017	Hgf	Hepatocyte growth factor
D08	Rn.6282	NM_178866	Igf1	Insulin-like growth factor 1
D09	Rn.118681	NM_031511	Igf2	Insulin-like growth factor 2
D10	Rn.198483	NM_133519	Il11	Interleukin 11
D11	Rn.207199	NM_053390	Il12a	Interleukin 12a
D12	Rn.11118	NM_019165	Il18	Interleukin 18
E01	Rn.12300	NM_017019	Il1a	Interleukin 1 alpha
E02	Rn.9869	NM_031512	Il1b	Interleukin 1 beta
E03	Rn.9871	NM_053836	Il2	Interleukin 2
E04	Rn.10652	NM_031513	Il3	Interleukin 3
E05	Rn.108255	NM_201270	Il4	Interleukin 4
E06	Rn.9873	NM_012589	Il6	Interleukin 6
E07	Rn.10793	NM_013110	Il7	Interleukin 7
E08	Rn.8831	NM_012590	Inha	Inhibin alpha
E09	Rn.9874	NM_017128	Inhba	Inhibin beta-A

Position	UniGene	GenBank	Symbol	Description
E10	Rn.35074	NM_080771	Inhbb	Inhibin beta-B
E11	Rn.44216	NM_021843	Kitlg	KIT ligand
E12	Rn.218667	NM_001109080	Lefty1	Left right determination factor 1
F01	Rn.44444	NM_013076	Lep	Leptin
F02	Rn.44379	NM_022196	Lif	Leukemia inhibitory factor
F03	Rn.17447	NM_030859	Mdk	Midkine
F04	Rn.44460	NM_019151	Mstn	Myostatin
F05	Rn.22168	XM_227525	Ngf	Nerve growth factor (beta polypeptide)
F06	Rn.218528	NM_001106394	Nodal	Nodal homolog (mouse)
F07	Rn.37438	NM_031588	Nrg1	Neuregulin 1
F08	Rn.9715	NM_031073	Ntf3	Neurotrophin 3
F09	Rn.44225	NM_013184	Ntf4	Neurotrophin 4
F10	Rn.10999	NM_012801	Pdgfa	Platelet-derived growth factor alpha polypeptide
F11	Rn.6960	NM_053595	Pgf	Placental growth factor
F12	Rn.5876	NM_019124	Rabep1	Rabaptin, RAB GTPase binding effector protein 1
G01	Rn.3233	NM_053485	S100a6	S100 calcium binding protein A6
G02	Rn.8871	NM_012881	Spp1	Secreted phosphoprotein 1
G03	N/A	XM_001056317	Tdgf1	Teratocarcinoma-derived growth factor 1
G04	Rn.10715	NM_057129	Tff1	Trefoil factor 1
G05	Rn.9952	NM_012671	Tgfa	Transforming growth factor alpha
G06	Rn.40136	NM_021578	Tgfb1	Transforming growth factor, beta 1
G07	Rn.24539	NM_031131	Tgfb2	Transforming growth factor, beta 2
G08	Rn.7018	NM_013174	Tgfb3	Transforming growth factor, beta 3
G09	Rn.1923	NM_031836	Vegfa	Vascular endothelial growth factor A
G10	Rn.198550	NM_053549	Vegfb	Vascular endothelial growth factor B
G11	Rn.6913	NM_053653	Vegfc	Vascular endothelial growth factor C
G12	Rn.1772	NM_001169120	Zfp91	Zinc finger protein 91
H01	Rn.94978	NM_031144	Actb	Actin, beta
H02	Rn.1868	NM_012512	B2m	Beta-2 microglobulin
H03	Rn.47	NM_012583	Hprt1	Hypoxanthine phosphoribosyltransferase 1
H04	Rn.107896	NM_017025	Ldha	Lactate dehydrogenase A
H05	Rn.973	NM_001007604	Rplp1	Ribosomal protein, large, P1
H06	N/A	U26919	RGDC	Rat Genomic DNA Contamination
H07	N/A	SA_00104	RTC	Reverse Transcription Control
H08	N/A	SA_00104	RTC	Reverse Transcription Control
H09	N/A	SA_00104	RTC	Reverse Transcription Control
H10	N/A	SA_00103	PPC	Positive PCR Control
H11	N/A	SA_00103	PPC	Positive PCR Control
H12	N/A	SA_00103	PPC	Positive PCR Control

## Related products

For optimal performance, RT<sup>2</sup> Profiler PCR Arrays should be used together with the RT<sup>2</sup> First Strand Kit for cDNA synthesis and RT<sup>2</sup> SYBR<sup>®</sup> Green qPCR Mastermixes for PCR.

Product	Contents	Cat. no.
RT <sup>2</sup> First Strand Kit (12)	Enzymes and reagents for cDNA synthesis	330401
RT <sup>2</sup> SYBR Green ROX <sup>™</sup> FAST Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the Rotor-Gene Q and other Rotor-Gene cyclers	330620

\* Larger kit sizes available; please inquire.

RT<sup>2</sup> Profiler PCR Array products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at [www.qiagen.com](http://www.qiagen.com) or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN<sup>®</sup>, Rotor-Gene<sup>®</sup>, Rotor-Disc<sup>™</sup> (QIAGEN Group); ROX<sup>™</sup> (Applied Biosystems or its subsidiaries); SYBR<sup>®</sup> (Molecular Probes, Inc.).

1067688 03/2011 © 2011 QIAGEN, all rights reserved.

[www.qiagen.com](http://www.qiagen.com)

Australia ■ 1-800-243-800

Austria ■ 0800/281010

Belgium ■ 0800-79612

Brazil ■ 0800-557779

Canada ■ 800-572-9613

China ■ 8621-3865-3865

Denmark ■ 80-885945

Finland ■ 0800-914416

France ■ 01-60-920-930

Germany ■ 02103-29-12000

Hong Kong ■ 800 933 965

Ireland ■ 1800 555 049

Italy ■ 800-787980

Japan ■ 03-6890-7300

Korea (South) ■ 080-000-7145

Luxembourg ■ 8002 2076

Mexico ■ 01-800-7742-436

The Netherlands ■ 0800 0229592

Norway ■ 800-18859

Singapore ■ 1800-742-4368

Spain ■ 91-630-7050

Sweden ■ 020-790282

Switzerland ■ 055-254-22-11

UK ■ 01293-422-911

USA ■ 800-426-8157



Sample & Assay Technologies