

RT² Profiler PCR Array (Rotor-Gene[®] Format)

Rat Cancer Inflammation & Immunity Crosstalk

Cat. no. 330231 PARN-181ZR

For pathway expression analysis

Format	For use with the following real-time cyclers
RT ² Profiler PCR Array, Format R	Rotor-Gene Q, other Rotor-Gene cyclers

Description

The Rat Cancer Inflammation & Immunity Crosstalk RT² Profiler PCR Array profiles the expression of 84 key genes involved in mediating communication between tumor cells and the cellular mediators of inflammation and immunity. In addition to epithelial and stromal compartments, the tumor microenvironment contains several cell types of the innate and adaptive immune systems including B and T lymphocytes, dendritic cells, and macrophages. In response to tumor-associated antigens presented via MHC Class I molecules, or to abnormal molecular patterns recognized by Toll-like receptors, the immune system eliminates target cells using a variety of effector enzymes and the engagement of pro-apoptotic signals including TRAIL and FAS ligand. If normal homeostasis is not resolved quickly, a state of chronic inflammation can ensue, including locally increased levels of reactive oxygen and nitrogen species that promote genomic instability. Immune cells produce a variety of cytokines that coordinate the inflammatory response, which is fueled by positive feedback loops commonly involving the STAT and NFκB signaling pathways in tumor cells. The resulting upregulation of antiapoptotic and immunosuppressive factors enables transformed cells to proliferate unchecked by the immune system. During cancer progression, the repertoire of chemokines, cytokines, and growth factors that orchestrates normal immune responses can be commandeered to create an immunosuppressive state that facilitates invasion and metastasis. The genes profiled with this array include mediators and effectors of the cross-talk between tumors and the immune system that influences the course of cancer progression. A set of controls present on each array enables data analysis using the Delta-Delta CT method of relative quantification as well as assessment of reverse transcription performance, genomic DNA contamination, and PCR performance. Using real-time PCR, researchers can easily and reliably analyze the expression of a focused panel of genes involved in cancer

inflammation and immune crosstalk with this array.

The RT² Profiler PCR Arrays are intended for molecular biology applications. This product is not intended for the diagnosis, prevention, or treatment of a disease.

For further details, consult the *RT² Profiler PCR Array Handbook*.

Shipping and storage

RT² Profiler PCR Arrays in the Rotor-Gene format are shipped at ambient temperature, on dry ice, or blue ice packs depending on destination and accompanying products.

For long term storage, keep plates at –20°C.

Note: Ensure that you have the correct RT² Profiler PCR Array format for your real-time cycler (see table above).

Note: Open the package and store the products appropriately immediately on receipt.



Array layout

The 96 real-time assays in the Rotor-Gene format are located in wells 1–96 of the Rotor-Disc™ (plate A1–A12=Rotor-Disc 1–12, plate B1–B12=Rotor-Disc 13–24, etc.). To maintain data analysis compatibility, wells 97–100 do not contain real-time assays but will contain master mix to account for weight balance.

Gene table: RT² Profiler PCR Array

Position	UniGene	GenBank	Symbol	Description
A01	Rn.118188	NM_001100779	Aicda	Activation-induced cytidine deaminase
A02	Rn.9996	NM_016993	Bcl2	B-cell CLL/lymphoma 2
A03	Rn.10323	NM_031535	Bcl2l1	Bcl2-like 1
A04	Rn.4772	NM_031530	Ccl2	Chemokine (C-C motif) ligand 2
A05	Rn.10722	NM_019233	Ccl20	Chemokine (C-C motif) ligand 20
A06	Rn.39658	NM_001008513	Ccl21	Chemokine (C-C motif) ligand 21
A07	Rn.48727	NM_057203	Ccl22	Chemokine (C-C motif) ligand 22
A08	Rn.195764	NM_053700	Ccl28	Chemokine (C-C motif) ligand 28
A09	Rn.37880	NM_053858	Ccl4	Chemokine (C-C motif) ligand 4
A10	Rn.8019	NM_031116	Ccl5	Chemokine (C-C motif) ligand 5
A11	Rn.34673	NM_020542	Ccr1	Chemokine (C-C motif) receptor 1
A12	Rn.49607	NM_001108836	Ccr10	Chemokine (C-C motif) receptor 10
B01	Rn.211983	NM_021866	Ccr2	Chemokine (C-C motif) receptor 2
B02	Rn.81076	NM_133532	Ccr4	Chemokine (C-C motif) receptor 4
B03	Rn.10736	NM_053960	Ccr5	Chemokine (C-C motif) receptor 5
B04	Rn.229736	NM_199489	Ccr7	Chemokine (C-C motif) receptor 7
B05	Rn.91140	NM_172329	Ccr9	Chemokine (C-C motif) receptor 9
B06	Rn.228198	NM_001191954	Cd274	CD274 molecule
B07	Rn.83632	NM_023981	Csf1	Colony stimulating factor 1 (macrophage)
B08	Rn.44285	NM_053852	Csf2	Colony stimulating factor 2 (granulocyte-macrophage)
B09	Rn.53973	NM_017104	Csf3	Colony stimulating factor 3 (granulocyte)
B10	Rn.10259	NM_031674	Ctla4	Cytotoxic T-lymphocyte-associated protein 4
B11	Rn.10907	NM_030845	Cxcl1	Chemokine (C-X-C motif) ligand 1 (melanoma growth stimulating activity, alpha)
B12	Rn.10584	NM_139089	Cxcl10	Chemokine (C-X-C motif) ligand 10
C01	Rn.13664	NM_182952	Cxcl11	Chemokine (C-X-C motif) ligand 11
C02	Rn.54439	NM_022177	Cxcl12	Chemokine (C-X-C motif) ligand 12 (stromal cell-derived factor 1)
C03	Rn.10230	NM_053647	Cxcl2	Chemokine (C-X-C motif) ligand 2
C04	Rn.44449	NM_022214	Cxcl5	Chemokine (C-X-C motif) ligand 5
C05	Rn.7391	NM_145672	Cxcl9	Chemokine (C-X-C motif) ligand 9
C06	Rn.138115	NM_019310	Cxcr1	Interleukin 8 receptor, alpha
C07	Rn.90347	NM_017183	Cxcr2	Chemokine (C-X-C motif) receptor 2
C08	Rn.24787	NM_053415	Cxcr3	Chemokine (C-X-C motif) receptor 3
C09	Rn.44431	NM_022205	Cxcr4	Chemokine (C-X-C motif) receptor 4
C10	Rn.127121	NM_053303	Cxcr5	Chemokine (C-X-C motif) receptor 5
C11	Rn.12959	NM_053352	Cxcr7	Chemokine (C-X-C motif) receptor 7
C12	Rn.6075	NM_012842	Egf	Epidermal growth factor
D01	Rn.37227	NM_031507	Egfr	Epidermal growth factor receptor
D02	Rn.9725	NM_012908	Faslg	Fas ligand (TNF superfamily, member 6)
D03	Rn.177272	NM_001108250	Foxp3	Forkhead box P3
D04	N/A	XM_006224278	Gbp1	Guanylate binding protein 1, interferon-inducible
D05	Rn.22283	NM_153468	Gzma	Granzyme A
D06	Rn.132115	NM_138517	Gzmb	Granzyme B
D07	Rn.10852	NM_024359	Hif1a	Hypoxia-inducible factor 1, alpha subunit (basic helix-loop-helix transcription factor)
D08	Rn.64513	NM_023973	Ido1	Indoleamine 2,3-dioxygenase 1
D09	Rn.10795	NM_138880	Igng	Interferon gamma
D10	Rn.201887	NM_178866	Igf1	Insulin-like growth factor 1
D11	Rn.9868	NM_012854	Il10	Interleukin 10
D12	Rn.207199	NM_053390	Il12a	Interleukin 12a
E01	Rn.48686	NM_022611	Il12b	Interleukin 12b
E02	Rn.9921	NM_053828	Il13	Interleukin 13
E03	Rn.2490	NM_013129	Il15	Interleukin 15
E04	Rn.218513	NM_001106897	Il17a	Interleukin 17A
E05	Rn.12300	NM_017019	Il1a	Interleukin 1 alpha
E06	Rn.9869	NM_031512	Il1b	Interleukin 1 beta
E07	Rn.9758	NM_013123	Il1r1	Interleukin 1 receptor, type 1
E08	Rn.9871	NM_053836	Il2	Interleukin 2

Position	UniGene	GenBank	Symbol	Description
E09	Rn.195879	NM_001191988	Il22	Interleukin 22
E10	Rn.81073	NM_130410	Il23a	Interleukin 23, alpha subunit p19
E11	Rn.108255	NM_201270	Il4	Interleukin 4
E12	Rn.9873	NM_012589	Il6	Interleukin 6
F01	Rn.6396	NM_012591	Irf1	Interferon regulatory factor 1
F02	Rn.44216	NM_021843	Kitlg	KIT ligand
F03	Rn.215725	NM_001017468	Micb	MHC class I polypeptide-related sequence B
F04	Rn.2661	NM_031051	Mif	Macrophage migration inhibitory factor
F05	Rn.12072	NM_012603	Myc	Myelocytomatosis oncogene
F06	Rn.37341	NM_198130	Myd88	Myeloid differentiation primary response gene 88
F07	Rn.2411	NM_001276711	Nfkb1	Nuclear factor of kappa light polypeptide gene enhancer in B-cells 1
F08	Rn.10400	NM_012611	Nos2	Nitric oxide synthase 2, inducible
F09	Rn.105023	NM_001106927	Pdcd1	Programmed cell death 1
F10	Rn.217585	NM_017232	Ptgs2	Prostaglandin-endoperoxide synthase 2
F11	Rn.233151	NM_001008827	RT1-A1	RT1 class Ia, locus A1
F12	Rn.234331	NM_001008854	RT1-N2	RT1 class Ib, locus N2
G01	Rn.8871	NM_012881	Spp1	Secreted phosphoprotein 1
G02	Rn.33229	NM_032612	Stat1	Signal transducer and activator of transcription 1
G03	Rn.10247	NM_012747	Stat3	Signal transducer and activator of transcription 3
G04	Rn.40136	NM_021578	Tgfb1	Transforming growth factor, beta 1
G05	Rn.46387	NM_198769	Tlr2	Toll-like receptor 2
G06	Rn.15273	NM_198791	Tlr3	Toll-like receptor 3
G07	Rn.14534	NM_019178	Tlr4	Toll-like receptor 4
G08	Rn.14819	NM_001097582	Tlr7	Toll-like receptor 7
G09	Rn.2275	NM_012675	Tnf	Tumor necrosis factor (TNF superfamily, member 2)
G10	Rn.83627	NM_145681	Tnfsf10	Tumor necrosis factor (ligand) superfamily, member 10
G11	Rn.54443	NM_030989	Tp53	Tumor protein p53
G12	Rn.1923	NM_031836	Vegfa	Vascular endothelial growth factor A
H01	Rn.94978	NM_031144	Actb	Actin, beta
H02	Rn.1868	NM_012512	B2m	Beta-2 microglobulin
H03	Rn.47	NM_012583	Hprt1	Hypoxanthine phosphoribosyltransferase 1
H04	Rn.107896	NM_017025	Ldha	Lactate dehydrogenase A
H05	Rn.973	NM_001007604	Rplp1	Ribosomal protein, large, P1
H06	N/A	U26919	RGDC	Rat Genomic DNA Contamination
H07	N/A	SA_00104	RTC	Reverse Transcription Control
H08	N/A	SA_00104	RTC	Reverse Transcription Control
H09	N/A	SA_00104	RTC	Reverse Transcription Control
H10	N/A	SA_00103	PPC	Positive PCR Control
H11	N/A	SA_00103	PPC	Positive PCR Control
H12	N/A	SA_00103	PPC	Positive PCR Control

Related products

For optimal performance, RT² Profiler PCR Arrays should be used together with the RT² First Strand Kit for cDNA synthesis and RT² SYBR[®] Green qPCR Mastermixes for PCR.

Product	Contents	Cat. no.
RT ² First Strand Kit (12)	Enzymes and reagents for cDNA synthesis	330401
RT ² SYBR Green ROX™ FAST Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the Rotor-Gene Q and other Rotor-Gene cyclers	330620

* Larger kit sizes available; please inquire.

RT² Profiler PCR Array products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN[®], Rotor-Gene[®], Rotor-Disc[™] (QIAGEN Group); ROX[™] (Applied Biosystems Corporation or its subsidiaries); SYBR[®] (Molecular Probes, Inc.).

1067688 03/2011 © 2011 QIAGEN, all rights reserved.

www.qiagen.com

Canada ■ 800-572-9613

China ■ 8621-3865-3865

Denmark ■ 80-885945

Finland ■ 0800-914416

France ■ 01-60-920-930

Germany ■ 02103-29-12000

Hong Kong ■ 800 933 965

Ireland ■ 1800 555 049

Italy ■ 800-787980

Japan ■ 03-6890-7300

Korea (South) ■ 080-000-7145

Luxembourg ■ 8002 2076

Mexico ■ 01-800-7742-436

The Netherlands ■ 0800 0229592

Norway ■ 800-18859

Singapore ■ 1800-742-4368

Spain ■ 91-630-7050

Sweden ■ 020-790282

Switzerland ■ 055-254-22-11

UK ■ 01293-422-911

USA ■ 800-426-8157

Australia ■ 1-800-243-800

Austria ■ 0800/281010

Belgium ■ 0800-79612

Brazil ■ 0800-557779



Sample & Assay Technologies