

RT² Profiler PCR Array (96-Well Format and 384-Well [4 x 96] Format)

Rat Dendritic & Antigen Presenting Cell

Cat. no. 330231 PARN-406ZA

For pathway expression analysis

Format	For use with the following real-time cyclers
RT ² Profiler PCR Array, Format A	Applied Biosystems [®] models 5700, 7000, 7300, 7500, 7700, 7900HT, ViiA™ 7 (96-well block); Bio-Rad [®] models iCycler [®] , iQ™ 5, MyiQ™, MyiQ2; Bio-Rad/MJ Research Chromo4™; Eppendorf [®] Mastercycler [®] ep realplex models 2, 2s, 4, 4s; Stratagene [®] models Mx3005P [®] , Mx3000P [®] ; Takara TP-800
RT ² Profiler PCR Array, Format C	Applied Biosystems models 7500 (Fast block), 7900HT (Fast block), StepOnePlus™, ViiA 7 (Fast block)
RT ² Profiler PCR Array, Format D	Bio-Rad CFX96™; Bio-Rad/MJ Research models DNA Engine Opticon [®] , DNA Engine Opticon 2; Stratagene Mx4000 [®]
RT ² Profiler PCR Array, Format E	Applied Biosystems models 7900HT (384-well block), ViiA 7 (384-well block); Bio-Rad CFX384™
RT ² Profiler PCR Array, Format F	Roche [®] LightCycler [®] 480 (96-well block)
RT ² Profiler PCR Array, Format G	Roche LightCycler 480 (384-well block)
RT ² Profiler PCR Array, Format H	Fluidigm [®] BioMark™



Sample & Assay Technologies

Description

The Rat Dendritic and Antigen Presenting Cell RT² Profiler PCR Array profiles the expression of 84 genes focused on dendritic cell activation and maturation. Genes important for dendritic cell activation and maturation such as cytokines, chemokines and their receptors are included on this array along with other related cell surface receptors and signal transduction molecules. Genes involved in antigen uptake, processing, and presentation are also represented on this array. In addition to being functionally defined, many of these genes on the array are highly expressed in mature dendritic cells or show significant changes in expression during cell differentiation. Using real-time PCR, you can easily and reliably analyze expression of a focused panel of genes related to dendritic and antigen presenting cells with this array.

For further details, consult the *RT² Profiler PCR Array Handbook*.

Shipping and storage

RT² Profiler PCR Arrays in formats A, C, D, E, F, and G are shipped at ambient temperature, on dry ice, or blue ice packs depending on destination and accompanying products. RT² Profiler PCR Arrays in format H are shipped on dry ice or blue ice packs.

For long term storage, keep plates at -20°C .

Note: Ensure that you have the correct RT² Profiler PCR Array format for your real-time cycler (see table above).

Note: Open the package and store the products appropriately immediately on receipt.

Array layout (96-well)

For 384-well 4 x 96 PCR arrays, genes are present in a staggered format. Refer to the *RT² Profiler PCR Array Handbook* for layout.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Ccl11	Ccl12	Ccl17	Ccl19	Ccl2	Ccl20	Cd3	Ccl4	Ccl5	Ccl7	Ccr1	Ccr2
B	Ccr3	Ccr5	Cd1d1	Cd2	Cd28	Cd36	Cd4	Cd40	Cd40lg	Cd44	Cd74	Cd80
C	Cd86	Cd8a	Cdc42	Cdkn1a	Cebpa	Clec4b2	Clec7a	Cmklr1	Csf1r	Csf2	Cxcl1	Cxcl10
D	Cxcl12	Cxcl2	Cxcr4	ErbB2	Fas	Fcer1a	Fcer2	Fcgr1a	Fcgrt	Flt3	Icam1	Icam2
E	Ifng	Il10	Il12a	Il12b	Il16	Il2	Il6	Il8ra	Irf7	Irf8	Ilgam	Ilgb2
F	Lrp1	Lyn	Mif	NfkB1	Ptprc	Rac1	Rag1	Rela	RT1-DMa	RT1-Ha	Stat3	Stat5a
G	Tap2	Tapbp	Tgfb1	Thbs1	Tlr1	Tlr2	Tlr7	Tlr9	Tnf	Tnfrsf11	Trpc4ap	Vcam1
H	Actb	B2m	Hprt1	Ldha	Rplp1	RGDC	RTC	RTC	RTC	PPC	PPC	PPC

Gene table: RT² Profiler PCR Array

Position	UniGene	GenBank	Symbol	Description
A01	Rn.10632	NM_019205	Ccl11	Chemokine (C-C motif) ligand 11
A02	Rn.137780	NM_001105822	Ccl12	Chemokine (C-C motif) ligand 12
A03	Rn.46679	NM_057151	Ccl17	Chemokine (C-C motif) ligand 17
A04	Rn.12445	NM_001108661	Ccl19	Chemokine (C-C motif) ligand 19
A05	Rn.4772	NM_031530	Ccl2	Chemokine (C-C motif) ligand 2
A06	Rn.10722	NM_019233	Ccl20	Chemokine (C-C motif) ligand 20
A07	Rn.10139	NM_013025	Ccl3	Chemokine (C-C motif) ligand 3
A08	Rn.37880	NM_053858	Ccl4	Chemokine (C-C motif) ligand 4
A09	Rn.8019	NM_031116	Ccl5	Chemokine (C-C motif) ligand 5
A10	Rn.26815	NM_001007612	Ccl7	Chemokine (C-C motif) ligand 7
A11	Rn.34673	NM_020542	Ccr1	Chemokine (C-C motif) receptor 1
A12	Rn.211983	NM_021866	Ccr2	Chemokine (C-C motif) receptor 2
B01	Rn.9890	NM_053958	Ccr3	Chemokine (C-C motif) receptor 3
B02	Rn.10736	NM_053960	Ccr5	Chemokine (C-C motif) receptor 5
B03	Rn.11120	NM_017079	Cd1d1	CD1d1 molecule
B04	Rn.10328	NM_012830	Cd2	Cd2 molecule
B05	Rn.10327	NM_013121	Cd28	Cd28 molecule
B06	Rn.102418	NM_031561	Cd36	CD36 molecule (thrombospondin receptor)
B07	Rn.10748	NM_012705	Cd4	Cd4 molecule
B08	Rn.25180	NM_134360	Cd40	CD40 molecule, TNF receptor superfamily member 5
B09	Rn.44218	NM_053353	Cd40lg	CD40 ligand
B10	Rn.1120	NM_012924	Cd44	Cd44 molecule
B11	Rn.33804	NM_013069	Cd74	Cd74 molecule, major histocompatibility complex, class II invariant chain
B12	Rn.10138	NM_012926	Cd80	Cd80 molecule
C01	Rn.6734	NM_020081	Cd86	CD86 molecule
C02	Rn.10306	NM_031538	Cd8a	CD8a molecule
C03	Rn.60067	NM_171994	Cdc42	Cell division cycle 42 (GTP binding protein)
C04	Rn.10089	NM_080782	Cdkn1a	Cyclin-dependent kinase inhibitor 1A
C05	Rn.204833	NM_012524	Cebpa	CCAAT/enhancer binding protein (C/EBP), alpha
C06	Rn.129075	NM_001005896	Clec4b2	C-type lectin domain family 4, member b2
C07	Rn.15743	XM_001067977	Clec7a	C-type lectin domain family 7, member a
C08	Rn.44465	NM_022218	Cmklr1	Chemokine-like receptor 1
C09	Rn.72599	NM_001029901	Csf1r	Colony stimulating factor 1 receptor
C10	Rn.44285	XM_340799	Csf2	Colony stimulating factor 2 (granulocyte-macrophage)
C11	Rn.10907	NM_030845	Cxcl1	Chemokine (C-X-C motif) ligand 1 (melanoma growth stimulating activity, alpha)
C12	Rn.10584	NM_139089	Cxcl10	Chemokine (C-X-C motif) ligand 10
D01	Rn.54439	NM_022177	Cxcl12	Chemokine (C-X-C motif) ligand 12 (stromal cell-derived factor 1)
D02	Rn.10230	NM_053647	Cxcl2	Chemokine (C-X-C motif) ligand 2
D03	Rn.44431	NM_022205	Cxcr4	Chemokine (C-X-C motif) receptor 4
D04	Rn.93966	NM_017003	ErbB2	V-erb-b2 erythroblastic leukemia viral oncogene homolog 2, neuro/glioblastoma derived oncogene homolog (avian)
D05	Rn.162521	NM_139194	Fas	Fas (TNF receptor superfamily, member 6)
D06	Rn.9677	NM_012724	Fcer1a	Fc fragment of IgE, high affinity I, receptor for; alpha polypeptide
D07	Rn.10326	NM_133550	Fcer2	Fc fragment of IgE, low affinity II, receptor for (CD23)
D08	Rn.16643	NM_001100836	Fcgr1a	Fc fragment of IgG, high affinity Ia, receptor (CD64)

Position	UniGene	GenBank	Symbol	Description
D09	Rn.15845	NM_033351	Fcgrt	Fc fragment of IgG, receptor, transporter, alpha
D10	Rn.6774	NM_001100822	FR3	Fms-related tyrosine kinase 3
D11	Rn.12	NM_012967	Icam1	Intercellular adhesion molecule 1
D12	Rn.162206	NM_001007725	Icam2	Intercellular adhesion molecule 2
E01	Rn.10795	NM_138880	Iifng	Interferon gamma
E02	Rn.9868	NM_012854	Il10	Interleukin 10
E03	Rn.207199	NM_053390	Il12a	Interleukin 12a
E04	Rn.48686	NM_022611	Il12b	Interleukin 12b
E05	Rn.104665	NM_001105749	Il16	Interleukin 16
E06	Rn.9871	NM_053836	Il2	Interleukin 2
E07	Rn.9873	NM_012589	Il6	Interleukin 6
E08	Rn.138115	NM_019310	Il8ra	Interleukin 8 receptor, alpha
E09	Rn.101159	NM_001033691	Irf7	Interferon regulatory factor 7
E10	Rn.3765	NM_001008722	Irf8	Interferon regulatory factor 8
E11	Rn.54465	NM_012711	Ilgam	Integrin, alpha M
E12	Rn.42962	NM_001037780	Ilgb2	Integrin, beta 2
F01	Rn.22436	XM_243524	Lrp1	Low density lipoprotein-related protein 1 (alpha-2-macroglobulin receptor)
F02	Rn.4338	NM_030857	Lyn	V-yes-1 Yamaguchi sarcoma viral related oncogene homolog
F03	Rn.2661	NM_031051	Mif	Macrophage migration inhibitory factor
F04	Rn.2411	XM_342346	Nfkb1	Nuclear factor of kappa light polypeptide gene enhancer in B-cells 1
F05	Rn.90166	NM_001109887	Ptprc	Protein tyrosine phosphatase, receptor type, C
F06	Rn.29157	NM_134366	Rac1	Ras-related C3 botulinum toxin substrate 1
F07	Rn.92344	NM_053468	Rag1	Recombination activating gene 1
F08	Rn.19480	NM_199267	Rela	V-rel reticuloendotheliosis viral oncogene homolog A (avian)
F09	Rn.35809	NM_198741	RT1-DMa	RT1 class II, locus DMA
F10	Rn.144557	NM_001008848	RT1-Ha	RT1 class II, locus Ha
F11	Rn.10247	NM_012747	Stat3	Signal transducer and activator of transcription 3
F12	Rn.154399	NM_017064	Stat5a	Signal transducer and activator of transcription 5A
G01	Rn.202962	NM_032056	Tap2	Transporter 2, ATP-binding cassette, sub-family B (MDR/TAP)
G02	Rn.51023	NM_033098	Tapbp	TAP binding protein
G03	Rn.40136	NM_021578	Tgfb1	Transforming growth factor, beta 1
G04	Rn.185771	NM_001013062	Thbs1	Thrombospondin 1
G05	Rn.107212	NM_001172120	Tlr1	Toll-like receptor 1
G06	Rn.46387	NM_198769	Tlr2	Toll-like receptor 2
G07	Rn.219862	NM_001097582	Tlr7	Toll-like receptor 7
G08	Rn.92495	NM_198131	Tlr9	Toll-like receptor 9
G09	Rn.2275	NM_012675	Tnf	Tumor necrosis factor (TNF superfamily, member 2)
G10	Rn.64517	NM_057149	Tnfsf11	Tumor necrosis factor (ligand) superfamily, member 11
G11	Rn.218020	NM_001100748	Trpc4ap	Transient receptor potential cation channel, subfamily C, member 4 associated protein
G12	Rn.11267	NM_012889	Vcam1	Vascular cell adhesion molecule 1
H01	Rn.94978	NM_031144	Actb	Actin, beta
H02	Rn.1868	NM_012512	B2m	Beta-2 microglobulin
H03	Rn.47	NM_012583	Hprt1	Hypoxanthine phosphoribosyltransferase 1
H04	Rn.107896	NM_017025	Ldha	Lactate dehydrogenase A
H05	Rn.973	NM_001007604	Rplp1	Ribosomal protein, large, P1
H06	N/A	U26919	RGDC	Rat Genomic DNA Contamination
H07	N/A	SA_00104	RTC	Reverse Transcription Control
H08	N/A	SA_00104	RTC	Reverse Transcription Control
H09	N/A	SA_00104	RTC	Reverse Transcription Control
H10	N/A	SA_00103	PPC	Positive PCR Control
H11	N/A	SA_00103	PPC	Positive PCR Control
H12	N/A	SA_00103	PPC	Positive PCR Control

Related products

For optimal performance, RT² Profiler PCR Arrays should be used together with the RT² First Strand Kit for cDNA synthesis and RT² SYBR[®] Green qPCR Mastermixes for PCR.

Product	Contents	Cat. no.
RT ² First Strand Kit (12)	Enzymes and reagents for cDNA synthesis	330401
RT ² SYBR Green qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with real-time cyclers that do not require a reference dye, including: Bio-Rad models CFX96, CFX384, DNA Engine Opticon 2; Bio-Rad/MJ Research Chromo4; Roche LightCycler 480 (96-well and 384-well); all other cyclers	330500
RT ² SYBR Green ROX™ qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the following real-time cyclers: Applied Biosystems models 5700, 7000, 7300, 7500 [Standard and FAST], 7700, 7900HT 96-well block [Standard and FAST] and 384-well block, StepOnePlus; Eppendorf Mastercycler ep realplex models 2, 2S, 4, 4S; Stratagene models Mx3000P, Mx3005P, Mx4000; Takara TP-800	330520
RT ² SYBR Green Fluor qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the following real-time cyclers: Bio-Rad models iCycler, iQ5, MyiQ, MyiQ2	330510

* Larger kit sizes available; please inquire.

RT² Profiler PCR Array products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN® (QIAGEN Group); Applied Biosystems®, ViiA™, StepOnePlus™, ROX™ (Applied Biosystems Corporation or its subsidiaries); Bio-Rad®, iCycler®, iQ™, MyiQ™, Chromo4™, CFX96™, DNA Engine Opticon®, CFX384™ (Bio-Rad Laboratories, Inc.); Stratagene®, Mx3005P®, Mx3000P®, Mx4000® (Stratagene); Eppendorf®, Mastercycler® (Eppendorf AG); Roche®, LightCycler® (Roche Group); Fluidigm® BioMark™ (Fluidigm Corporation); SYBR® (Molecular Probes, Inc.).

1066029 03/2011 © 2011 QIAGEN, all rights reserved.

www.qiagen.com

Canada ■ 800-572-9613

Ireland ■ 1800 555 049

Norway ■ 800-18859

China ■ 8621-3865-3865

Italy ■ 800-787980

Singapore ■ 1800-742-4368

Denmark ■ 80-885945

Japan ■ 03-6890-7300

Spain ■ 91-630-7050

Australia ■ 1-800-243-800

Finland ■ 0800-914416

Korea (South) ■ 080-000-7145

Sweden ■ 020-790282

Austria ■ 0800/281010

France ■ 01-60-920-930

Luxembourg ■ 8002 2076

Switzerland ■ 055-254-22-11

Belgium ■ 0800-79612

Germany ■ 02103-29-12000

Mexico ■ 01-800-7742-436

UK ■ 01293-422-911

Brazil ■ 0800-557779

Hong Kong ■ 800 933 965

The Netherlands ■ 0800 0229592

USA ■ 800-426-8157



Sample & Assay Technologies