

RT² Profiler PCR Array (96-Well Format and 384-Well [4 x 96] Format)

Rat Cancer Drug Targets

Cat. no. 330231 PARN-507ZA

For pathway expression analysis

Format	For use with the following real-time cyclers
RT ² Profiler PCR Array, Format A	Applied Biosystems® models 5700, 7000, 7300, 7500, 7700, 7900HT, ViiA™ 7 (96-well block); Bio-Rad® models iCycler®, iQ™ 5, MyiQ™, MyiQ2; Bio-Rad/MJ Research Chromo4™; Eppendorf® Mastercycler® ep realplex models 2, 2s, 4, 4s; Stratagene® models Mx3005P®, Mx3000P®; Takara TP-800
RT ² Profiler PCR Array, Format C	Applied Biosystems models 7500 (Fast block), 7900HT (Fast block), StepOnePlus™, ViiA 7 (Fast block)
RT ² Profiler PCR Array, Format D	Bio-Rad CFX96™; Bio-Rad/MJ Research models DNA Engine Opticon®, DNA Engine Opticon 2; Stratagene Mx4000®
RT ² Profiler PCR Array, Format E	Applied Biosystems models 7900HT (384-well block), ViiA 7 (384-well block); Bio-Rad CFX384™
RT ² Profiler PCR Array, Format F	Roche® LightCycler® 480 (96-well block)
RT ² Profiler PCR Array, Format G	Roche LightCycler 480 (384-well block)
RT ² Profiler PCR Array, Format H	Fluidigm® BioMark™



Sample & Assay Technologies

Description

The Rat Cancer Drug Targets RT² Profiler PCR Array profiles the expression of 84 actively sought targets for anticancer therapeutics and drug development. Cancer is a heterogeneous disease with a variety of survival mechanisms resulting from accumulated mutations that alter gene expression. In one of its key roles, cancer research continually identifies novel dysregulated carcinogenesis-related genes elucidating new mechanisms of cancer progression or treatment evasion, and potentially leading to new avenues for drug development. Further research into the expression of these genes may identify how and when they are dysregulated and potentially discover the underlying mechanism(s) behind cancer growth and progression. This array includes genes dysregulated during carcinogenesis, including those involved in key cellular growth pathways such as apoptosis, DNA damage repair, epigenetics, and growth factor and other signaling pathways. Using real-time PCR, your research study can easily and reliably analyze the expression of a focused panel of genes involved in oncogenesis with this array.

For further details, consult the *RT² Profiler PCR Array Handbook*.

Shipping and storage

RT² Profiler PCR Arrays in formats A, C, D, E, F, and G are shipped at ambient temperature, on dry ice, or blue ice packs depending on destination and accompanying products. RT² Profiler PCR Arrays in format H are shipped on dry ice or blue ice packs.

For long term storage, keep plates at –20°C.

Note: Ensure that you have the correct RT² Profiler PCR Array format for your real-time cycler (see table above).

Note: Open the package and store the products appropriately immediately on receipt.

Array layout (96-well)

For 384-well 4 x 96 PCR arrays, genes are present in a staggered format. Refer to the *RT² Profiler PCR Array Handbook* for layout.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Abcc1	Akt1	Akt2	Aif2	Aurka	Aurkb	Aurkc	Bcl2	Birc5	Cdc25a	Cdk1	Cdk2
B	Cdk4	Cdk5	Cdk7	Cdk9	Ctsb	Ctsd	Ctsl1	Ctss	Egfr	ErbB2	ErbB3	ErbB4
C	Esr1	Esr2	Figf	Flt1	Flt4	Grb2	Gstp1	Hdac1	Hdac11	Hdac2	Hdac3	Hdac4
D	Hdac6	Hdac7	Hdac8	Hif1a	Hras	Hsp90aa1	Hsp90b1	Igf1	Igf1r	Igf2	Ir5	Kdr
E	Kit	Kras	Mdm2	Mdm4	Mtor	Nfkb1	Nras	Ntn3	Parp1	Parp2	Pdgfra	Pdgfrb
F	Pgr	Pik3c2a	Pik3c3	Pik3ca	Plk1	Plk2	Plk3	Plk4	Prkca	Prkcb	Prkcd	Prkce
G	Pigs2	RGD1560888	Rhoa	Rhob	Tert	Tnks	Top2a	Top2b	Txn1	Txnrd1	Vegfb	Vegfc
H	Actb	B2m	Hprt1	Ldha	Rplp1	RGDC	RTC	RTC	RTC	PPC	PPC	PPC

Gene table: RT² Profiler PCR Array

Position	UniGene	GenBank	Symbol	Description
A01	Rn.10495	NM_022281	Abcc1	ATP-binding cassette, subfamily C (CFTR/MRP), member 1
A02	Rn.11422	NM_033230	Akt1	V-akt murine thymoma viral oncogene homolog 1
A03	Rn.87066	NM_017093	Akt2	V-akt murine thymoma viral oncogene homolog 2
A04	Rn.9825	NM_031018	Aif2	Activating transcription factor 2
A05	Rn.161874	NM_153296	Aurka	Aurora kinase A
A06	Rn.10865	NM_053749	Aurkb	Aurora kinase B
A07	Rn.41092	NM_001106221	Aurkc	Aurora kinase C
A08	Rn.9996	NM_016993	Bcl2	B-cell CLL/lymphoma 2
A09	Rn.54471	NM_022274	Birc5	Baculoviral IAP repeat-containing 5
A10	Rn.11390	NM_133571	Cdc25a	Cell division cycle 25 homolog A (S. pombe)
A11	Rn.6934	NM_019296	Cdk1	Cyclin-dependent kinase 1
A12	Rn.104460	NM_199501	Cdk2	Cyclin dependent kinase 2
B01	Rn.6115	NM_053593	Cdk4	Cyclin-dependent kinase 4
B02	Rn.10749	NM_080885	Cdk5	Cyclin-dependent kinase 5
B03	Rn.98896	XM_215467	Cdk7	Cyclin-dependent kinase 7
B04	Rn.98228	NM_001007743	Cdk9	Cyclin-dependent kinase 9
B05	Rn.100909	NM_022597	Ctsb	Cathepsin B
B06	Rn.11085	NM_134334	Ctsd	Cathepsin D
B07	Rn.1294	NM_013156	Ctsl1	Cathepsin L1
B08	Rn.11347	NM_017320	Ctss	Cathepsin S
B09	Rn.37227	NM_031507	Egfr	Epidermal growth factor receptor
B10	Rn.93966	NM_017003	ErbB2	V-erb-b2 erythroblastic leukemia viral oncogene homolog 2, neuro/glioblastoma derived oncogene homolog (avian)
B11	Rn.10228	NM_017218	ErbB3	V-erb-b2 erythroblastic leukemia viral oncogene homolog 3 (avian)
B12	Rn.163078	NM_021687	ErbB4	V-erb-a erythroblastic leukemia viral oncogene homolog 4 (avian)
C01	Rn.10595	NM_012689	Esr1	Estrogen receptor 1
C02	Rn.37460	NM_012754	Esr2	Estrogen receptor 2 (ER beta)
C03	Rn.10796	NM_031761	Figf	C-fos induced growth factor
C04	Rn.10239	NM_019306	Flt1	Fms-related tyrosine kinase 1
C05	Rn.81043	NM_053652	Flt4	Fms-related tyrosine kinase 4
C06	Rn.3360	NM_030846	Grb2	Growth factor receptor bound protein 2
C07	Rn.87063	NM_012577	Gstp1	Glutathione S-transferase pi 1
C08	Rn.1863	NM_001025409	Hdac1	Histone deacetylase 1
C09	Rn.28065	NM_001106610	Hdac11	Histone deacetylase 11
C10	Rn.1797	NM_053447	Hdac2	Histone deacetylase 2
C11	Rn.17284	NM_053448	Hdac3	Histone deacetylase 3
C12	Rn.23483	XM_343629	Hdac4	Histone deacetylase 4
D01	Rn.13453	XM_228753	Hdac6	Histone deacetylase 6
D02	Rn.203327	XM_345868	Hdac7	Histone deacetylase 7
D03	Rn.208476	NM_001126373	Hdac8	Histone deacetylase 8
D04	Rn.10852	NM_024359	Hif1a	Hypoxia-inducible factor 1, alpha subunit (basic helix-loop-helix transcription factor)
D05	Rn.102180	NM_001098241	Hras	Harvey rat sarcoma virus oncogene
D06	Rn.119867	NM_175761	Hsp90aa1	Heat shock protein 90, alpha (cytosolic), class A member 1
D07	Rn.101146	NM_001012197	Hsp90b1	Heat shock protein 90, beta, member 1

Position	UniGene	GenBank	Symbol	Description
D08	Rn.6282	NM_178866	Igf1	Insulin-like growth factor 1
D09	Rn.10957	NM_052807	Igf1r	Insulin-like growth factor 1 receptor
D10	Rn.118681	NM_031511	Igf2	Insulin-like growth factor 2
D11	Rn.203787	NM_001106586	Irf5	Interferon regulatory factor 5
D12	Rn.88869	NM_013062	Kdr	Kinase insert domain receptor
E01	Rn.54004	NM_022264	Kit	V-kit Hardy-Zuckerman 4 feline sarcoma viral oncogene homolog
E02	Rn.24554	NM_031515	Kras	V-Ki-ras2 Kirsten rat sarcoma viral oncogene homolog
E03	Rn.91829	NM_001108099	Mdm2	Mdm2 p53 binding protein homolog (mouse)
E04	Rn.90967	NM_001012026	Mdm4	Mdm4 p53 binding protein homolog (mouse)
E05	Rn.11008	NM_019906	Mtor	Mechanistic target of rapamycin (serine/threonine kinase)
E06	Rn.2411	XM_342346	Nfkb1	Nuclear factor of kappa light polypeptide gene enhancer in B-cells 1
E07	Rn.217722	NM_080766	Nras	Neuroblastoma ras oncogene
E08	Rn.95949	NM_053732	Ntn3	Netrin 3
E09	Rn.11327	NM_013063	Parp1	Poly (ADP-ribose) polymerase 1
E10	Rn.22730	NM_001106030	Parp2	Poly (ADP-ribose) polymerase 2
E11	Rn.55127	NM_012802	Pdgfra	Platelet derived growth factor receptor, alpha polypeptide
E12	Rn.98311	NM_031525	Pdgfrb	Platelet derived growth factor receptor, beta polypeptide
F01	Rn.10303	NM_022847	Pgr	Progesterone receptor
F02	Rn.60820	NM_001108500	Pik3c2a	Phosphoinositide-3-kinase, class 2, alpha polypeptide
F03	Rn.30010	NM_022958	Pik3c3	Phosphoinositide-3-kinase, class 3
F04	Rn.44193	NM_133399	Pik3ca	Phosphoinositide-3-kinase, catalytic, alpha polypeptide
F05	Rn.11034	NM_017100	Plk1	Polo-like kinase 1 (Drosophila)
F06	Rn.12100	NM_031821	Plk2	Polo-like kinase 2 (Drosophila)
F07	Rn.9284	NM_022187	Plk3	Polo-like kinase 3 (Drosophila)
F08	Rn.14848	NM_001107669	Plk4	Polo-like kinase 4 (Drosophila)
F09	Rn.207908	NM_001105713	Prkca	Protein kinase C, alpha
F10	Rn.91118	NM_012713	Prkcb	Protein kinase C, beta
F11	Rn.98279	NM_133307	Prkcd	Protein kinase C, delta
F12	Rn.216481	NM_017171	Prkce	Protein kinase C, epsilon
G01	Rn.44369	NM_017232	Ptgs2	Prostaglandin-endoperoxide synthase 2
G02	Rn.1331	NM_001109061	RGD1560888	Similar to Cell division protein kinase 8 (Protein kinase K35)
G03	Rn.107401	NM_057132	Rhoa	Ras homolog gene family, member A
G04	Rn.2042	NM_022542	Rhob	Ras homolog gene family, member B
G05	Rn.48802	NM_053423	Tert	Telomerase reverse transcriptase
G06	Rn.213840	NM_001106084	Tnks	Tankyrase, TRF1-interacting ankyrin-related ADP-ribose polymerase
G07	Rn.90996	NM_022183	Top2a	Topoisomerase (DNA) II alpha
G08	Rn.54443	NM_030989	Tp53	Tumor protein p53
G09	Rn.29777	NM_053800	Txn1	Thioredoxin 1
G10	Rn.67581	NM_031614	Txnrd1	Thioredoxin reductase 1
G11	Rn.198550	NM_053549	Vegfb	Vascular endothelial growth factor B
G12	Rn.6913	NM_053653	Vegfc	Vascular endothelial growth factor C
H01	Rn.94978	NM_031144	Actb	Actin, beta
H02	Rn.1868	NM_012512	B2m	Beta-2 microglobulin
H03	Rn.47	NM_012583	Hprt1	Hypoxanthine phosphoribosyltransferase 1
H04	Rn.107896	NM_017025	Ldha	Lactate dehydrogenase A
H05	Rn.973	NM_001007604	Rplp1	Ribosomal protein, large, P1
H06	N/A	U26919	RGDC	Rat Genomic DNA Contamination
H07	N/A	SA_00104	RTC	Reverse Transcription Control
H08	N/A	SA_00104	RTC	Reverse Transcription Control
H09	N/A	SA_00104	RTC	Reverse Transcription Control
H10	N/A	SA_00103	PPC	Positive PCR Control
H11	N/A	SA_00103	PPC	Positive PCR Control
H12	N/A	SA_00103	PPC	Positive PCR Control

Related products

For optimal performance, RT² Profiler PCR Arrays should be used together with the RT² First Strand Kit for cDNA synthesis and RT² SYBR[®] Green qPCR Mastermixes for PCR.

Product	Contents	Cat. no.
RT ² First Strand Kit (12)	Enzymes and reagents for cDNA synthesis	330401
RT ² SYBR Green qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with real-time cyclers that do not require a reference dye, including: Bio-Rad models CFX96, CFX384, DNA Engine Opticon 2; Bio-Rad/MJ Research Chromo4; Roche LightCycler 480 (96-well and 384-well); all other cyclers	330500
RT ² SYBR Green ROX [™] qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the following real-time cyclers: Applied Biosystems models 5700, 7000, 7300, 7500 [Standard and FAST], 7700, 7900HT 96-well block [Standard and FAST] and 384-well block, StepOnePlus; Eppendorf Mastercycler ep realplex models 2, 2S, 4, 4S; Stratagene models Mx3000P, Mx3005P, Mx4000; Takara TP-800	330520
RT ² SYBR Green Fluor qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the following real-time cyclers: Bio-Rad models iCycler, iQ5, MyiQ, MyiQ2	330510

* Larger kit sizes available; please inquire.

RT² Profiler PCR Array products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN® (QIAGEN Group); Applied Biosystems®, ViiA™, StepOnePlus™, ROX™ (Applied Biosystems Corporation or its subsidiaries); Bio-Rad®, iCycler®, iQ™, MyiQ™, Chromo4™, CFX96™, DNA Engine Opticon®, CFX384™ (Bio-Rad Laboratories, Inc.); Stratagene®, Mx3005P®, Mx3000P®, Mx4000® (Stratagene); Eppendorf®, Mastercycler® (Eppendorf AG); Roche®, LightCycler® (Roche Group); Fluidigm® BioMark™ (Fluidigm Corporation); SYBR® (Molecular Probes, Inc.).

1066029 03/2011 © 2011 QIAGEN, all rights reserved.

www.qiagen.com

Canada ■ 800-572-9613

China ■ 8621-3865-3865

Denmark ■ 80-885945

Finland ■ 0800-914416

France ■ 01-60-920-930

Germany ■ 02103-29-12000

Hong Kong ■ 800 933 965

Ireland ■ 1800 555 049

Italy ■ 800-787980

Japan ■ 03-6890-7300

Korea (South) ■ 080-000-7145

Luxembourg ■ 8002 2076

Mexico ■ 01-800-7742-436

The Netherlands ■ 0800 0229592

Norway ■ 800-18859

Singapore ■ 1800-742-4368

Spain ■ 91-630-7050

Sweden ■ 020-790282

Switzerland ■ 055-254-22-11

UK ■ 01293-422-911

USA ■ 800-426-8157

Australia ■ 1-800-243-800

Austria ■ 0800/281010

Belgium ■ 0800-79612

Brazil ■ 0800-557779



Sample & Assay Technologies