

RT² Profiler PCR Array (Rotor-Gene® Format)

Rat miR-210 Targets

Cat. no. 330231 PARN-6009ZR

For pathway expression analysis

Format	For use with the following real-time cyclers
RT ² Profiler PCR Array, Format R	Rotor-Gene Q, other Rotor-Gene cyclers

Description

The Rat miR-210 Targets RT² Profiler PCR Array profiles the expression of 84 rno-miR-210-3p target genes. This panel of 84 genes includes currently known experimentally verified plus bioinformatically predicted target genes regulated by rno-miR-210-3p. miRNA target gene expression analysis provides further insight into the function of these specific miRNAs. Using real-time PCR, research studies can easily and reliably analyze the expression of a focused panel of genes likely to be regulated by miR-210 with this array.

For further details, consult the *RT² Profiler PCR Array Handbook*.

Shipping and storage

RT² Profiler PCR Arrays in the Rotor-Gene format are shipped at ambient temperature, on dry ice, or blue ice packs depending on destination and accompanying products.

For long term storage, keep plates at –20°C.

Note: Ensure that you have the correct RT² Profiler PCR Array format for your real-time cycler (see table above).

Note: Open the package and store the products appropriately immediately on receipt.



Array layout

The 96 real-time assays in the Rotor-Gene format are located in wells 1–96 of the Rotor-Disc™ (plate A1–A12=Rotor-Disc 1–12, plate B1–B12=Rotor-Disc 13–24, etc.). To maintain data analysis compatibility, wells 97–100 do not contain real-time assays but will contain master mix to account for weight balance.

Gene table: RT² Profiler PCR Array

Position	UniGene	GenBank	Symbol	Description
A01	Rn.214018	NM_199230	Acvr1b	Activin A receptor, type IB
A02	Rn.107051	NM_198761	Adams5	ADAM metalloproteinase with thrombospondin type 1 motif, 5
A03	Rn.214145	NM_001107239	Adcy1	Adenylate cyclase 1 (brain)
A04	Rn.216683	NM_001106980	Adcy9	Adenylate cyclase 9
A05	Rn.138718	NM_001013977	Aifm3	Apoptosis-inducing factor, mitochondrion-associated 3
A06	Rn.10070	XM_214478	Aldh5a1	Aldehyde dehydrogenase 5 family, member A1
A07	Rn.42835	NM_001024980	Ankrd34a	Ankyrin repeat domain 34A
A08	Rn.29180	NM_001033703	Apobec3f	Apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 3F
A09	Rn.17507	NM_001106425	Armcl	Armadillo repeat containing 1
A10	Rn.162765	NM_001012097	Atg7	ATG7 autophagy related 7 homolog (S. cerevisiae)
A11	Rn.107246	NM_001108608	B4galt5	UDP-Gal:betaGlcNAc beta 1,4-galactosyltransferase, polypeptide 5
A12	Rn.11266	NM_012513	Bdnf	Brain-derived neurotrophic factor
B01	Rn.198715	NM_001107921	Casp8ap2	Caspase 8 associated protein 2
B02	Rn.10771	NM_078621	Ccbp2	Chemokine binding protein 2
B03	Rn.44186	NM_001106235	Ccdc97	Coiled-coil domain containing 97
B04	Rn.11166	NM_032083	Chn1	Chimerin (chimaerin) 1
B05	Rn.44611	NM_012528	Chrn1	Cholinergic receptor, nicotinic, beta 1 (muscle)
B06	Rn.8477	XM_217183	Coro2b	Coronin, actin binding protein, 2B
B07	Rn.162283	NM_001108361	Cpeb2	Cytoplasmic polyadenylation element binding protein 2
B08	Rn.7398	NM_022297	Ddah1	Dimethylarginine dimethylaminohydrolase 1
B09	Rn.25178	NM_001106408	Dimt1l	DIM1 dimethyladenosine transferase 1-like (S. cerevisiae)
B10	Rn.8092	NM_001106050	Dnajc15	DnaJ (Hsp40) homolog, subfamily C, member 15
B11	Rn.138110	NM_012768	Drd5	Dopamine receptor D5
B12	Rn.99931	NM_001107402	Dsc3	Desmocollin 3
C01	Rn.20041	NM_153303	Dyt1	Dystonia 1
C02	Rn.73967	NM_001137626	E2f3	E2F transcription factor 3
C03	Rn.8427	NM_053599	Efna1	Ephrin A1
C04	Rn.214210	XM_574979	Efna3	Ephrin A3
C05	Rn.198838	NM_001012744	Enpp5	Ectonucleotide pyrophosphatase/phosphodiesterase 5
C06	Rn.4177	NM_001134467	Fam116a	Family with sequence similarity 116, member A
C07	Rn.17259	NM_001106566	Fam73b	Family with sequence similarity 73, member B
C08	Rn.44953	NM_199114	Fgfr1l	Fibroblast growth factor receptor-like 1
C09	Rn.10470	NM_012742	Foxa1	Forkhead box A1
C10	Rn.163036	XM_343498	Gpd1l	Glycerol-3-phosphate dehydrogenase 1-like
C11	Rn.19721	NM_053906	Gsr	Glutathione reductase
C12	Rn.204425	XM_238042	Hhip	Hedgehog-interacting protein
D01	Rn.161814	NM_022528	Hif3a	Hypoxia inducible factor 3, alpha subunit
D02	Rn.9780	NM_013075	Hoxa1	Homeo box A1
D03	Rn.85723	NM_001109233	Hoxa9	Homeobox A9
D04	Rn.87132	NM_022938	Htr7	5-hydroxytryptamine (serotonin) receptor 7
D05	Rn.121	NM_001105936	Iscu	Iron-sulfur cluster scaffold homolog (E. coli)
D06	Rn.100904	XM_001069861	Kcmf1	Potassium channel modulatory factor 1
D07	Rn.27301	NM_001108831	Kctd11	Potassium channel tetramerisation domain containing 11
D08	Rn.61602	NM_001107281	Klf12	Kruppel-like factor 12
D09	Rn.149608	XM_001070337	Lair1	Leukocyte-associated immunoglobulin-like receptor 1
D10	Rn.40177	NM_012857	Lamp1	Lysosomal-associated membrane protein 1
D11	Rn.76594	XM_001079869	LOC691853	Similar to COX10 homolog, cytochrome c oxidase assembly protein, heme A: farnesyltransferase
D12	Rn.162435	NM_001014035	Mef2a	Myocyte enhancer factor 2a
E01	Rn.8267	NM_206950	Mid1ip1	MID1 interacting protein 1 (gastrulation specific G12 homolog (zebrafish))
E02	Rn.30696	NM_001105807	Mnt	Max binding protein
E03	Rn.11283	NM_031521	Ncam1	Neural cell adhesion molecule 1
E04	Rn.793	NM_001127684	Ndufa4	NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 4
E05	Rn.64485	NM_022302	Necab1	N-terminal EF-hand calcium binding protein 1
E06	Rn.10724	NM_019326	Neurod2	Neurogenic differentiation 2
E07	Rn.2411	XM_342346	Nfkb1	Nuclear factor of kappa light polypeptide gene enhancer in B-cells 1
E08	Rn.54707	NM_153735	Nptx1	Neuronal pentraxin 1

Position	UniGene	GenBank	Symbol	Description
E09	Rn.10055	NM_147210	Nr1d2	Nuclear receptor subfamily 1, group D, member 2
E10	Rn.7549	NM_001107141	Pptc7	PTC7 protein phosphatase homolog (S. cerevisiae)
E11	Rn.11317	NM_012637	Ptpn1	Protein tyrosine phosphatase, non-receptor type 1
E12	Rn.3788	NM_031091	Rab3b	RAB3B, member RAS oncogene family
F01	Rn.8154	NM_001106617	Rad52	RAD52 homolog (S. cerevisiae)
F02	Rn.24428	NM_175578	Rcan2	Regulator of calcineurin 2
F03	Rn.22326	NM_001107524	Rgma	RGM domain family, member A
F04	Rn.55126	NM_053865	Rtn1	Reticulon 1
F05	Rn.220523	NM_001130570	Scr1	Scratch homolog 1, zinc finger protein (Drosophila)
F06	Rn.32183	NM_001100882	Sema3g	Sema domain, immunoglobulin domain (Ig), short basic domain, secreted, (semaphorin) 3G
F07	Rn.23507	NM_001107002	Sept8	Septin 8
F08	Rn.96303	NM_001108638	Slc25a26	Solute carrier family 25, member 26
F09	Rn.203043	XM_001059906	Smg5	Smg-5 homolog, nonsense mediated mRNA decay factor (C. elegans)
F10	Rn.7582	NM_001047094	Spred2	Sprouty-related, EVH1 domain containing 2
F11	Rn.9908	NM_181092	Syngap1	Synaptic Ras GTPase activating protein 1 homolog (rat)
F12	Rn.105849	XM_001054844	Tcf7l2	Transcription factor 7-like 2 (T-cell specific, HMG-box)
G01	Rn.208326	XM_001055908	Thsd7a	Thrombospondin, type I, domain containing 7A
G02	Rn.214112	NM_001191804	Tmem151b	Transmembrane protein 151B
G03	Rn.138388	NM_001108461	Tppp	Tubulin polymerization promoting protein
G04	Rn.4272	NM_172032	Txndc9	Thioredoxin domain containing 9
G05	Rn.2758	NM_001008767	Txnip	Thioredoxin interacting protein
G06	Rn.19478	NM_001130582	Xpnpep3	X-prolyl aminopeptidase (aminopeptidase P) 3, putative
G07	Rn.12321	NM_001106129	Zadh2	Zinc binding alcohol dehydrogenase, domain containing 2
G08	Rn.136939	NM_001012175	Zfand3	Zinc finger, AN1-type domain 3
G09	Rn.64671	NM_031615	Zfp148	Zinc finger protein 148
G10	Rn.91545	NM_001107227	Zfyve28	Zinc finger, FYVE domain containing 28
G11	Rn.2925	NM_022548	Zmat3	Zinc finger, matrin type 3
G12	Rn.1712	NM_001108393	Zmiz1	Zinc finger, MIZ-type containing 1
H01	Rn.94978	NM_031144	Actb	Actin, beta
H02	Rn.1868	NM_012512	B2m	Beta-2 microglobulin
H03	Rn.47	NM_012583	Hprt1	Hypoxanthine phosphoribosyltransferase 1
H04	Rn.107896	NM_017025	Ldha	Lactate dehydrogenase A
H05	Rn.973	NM_001007604	Rplp1	Ribosomal protein, large, P1
H06	N/A	U26919	RGDC	Rat Genomic DNA Contamination
H07	N/A	SA_00104	RTC	Reverse Transcription Control
H08	N/A	SA_00104	RTC	Reverse Transcription Control
H09	N/A	SA_00104	RTC	Reverse Transcription Control
H10	N/A	SA_00103	PPC	Positive PCR Control
H11	N/A	SA_00103	PPC	Positive PCR Control
H12	N/A	SA_00103	PPC	Positive PCR Control

Related products

For optimal performance, RT² Profiler PCR Arrays should be used together with the RT² First Strand Kit for cDNA synthesis and RT² SYBR[®] Green qPCR Mastermixes for PCR.

Product	Contents	Cat. no.
RT ² First Strand Kit (12)	Enzymes and reagents for cDNA synthesis	330401
RT ² SYBR Green ROX [™] FAST Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the Rotor-Gene Q and other Rotor-Gene cyclers	330620

* Larger kit sizes available; please inquire.

RT² Profiler PCR Array products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN[®], Rotor-Gene[®], Rotor-Disc[™] (QIAGEN Group); ROX[™] (Applied Biosystems or its subsidiaries); SYBR[®] (Molecular Probes, Inc.).

1067688 03/2011 © 2011 QIAGEN, all rights reserved.

www.qiagen.com

Australia ■ 1-800-243-800

Austria ■ 0800/281010

Belgium ■ 0800-79612

Brazil ■ 0800-557779

Canada ■ 800-572-9613

China ■ 8621-3865-3865

Denmark ■ 80-885945

Finland ■ 0800-914416

France ■ 01-60-920-930

Germany ■ 02103-29-12000

Hong Kong ■ 800 933 965

Ireland ■ 1800 555 049

Italy ■ 800-787980

Japan ■ 03-6890-7300

Korea (South) ■ 080-000-7145

Luxembourg ■ 8002 2076

Mexico ■ 01-800-7742-436

The Netherlands ■ 0800 0229592

Norway ■ 800-18859

Singapore ■ 1800-742-4368

Spain ■ 91-630-7050

Sweden ■ 020-790282

Switzerland ■ 055-254-22-11

UK ■ 01293-422-911

USA ■ 800-426-8157



Sample & Assay Technologies