

RT² Profiler PCR Array (Rotor-Gene[®] Format)

Rat miR-29 Targets

Cat. no. 330231 PARN-6012ZR

For pathway expression analysis

Format	For use with the following real-time cyclers
RT ² Profiler PCR Array, Format R	Rotor-Gene Q, other Rotor-Gene cyclers

Description

The Rat miR-29 Targets RT² Profiler PCR Array profiles the expression of 84 rno-miR-29a-3p target genes. This panel of 84 genes includes currently known experimentally verified plus bioinformatically predicted target genes regulated by rno-miR-29a-3p. This array also includes target genes regulated by other miRNAs that have the same seed sequence as rno-miR-29a-3p, including rno-miR-29b-3p and rno-miR-29c-3p. miRNA target gene expression analysis provides further insight into the function of these specific miRNAs. A set of controls present on each array enables data analysis using the $\Delta\Delta CT$ method of relative quantification as well as assessment of reverse transcription performance, genomic DNA contamination, and PCR performance. Using real-time PCR, research studies can easily and reliably analyze the expression of a focused panel of genes likely to be regulated by miR-29 with this array.

For further details, consult the *RT² Profiler PCR Array Handbook*.

Shipping and storage

RT² Profiler PCR Arrays in the Rotor-Gene format are shipped at ambient temperature, on dry ice, or blue ice packs depending on destination and accompanying products.

For long term storage, keep plates at -20°C .

Note: Ensure that you have the correct RT² Profiler PCR Array format for your real-time cycler (see table above).

Note: Open the package and store the products appropriately immediately on receipt.



Array layout

The 96 real-time assays in the Rotor-Gene format are located in wells 1–96 of the Rotor-Disc™ (plate A1–A12=Rotor-Disc 1–12, plate B1–B12=Rotor-Disc 13–24, etc.). To maintain data analysis compatibility, wells 97–100 do not contain real-time assays but will contain master mix to account for weight balance.

Gene table: RT² Profiler PCR Array

Position	UniGene	GenBank	Symbol	Description
A01	Rn.161783	NM_031571	Acvr2a	Activin A receptor, type IIA
A02	Rn.161136	XM_001054670	Adam12	ADAM metalloproteinase domain 12
A03	Rn.162758	NM_001107877	Adams9	A disintegrin-like and metalloprotease (reprolysin type) with thrombospondin type 1 motif, 9
A04	Rn.60	NM_013218	Ak3	Adenylate kinase 3
A05	Rn.207201	NM_019204	Bace1	Beta-site APP cleaving enzyme 1
A06	Rn.14598	NM_053812	Bak1	BCL2-antagonist/killer 1
A07	Rn.25176	NM_173837	Bbc3	Bcl-2 binding component 3
A08	Rn.9996	NM_016993	Bcl2	B-cell CLL/lymphoma 2
A09	Rn.82709	NM_022612	Bcl2l1	BCL2-like 11 (apoptosis facilitator)
A10	Rn.213264	NM_139258	Bmf	Bcl2 modifying factor
A11	Rn.76129	NM_182824	Cd276	Cd276 molecule
A12	Rn.60067	NM_171994	Cdc42	Cell division cycle 42 (GTP binding protein)
B01	Rn.162731	XM_342638	Cdk6	Cyclin-dependent kinase 6
B02	Rn.31832	XM_216399	Col15a1	Collagen, type XV, alpha 1
B03	Rn.2953	NM_053304	Col1a1	Collagen, type I, alpha 1
B04	Rn.107239	NM_053356	Col1a2	Collagen, type I, alpha 2
B05	Rn.10124	NM_012929	Col2a1	Collagen, type II, alpha 1
B06	Rn.3247	NM_032085	Col3a1	Collagen, type III, alpha 1
B07	Rn.53801	NM_001135009	Col4a1	Collagen, type IV, alpha 1
B08	Rn.2237	XM_225043	Col4a2	Collagen, type IV, alpha 2
B09	Rn.2875	NM_053488	Col5a2	Collagen, type V, alpha 2
B10	Rn.38654	NM_021760	Col5a3	Collagen, type V, alpha 3
B11	Rn.24307	NM_001106858	Col7a1	Procollagen, type VII, alpha 1
B12	Rn.36076	XM_001075213	Ctnnbip1	Catenin, beta-interacting protein 1
C01	Rn.17832	XM_347258	Dgkd	Diacylglycerol kinase, delta
C02	Rn.205881	XM_001068155	Dicer1	Dicer 1, ribonuclease type III
C03	Rn.14603	NM_001015021	Dnajb11	DnaJ (Hsp40) homolog, subfamily B, member 11
C04	Rn.6955	NM_053354	Dnmt1	DNA (cytosine-5-)-methyltransferase 1
C05	Rn.92659	NM_001003958	Dnmt3a	DNA (cytosine-5-)-methyltransferase 3 alpha
C06	Rn.117353	NM_001003959	Dnmt3b	DNA (cytosine-5-)-methyltransferase 3 beta
C07	Rn.136933	NM_001012089	Dusp2	Dual specificity phosphatase 2
C08	Rn.107459	NM_001012181	Elf2	E74-like factor 2
C09	Rn.54384	NM_012722	Eln	Elastin
C10	Rn.208436	XM_001061749	Eomes	Eomesodermin homolog (Xenopus laevis)
C11	Rn.12759	NM_031825	Fbn1	Fibrillin 1
C12	Rn.219320	NM_001108157	Fem1b	Fem-1 homolog b (C. elegans)
D01	Rn.98846	NM_001008724	Fga	Fibrinogen alpha chain
D02	Rn.11416	NM_020071	Fgb	Fibrinogen beta chain
D03	Rn.1702	NM_012559	Fgg	Fibrinogen gamma chain
D04	Rn.63642	XM_001060923	Foxj2	Forkhead box J2
D05	Rn.2204	NM_017073	Glul	Glutamate-ammonia ligase (glutamine synthetase)
D06	Rn.5820	NM_017113	Grn	Granulin
D07	Rn.23483	XM_343629	Hdac4	Histone deacetylase 4
D08	Rn.89639	NM_057130	Hrk	Harakiri, BCL2 interacting protein (contains only BH3 domain)
D09	Rn.3370	NM_001030026	Ifi30	Interferon gamma inducible protein 30
D10	Rn.10132	NM_022863	Ireb2	Iron responsive element binding protein 2
D11	Rn.64478	NM_001108156	Itga11	Integrin, alpha 11
D12	Rn.7145	NM_053966	Lamc1	Laminin, gamma 1
E01	Rn.3834	NM_012598	Lpl	Lipoprotein lipase
E02	Rn.23232	NM_130749	Mark3	MAP/microtubule affinity-regulating kinase 3
E03	Rn.129914	NM_021846	Mcl1	Myeloid cell leukemia sequence 1
E04	Rn.59431	NM_001009617	Mest	Mesoderm specific transcript homolog (mouse)
E05	Rn.165433	NM_001106168	Mmp15	Matrix metalloproteinase 15
E06	Rn.3117	NM_031757	Mmp24	Matrix metalloproteinase 24
E07	Rn.81116	NM_001013096	Mycn	V-myc myelocytomatosis viral related oncogene, neuroblastoma derived (avian)
E08	Rn.91711	NM_001191782	Nav3	Neuron navigator 3

Position	UniGene	GenBank	Symbol	Description
E09	Rn.105658	XM_213954	Nid1	Nidogen 1
E10	Rn.221910	NM_178096	Nrep	Neuronal regeneration related protein
E11	Rn.22308	NM_053940	Pcdha12	Protocadherin alpha 12
E12	Rn.10599	NM_013005	Pik3r1	Phosphoinositide-3-kinase, regulatory subunit 1 (alpha)
F01	Rn.1476	NM_017037	Pmp22	Peripheral myelin protein 22
F02	Rn.15540	NM_001105825	Ppm1d	Protein phosphatase 1D magnesium-dependent, delta isoform
F03	Rn.30183	NM_001108062	Ppp1r13b	Protein phosphatase 1, regulatory (inhibitor) subunit 13B
F04	Rn.22158	NM_031606	Pten	Phosphatase and tensin homolog
F05	Rn.40541	XM_001053103	Pxdn	Peroxidasin homolog (Drosophila)
F06	Rn.161854	NM_001015012	Rab30	RAB30, member RAS oncogene family
F07	Rn.50631	XM_233485	Rlf	Rearranged L-myc fusion
F08	Rn.8937	NM_013191	S100b	S100 calcium binding protein B
F09	Rn.95177	NM_001007732	Serpinb9	Serpin peptidase inhibitor, clade B (ovalbumin), member 9
F10	Rn.18006	NM_001134514	Sestd1	SEC14 and spectrin domains 1
F11	Rn.54645	NM_001025271	Sfpq	Splicing factor proline/glutamine rich (polypyrimidine tract binding protein associated)
F12	Rn.44609	NM_012655	Sp1	Sp1 transcription factor
G01	Rn.98989	NM_012656	Sparc	Secreted protein, acidic, cysteine-rich (osteonectin)
G02	Rn.221903	NM_001106427	Spry1	Sprouty homolog 1, antagonist of FGF signaling (Drosophila)
G03	Rn.162118	NM_001025738	Srsf10	Serine/arginine-rich splicing factor 10
G04	Rn.144930	NM_001107043	Tbx21	T-box 21
G05	Rn.218634	NM_001109601	Tcl1a	T-cell leukemia/lymphoma 1A
G06	Rn.98685	NM_053729	Tdg	Thymine-DNA glycosylase
G07	Rn.219960	NM_001107643	Tet1	Tet oncogene 1
G08	Rn.16142	NM_201420	Tfap2c	Transcription factor AP-2 gamma
G09	Rn.7018	NM_013174	Tgfb3	Transforming growth factor, beta 3
G10	Rn.1923	NM_031836	Vegfa	Vascular endothelial growth factor A
G11	Rn.82737	NM_133290	Zfp36	Zinc finger protein 36
G12	Rn.6142	NM_017172	Zfp361l	Zinc finger protein 36, C3H type-like 1
H01	Rn.94978	NM_031144	Actb	Actin, beta
H02	Rn.1868	NM_012512	B2m	Beta-2 microglobulin
H03	Rn.47	NM_012583	Hprt1	Hypoxanthine phosphoribosyltransferase 1
H04	Rn.107896	NM_017025	Ldha	Lactate dehydrogenase A
H05	Rn.973	NM_001007604	Rplp1	Ribosomal protein, large, P1
H06	N/A	U26919	RGDC	Rat Genomic DNA Contamination
H07	N/A	SA_00104	RTC	Reverse Transcription Control
H08	N/A	SA_00104	RTC	Reverse Transcription Control
H09	N/A	SA_00104	RTC	Reverse Transcription Control
H10	N/A	SA_00103	PPC	Positive PCR Control
H11	N/A	SA_00103	PPC	Positive PCR Control
H12	N/A	SA_00103	PPC	Positive PCR Control

Related products

For optimal performance, RT² Profiler PCR Arrays should be used together with the RT² First Strand Kit for cDNA synthesis and RT² SYBR[®] Green qPCR Mastermixes for PCR.

Product	Contents	Cat. no.
RT ² First Strand Kit (12)	Enzymes and reagents for cDNA synthesis	330401
RT ² SYBR Green ROX [™] FAST Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the Rotor-Gene Q and other Rotor-Gene cyclers	330620

* Larger kit sizes available; please inquire.

RT² Profiler PCR Array products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN[®], Rotor-Gene[®], Rotor-Disc[™] (QIAGEN Group); ROX[™] (Applied Biosystems or its subsidiaries); SYBR[®] (Molecular Probes, Inc.).

1067688 03/2011 © 2011 QIAGEN, all rights reserved.

www.qiagen.com

Australia ■ 1-800-243-800

Austria ■ 0800/281010

Belgium ■ 0800-79612

Brazil ■ 0800-557779

Canada ■ 800-572-9613

China ■ 8621-3865-3865

Denmark ■ 80-885945

Finland ■ 0800-914416

France ■ 01-60-920-930

Germany ■ 02103-29-12000

Hong Kong ■ 800 933 965

Ireland ■ 1800 555 049

Italy ■ 800-787980

Japan ■ 03-6890-7300

Korea (South) ■ 080-000-7145

Luxembourg ■ 8002 2076

Mexico ■ 01-800-7742-436

The Netherlands ■ 0800 0229592

Norway ■ 800-18859

Singapore ■ 1800-742-4368

Spain ■ 91-630-7050

Sweden ■ 020-790282

Switzerland ■ 055-254-22-11

UK ■ 01293-422-911

USA ■ 800-426-8157



Sample & Assay Technologies