

RT² Profiler PCR Array (Rotor-Gene® Format)

Pig Drug Metabolism: Phase I Enzymes

Cat. no. 330231 PASS-068ZR

For pathway expression analysis

Format	For use with the following real-time cyclers
RT ² Profiler PCR Array, Format R	Rotor-Gene Q, other Rotor-Gene cyclers

Description

The Pig Drug Metabolism: Phase I Enzymes RT² Profiler™ PCR Array contains 84 genes involved in phase I drug metabolism. Phase I drug metabolism enzymes make compounds more hydrophilic and add functional groups necessary for the completion of Phase II drug metabolism. This array represents genes involved in Phase I drug metabolism reactions including oxidation, reduction, hydrolysis, cyclization, and decyclization. Members of the Cytochrome P450 enzyme family that play a key role in mediating phase I drug metabolism reactions are also included on this array. Using real-time PCR, you can easily and reliably analyze expression of a focused panel of genes related to drug phase I metabolism with this array.

For further details, consult the *RT² Profiler PCR Array Handbook*.



Shipping and storage

RT² Profiler PCR Arrays in the Rotor-Gene format are shipped at ambient temperature, on dry ice, or blue ice packs depending on destination and accompanying products.

For long term storage, keep plates at –20°C.

Note: Ensure that you have the correct RT² Profiler PCR Array format for your real-time cyclers (see table above).

Note: Open the package and store the products appropriately immediately on receipt.

Array layout

The 96 real-time assays in the Rotor-Gene format are located in wells 1–96 of the Rotor-Disc™ (plate A1–A12=Rotor-Disc 1–12, plate B1–B12=Rotor-Disc 13–24, etc.). To maintain data analysis compatibility, wells 97–100 do not contain real-time assays but will contain master mix to account for weight balance.

Gene table: RT² Profiler PCR Array

Position	UniGene	GenBank	Symbol	Description
1	Ssc.14069	NM_001243939	ADH1C	Alcohol dehydrogenase 1C (class I), gamma polypeptide
2	N/A	NM_001315670	ADH4	Alcohol dehydrogenase 4
3	Ssc.4764	NM_001244833	ADH5	Alcohol dehydrogenase 5 (class III), chi polypeptide
4	N/A	XM_013998428	ADH7	Alcohol dehydrogenase class 4 mu/sigma chain-like
5	Ssc.6928	XM_013993622	ALDH1A1	Retinal dehydrogenase 1-like
6	Ssc.87537	XM_005659561	ALDH1A2	Aldehyde dehydrogenase 1 family member A2
7	N/A	XM_013993236	ALDH1A3	Aldehyde dehydrogenase 1 family member A3
8	Ssc.84202	XM_003353586	ALDH1B1	Aldehyde dehydrogenase 1 family, member B1
9	Ssc.11147	NM_001044611	ALDH2	Aldehyde dehydrogenase 2 family (mitochondrial)
10	N/A	XM_013981361	ALDH3A1	Aldehyde dehydrogenase 3 family, member A1
11	Ssc.5884	XM_003132032	ALDH3A2	Fatty aldehyde dehydrogenase-like
12	Ssc.79550	NM_001315593	ALDH3B1	Aldehyde dehydrogenase 3 family member B1
13	Ssc.54494	NM_001252213	ALDH4A1	Aldehyde dehydrogenase 4 family member A1
14	Ssc.27232	NM_001244467	ALDH5A1	Aldehyde dehydrogenase 5 family, member A1
15	Ssc.32192	XM_005656399	ALDH6A1	Aldehyde dehydrogenase 6 family, member A1
16	Ssc.11751	XM_013995069	ALDH7A1	Aldehyde dehydrogenase 7 family member A1
17	N/A	XM_013988074	ALDH8A1	Aldehyde dehydrogenase 8 family member A1
18	Ssc.3723	XM_005663148	ALDH9A1	Aldehyde dehydrogenase 9 family, member A1
19	N/A	XM_003353694	CEL	Carboxyl ester lipase (bile salt-stimulated lipase)
20	Ssc.8776	NM_214427	CYP11A1	Cytochrome P450, family 11, subfamily A, polypeptide 1
21	Ssc.51528	NM_214428	CYP17A1	Cytochrome P450 17A1
22	Ssc.56422	NM_214429	CYP19A1	Cytochrome P450 19A1
23	Ssc.82444	NM_214412	CYP1A1	Cytochrome P450 1A1
24	Ssc.84322	NM_001159614	CYP1A2	Cytochrome P450, family 1, subfamily A, polypeptide 2
25	Ssc.50259	XM_013992312	CYP20A1	Cytochrome P450, family 20, subfamily A, polypeptide 1
26	Ssc.52818	NM_214433	CYP21A2	Cytochrome P450, family 21, subfamily A, polypeptide 2
27	Ssc.268	NM_214075	CYP24A1	Cytochrome P450, family 24, subfamily A, polypeptide 1
28	Ssc.97040	XM_003483530	CYP26A1	Cytochrome P450, family 26, subfamily A, polypeptide 1
29	Ssc.3804	NM_001243304	CYP27A1	Cytochrome P450, family 27, subfamily A, polypeptide 1
30	Ssc.16	NM_213995	CYP27B1	25-hydroxyvitamin D3 1alpha-hydroxylase
31	Ssc.15742	NM_214417	CYP2A19	Cytochrome P450 2A19
32	Ssc.11267	NM_214413	CYP2B22	Cytochrome P450 2B22
33	Ssc.4617	NM_214394	CYP2D25	Vitamin D3 25-Hydroxylase
34	Ssc.15741	NM_214421	CYP2E1	Cytochrome P450, family 2, subfamily E, polypeptide 1
35	Ssc.80416	NM_001244633	CYP2J34	Cytochrome P450, family 2, subfamily J, polypeptide 34
36	Ssc.15486	XM_003480731	CYP2R1	Cytochrome P450, family 2, subfamily R, polypeptide 1
37	Ssc.81969	NM_001195509	CYP3A22	Cytochrome P450, subfamily IIIA, polypeptide 22
38	Ssc.204	NM_214423	CYP3A29	Cytochrome P450 3A29

Position	UniGene	GenBank	Symbol	Description
39	Ssc.203	NM_214422	CYP3A39	Cytochrome P450 3A39
40	Ssc.203	NM_001134824	CYP3A46	Cytochrome P450 3A46
41	Ssc.15895	NM_214425	CYP4A21	Cytochrome P450, family 4, subfamily A, polypeptide 21
42	Ssc.51897	XM_003123433	CYP4F2	Cytochrome P450, family 4, subfamily F, polypeptide 2
43	Ssc.51897	NM_001244636	CYP4F55	Cytochrome P450, family 4, subfamily F, polypeptide 55
44	Ssc.28948	XM_013984059	CYP4V2	Cytochrome P450, family 4, subfamily v, polypeptide 2
45	Ssc.30521	NM_001005352	CYP7A1	Cytochrome P450, family 7, subfamily A, polypeptide 1
46	Ssc.75265	XM_013996723	CYP7B1	Cytochrome P450, family 7, subfamily B, polypeptide 1
47	Ssc.42789	NM_214426	CYP8B1	Cytochrome P-450 8B1
48	Ssc.153	NM_214044	DPYD	Dihydropyrimidine dehydrogenase
49	Ssc.217	NM_214060	ESD	Esterase D
50	Ssc.229	NM_214064	FMO1	Flavin containing monooxygenase 1
51	N/A	XM_003130104	FMO2	Dimethylaniline monooxygenase [N-oxide-forming] 2-like
52	Ssc.45970	XM_003130094	FMO3	Dimethylaniline monooxygenase [N-oxide-forming] 3-like
53	Ssc.98212	XM_005656647	FMO4	Dimethylaniline monooxygenase [N-oxide-forming] 4
54	Ssc.96266	XM_001928594	FMO5	Flavin containing monooxygenase 5
55	Ssc.53099	NM_001143709	GZMA	Granzyme A (granzyme 1, cytotoxic T-lymphocyte-associated serine esterase 3)
56	Ssc.4939	NM_001143710	GZMB	Granzyme B (granzyme 2, cytotoxic T-lymphocyte-associated serine esterase 1)
57	Ssc.1301	XM_003135126	HSD17B10	3-hydroxyacyl-CoA dehydrogenase type-2-like
58	N/A	XM_003125033	LOC100521659	Cytochrome P450 26B1-like
59	N/A	XM_013987674	LOC100522669	Cytochrome P450 2W1-like
60	N/A	XM_005655855	LOC100523190	Cytochrome P450 2F1-like
61	Ssc.93967	XM_003128017	LOC100523909	Cytochrome P450 4B1-like
62	N/A	XM_003354125	LOC100620265	Cytochrome P450 4F22-like
63	N/A	XM_003361542	LOC100623616	Arylacetamide deacetylase-like
64	N/A	XM_013983577	LOC100624257	Cytochrome P450 2C42-like
65	Ssc.44467	XM_013991464	LOC100625479	Cytochrome P450 1B1-like
66	N/A	XM_013998513	LOC100737508	Cytochrome P450 2G1
67	N/A	XM_003483531	LOC100739552	Cytochrome P450 26C1
68	N/A	XM_013998511	LOC102165015	Cytochrome P450 2B11-like
69	N/A	XM_013998507	LOC106510546	Cytochrome P450 2S1
70	N/A	XM_021074021	LOC110256857	Cytochrome P450 2C9-like
71	N/A	XM_021083598	LOC110259328	Phylloquinone omega-hydroxylase CYP4F2-like
72	N/A	XM_021083606	LOC110259329	Docosahexaenoic acid omega-hydroxylase CYP4F3-like
73	N/A	XM_021086162	LOC110259856	Cytochrome P450 3A29-like
74	N/A	XM_021088772	LOC110260194	Cytochrome P450 11B1, mitochondrial
75	N/A	XM_021097627	LOC110261409	Cytochrome P450 4A25-like
76	N/A	XM_021100552	LOC110261964	Alcohol dehydrogenase 6-like
77	N/A	XM_021101730	LOC110262184	Alcohol dehydrogenase 1-like
78	Ssc.7297	NM_001001640	MAOA	Monoamine oxidase A

Position	UniGene	GenBank	Symbol	Description
79	Ssc.7297	NM_001001864	MAOB	Monoamine oxidase B
80	Ssc.1986	XM_001926129	PTGS1	Prostaglandin-endoperoxide synthase 1 (prostaglandin G/H synthase and cyclooxygenase)
81	Ssc.23994	NM_214321	PTGS2	Prostaglandin G/H synthase-2
82	Ssc.2073	NM_213763	UCHL1	Ubiquitin carboxyl-terminal esterase L1 (ubiquitin thiolesterase)
83	Ssc.13668	NM_001077227	UCHL3	Ubiquitin carboxyl-terminal esterase L3 (ubiquitin thiolesterase)
84	Ssc.30588	NM_001285974	XDH	Xanthine dehydrogenase
85	Ssc.10316	XM_003124280	ACTB	Actin, beta
86	Ssc.73773	NM_213978	B2M	Beta-2-microglobulin
87	Ssc.16135	NM_001206359	GAPDH	Glyceraldehyde-3-phosphate dehydrogenase
88	Ssc.4158	NM_001032376	HPRT1	Hypoxanthine phosphoribosyltransferase 1
89	Ssc.27927	NM_001244068	RPL13A	Ribosomal protein L13a
90	N/A	SA_00133	SGDC	Pig Genomic DNA Contamination
91	N/A	SA_00104	RTC	Reverse Transcription Control
92	N/A	SA_00104	RTC	Reverse Transcription Control
93	N/A	SA_00104	RTC	Reverse Transcription Control
94	N/A	SA_00103	PPC	Positive PCR Control
95	N/A	SA_00103	PPC	Positive PCR Control
96	N/A	SA_00103	PPC	Positive PCR Control

Related products

For optimal performance, RT² Profiler PCR Arrays should be used together with the RT² First Strand Kit for cDNA synthesis and RT² SYBR[®] Green qPCR Mastermixes for PCR.

Product	Contents	Cat. no.
RT ² First Strand Kit (12)	Enzymes and reagents for cDNA synthesis	330401
RT ² SYBR Green ROX [™] FAST Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the Rotor-Gene Q and other Rotor-Gene cyclers	330620

* Larger kit sizes available; please inquire.

RT² Profiler PCR Array products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN[®], Rotor-Gene[®], Rotor-Disc[™] (QIAGEN Group); ROX[™] (Applied Biosystems Corporation or its subsidiaries); SYBR[®] (Molecular Probes, Inc.).

1067688 03/2011 © 2011 QIAGEN, all rights reserved.

www.qiagen.com

Canada ■ 800-572-9613

Ireland ■ 1800 555 049

Norway ■ 800-18859

China ■ 8621-3865-3865

Italy ■ 800-787980

Singapore ■ 1800-742-4368

Denmark ■ 80-885945

Japan ■ 03-6890-7300

Spain ■ 91-630-7050

Australia ■ 1-800-243-800

Finland ■ 0800-914416

Korea (South) ■ 080-000-7145

Sweden ■ 020-790282

Austria ■ 0800/281010

France ■ 01-60-920-930

Luxembourg ■ 8002 2076

Switzerland ■ 055-254-22-11

Belgium ■ 0800-79612

Germany ■ 02103-29-12000

Mexico ■ 01-800-7742-436

UK ■ 01293-422-911

Brazil ■ 0800-557779

Hong Kong ■ 800 933 965

The Netherlands ■ 0800 0229592

USA ■ 800-426-8157



Sample & Assay Technologies