

QIAseq Targeted DNA Pro Panels

Human Lung Cancer Research

Cat. nos. 333651 PHS-005Z-12 and 333655 PHS-005Z-96

For ultrasensitive targeted next-generation sequencing (NGS) of DNA

The QIAseq Targeted DNA Pro Human Lung Cancer Research Panel permits the targeted next generation sequencing analysis of the whole coding region of 61 oncogenes and tumor suppressor genes relevant to lung cancer, covering 207,663 bases with 2,571 primers. The panel also targets actionable and interpretable cancer-related somatic variants relevant for tumor classification of adenocarcinomas, non-small cell lung cancers (NSCLC), small cell lung cancers (SCLC), squamous cell carcinomas (SCC), and other lung cancers. Targeted actionable and interpretable genes and variants were collected from 1) clinical practice guidelines recommended by multiple organizations, including the American College of Medical Genetics (ACMG), the American Society Of Hematology (ASH), the Association for Molecular Pathology (AMP), the American Society of Clinical Oncology (ASCO), the College of American Pathologists (CAP), the European Society for Medical Oncology (ESMO), and the National Comprehensive Cancer Network (NCCN); 2) the NIH NCBI ClinVar and Online Mendelian Inheritance in Man (OMIM) resources that compile variants with known associations to cancers cited by peer-reviewed publications, 3) clinical trial documents from resources such as Cancer.gov, ClinicalTrials.gov, and the UMIN Clinical Trials Registry; 4) requests and suggestions from QIAGEN KOLs and SMEs; and 5) other peer-reviewed publications. As a result, this primer pool gives researchers a comprehensive insight into lung cancer genomics with the detection of SNP/indels in a focused panel of relevant genes using digital DNA sequencing.

For further details, consult the *QIAseq Targeted DNA Pro Panels Handbook*.

Shipping and storage

The QIAseq Targeted DNA Pro Kits (except QIAseq Beads and QIAseq Bead Binding Buffer) are shipped on dry ice and should be stored at -30°C to -15°C in a constant-temperature freezer upon arrival. The QIAseq Beads and QIAseq Bead Binding Buffer are shipped on cold packs and should be stored at $2-8^{\circ}\text{C}$ upon arrival. The QIAseq Targeted DNA Pro Index Kits are shipped on dry ice and should be stored at -30°C to -15°C upon arrival. When stored correctly, the QIAseq Targeted DNA Pro Kits are good until the expiration date printed on the kit label.

Note: Open the package and store the products appropriately immediately upon receipt.



Panel Contents

QIAseq Targeted DNA Pro Catalog Number	333651	333655
No. of samples	12	96
One pool of region-specific primers	90 μ L	700 μ L
FFPE repair reagent	15 μ L	80 μ L
FFPE repair enzyme	17 μ L	120 μ L
FX Buffer, 10x	23 μ L	170 μ L
5X WGS FX Mix	42 μ L	330 μ L
FG Solution	170 μ L	170 μ L
UPH Ligation Buffer, 2.5x	288 μ L	1152 μ L
DNA Ligase	75 μ L	330 μ L
Ligation Cleanup Reagent	30 μ L	240 μ L
TEPCR Buffer, 5x	3 x 60 μ L	3 x 500 μ L
TEPCR modifier	30 μ L	230 μ L
TEPCR Cleanup Reagent	80 μ L	600 μ L
UPCR Buffer, 5x	500 μ L	4 x 500 μ L
QN Taq Polymerase	2 x 60 μ L	2 x 450 μ L
Nuclease-free water	1.5 mL	10 mL
One bottle containing QIAseq Bead Binding Buffer	10.2 mL	10.2 mL
One bottle containing QIAseq Beads	10 mL	10 mL

Gene Table

Gene Symbol	Gene Symbol	Gene Symbol	Gene Symbol
ABL1	AKT1	ALK	ATM
BAP1	BARD1	BRAF	BRCA1
BRCA2	BRIP1	CCND1	CD274
CDK12	CDK4	CDKN2A	CHEK1
CHEK2	CTNNB1	DDR2	EGFR
ERBB2	FANCL	FGFR1	FGFR3
HRAS	IDH1	IRF1	KIT
KRAS	MAP2K1	MAP3K8	MET
MLH1	MLH3	MSH2	MSH3
MSH6	MYC	MYCN	NF2
NRAS	NTRK1	PALB2	PDGFRA
PIK3CA	PMS1	PMS2	PPP2R1B
PRKN	PTEN	RAD51B	RAD51C
RAD51D	RAD54L	RB1	RET
ROS1	SLC22A18	STK11	TP53
TSC2			

Functional Groupings

Adenocarcinomas: ABL1, AKT1, ALK, ATM, BRAF, CCND1, CD274, CDK4, CDKN2A, CTNNB1, DDR2, EGFR, ERBB2, FGFR1, FGFR3, HRAS, KRAS, MAP2K1, MAP3K8, MET, MYC, MYCN, NRAS, NTRK1, PIK3CA, PPP2R1B, PRKN, PTEN, RET, ROS1, SLC22A18, STK11, TP53, TSC2.

Malignant Pleural Mesotheliomas: BAP1, NF2.

Non-Small Cell Lung Cancers (NSCLC): AKT1, ALK, ATM, BARD1, BRAF, BRCA1, BRCA2, BRIP1, CCND1, CDK12, CDK4, CDKN2A, CHEK1, CHEK2, CTNNB1, EGFR, ERBB2, FANCL, FGFR3, HRAS, IDH1, IRF1, KIT, KRAS, MAP2K1, MET, MLH1, MLH3, MSH2, MSH3, MSH6, MYC, NRAS, PALB2, PDGFRA, PIK3CA, PMS1, PMS2, PTEN, RAD51B, RAD51C, RAD51D, RAD54L, RET, ROS1, TP53.

Small Cell Lung Cancers (SCLC): CD274, HRAS, MYC, MYCN, NTRK1, PIK3CA, PPP2R1B, PTEN, RB1, TP53.

Squamous Cell Carcinomas (SCC): BRAF, CD274, CDKN2A, DDR2, EGFR, FGFR1, FGFR3, HRAS, MAP3K8, MYC, MYCN, NTRK1, PIK3CA, PPP2R1B, PTEN, SLC22A18, TP53, TSC2.

Ordering Information

Product	Contents	Cat. no.
QIAseq Targeted DNA Pro (12)	All reagents (except indexes) for targeted DNA sequencing; fixed panel for 12 samples; less than 200 genes	333651
QIAseq Targeted DNA Pro (96)	All reagents (except indexes) for targeted DNA sequencing; fixed panel for 96 samples; less than 200 genes	333655
QIAseq Targeted DNA Pro HC (12)	All reagents (except indexes) for targeted DNA sequencing; fixed panel for 12 samples; more than 200 genes	333661
QIAseq Targeted DNA Pro HC (96)	All reagents (except indexes) for targeted DNA sequencing; fixed panel for 96 samples; more than 200 genes	333665
QIAseq Targeted DNA Pro Custom (96)	All reagents (except indexes) for targeted DNA sequencing; custom panel for 96 samples	333675
QIAseq Targeted DNA Pro Booster (96)	Pool of primers used in combination with either catalogued or custom panels	333685

QIAseq Targeted DNA Pro Unique Dual Indices

QIAseq Targeted DNA Pro UDI Set A (96)	Box containing unique molecularly indexed adapters and primers, enough for a total of 96 samples, for indexing up to 96 samples for targeted panel sequencing on Illumina platforms; Set A (of A, B, C, D, E, F, G, and H) required for multiplexing 768 samples in one run	333455
QIAseq Targeted DNA Pro UDI Set B (96)	Box containing unique molecularly indexed adapters and primers, enough for a total of 96 samples, for indexing up to 96 samples for targeted panel sequencing on Illumina platforms; Set B (of A, B, C, D, E, F, G, and H) required for multiplexing 768 samples in one run	333465
QIAseq Targeted DNA Pro UDI Set C (96)	Box containing unique molecularly indexed adapters and primers, enough for a total of 96 samples, for indexing up to 96 samples for targeted panel sequencing on Illumina platforms; Set C (of A, B, C, D, E, F, G, and H) required for multiplexing 768 samples in one run	333475

Product	Contents	Cat. no.
QIAseq Targeted DNA Pro UDI Set D (96)	Box containing unique molecularly indexed adapters and primers, enough for a total of 96 samples, for indexing up to 96 samples for targeted panel sequencing on Illumina platforms; Set D (of A, B, C, D, E, F, G, and H) required for multiplexing 768 samples in one run	333485
QIAseq Targeted DNA Pro UDI Set E (96)	Box containing unique molecularly indexed adapters and primers, enough for a total of 96 samples, for indexing up to 96 samples for targeted panel sequencing on Illumina platforms; Set E (of A, B, C, D, E, F, G, and H) required for multiplexing 768 samples in one run	334015
QIAseq Targeted DNA Pro UDI Set F (96)	Box containing unique molecularly indexed adapters and primers, enough for a total of 96 samples, for indexing up to 96 samples for targeted panel sequencing on Illumina platforms; Set F (of A, B, C, D, E, F, G, and H) required for multiplexing 768 samples in one run	334025
QIAseq Targeted DNA Pro UDI Set G (96)	Box containing unique molecularly indexed adapters and primers, enough for a total of 96 samples, for indexing up to 96 samples for targeted panel sequencing on Illumina platforms; Set G (of A, B, C, D, E, F, G, and H) required for multiplexing 768 samples in one run	334035
QIAseq Targeted DNA Pro UDI Set H (96)	Box containing unique molecularly indexed adapters and primers, enough for a total of 96 samples, for indexing up to 96 samples for targeted panel sequencing on Illumina platforms; Set H (of A, B, C, D, E, F, G, and H) required for multiplexing 768 samples in one run	334045
QIAseq Targeted DNA Pro UDI (12)	Box containing unique molecularly-indexed adapters and primers, enough for a total of 12 samples, for indexing up to 12 samples for targeted panel sequencing on Illumina platforms	333441
QIAseq Targeted DNA Pro Indices for Ion Torrent platforms		
QIAseq Targeted DNA Pro Index L (12)	Box containing molecularly-indexed adapters and primers, enough for a total of 12 samples – for indexing up to 12 samples for targeted panel sequencing on Ion Torrent platforms	333491

Product	Contents	Cat. no.
QIAseq Targeted DNA Pro Index L (24)	Box containing molecularly-indexed adapters and primers, enough for a total of 24 samples – for indexing up to 24 samples for targeted panel sequencing on Ion Torrent platforms	333492
Related products		
QIAseq Library Quant Assay Kit	Reagents for quantification of libraries prepared for Illumina or Ion Torrent platforms; assay format	333314
QIAseq Universal Normalizer Kit (96)	For 96 reactions: Reagents for the normalization of libraries using indexed library amplification primers.	180615
QIAseq Universal Normalizer Kit (24)	For 24 reactions: Reagents for the normalization of libraries using indexed library amplification primers.	180613
QIAamp DNA Mini Kit (50)	For 50 DNA preps: 50 QIAamp Mini Spin Columns, QIAGEN Proteinase K, Collection Tubes (2 mL), reagents and buffers	51304
QIAamp Circulating Nucleic Acid Kit (50)	For 50 DNA preps: QIAamp Mini Columns, Tube Extenders (20 mL), QIAGEN Proteinase K, Carrier RNA, Buffers, VacConnectors, and Collection Tubes (1.5 mL and 2 mL)	55114
QIAamp DNA FFPE Advanced UNG Kit (50)	For 50 preps: Uracil-N-glycosylase, QIAamp UCP MinElute columns, collection tubes, Deparaffinization Solution, Proteinase K, RNase A, RNase-free water and buffers	56704

QIAseq Targeted DNA Pro Panels are intended for molecular biology applications. These products are not intended for the diagnosis, prevention or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN® (QIAGEN Group)

1100784 02/2016 © 2016 QIAGEN, all rights reserved.

www.qiagen.com

Australia ■ 1-800-243-800

Austria ■ 0800/281010

Belgium ■ 0800-79612

Brazil ■ 0800-557779

Canada ■ 800-572-9613

China ■ 8621-3865-3865

Denmark ■ 80-885945

Finland ■ 0800-914416

France ■ 01-60-920-930

Germany ■ 02103-29-12000

Hong Kong ■ 800 933 965

Ireland ■ 1800 555 049

Italy ■ 800-787980

Japan ■ 03-6890-7300

Korea (South) ■ 080-000-7145

Luxembourg ■ 8002 2076

Mexico ■ 01-800-7742-436

The Netherlands ■ 0800 0229592

Norway ■ 800-18859

Singapore ■ 1800-742-4368

Spain ■ 91-630-7050

Sweden ■ 020-790282

Switzerland ■ 055-254-22-11

UK ■ 01293-422-911

USA ■ 800-426-8157



Sample to Insight