


Validation Data Sheet

Quant Nova LNA PCR Assay	IL17C
Assay Catalog No.	SBH0371753
dPCR Instrument	QIAcuity digital PCR instrument
dPCR Master Mix	QIAcuity EG PCR Master Mix
Experimental Sample(s)	gBlock

Gene Expression Analysis: Concentration values [copies/ μ l] for different gBlock sample amounts (0 copies, 15600 copies, 23400 copies, 31200 copies and 46800 copies). Technical triplicates are shown.



1D Scatter Plot: Single-well data (representative well of technical triplicates) for the different gBlock input amounts (FAM channel for EvaGreen detection; gBlock amounts in the 5 wells: 0 copies/ μ l, 1300 copies/ μ l, 1950 copies/ μ l, 2600 copies/ μ l and 3900 copies/ μ l).

