

Validation Data Sheet

Quant Nova LNA PCR Assay	IL3
Assay Catalog No.	SBH0584080
dPCR Instrument	QIAcuity digital PCR instrument
dPCR Master Mix	QIAcuity EG PCR Master Mix
Experimental Sample(s)	gBlock

Gene Expression Analysis: Concentration values [copies/ μ l] for different gBlock sample amounts (0 copies, 7560 copies, 11340 copies, 15120 copies and 22680 copies). Technical triplicates are shown.



1D Scatter Plot: Single-well data (representative well of technical triplicates) for the different gBlock input amounts (FAM channel for EvaGreen detection; gBlock amounts in the 5 wells: 0 copies/ μ l, 630 copies/ μ l, 945 copies/ μ l, 1260 copies/ μ l and 1890 copies/ μ l).

