


Validation Data Sheet

QuantiNova LNA PCR Assay	GSC
Assay Catalog No.	SBH1220040
dPCR Instrument	QIAcuity digital PCR instrument
dPCR Master Mix	QIAcuity EG PCR Master Mix
Experimental Sample(s)	gBlock

Gene Expression Analysis: Concentration values [copies/ μ l] for different gBlock sample amounts (0 copies, 6000 copies, 9000 copies, 12000 copies and 18000 copies). Technical triplicates are shown.



1D Scatter Plot: Single-well data (representative well of technical triplicates) for the different gBlock input amounts (FAM channel for EvaGreen detection; gBlock amounts in the 5 wells: 0 copies/ μ l, 500 copies/ μ l, 750 copies/ μ l, 1000 copies/ μ l and 1500 copies/ μ l).

