


# Validation Data Sheet

QuantiNova LNA PCR Assay	SNAI3
Assay Catalog No.	SBH1220411
dPCR Instrument	QIAcuity digital PCR instrument
dPCR Master Mix	QIAcuity EG PCR Master Mix
Experimental Sample(s)	gBlock

**Gene Expression Analysis:** Concentration values [copies/ $\mu$ l] for different gBlock sample amounts (0 copies, 2040 copies, 3060 copies, 4080 copies and 6120 copies). Technical triplicates are shown.



**1D Scatter Plot:** Single-well data (representative well of technical triplicates) for the different gBlock input amounts (FAM channel for EvaGreen detection; gBlock amounts in the 5 wells: 0 copies/ $\mu$ l, 170 copies/ $\mu$ l, 255 copies/ $\mu$ l, 340 copies/ $\mu$ l and 510 copies/ $\mu$ l).

