

QuantiNova® LNA® PCR Focus Panels (Rotor-Gene® Format)

Human Epithelial to Mesenchymal Transition (EMT)

Cat. no. 249950 SBHS-090ZR

For study focus gene expression analysis

Shipping and storage

QuantiNova LNA PCR Focus Panels are shipped at ambient temperature. Immediately upon receipt, they should be stored at 2–8°C for short term storage or at –30°C to –15°C for long time storage. Under these conditions, all components are stable for at least 12 months.

Note: Open the package and store the products appropriately immediately upon receipt.

For optimal performance, QuantiNova LNA PCR Focus Panels should be used together with the QuantiNova Reverse Transcription Kit for cDNA synthesis and the QuantiNova SYBR® Green PCR Kit (Mastermix) for PCR.

Panel layout (Rotor-Gene): QuantiNova LNA PCR Focus Panel

The 96 real-time assays in the Rotor-Gene format are located in wells 1–96 of the Rotor-Disc® (plate A1–A12=Rotor-Disc 1–12, plate B1–B12=Rotor-Disc 13–24, etc.). To maintain data analysis compatibility, wells 97–100 do not contain real-time assays but will contain master mix to account for weight balance. Refer to the QuantiNova LNA PCR System Handbook at www.qiagen.com for further details.

	1	2	3	4	5	6	7	8	9	10	11	12
A	AHNK	AKT1	BMP1	BMP2	BMP7	CALD1	CAMK2N1	CAV2	CDH1	CDH2	COL1A2	COL3A1
B	COL5A2	CTNNB1	DSC2	DSP	EGFR	ERBB3	ESR1	F11R	FGFBP1	FN1	FOXC2	FZD7
C	GNG11	GSC	GSK3B	IGFBP4	IL1RN	ILK	ITGA5	ITGAV	ITGB1	JAG1	KRT14	KRT19
D	KRT7	MAP1B	MMP2	MMP3	MMP9	MSN	MST1R	NODAL	NOTCH1	NUDT13	OCLN	PDGFRB
E	PLEK2	DES1	PTK2	PTP4A1	RAC1	RGS2	SERPINE1	GEMIN2	SMAD2	SNAI1	SNAI2	SNAI3
F	SOX10	SPARC	SPP1	STAT3	STEAP1	TCF3	TCF4	TFPI2	TGFB1	TGFB2	TGFB3	TIMP1
G	TMEFF1	TMEM132A	TSPAN13	TWIST1	VCAN	VIM	VPS13A	WNT11	WNT5A	WNT5B	ZEB1	ZEB2
H	ACTB	B2M	GAPDH	HPRT1	RPLP0	HGDC	QIC	QIC	QIC	PPC	PPC	PPC

Gene table: QuantiNova LNA PCR Focus Panel

Position	Assay	Name	Symbol	Ensembl ID	Description
A01	SBH0115423	ENST00000257247.11	AHNAK	ENSG00000124942	AHNAK nucleoprotein Source HGNC Symbol Acc HGNC 347
A02	SBH0095396	ENST00000555528.5	AKT1	ENSG00000142208	AKT serine/threonine kinase 1 Source HGNC Symbol Acc HGNC 391
A03	SBH1219801	ENST00000354870.5	BMP1	ENSG00000168487	bone morphogenetic protein 1 Source HGNC Symbol Acc HGNC 1067
A04	SBH1219802	ENST00000378827.5	BMP2	ENSG00000125845	bone morphogenetic protein 2 Source HGNC Symbol Acc HGNC 1069
A05	SBH1219806	ENST00000450594.6	BMP7	ENSG00000101144	bone morphogenetic protein 7 Source HGNC Symbol Acc HGNC 1074
A06	SBH0414365	ENST00000361901.6	CALD1	ENSG00000122786	caldesmon 1 Source HGNC Symbol Acc HGNC 1441
A07	SBH1219821	ENST00000375078.4	CAMK2N1	ENSG00000162545	calcium/calmodulin dependent protein kinase II inhibitor 1 Source HGNC Symbol Acc HGNC 24190
A08	SBH0607858	ENST00000484871.5	CAV2	ENSG00000105971	caveolin 2 Source HGNC Symbol Acc HGNC 1528
A09	SBH1219869	ENST00000261769.10	CDH1	ENSG00000039068	cadherin 1 Source HGNC Symbol Acc HGNC 1748
A10	SBH1219870	ENST00000269141.8	CDH2	ENSG00000170558	cadherin 2 Source HGNC Symbol Acc HGNC 1759
A11	SBH0096733	ENST00000297268.10	COL1A2	ENSG00000164692	collagen type I alpha 2 chain Source HGNC Symbol Acc HGNC 2198
A12	SBH0521348	ENST00000304636.7	COL3A1	ENSG00000168542	collagen type III alpha 1 chain Source HGNC Symbol Acc HGNC 2201
B01	SBH1219899	ENST00000649066.1	COL5A2	ENSG00000204262	collagen type V alpha 2 chain Source HGNC Symbol Acc HGNC 2210
B02	SBH0588482	ENST00000396183.7	CTNNB1	ENSG00000168036	catenin beta 1 Source HGNC Symbol Acc HGNC 2514
B03	SBH1219959	ENST00000280904.10	DSC2	ENSG00000134755	desmocollin 2 Source HGNC Symbol Acc HGNC 3036
B04	SBH1219960	ENST00000379802.8	DSP	ENSG00000096696	desmoplakin Source HGNC Symbol Acc HGNC 3052
B05	SBH1219970	ENST00000454757.6	EGFR	ENSG00000146648	epidermal growth factor receptor Source HGNC Symbol Acc HGNC 3236
B06	SBH1219981	ENST00000549832.1	ERBB3	ENSG00000065361	erb-b2 receptor tyrosine kinase 3 Source HGNC Symbol Acc HGNC 3431
B07	SBH0125383	ENST00000206249.7	ESR1	ENSG00000091831	estrogen receptor 1 Source HGNC Symbol Acc HGNC 3467
B08	SBH1219988	ENST00000368026.11	F11R	ENSG00000158769	F11 receptor Source HGNC Symbol Acc HGNC 14685
B09	SBH0619645	ENST00000382333.1	FGFBP1	ENSG00000137440	fibroblast growth factor binding protein 1 Source HGNC Symbol Acc HGNC 19695
B10	SBH1220003	ENST00000354785.9	FN1	ENSG00000115414	fibronectin 1 Source HGNC Symbol Acc HGNC 3778
B11	SBH1220006	ENST00000320354.6	FOXC2	ENSG00000176692	forkhead box C2 Source HGNC Symbol Acc HGNC 3801
B12	SBH0357958	ENST00000286201.2	FZD7	ENSG00000155760	frizzled class receptor 7 Source HGNC Symbol Acc HGNC 4045
C01	SBH0478483	ENST00000248564.5	GNG11	ENSG00000127920	G protein subunit gamma 11 Source HGNC Symbol Acc HGNC 4403
C02	SBH1220040	ENST00000238558.4	GSC	ENSG00000133937	goosecoid homeobox Source HGNC Symbol Acc HGNC 4612
C03	SBH0579883	ENST00000316626.5	GSK3B	ENSG00000082701	glycogen synthase kinase 3 beta Source HGNC Symbol Acc HGNC 4617
C04	SBH1220093	ENST00000269593.5	IGFBP4	ENSG00000141753	insulin like growth factor binding protein 4 Source HGNC Symbol Acc HGNC 5473
C05	SBH0473919	ENST00000354115.6	IL1RN	ENSG00000136689	interleukin 1 receptor antagonist Source HGNC Symbol Acc HGNC 6000
C06	SBH0381135	ENST00000299421.8	ILK	ENSG00000166333	integrin linked kinase Source HGNC Symbol Acc HGNC 6040
C07	SBH1220133	ENST00000293379.9	ITGA5	ENSG00000161638	integrin subunit alpha 5 Source HGNC Symbol Acc HGNC 6141
C08	SBH0064907	ENST00000460641.1	ITGAV	ENSG00000138448	integrin subunit alpha V Source HGNC Symbol Acc HGNC 6150
C09	SBH1220136	ENST00000302278.8	ITGB1	ENSG00000150093	integrin subunit beta 1 Source HGNC Symbol Acc HGNC 6153
C10	SBH0407654	ENST00000254958.10	JAG1	ENSG00000101384	jagged 1 Source HGNC Symbol Acc HGNC 6188
		ENST00000167		ENSG000000	

Position	Assay	Name	Symbol	Ensembl ID	Description
C11	SBH1220158	586.7	KRT14	186847	keratin 14 Source HGNC Symbol Acc HGNC 6416
C12	SBH1220159	ENST00000361566.7	KRT19	ENSG00000171345	keratin 19 Source HGNC Symbol Acc HGNC 6436
D01	SBH1220160	ENST00000331817.6	KRT7	ENSG00000135480	keratin 7 Source HGNC Symbol Acc HGNC 6445
D02	SBH1220187	ENST00000296755.12	MAP1B	ENSG00000131711	microtubule associated protein 1B Source HGNC Symbol Acc HGNC 6836
D03	SBH1220222	ENST00000570308.5	MMP2	ENSG00000087245	matrix metalloproteinase 2 Source HGNC Symbol Acc HGNC 7166
D04	SBH1220223	ENST00000299855.10	MMP3	ENSG00000149968	matrix metalloproteinase 3 Source HGNC Symbol Acc HGNC 7173
D05	SBH0471278	ENST00000372330.3	MMP9	ENSG00000100985	matrix metalloproteinase 9 Source HGNC Symbol Acc HGNC 7176
D06	SBH1220230	ENST00000360270.7	MSN	ENSG00000147065	moesin Source HGNC Symbol Acc HGNC 7373
D07	SBH0127162	ENST00000621387.4	MST1R	ENSG00000164078	macrophage stimulating 1 receptor Source HGNC Symbol Acc HGNC 7381
D08	SBH0463463	ENST00000287139.7	NODAL	ENSG00000156574	nodal growth differentiation factor Source HGNC Symbol Acc HGNC 7865
D09	SBH0615258	ENST00000277541.7	NOTCH1	ENSG00000148400	notch 1 Source HGNC Symbol Acc HGNC 7881
D10	SBH0402961	ENST00000357321.9	NUDT13	ENSG00000166321	nudix hydrolase 13 Source HGNC Symbol Acc HGNC 18827
D11	SBH1220284	ENST00000355237.6	OCLN	ENSG00000197822	occludin Source HGNC Symbol Acc HGNC 8104
D12	SBH1220293	ENST00000261799.9	PDGFRB	ENSG00000113721	platelet derived growth factor receptor beta Source HGNC Symbol Acc HGNC 8804
E01	SBH1220316	ENST00000216446.9	PLEK2	ENSG00000100558	pleckstrin 2 Source HGNC Symbol Acc HGNC 19238
E02	SBH1219945	ENST00000263256.7	DESI1	ENSG00000100418	desumoylating isopeptidase 1 Source HGNC Symbol Acc HGNC 24577
E03	SBH1220345	ENST00000523539.5	PTK2	ENSG00000169398	protein tyrosine kinase 2 Source HGNC Symbol Acc HGNC 9611
E04	SBH0121737	ENST00000626021.2	PTP4A1	ENSG00000112245	protein tyrosine phosphatase type IVA, member 1 Source HGNC Symbol Acc HGNC 9634
E05	SBH1220352	ENST00000356142.4	RAC1	ENSG00000136238	Rac family small GTPase 1 Source HGNC Symbol Acc HGNC 9801
E06	SBH1220366	ENST00000235382.7	RGS2	ENSG00000116741	regulator of G protein signaling 2 Source HGNC Symbol Acc HGNC 9998
E07	SBH1220389	ENST00000223095.4	SERPINE1	ENSG00000106366	serpin family E member 1 Source HGNC Symbol Acc HGNC 8583
E08	SBH0401216	ENST00000308317.10	GEMIN2	ENSG00000092208	gem nuclear organelle associated protein 2 Source HGNC Symbol Acc HGNC 10884
E09	SBH1220405	ENST00000262160.11	SMAD2	ENSG00000175387	SMAD family member 2 Source HGNC Symbol Acc HGNC 6768
E10	SBH1220409	ENST00000244050.3	SNAI1	ENSG00000124216	snail family transcriptional repressor 1 Source HGNC Symbol Acc HGNC 11128
E11	SBH1220410	ENST00000642303.1	SNAI2	ENSG00000019549	snail family transcriptional repressor 2 Source HGNC Symbol Acc HGNC 11094
E12	SBH1220411	ENST00000332281.6	SNAI3	ENSG00000185669	snail family transcriptional repressor 3 Source NCBI gene Acc 333929
F01	SBH1220416	ENST00000360880.6	SOX10	ENSG00000100146	SRY-box 10 Source HGNC Symbol Acc HGNC 11190
F02	SBH1220420	ENST00000521569.1	SPARC	ENSG00000113140	secreted protein acidic and cysteine rich Source HGNC Symbol Acc HGNC 11219
F03	SBH0180162	ENST00000237623.11	SPP1	ENSG00000118785	secreted phosphoprotein 1 Source HGNC Symbol Acc HGNC 11255
F04	SBH0341614	ENST00000404395.3	STAT3	ENSG00000168610	signal transducer and activator of transcription 3 Source HGNC Symbol Acc HGNC 11364
F05	SBH1220425	ENST00000297205.7	STEAP1	ENSG00000164647	STEAP family member 1 Source HGNC Symbol Acc HGNC 11378
F06	SBH0362567	ENST00000395423.7	TCF3	ENSG00000071564	transcription factor 3 Source HGNC Symbol Acc HGNC 11633
F07	SBH0089633	ENST00000356073.8	TCF4	ENSG00000196628	transcription factor 4 Source HGNC Symbol Acc HGNC 11634
F08	SBH0102008	ENST00000649730.1	TFPI2	ENSG00000105825	tissue factor pathway inhibitor 2 Source HGNC Symbol Acc HGNC 11761
F09	SBH1220443	ENST00000598758.5	TGFB1	ENSG00000105329	transforming growth factor beta 1 Source NCBI gene Acc 7040
F10	SBH1220444	ENST00000366930.9	TGFB2	ENSG00000092969	transforming growth factor beta 2 Source HGNC Symbol Acc HGNC 11768

Position	Assay	Name	Symbol	Ensembl ID	Description
F11	SBH0179529	ENST00000238682.7	TGFB3	ENSG00000119699	transforming growth factor beta 3 Source HGNC Symbol Acc HGNC 11769
F12	SBH1220454	ENST00000218388.9	TIMP1	ENSG00000102265	TIMP metalloproteinase inhibitor 1 Source HGNC Symbol Acc HGNC 11820
G01	SBH1220467	ENST00000374879.5	TMEFF1	ENSG00000241697	transmembrane protein with EGF like and two follistatin like domains 1 Source HGNC Symbol Acc HGNC 11866
G02	SBH0430856	ENST00000453848.6	TMEM132A	ENSG00000006118	transmembrane protein 132A Source HGNC Symbol Acc HGNC 31092
G03	SBH0198926	ENST00000466195.1	TSPAN13	ENSG00000106537	tetraspanin 13 Source HGNC Symbol Acc HGNC 21643
G04	SBH1220496	ENST00000242261.6	TWIST1	ENSG00000122691	twist family bHLH transcription factor 1 Source HGNC Symbol Acc HGNC 12428
G05	SBH1220516	ENST00000343200.9	VCAN	ENSG00000038427	versican Source HGNC Symbol Acc HGNC 2464
G06	SBH1220518	ENST00000544301.7	VIM	ENSG00000026025	vimentin Source HGNC Symbol Acc HGNC 12692
G07	SBH1220520	ENST00000643348.1	VPS13A	ENSG00000197969	vacuolar protein sorting 13 homolog A Source HGNC Symbol Acc HGNC 1908
G08	SBH1220528	ENST00000322563.8	WNT11	ENSG00000085741	Wnt family member 11 Source HGNC Symbol Acc HGNC 12776
G09	SBH0548767	ENST00000264634.8	WNT5A	ENSG00000114251	Wnt family member 5A Source HGNC Symbol Acc HGNC 12784
G10	SBH1220534	ENST00000397196.7	WNT5B	ENSG00000111186	Wnt family member 5B Source HGNC Symbol Acc HGNC 16265
G11	SBH0357652	ENST00000320985.14	ZEB1	ENSG00000148516	zinc finger E-box binding homeobox 1 Source HGNC Symbol Acc HGNC 11642
G12	SBH1220542	ENST00000558170.6	ZEB2	ENSG00000169554	zinc finger E-box binding homeobox 2 Source HGNC Symbol Acc HGNC 14881
H01	SBH1220543	ENST00000646664.1	ACTB	ENSG00000075624	actin beta Source HGNC Symbol Acc HGNC 132
H02	SBH1220550	ENST00000558401.6	B2M	ENSG00000166710	beta-2-microglobulin Source HGNC Symbol Acc HGNC 914
H03	SBH1220545	ENST00000396861.5	GAPDH	ENSG00000111640	glyceraldehyde-3-phosphate dehydrogenase Source HGNC Symbol Acc HGNC 4141
H04	SBH1220546	ENST00000298556.8	HPRT1	ENSG00000165704	hypoxanthine phosphoribosyltransferase 1 Source HGNC Symbol Acc HGNC 5157
H05	SBH1220553	ENST00000546989.5	RPLP0	ENSG00000089157	ribosomal protein lateral stalk subunit P0 Source HGNC Symbol Acc HGNC 10371
H06	SBH1218553	Sybr_HGDC	HGDC	Sybr_HGDC	Human Genomic DNA Contamination
H07	SBH1218551	Sybr_QIC	QIC	Sybr_QIC	QuantiNova Internal Control
H08	SBH1218551	Sybr_QIC	QIC	Sybr_QIC	QuantiNova Internal Control
H09	SBH1218551	Sybr_QIC	QIC	Sybr_QIC	QuantiNova Internal Control
H10	SBH1218550	Sybr_PPC	PPC	Sybr_PPC	Positive PCR Control
H11	SBH1218550	Sybr_PPC	PPC	Sybr_PPC	Positive PCR Control
H12	SBH1218550	Sybr_PPC	PPC	Sybr_PPC	Positive PCR Control

Related products

Product	Contents	Cat. no.
QuantiNova LNA PCR QC Panel	These panels are designed to assess the quality of RNA samples before characterization using QuantiNova LNA PCR Focus Panels; available in 96-well, 384-well, and Rotor-Disc 100 formats	249940
QuantiNova Reverse Transcription Kit (10)*	For 10 x 20 µl reactions: 20 µl 8x gDNA Removal Mix, 10 µl Reverse Transcription Enzyme, 40 µl Reverse Transcription Mix (containing RT primers), 20 µl Internal Control RNA, 1.9 ml RNase-Free Water	205410
QuantiNova SYBR Green RT-PCR Kit (100)*	For 100 x 20 µl reactions: 1 ml QuantiNova SYBR Green RT-PCR Master Mix, 20 µl QuantiNova SYBR Green RT Mix, 20 µl Internal Control RNA, 500 µl Yellow Template Dilution Buffer, 250 µl ROX Reference Dye, 1.9 µl RNase-Free Water	208152
QuantiNova SYBR Green PCR Kit (100)*	For 100 x 20 µl reactions: 1 ml 2x QuantiNova SYBR Green PCR Master Mix, 500 µl QuantiNova Yellow Template Dilution Buffer, 250 µl QN ROX Reference Dye, 1.9 ml Water	208052

*Larger kit sizes available.

The QuantiNova LNA PCR Focus Panels are intended for molecular biology applications. These products are not intended for the diagnosis, prevention or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, LNA®, QuantiNova®, Sample to Insight® (QIAGEN Group); SYBR® (Life Technologies Corp.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

09/2019 © 2019 QIAGEN, all rights reserved.