

QuantiNova® LNA® PCR Focus Panels (Rotor-Gene® Format)

Human Fibrosis

Cat. no. 249950 SBHS-120ZR

For study focus gene expression analysis

Shipping and storage

QuantiNova LNA PCR Focus Panels are shipped at ambient temperature. Immediately upon receipt, they should be stored at 2–8°C for short term storage or at –30°C to –15°C for long time storage. Under these conditions, all components are stable for at least 12 months.

Note: Open the package and store the products appropriately immediately upon receipt.

For optimal performance, QuantiNova LNA PCR Focus Panels should be used together with the QuantiNova Reverse Transcription Kit for cDNA synthesis and the QuantiNova SYBR® Green PCR Kit (Mastermix) for PCR.

Panel layout (Rotor-Gene): QuantiNova LNA PCR Focus Panel

The 96 real-time assays in the Rotor-Gene format are located in wells 1–96 of the Rotor-Disc® (plate A1–A12=Rotor-Disc 1–12, plate B1–B12=Rotor-Disc 13–24, etc.). To maintain data analysis compatibility, wells 97–100 do not contain real-time assays but will contain master mix to account for weight balance. Refer to the QuantiNova LNA PCR System Handbook at www.qiagen.com for further details.

	1	2	3	4	5	6	7	8	9	10	11	12
A	ACTA2	AGT	AKT1	BCL2	BMP7	CAV1	CCL11	CCL2	CCL3	CCR2	CEBPB	COL1A2
B	COL3A1	CCN2	CXCR4	DCN	EDN1	EGF	ENG	FASLG	GREM1	HGF	IFNG	IL10
C	IL13	IL13RA2	IL1A	IL1B	IL4	IL5	ILK	INHBE	ITGA1	ITGA2	ITGA3	ITGAV
D	ITGB1	ITGB3	ITGB5	ITGB6	ITGB8	JUN	LOX	LTBP1	MMP1	MMP13	MMP14	MMP2
E	MMP3	MMP8	MMP9	MYC	NFKB1	PDGFA	PDGFB	PLAT	PLAU	PLG	SERPINA1	SERPINE1
F	SERPINH1	SMAD2	SMAD3	SMAD4	SMAD6	SMAD7	SNAI1	SP1	STAT1	STAT6	TGFB1	TGFB2
G	TGFB3	TGFBRI	TGFBRI2	TGFI1	THBS1	THBS2	TIMP1	TIMP2	TIMP3	TIMP4	TNF	VEGFA
H	ACTB	B2M	GAPDH	HPRT1	RPLP0	HGDC	QIC	QIC	QIC	PPC	PPC	PPC

Gene table: QuantiNova LNA PCR Focus Panel

Position	Assay	Name	Symbol	Ensembl ID	Description
A01	SBH0247260	ENST00000224784.10	ACTA2	ENSG00000107796	actin, alpha 2, smooth muscle, aorta Source HGNC Symbol Acc HGNC 130
A02	SBH1219729	ENST00000366667.5	AGT	ENSG00000135744	angiotensinogen Source HGNC Symbol Acc HGNC 333
A03	SBH0095396	ENST00000555528.5	AKT1	ENSG00000142208	AKT serine/threonine kinase 1 Source HGNC Symbol Acc HGNC 391
A04	SBH1219786	ENST00000398117.1	BCL2	ENSG00000171791	BCL2, apoptosis regulator Source HGNC Symbol Acc HGNC 990
A05	SBH1219806	ENST00000450594.6	BMP7	ENSG00000101144	bone morphogenetic protein 7 Source HGNC Symbol Acc HGNC 1074
A06	SBH0105254	ENST00000451122.5	CAV1	ENSG00000105974	caveolin 1 Source HGNC Symbol Acc HGNC 1527
A07	SBH0204041	ENST00000305869.3	CCL11	ENSG00000172156	C-C motif chemokine ligand 11 Source HGNC Symbol Acc HGNC 10610
A08	SBH0228134	ENST00000225831.4	CCL2	ENSG00000108691	C-C motif chemokine ligand 2 Source HGNC Symbol Acc HGNC 10618
A09	SBH1219838	ENST00000613922.2	CCL3	ENSG00000277632	C-C motif chemokine ligand 3 Source HGNC Symbol Acc HGNC 10627
A10	SBH0387563	ENST00000445132.2	CCR2	ENSG00000121807	C-C motif chemokine receptor 2 Source HGNC Symbol Acc HGNC 1603
A11	SBH0569983	ENST00000303004.4	CEBPB	ENSG00000172216	CCAAT enhancer binding protein beta Source HGNC Symbol Acc HGNC 1834
A12	SBH0096733	ENST00000297268.10	COL1A2	ENSG00000164692	collagen type I alpha 2 chain Source HGNC Symbol Acc HGNC 2198
B01	SBH0521348	ENST00000304636.7	COL3A1	ENSG00000168542	collagen type III alpha 1 chain Source HGNC Symbol Acc HGNC 2201
B02	SBH1219917	ENST00000367976.4	CCN2	ENSG00000118523	cellular communication network factor 2 Source HGNC Symbol Acc HGNC 2500
B03	SBH0591410	ENST00000241393.3	CXCR4	ENSG00000121966	C-X-C motif chemokine receptor 4 Source HGNC Symbol Acc HGNC 2561
B04	SBH1219942	ENST00000420120.6	DCN	ENSG00000011465	decorin Source HGNC Symbol Acc HGNC 2705
B05	SBH1219968	ENST00000379375.6	EDN1	ENSG00000078401	endothelin 1 Source HGNC Symbol Acc HGNC 3176
B06	SBH0321686	ENST00000265171.9	EGF	ENSG00000138798	epidermal growth factor Source HGNC Symbol Acc HGNC 3229
B07	SBH1219975	ENST00000480266.5	ENG	ENSG00000106991	endoglin Source HGNC Symbol Acc HGNC 3349
B08	SBH1219995	ENST00000367721.3	FASLG	ENSG00000117560	Fas ligand Source HGNC Symbol Acc HGNC 11936
B09	SBH1220039	ENST00000651154.1	GREM1	ENSG00000166923	gremlin 1, DAN family BMP antagonist Source HGNC Symbol Acc HGNC 2001
B10	SBH1220058	ENST00000457544.7	HGF	ENSG00000019991	hepatocyte growth factor Source HGNC Symbol Acc HGNC 4893
B11	SBH1220090	ENST00000229135.4	IFNG	ENSG00000111537	interferon gamma Source HGNC Symbol Acc HGNC 5438
B12	SBH1220095	ENST00000423557.1	IL10	ENSG00000136634	interleukin 10 Source HGNC Symbol Acc HGNC 5962
C01	SBH0375568	ENST00000304506.7	IL13	ENSG00000169194	interleukin 13 Source HGNC Symbol Acc HGNC 5973
C02	SBH1220100	ENST00000243213.2	IL13RA2	ENSG00000123496	interleukin 13 receptor subunit alpha 2 Source HGNC Symbol Acc HGNC 5975
C03	SBH0663647	ENST00000263339.3	IL1A	ENSG00000115008	interleukin 1 alpha Source HGNC Symbol Acc HGNC 5991
C04	SBH0079231	ENST00000263341.6	IL1B	ENSG00000125538	interleukin 1 beta Source HGNC Symbol Acc HGNC 5992
C05	SBH1220109	ENST00000350025.2	IL4	ENSG00000113520	interleukin 4 Source HGNC Symbol Acc HGNC 6014
C06	SBH1220110	ENST00000231454.6	IL5	ENSG00000113525	interleukin 5 Source HGNC Symbol Acc HGNC 6016
C07	SBH0381135	ENST00000299421.8	ILK	ENSG00000166333	integrin linked kinase Source HGNC Symbol Acc HGNC 6040
C08	SBH1220118	ENST00000266646.3	INHBE	ENSG00000139269	inhibin subunit beta E Source HGNC Symbol Acc HGNC 24029
C09	SBH1220129	ENST00000282588.6	ITGA1	ENSG00000213949	integrin subunit alpha 1 Source HGNC Symbol Acc HGNC 6134
C10	SBH1220130	ENST00000296585.10	ITGA2	ENSG00000164171	integrin subunit alpha 2 Source HGNC Symbol Acc HGNC 6137
		ENST000000007		ENSG000000	

Position	Assay	Name	Symbol	Ensembl ID	Description
C11	SBH1220131	722.11	ITGA3	005884	integrin subunit alpha 3 Source HGNC Symbol Acc HGNC 6139
C12	SBH0064907	ENST00000460641.1	ITGAV	ENSG00000138448	integrin subunit alpha V Source HGNC Symbol Acc HGNC 6150
D01	SBH1220136	ENST00000302278.8	ITGB1	ENSG00000150093	integrin subunit beta 1 Source HGNC Symbol Acc HGNC 6153
D02	SBH1220137	ENST00000559488.5	ITGB3	ENSG00000259207	integrin subunit beta 3 Source HGNC Symbol Acc HGNC 6156
D03	SBH1220139	ENST00000608657.5	ITGB5	ENSG00000082781	integrin subunit beta 5 Source HGNC Symbol Acc HGNC 6160
D04	SBH1220140	ENST00000409872.1	ITGB6	ENSG00000115221	integrin subunit beta 6 Source HGNC Symbol Acc HGNC 6161
D05	SBH1220141	ENST00000222573.5	ITGB8	ENSG00000105855	integrin subunit beta 8 Source HGNC Symbol Acc HGNC 6163
D06	SBH0613340	ENST00000371222.3	JUN	ENSG00000177606	Jun proto-oncogene, AP-1 transcription factor subunit Source HGNC Symbol Acc HGNC 6204
D07	SBH1220174	ENST00000639739.1	LOX	ENSG00000113083	lysyl oxidase Source HGNC Symbol Acc HGNC 6664
D08	SBH1220179	ENST00000407925.5	LTBP1	ENSG00000049323	latent transforming growth factor beta binding protein 1 Source HGNC Symbol Acc HGNC 6714
D09	SBH1220215	ENST00000315274.7	MMP1	ENSG00000196611	matrix metalloproteinase 1 Source HGNC Symbol Acc HGNC 7155
D10	SBH1220218	ENST00000615555.4	MMP13	ENSG00000137745	matrix metalloproteinase 13 Source HGNC Symbol Acc HGNC 7159
D11	SBH1220219	ENST00000311852.11	MMP14	ENSG00000157227	matrix metalloproteinase 14 Source HGNC Symbol Acc HGNC 7160
D12	SBH1220222	ENST00000570308.5	MMP2	ENSG00000087245	matrix metalloproteinase 2 Source HGNC Symbol Acc HGNC 7166
E01	SBH1220223	ENST00000299855.10	MMP3	ENSG00000149968	matrix metalloproteinase 3 Source HGNC Symbol Acc HGNC 7173
E02	SBH1220225	ENST00000236826.8	MMP8	ENSG00000118113	matrix metalloproteinase 8 Source HGNC Symbol Acc HGNC 7175
E03	SBH0471278	ENST00000372330.3	MMP9	ENSG00000100985	matrix metalloproteinase 9 Source HGNC Symbol Acc HGNC 7176
E04	SBH0426145	ENST00000524013.1	MYC	ENSG00000136997	MYC proto-oncogene, bHLH transcription factor Source HGNC Symbol Acc HGNC 7553
E05	SBH1220264	ENST00000651197.1	NFKB1	ENSG00000109320	nuclear factor kappa B subunit 1 Source HGNC Symbol Acc HGNC 7794
E06	SBH0498934	ENST00000354513.9	PDGFA	ENSG00000197461	platelet derived growth factor subunit A Source HGNC Symbol Acc HGNC 8799
E07	SBH0091370	ENST00000331163.10	PDGFB	ENSG00000100311	platelet derived growth factor subunit B Source HGNC Symbol Acc HGNC 8800
E08	SBH0432338	ENST00000220809.8	PLAT	ENSG00000104368	plasminogen activator, tissue type Source HGNC Symbol Acc HGNC 9051
E09	SBH1220315	ENST00000446342.5	PLAU	ENSG00000122861	plasminogen activator, urokinase Source HGNC Symbol Acc HGNC 9052
E10	SBH0191454	ENST00000308192.13	PLG	ENSG00000122194	plasminogen Source HGNC Symbol Acc HGNC 9071
E11	SBH0149659	ENST00000393087.8	SERPINA1	ENSG00000197249	serpin family A member 1 Source HGNC Symbol Acc HGNC 8941
E12	SBH1220389	ENST00000223095.4	SERPINE1	ENSG00000106366	serpin family E member 1 Source HGNC Symbol Acc HGNC 8583
F01	SBH0596911	ENST00000358171.7	SERPINH1	ENSG00000149257	serpin family H member 1 Source HGNC Symbol Acc HGNC 1546
F02	SBH1220405	ENST00000262160.11	SMAD2	ENSG00000175387	SMAD family member 2 Source HGNC Symbol Acc HGNC 6768
F03	SBH0216540	ENST00000558428.5	SMAD3	ENSG00000166949	SMAD family member 3 Source HGNC Symbol Acc HGNC 6769
F04	SBH1220406	ENST00000588745.5	SMAD4	ENSG00000141646	SMAD family member 4 Source HGNC Symbol Acc HGNC 6770
F05	SBH1220408	ENST00000558937.1	SMAD6	ENSG00000137834	SMAD family member 6 Source HGNC Symbol Acc HGNC 6772
F06	SBH0167278	ENST00000262158.7	SMAD7	ENSG00000101665	SMAD family member 7 Source HGNC Symbol Acc HGNC 6773
F07	SBH1220409	ENST00000244050.3	SNAI1	ENSG00000124216	snail family transcriptional repressor 1 Source HGNC Symbol Acc HGNC 11128
F08	SBH1220419	ENST00000426431.2	SP1	ENSG00000185591	Sp1 transcription factor Source HGNC Symbol Acc HGNC 11205
F09	SBH0333289	ENST00000361099.7	STAT1	ENSG00000115415	signal transducer and activator of transcription 1 Source HGNC Symbol Acc HGNC 11362
F10	SBH1220424	ENST00000300134.8	STAT6	ENSG00000166888	signal transducer and activator of transcription 6 Source HGNC Symbol Acc HGNC 11368

Position	Assay	Name	Symbol	Ensembl ID	Description
F11	SBH1220443	ENST00000598758.5	TGFB1	ENSG00000105329	transforming growth factor beta 1 Source NCBI gene Acc 7040
F12	SBH1220444	ENST00000366930.9	TGFB2	ENSG00000092969	transforming growth factor beta 2 Source HGNC Symbol Acc HGNC 11768
G01	SBH0179529	ENST00000238682.7	TGFB3	ENSG00000119699	transforming growth factor beta 3 Source HGNC Symbol Acc HGNC 11769
G02	SBH1220446	ENST00000374994.9	TGFBR1	ENSG00000106799	transforming growth factor beta receptor 1 Source HGNC Symbol Acc HGNC 11772
G03	SBH0598842	ENST00000295754.9	TGFBR2	ENSG00000163513	transforming growth factor beta receptor 2 Source HGNC Symbol Acc HGNC 11773
G04	SBH0261061	ENST00000334820.9	TGIF1	ENSG00000177426	TGFB induced factor homeobox 1 Source HGNC Symbol Acc HGNC 11776
G05	SBH1220450	ENST00000260356.5	THBS1	ENSG00000137801	thrombospondin 1 Source HGNC Symbol Acc HGNC 11785
G06	SBH0457929	ENST00000366787.7	THBS2	ENSG00000186340	thrombospondin 2 Source HGNC Symbol Acc HGNC 11786
G07	SBH1220454	ENST00000218388.9	TIMP1	ENSG00000102265	TIMP metalloproteinase inhibitor 1 Source HGNC Symbol Acc HGNC 11820
G08	SBH0450624	ENST00000262768.11	TIMP2	ENSG00000035862	TIMP metalloproteinase inhibitor 2 Source HGNC Symbol Acc HGNC 11821
G09	SBH1220455	ENST00000266085.6	TIMP3	ENSG00000100234	TIMP metalloproteinase inhibitor 3 Source HGNC Symbol Acc HGNC 11822
G10	SBH1220456	ENST00000287814.5	TIMP4	ENSG00000157150	TIMP metalloproteinase inhibitor 4 Source HGNC Symbol Acc HGNC 11823
G11	SBH1220471	ENST00000449264.3	TNF	ENSG00000232810	tumor necrosis factor Source HGNC Symbol Acc HGNC 11892
G12	SBH0420322	ENST00000425836.6	VEGFA	ENSG00000112715	vascular endothelial growth factor A Source HGNC Symbol Acc HGNC 12680
H01	SBH1220543	ENST00000646664.1	ACTB	ENSG00000075624	actin beta Source HGNC Symbol Acc HGNC 132
H02	SBH1220550	ENST00000558401.6	B2M	ENSG00000166710	beta-2-microglobulin Source HGNC Symbol Acc HGNC 914
H03	SBH1220545	ENST00000396861.5	GAPDH	ENSG00000111640	glyceraldehyde-3-phosphate dehydrogenase Source HGNC Symbol Acc HGNC 4141
H04	SBH1220546	ENST00000298556.8	HPRT1	ENSG00000165704	hypoxanthine phosphoribosyltransferase 1 Source HGNC Symbol Acc HGNC 5157
H05	SBH1220553	ENST00000546989.5	RPLP0	ENSG00000089157	ribosomal protein lateral stalk subunit P0 Source HGNC Symbol Acc HGNC 10371
H06	SBH1218553	Sybr_HGDC	HGDC	Sybr_HGDC	Human Genomic DNA Contamination
H07	SBH1218551	Sybr_QIC	QIC	Sybr_QIC	QuantiNova Internal Control
H08	SBH1218551	Sybr_QIC	QIC	Sybr_QIC	QuantiNova Internal Control
H09	SBH1218551	Sybr_QIC	QIC	Sybr_QIC	QuantiNova Internal Control
H10	SBH1218550	Sybr_PPC	PPC	Sybr_PPC	Positive PCR Control
H11	SBH1218550	Sybr_PPC	PPC	Sybr_PPC	Positive PCR Control
H12	SBH1218550	Sybr_PPC	PPC	Sybr_PPC	Positive PCR Control



Related products

Product	Contents	Cat. no.
QuantiNova LNA PCR QC Panel	These panels are designed to assess the quality of RNA samples before characterization using QuantiNova LNA PCR Focus Panels; available in 96-well, 384-well, and Rotor-Disc 100 formats	249940
QuantiNova Reverse Transcription Kit (10)*	For 10 x 20 μ l reactions: 20 μ l 8x gDNA Removal Mix, 10 μ l Reverse Transcription Enzyme, 40 μ l Reverse Transcription Mix (containing RT primers), 20 μ l Internal Control RNA, 1.9 ml RNase-Free Water	205410
QuantiNova SYBR Green RT-PCR Kit (100)*	For 100 x 20 μ l reactions: 1 ml QuantiNova SYBR Green RT-PCR Master Mix, 20 μ l QuantiNova SYBR Green RT Mix, 20 μ l Internal Control RNA, 500 μ l Yellow Template Dilution Buffer, 250 μ l ROX Reference Dye, 1.9 μ l RNase-Free Water	208152
QuantiNova SYBR Green PCR Kit (100)*	For 100 x 20 μ l reactions: 1 ml 2x QuantiNova SYBR Green PCR Master Mix, 500 μ l QuantiNova Yellow Template Dilution Buffer, 250 μ l QN ROX Reference Dye, 1.9 ml Water	208052

*Larger kit sizes available.

The QuantiNova LNA PCR Focus Panels are intended for molecular biology applications. These products are not intended for the diagnosis, prevention or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, LNA®, QuantiNova®, Sample to Insight® (QIAGEN Group); SYBR® (Life Technologies Corp.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

09/2019 © 2019 QIAGEN, all rights reserved.