

QuantiNova® LNA® PCR Focus Panels (Rotor-Gene® Format)

Rat Mitochondria

Cat. no. 249950 SBRN-087ZR

For study focus gene expression analysis

Shipping and storage

QuantiNova LNA PCR Focus Panels are shipped at ambient temperature. Immediately upon receipt, they should be stored at 2–8°C for short term storage or at –30°C to –15°C for long time storage. Under these conditions, all components are stable for at least 12 months.

Note: Open the package and store the products appropriately immediately upon receipt.

For optimal performance, QuantiNova LNA PCR Focus Panels should be used together with the QuantiNova Reverse Transcription Kit for cDNA synthesis and the QuantiNova SYBR® Green PCR Kit (Mastermix) for PCR.

Panel layout (Rotor-Gene): QuantiNova LNA PCR Focus Panel

The 96 real-time assays in the Rotor-Gene format are located in wells 1–96 of the Rotor-Disc® (plate A1–A12=Rotor-Disc 1–12, plate B1–B12=Rotor-Disc 13–24, etc.). To maintain data analysis compatibility, wells 97–100 do not contain real-time assays but will contain master mix to account for weight balance. Refer to the QuantiNova LNA PCR System Handbook at www.qiagen.com for further details.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Aifm2	Aip	Akt1	Bak1	Bmf	Bcl2	Bcl2l1	Bid	Bnip3	Cav2	Cdkn2a	Cln8
B	Cox18	Cpt1b	Cpt2	Dnm1l	Fis1	Timm10b	Glc	Gclm	Gpx1	Grpel1	Hsp90aa1	Hspd1
C	Immp1l	Cox10	Mfn1	LOC1009114 85	Mipep	Msto1	Mtr2	Nefl	Opa1	Pmaip1	Ppargc1a	Rhot2
D	Rhot1	Sfn	Sh3glb1	Slc25a10	Slc25a12	Slc25a13	Slc25a14	Slc25a15	Slc25a16	Slc25a17	Slc25a19	Slc25a2
E	Slc25a20	Slc25a21	LOC1009114 40	Slc25a23	Slc25a24	Slc25a25	Slc25a27	Slc25a3	Slc25a30	Slc25a36	Slc25a37	Slc25a4
F	Slc25a5	Sod1	Sod2	Stard3	Taz	Timm10	Timm17a	LOC1009111 30	Timm22	Timm44	Timm8a1	Timm8b
G	Timm9	Tomm22	Tomm34	Tomm40	Tomm40l	Tomm70	Tp53	Tspo	Ucp1	Ucp2	Ucp3	Uxt
H	Actb	B2m	Hprt1	Ldha	Rplp1	RGDC	QIC	QIC	QIC	PPC	PPC	PPC

Gene table: QuantiNova LNA PCR Focus Panel

Position	Assay	Name	Symbol	Ensembl ID	Description
A01	SBR1127028	ENSRNOT00000082463.1	Aifm2	ENSRNOG0000059445	apoptosis inducing factor, mitochondria associated 2 Source RGD Symbol Acc 1304964
A02	SBR1200763	ENSRNOT00000091542.1	Aip	ENSRNOG0000022289	aryl-hydrocarbon receptor-interacting protein Source RGD Symbol Acc 628619
A03	SBR1163078	ENSRNOT00000031164.3	Akt1	ENSRNOG0000028629	AKT serine/threonine kinase 1 Source RGD Symbol Acc 2081
A04	SBR1162844	ENSRNOT00000000576.5	Bak1	ENSRNOG0000000485	BCL2-antagonist/killer 1 Source RGD Symbol Acc 621635
A05	SBR1196608	ENSRNOT00000010381.4	Bmf	ENSRNOG0000007529	Bcl2 modifying factor Source RGD Symbol Acc 628658
A06	SBR1170984	ENSRNOT00000003768.2	Bcl2	ENSRNOG0000002791	BCL2, apoptosis regulator Source RGD Symbol Acc 2199
A07	SBR1205791	ENSRNOT00000010762.7	Bcl2l1	ENSRNOG0000007946	Bcl2-like 1 Source RGD Symbol Acc 2200
A08	SBR1147619	ENSRNOT00000016776.2	Bid	ENSRNOG0000012439	BH3 interacting domain death agonist Source RGD Symbol Acc 620160
A09	SBR1178674	ENSRNOT00000023477.6	Bnip3	ENSRNOG0000017243	BCL2 interacting protein 3 Source RGD Symbol Acc 620800
A10	SBR1201369	ENSRNOT00000080271.1	Cav2	ENSRNOG0000057713	caveolin 2 Source RGD Symbol Acc 620348
A11	SBR1138605	ENSRNOT00000084293.1	Cdkn2a	ENSRNOG0000059837	cyclin-dependent kinase inhibitor 2A Source RGD Symbol Acc 2323
A12	SBR1191088	ENSRNOT00000016756.4	Cln8	ENSRNOG0000012565	CLN8, transmembrane ER and ERGIC protein Source RGD Symbol Acc 1359518
B01	SBR1193638	ENSRNOT00000004074.5	Cox18	ENSRNOG0000003052	cytochrome c oxidase assembly factor COX18 Source RGD Symbol Acc 1559547
B02	SBR1162370	ENSRNOT00000013985.3	Cpt1b	ENSRNOG0000010438	carnitine palmitoyltransferase 1B Source RGD Symbol Acc 2397
B03	SBR1129520	ENSRNOT00000016954.3	Cpt2	ENSRNOG0000012443	carnitine palmitoyltransferase 2 Source RGD Symbol Acc 2398
B04	SBR1187963	ENSRNOT00000002482.7	Dnm1l	ENSRNOG0000001813	dynamin 1-like Source RGD Symbol Acc 620416
B05	SBR1143750	ENSRNOT00000083471.1	Fis1	ENSRNOG0000001420	fission, mitochondrial 1 Source RGD Symbol Acc 1306668
B06	SBR1105737	ENSRNOT00000076961.2	Timm10b	ENSRNOG0000050846	translocase of inner mitochondrial membrane 10B Source RGD Symbol Acc 71097
B07	SBR1095943	ENSRNOT00000033196.3	Gclc	ENSRNOG0000006302	glutamate-cysteine ligase, catalytic subunit Source RGD Symbol Acc 619868
B08	SBR1197356	ENSRNOT00000018343.4	Gclm	ENSRNOG0000013409	glutamate cysteine ligase, modifier subunit Source RGD Symbol Acc 619871
B09	SBR1185178	ENSRNOT00000073247.2	Gpx1	ENSRNOG0000048812	glutathione peroxidase 1 Source RGD Symbol Acc 2729
B10	SBR1196027	ENSRNOT00000082180.1	Grpel1	ENSRNOG0000006593	GrpE-like 1, mitochondrial Source RGD Symbol Acc 70947
B11	SBR1113341	ENSRNOT00000086310.1	Hsp90aa1	ENSRNOG0000059714	heat shock protein 90 alpha family class A member 1 Source RGD Symbol Acc 631409
B12	SBR1128872	ENSRNOT00000066589.2	Hspd1	ENSRNOG0000014525	heat shock protein family D member 1 Source RGD Symbol Acc 621314
C01	SBR1128331	ENSRNOT00000006406.6	Immp1l	ENSRNOG0000004829	inner mitochondrial membrane peptidase subunit 1 Source RGD Symbol Acc 1587441
C02	SBR1182078	ENSRNOT00000034261.5	Cox10	ENSRNOG0000024972	cytochrome c oxidase assembly factor heme A farnesyltransferase COX10 Source RGD Symbol Acc 1594623
C03	SBR1146003	ENSRNOT00000077170.1	Mfn1	ENSRNOG0000011057	mitofusin 1 Source RGD Symbol Acc 621460
C04	SBR1151248	ENSRNOT00000091132.1	LOC100911485	ENSRNOG0000046424	mitofusin 2 Source RGD Symbol Acc 628843
C05	SBR1110971	ENSRNOT00000085785.1	Mipep	ENSRNOG0000013876	mitochondrial intermediate peptidase Source RGD Symbol Acc 621680
C06	SBR1165108	ENSRNOT00000027579.4	Msto1	ENSRNOG0000020357	misato mitochondrial distribution and morphology regulator 1 Source RGD Symbol Acc 1306110
C07	SBR1184774	ENSRNOT00000002134.5	Mtx2	ENSRNOG0000001559	metaxin 2 Source RGD Symbol Acc 1306473
C08	SBR1148387	ENSRNOT00000018599.3	Nefl	ENSRNOG0000013658	neurofilament light Source RGD Symbol Acc 621458
C09	SBR1161494	ENSRNOT00000002343.7	Opa1	ENSRNOG0000001717	OPA1, mitochondrial dynamin like GTPase Source RGD Symbol Acc 708423
C10	SBR1112809	ENSRNOT00000025362.3	Pmaip1	ENSRNOG0000018770	phorbol-12-myristate-13-acetate-induced protein 1 Source RGD Symbol Acc 1359266
		ENSRNOT000000		ENSRNOG00	

Position	Assay	Name	Symbol	Ensembl ID	Description
C11	SBR1191958	006071.5	Ppargc1a	000004473	PPARG coactivator 1 alpha Source RGD Symbol Acc 620925
C12	SBR1101823	ENSRNOT00000027083.5	Rhot2	ENSRNOG0000019930	ras homolog family member T2 Source RGD Symbol Acc 727970
D01	SBR1174428	ENSRNOT00000006218.5	Rhot1	ENSRNOG0000004093	ras homolog family member T1 Source RGD Symbol Acc 1307023
D02	SBR1202515	ENSRNOT00000043959.4	Sfn	ENSRNOG0000033153	stratifin Source RGD Symbol Acc 1304729
D03	SBR1101841	ENSRNOT00000017770.5	Sh3glb1	ENSRNOG0000012957	SH3 domain -containing GRB2-like endophilin B1 Source RGD Symbol Acc 1304859
D04	SBR1203003	ENSRNOT00000054963.2	Slc25a10	ENSRNOG0000036693	solute carrier family 25 member 10 Source RGD Symbol Acc 621430
D05	SBR1099399	ENSRNOT00000036728.6	Slc25a12	ENSRNOG0000022922	solute carrier family 25 member 12 Source RGD Symbol Acc 1305181
D06	SBR1122221	ENSRNOT00000013817.7	Slc25a13	ENSRNOG0000009957	solute carrier family 25 member 13 Source RGD Symbol Acc 1565889
D07	SBR1146752	ENSRNOT00000009637.5	Slc25a14	ENSRNOG0000006871	solute carrier family 25 member 14 Source RGD Symbol Acc 621433
D08	SBR1116441	ENSRNOT00000078971.1	Slc25a15	ENSRNOG0000011881	solute carrier family 25 member 15 Source RGD Symbol Acc 1311488
D09	SBR1205766	ENSRNOT00000000438.7	Slc25a16	ENSRNOG0000000387	solute carrier family 25 member 16 Source RGD Symbol Acc 1311311
D10	SBR1196542	ENSRNOT00000025758.5	Slc25a17	ENSRNOG0000018920	solute carrier family 25 member 17 Source RGD Symbol Acc 1311498
D11	SBR1186421	ENSRNOT00000005292.4	Slc25a19	ENSRNOG0000003918	solute carrier family 25 member 19 Source RGD Symbol Acc 1359554
D12	SBR1192004	ENSRNOT00000027128.3	Slc25a2	ENSRNOG0000027359	solute carrier family 25 member 2 Source RGD Symbol Acc 1310667
E01	SBR1142250	ENSRNOT00000027520.4	Slc25a20	ENSRNOG0000020288	solute carrier family 25 member 20 Source RGD Symbol Acc 621443
E02	SBR1129277	ENSRNOT00000034342.5	Slc25a21	ENSRNOG0000008931	solute carrier family 25 member 21 Source RGD Symbol Acc 621444
E03	SBR1195919	ENSRNOT00000025240.6	LOC100911440	ENSRNOG0000049944	solute carrier family 25 member 22 Source RGD Symbol Acc 1307826
E04	SBR1140021	ENSRNOT00000086575.1	Slc25a23	ENSRNOG0000047781	solute carrier family 25 member 23 Source RGD Symbol Acc 1588586
E05	SBR1152762	ENSRNOT00000027765.5	Slc25a24	ENSRNOG0000020470	solute carrier family 25 member 24 Source RGD Symbol Acc 1311982
E06	SBR1170427	ENSRNOT00000020206.5	Slc25a25	ENSRNOG0000014338	solute carrier family 25 member 25 Source RGD Symbol Acc 628666
E07	SBR1206973	ENSRNOT00000052068.3	Slc25a27	ENSRNOG0000010592	solute carrier family 25, member 27 Source RGD Symbol Acc 620787
E08	SBR1204776	ENSRNOT00000080824.1	Slc25a3	ENSRNOG0000008289	solute carrier family 25 member 3 Source RGD Symbol Acc 70986
E09	SBR1157856	ENSRNOT00000001389.7	Slc25a30	ENSRNOG0000001052	solute carrier family 25, member 30 Source RGD Symbol Acc 1359702
E10	SBR1131953	ENSRNOT00000037120.5	Slc25a36	ENSRNOG0000013802	solute carrier family 25 member 36 Source RGD Symbol Acc 1591928
E11	SBR1113752	ENSRNOT00000021636.4	Slc25a37	ENSRNOG0000015495	solute carrier family 25 member 37 Source RGD Symbol Acc 1359361
E12	SBR1180410	ENSRNOT00000014704.4	Slc25a4	ENSRNOG0000010830	solute carrier family 25 member 4 Source RGD Symbol Acc 620352
F01	SBR1173244	ENSRNOT00000015085.6	Slc25a5	ENSRNOG0000039980	solute carrier family 25 member 5 Source RGD Symbol Acc 620353
F02	SBR1136942	ENSRNOT00000002885.6	Sod1	ENSRNOG0000002115	superoxide dismutase 1 Source RGD Symbol Acc 3731
F03	SBR1204773	ENSRNOT00000025794.4	Sod2	ENSRNOG0000019048	superoxide dismutase 2 Source RGD Symbol Acc 3732
F04	SBR1105956	ENSRNOT00000064190.3	Stard3	ENSRNOG0000042044	StAR-related lipid transfer domain containing 3 Source RGD Symbol Acc 1549768
F05	SBR1121508	ENSRNOT00000089778.1	Taz	ENSRNOG0000053023	tafazzin Source RGD Symbol Acc 1588532
F06	SBR1162876	ENSRNOT00000010619.4	Timm10	ENSRNOG0000007883	translocase of inner mitochondrial membrane 10 Source RGD Symbol Acc 621741
F07	SBR1148602	ENSRNOT00000009398.4	Timm17a	ENSRNOG0000007040	translocase of inner mitochondrial membrane 17A Source RGD Symbol Acc 3862
F08	SBR1101808	ENSRNOT00000066373.3	LOC100911130	ENSRNOG0000048613	mitochondrial import inner membrane translocase subunit Tim17-B-like Source RGD Symbol Acc 6494590
F09	SBR1115684	ENSRNOT00000010778.7	Timm22	ENSRNOG0000007988	translocase of inner mitochondrial membrane 22 Source RGD Symbol Acc 68406
F10	SBR1174386	ENSRNOT00000001409.4	Timm44	ENSRNOG0000001058	translocase of inner mitochondrial membrane 44 Source RGD Symbol Acc 3864

Position	Assay	Name	Symbol	Ensembl ID	Description
F11	SBR1197680	ENSRNOT00000015043.7	Timm8a1	ENSRNOG00000011226	translocase of inner mitochondrial membrane 8A1 Source RGD Symbol Acc 621801
F12	SBR1166379	ENSRNOT00000013188.5	Timm8b	ENSRNOG0000009888	translocase of inner mitochondrial membrane 8 homolog B Source RGD Symbol Acc 619840
G01	SBR1098039	ENSRNOT00000088790.1	Timm9	ENSRNOG0000008222	translocase of inner mitochondrial membrane 9 Source RGD Symbol Acc 621656
G02	SBR1140216	ENSRNOT00000019323.5	Tomm22	ENSRNOG0000014058	translocase of outer mitochondrial membrane 22 Source RGD Symbol Acc 1303260
G03	SBR1095644	ENSRNOT00000040156.4	Tomm34	ENSRNOG0000029799	translocase of outer mitochondrial membrane 34 Source RGD Symbol Acc 1309029
G04	SBR1182057	ENSRNOT00000025281.2	Tomm40	ENSRNOG0000018556	translocase of outer mitochondrial membrane 40 Source RGD Symbol Acc 1303022
G05	SBR1153269	ENSRNOT00000004641.3	Tomm40l	ENSRNOG0000003398	translocase of outer mitochondrial membrane 40 like Source RGD Symbol Acc 1562006
G06	SBR1163402	ENSRNOT00000002238.6	Tomm70	ENSRNOG0000001640	translocase of outer mitochondrial membrane 70 Source RGD Symbol Acc 1303049
G07	SBR1172374	ENSRNOT00000046490.3	Tp53	ENSRNOG0000010756	tumor protein p53 Source RGD Symbol Acc 3889
G08	SBR1101780	ENSRNOT00000014089.4	Tspo	ENSRNOG0000010549	translocator protein Source RGD Symbol Acc 2228
G09	SBR1129529	ENSRNOT00000004900.4	Ucp1	ENSRNOG0000003580	uncoupling protein 1 Source RGD Symbol Acc 3931
G10	SBR1104067	ENSRNOT00000024156.5	Ucp2	ENSRNOG0000017854	uncoupling protein 2 Source RGD Symbol Acc 3932
G11	SBR1125659	ENSRNOT00000024005.5	Ucp3	ENSRNOG0000017716	uncoupling protein 3 Source RGD Symbol Acc 3933
G12	SBR1095475	ENSRNOT00000013476.7	Uxt	ENSRNOG0000009893	ubiquitously-expressed, prefoldin-like chaperone Source RGD Symbol Acc 1359326
H01	SBR1220567	ENSRNOT00000042459.4	Actb	ENSRNOG00000034254	actin, beta Source RGD Symbol Acc 628837
H02	SBR1220568	ENSRNOT00000023017.5	B2m	ENSRNOG0000017123	beta-2 microglobulin Source RGD Symbol Acc 2189
H03	SBR1225377	ENSRNOT00000065935.3	Hprt1	ENSRNOG0000048561	hypoxanthine phosphoribosyltransferase 1 Source RGD Symbol Acc 2826
H04	SBR1122313	ENSRNOT00000017468.2	Ldha	ENSRNOG0000013009	lactate dehydrogenase A Source RGD Symbol Acc 2996
H05	SBR1220572	ENSRNOT00000018820.5	Rplp1	ENSRNOG0000013874	ribosomal protein lateral stalk subunit P1 Source RGD Symbol Acc 621774
H06	SBR1218555	Sybr_RGDC	RGDC	Sybr_RGDC	Rat Genomic DNA Contamination
H07	SBH1218551	Sybr_QIC	QIC	Sybr_QIC	QuantiNova Internal Control
H08	SBH1218551	Sybr_QIC	QIC	Sybr_QIC	QuantiNova Internal Control
H09	SBH1218551	Sybr_QIC	QIC	Sybr_QIC	QuantiNova Internal Control
H10	SBH1218550	Sybr_PPC	PPC	Sybr_PPC	Positive PCR Control
H11	SBH1218550	Sybr_PPC	PPC	Sybr_PPC	Positive PCR Control
H12	SBH1218550	Sybr_PPC	PPC	Sybr_PPC	Positive PCR Control

Related products

Product	Contents	Cat. no.
QuantiNova LNA PCR QC Panel	These panels are designed to assess the quality of RNA samples before characterization using QuantiNova LNA PCR Focus Panels; available in 96-well, 384-well, and Rotor-Disc 100 formats	249940
QuantiNova Reverse Transcription Kit (10)*	For 10 x 20 μ l reactions: 20 μ l 8x gDNA Removal Mix, 10 μ l Reverse Transcription Enzyme, 40 μ l Reverse Transcription Mix (containing RT primers), 20 μ l Internal Control RNA, 1.9 ml RNase-Free Water	205410
QuantiNova SYBR Green RT-PCR Kit (100)*	For 100 x 20 μ l reactions: 1 ml QuantiNova SYBR Green RT-PCR Master Mix, 20 μ l QuantiNova SYBR Green RT Mix, 20 μ l Internal Control RNA, 500 μ l Yellow Template Dilution Buffer, 250 μ l ROX Reference Dye, 1.9 μ l RNase-Free Water	208152
QuantiNova SYBR Green PCR Kit (100)*	For 100 x 20 μ l reactions: 1 ml 2x QuantiNova SYBR Green PCR Master Mix, 500 μ l QuantiNova Yellow Template Dilution Buffer, 250 μ l QN ROX Reference Dye, 1.9 ml Water	208052

*Larger kit sizes available.

The QuantiNova LNA PCR Focus Panels are intended for molecular biology applications. These products are not intended for the diagnosis, prevention or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, LNA®, QuantiNova®, Sample to Insight® (QIAGEN Group); SYBR® (Life Technologies Corp.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

09/2019 © 2019 QIAGEN, all rights reserved.