

QuantiNova® LNA® Probe PCR Focus Panels (Rotor-Gene® Format)

Human Endothelial Cell Biology

Cat. no. 249955 UPHS-015ZR

For study focus gene expression analysis

Shipping and storage

QuantiNova LNA Probe PCR Focus Panels are shipped at room temperature. Immediately upon receipt, they should be stored protected from light at 2–8°C for short term storage or at –30°C to –15°C for long time storage. Under these conditions, all components are stable for at least 12 months.

Note: Open the package and store the products appropriately immediately upon receipt.

For optimal performance, QuantiNova LNA Probe PCR Focus Panels should be used together with the QuantiNova Reverse Transcription Kit for cDNA synthesis and the QuantiNova Probe PCR Kit (Mastermix) for PCR.

Panel layout (Rotor-Gene): QuantiNova LNA Probe PCR Focus Panel

The 96 real-time assays in the Rotor-Gene format are located in wells 1–96 of the Rotor-Disc® (plate A1–A12=Rotor-Disc 1–12, plate B1–B12=Rotor-Disc 13–24, etc.). To maintain data analysis compatibility, wells 97–100 do not contain real-time assays but will contain master mix to account for weight balance. Refer to the QuantiNova LNA Probe PCR Handbook at www.qiagen.com for further details.

	1	2	3	4	5	6	7	8	9	10	11	12
A	ACE	ADAM17	AGT	AGTR1	ALOX5	ANGPT1	ANXA5	APOE	BAX	BCL2	BCL2L1	CALCA
B	CASP1	CASP3	CAV1	CCL2	CCL5	CDH5	CFLAR	COL18A1	CX3CL1	EDN1	EDN2	EDNRA
C	ENG	F2R	F3	FAS	FASLG	FGF1	FGF2	FLT1	FN1	HIF1A	HMOX1	ICAM1
D	IL11	IL1B	IL3	IL6	IL7	ITGA5	ITGAV	ITGB1	ITGB3	KDR	KIT	KLK3
E	MMP1	MMP2	MMP9	NOS3	NPPB	NPR1	OCLN	PDGFRA	PECAM1	PF4	PGF	PLAT
F	PLAU	PLG	PROCR	PTGIS	PTGS2	PTK2	SELE	SELL	SELPLG	SERPINE1	SOD1	SPHK1
G	TEK	TFPI	TGFB1	THBD	THBS1	TIMP1	TNF	TNFSF10	TYMP	VCAM1	VEGFA	VWF
H	ACTB	B2M	GAPDH	HPRT1	RPLP0	HGDC	QIC	QIC	QIC	PPC	PPC	PPC

Gene table: QuantiNova LNA Probe PCR Focus Panel

Position	Assay	Name	Symbol	Ensembl ID	Description
A01	UPFH1132210	ENST00000428043.5	ACE	ENSG00000159640	angiotensin I converting enzyme Source HGNC Symbol Acc HGNC 2707
A02	UPFH1132216	ENST00000497134.6	ADAM17	ENSG00000151694	ADAM metallopeptidase domain 17 Source HGNC Symbol Acc HGNC 195
A03	UPFH1132773	ENST00000366667.5	AGT	ENSG00000135744	angiotensinogen Source HGNC Symbol Acc HGNC 333
A04	UPFH1132221	ENST00000418473.6	AGTR1	ENSG00000144891	angiotensin II receptor type 1 Source HGNC Symbol Acc HGNC 336
A05	UPFH1132227	ENST00000612635.4	ALOX5	ENSG00000012779	arachidonate 5-lipoxygenase Source HGNC Symbol Acc HGNC 435
A06	UPFH1132231	ENST00000517746.6	ANGPT1	ENSG00000154188	angiopoietin 1 Source HGNC Symbol Acc HGNC 484
A07	UPFH1132775	ENST00000515017.5	ANXA5	ENSG00000164111	annexin A5 Source HGNC Symbol Acc HGNC 543
A08	UPFH0539669	ENST00000425718.1	APOE	ENSG00000130203	apolipoprotein E Source HGNC Symbol Acc HGNC 613
A09	UPFH0540159	ENST00000293288.12	BAX	ENSG00000087088	BCL2 associated X, apoptosis regulator Source HGNC Symbol Acc HGNC 959
A10	UPFH1132900	ENST00000333681.5	BCL2	ENSG00000171791	BCL2, apoptosis regulator Source HGNC Symbol Acc HGNC 990
A11	UPFH1132271	ENST00000376062.6	BCL2L1	ENSG00000171552	BCL2 like 1 Source HGNC Symbol Acc HGNC 992
A12	UPFH0317713	ENST00000331587.8	CALCA	ENSG00000110680	calcitonin related polypeptide alpha Source HGNC Symbol Acc HGNC 1437
B01	UPFH0285144	ENST00000436863.7	CASP1	ENSG00000137752	caspase 1 Source HGNC Symbol Acc HGNC 1499
B02	UPFH1132892	ENST00000523916.5	CASP3	ENSG00000164305	caspase 3 Source HGNC Symbol Acc HGNC 1504
B03	UPFH0192343	ENST00000393468.1	CAV1	ENSG00000105974	caveolin 1 Source HGNC Symbol Acc HGNC 1527
B04	UPFH1132783	ENST00000225831.4	CCL2	ENSG00000108691	C-C motif chemokine ligand 2 Source HGNC Symbol Acc HGNC 10618
B05	UPFH1132786	ENST00000603197.6	CCL5	ENSG00000271503	C-C motif chemokine ligand 5 Source HGNC Symbol Acc HGNC 10632
B06	UPFH1132306	ENST00000563425.2	CDH5	ENSG00000179776	cadherin 5 Source HGNC Symbol Acc HGNC 1764
B07	UPFH1132312	ENST00000479953.6	CFLAR	ENSG00000003402	CASP8 and FADD like apoptosis regulator Source HGNC Symbol Acc HGNC 1876
B08	UPFH1132325	ENST00000359759.8	COL18A1	ENSG00000182871	collagen type XVIII alpha 1 chain Source HGNC Symbol Acc HGNC 2195
B09	UPFH1132348	ENST00000006053.7	CX3CL1	ENSG00000006210	C-X3-C motif chemokine ligand 1 Source HGNC Symbol Acc HGNC 10647
B10	UPFH1132801	ENST00000379375.6	EDN1	ENSG00000078401	endothelin 1 Source HGNC Symbol Acc HGNC 3176
B11	UPFH1132378	ENST00000372587.5	EDN2	ENSG00000127129	endothelin 2 Source HGNC Symbol Acc HGNC 3177
B12	UPFH1132379	ENST00000511804.5	EDNRA	ENSG00000151617	endothelin receptor type A Source HGNC Symbol Acc HGNC 3179
C01	UPFH0535657	ENST00000344849.4	ENG	ENSG00000106991	endoglin Source HGNC Symbol Acc HGNC 3349
C02	UPFH1132392	ENST00000319211.5	F2R	ENSG00000181104	coagulation factor II thrombin receptor Source HGNC Symbol Acc HGNC 3537
C03	UPFH1132393	ENST00000334047.12	F3	ENSG00000117525	coagulation factor III, tissue factor Source HGNC Symbol Acc HGNC 3541
C04	UPFH1132395	ENST00000357339.6	FAS	ENSG00000026103	Fas cell surface death receptor Source HGNC Symbol Acc HGNC 11920
C05	UPFH1132396	ENST00000367721.3	FASLG	ENSG00000117560	Fas ligand Source HGNC Symbol Acc HGNC 11936
C06	UPFH0087975	ENST00000612258.4	FGF1	ENSG00000113578	fibroblast growth factor 1 Source HGNC Symbol Acc HGNC 3665
C07	UPFH0613093	ENST00000264498.7	FGF2	ENSG00000138685	fibroblast growth factor 2 Source HGNC Symbol Acc HGNC 3676
C08	UPFH1132400	ENST00000541932.5	FLT1	ENSG00000102755	fms related tyrosine kinase 1 Source HGNC Symbol Acc HGNC 3763
C09	UPFH0605066	ENST00000336916.8	FN1	ENSG00000115414	fibronectin 1 Source HGNC Symbol Acc HGNC 3778
C10	UPFH1132447	ENST00000394997.5	HIF1A	ENSG00000100644	hypoxia inducible factor 1 subunit alpha Source HGNC Symbol Acc HGNC 4910
		ENST00000216		ENSG000000	

Position	Assay	Name	Symbol	Ensembl ID	Description
C11	UPFH0433285	117.9	HMOX1	100292	heme oxygenase 1 Source HGNC Symbol Acc HGNC 5013
C12	UPFH1132462	ENST00000264832.8	ICAM1	ENSG00000090339	intercellular adhesion molecule 1 Source HGNC Symbol Acc HGNC 5344
D01	UPFH1132477	ENST0000058513.1	IL11	ENSG00000095752	interleukin 11 Source HGNC Symbol Acc HGNC 5966
D02	UPFH0163764	ENST00000263341.6	IL1B	ENSG00000125538	interleukin 1 beta Source HGNC Symbol Acc HGNC 5992
D03	UPFH0282899	ENST00000296870.2	IL3	ENSG00000164399	interleukin 3 Source HGNC Symbol Acc HGNC 6011
D04	UPFH1172910	ENST00000258743.10	IL6	ENSG00000136244	interleukin 6 Source HGNC Symbol Acc HGNC 6018
D05	UPFH1132812	ENST00000263851.9	IL7	ENSG00000104432	interleukin 7 Source HGNC Symbol Acc HGNC 6023
D06	UPFH1132497	ENST00000293379.9	ITGA5	ENSG00000161638	integrin subunit alpha 5 Source HGNC Symbol Acc HGNC 6141
D07	UPFH1132816	ENST00000261023.8	ITGAV	ENSG00000138448	integrin subunit alpha V Source HGNC Symbol Acc HGNC 6150
D08	UPFH1132498	ENST00000423113.5	ITGB1	ENSG00000150093	integrin subunit beta 1 Source HGNC Symbol Acc HGNC 6153
D09	UPFH1132500	ENST00000559488.5	ITGB3	ENSG00000259207	integrin subunit beta 3 Source HGNC Symbol Acc HGNC 6156
D10	UPFH0596732	ENST00000263923.5	KDR	ENSG00000128052	kinase insert domain receptor Source HGNC Symbol Acc HGNC 6307
D11	UPFH0545239	ENST00000288135.5	KIT	ENSG00000157404	KIT proto-oncogene receptor tyrosine kinase Source HGNC Symbol Acc HGNC 6342
D12	UPFH1132511	ENST00000617027.4	KLK3	ENSG00000142515	kallikrein related peptidase 3 Source HGNC Symbol Acc HGNC 6364
E01	UPFH0322484	ENST00000315274.7	MMP1	ENSG00000196611	matrix metallopeptidase 1 Source HGNC Symbol Acc HGNC 7155
E02	UPFH1132551	ENST00000437642.6	MMP2	ENSG00000087245	matrix metallopeptidase 2 Source HGNC Symbol Acc HGNC 7166
E03	UPFH0367626	ENST00000372330.3	MMP9	ENSG00000100985	matrix metallopeptidase 9 Source HGNC Symbol Acc HGNC 7176
E04	UPFH1132897	ENST00000297494.8	NOS3	ENSG00000164867	nitric oxide synthase 3 Source HGNC Symbol Acc HGNC 7876
E05	UPFH1132597	ENST00000376468.4	NPPB	ENSG00000120937	natriuretic peptide B Source HGNC Symbol Acc HGNC 7940
E06	UPFH1132933	ENST00000368680.4	NPR1	ENSG00000169418	natriuretic peptide receptor 1 Source HGNC Symbol Acc HGNC 7943
E07	UPFH1132600	ENST00000355237.6	OCLN	ENSG00000197822	occludin Source HGNC Symbol Acc HGNC 8104
E08	UPFH1132609	ENST00000508170.5	PDGFRA	ENSG00000134853	platelet derived growth factor receptor alpha Source HGNC Symbol Acc HGNC 8803
E09	UPFH1132613	ENST00000563924.6	PECAM1	ENSG00000261371	platelet and endothelial cell adhesion molecule 1 Source HGNC Symbol Acc HGNC 8823
E10	UPFH0614943	ENST00000296029.3	PF4	ENSG00000163737	platelet factor 4 Source HGNC Symbol Acc HGNC 8861
E11	UPFH1132616	ENST00000555567.6	PGF	ENSG00000119630	placental growth factor Source HGNC Symbol Acc HGNC 8893
E12	UPFH1132830	ENST00000352041.7	PLAT	ENSG00000104368	plasminogen activator, tissue type Source HGNC Symbol Acc HGNC 9051
F01	UPFH1132831	ENST00000446342.5	PLAU	ENSG00000122861	plasminogen activator, urokinase Source HGNC Symbol Acc HGNC 9052
F02	UPFH0091404	ENST00000308192.13	PLG	ENSG00000122194	plasminogen Source HGNC Symbol Acc HGNC 9071
F03	UPFH1132638	ENST00000216968.5	PROCR	ENSG00000101000	protein C receptor Source HGNC Symbol Acc HGNC 9452
F04	UPFH1132641	ENST00000244043.5	PTGIS	ENSG00000124212	prostaglandin I2 synthase Source HGNC Symbol Acc HGNC 9603
F05	UPFH1132642	ENST00000367468.10	PTGS2	ENSG00000073756	prostaglandin-endoperoxide synthase 2 Source HGNC Symbol Acc HGNC 9605
F06	UPFH1132643	ENST00000522684.5	PTK2	ENSG00000169398	protein tyrosine kinase 2 Source HGNC Symbol Acc HGNC 9611
F07	UPFH1132838	ENST00000367777.5	SELE	ENSG00000007908	selectin E Source HGNC Symbol Acc HGNC 10718
F08	UPFH0444664	ENST00000236147.5	SELL	ENSG00000188404	selectin L Source HGNC Symbol Acc HGNC 10720
F09	UPFH0032702	ENST00000550948.1	SELPLG	ENSG00000110876	selectin P ligand Source HGNC Symbol Acc HGNC 10722
F10	UPFH0384736	ENST00000223095.4	SERPINE1	ENSG00000106366	serpin family E member 1 Source HGNC Symbol Acc HGNC 8583

Position	Assay	Name	Symbol	Ensembl ID	Description
F11	UPFH0220920	ENST00000270142.10	SOD1	ENSG00000142168	superoxide dismutase 1 Source HGNC Symbol Acc HGNC 11179
F12	UPFH0560569	ENST00000590959.5	SPHK1	ENSG00000176170	sphingosine kinase 1 Source NCBI gene Acc 8877
G01	UPFH1132711	ENST00000380036.9	TEK	ENSG00000120156	TEK receptor tyrosine kinase Source HGNC Symbol Acc HGNC 11724
G02	UPFH1132716	ENST00000392365.5	TFPI	ENSG00000003436	tissue factor pathway inhibitor Source HGNC Symbol Acc HGNC 11760
G03	UPFH0193430	ENST00000221930.5	TGFB1	ENSG00000105329	transforming growth factor beta 1 Source NCBI gene Acc 7040
G04	UPFH0356052	ENST00000377103.2	THBD	ENSG00000178726	thrombomodulin Source HGNC Symbol Acc HGNC 11784
G05	UPFH1132847	ENST00000260356.5	THBS1	ENSG00000137801	thrombospondin 1 Source HGNC Symbol Acc HGNC 11785
G06	UPFH1132725	ENST00000456754.6	TIMP1	ENSG00000102265	TIMP metalloproteinase inhibitor 1 Source HGNC Symbol Acc HGNC 11820
G07	UPFH1132978	ENST00000449264.3	TNF	ENSG00000232810	tumor necrosis factor Source HGNC Symbol Acc HGNC 11892
G08	UPFH1132733	ENST00000241261.7	TNFSF10	ENSG00000121858	TNF superfamily member 10 Source HGNC Symbol Acc HGNC 11925
G09	UPFH1132746	ENST00000252029.8	TYMP	ENSG00000025708	thymidine phosphorylase Source HGNC Symbol Acc HGNC 3148
G10	UPFH1132856	ENST00000294728.7	VCAM1	ENSG00000162692	vascular cell adhesion molecule 1 Source HGNC Symbol Acc HGNC 12663
G11	UPFH0281656	ENST00000425836.6	VEGFA	ENSG00000112715	vascular endothelial growth factor A Source HGNC Symbol Acc HGNC 12680
G12	UPFH0533670	ENST00000261405.9	VWF	ENSG00000110799	von Willebrand factor Source HGNC Symbol Acc HGNC 12726
H01	UPFH1132936	ENST00000646664.1	ACTB	ENSG00000075624	actin beta Source HGNC Symbol Acc HGNC 132
H02	UPFH1132937	ENST00000544417.5	B2M	ENSG00000166710	beta-2-microglobulin Source HGNC Symbol Acc HGNC 914
H03	UPFH1132938	ENST00000229239.10	GAPDH	ENSG00000111640	glyceraldehyde-3-phosphate dehydrogenase Source HGNC Symbol Acc HGNC 4141
H04	UPFH1132939	ENST00000298556.8	HPRT1	ENSG00000165704	hypoxanthine phosphoribosyltransferase 1 Source HGNC Symbol Acc HGNC 5157
H05	UPFH1132941	ENST00000392514.9	RPLP0	ENSG00000089157	ribosomal protein lateral stalk subunit P0 Source HGNC Symbol Acc HGNC 10371
H06	UPFH1126608	UPL_HGDC	HGDC	UPL_HGDC	Human Genomic DNA Contamination
H07	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H08	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H09	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H10	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control
H11	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control
H12	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control

Related products

Product	Contents	Cat. no.
QuantiNova LNA Probe PCR QC Panel	These panels are designed to assess the quality of RNA samples before characterization using QuantiNova LNA Probe PCR Focus Panels; available in 96-well, 384-well, and Rotor-Disc 100 formats	249945
QuantiNova Reverse Transcription Kit (10)*	For 10 x 20 μ l reactions: 20 μ l 8x gDNA Removal Mix, 10 μ l Reverse Transcription Enzyme, 40 μ l Reverse Transcription Mix (containing RT primers), 20 μ l Internal Control RNA, 1.9 ml RNase-Free Water	205410
QuantiNova Probe RT-PCR Kit (100)*	For 100 x 20 μ l reactions: 1 ml QuantiNova Probe RT-PCR Master Mix, 20 μ l QuantiNova Probe RT Mix, 20 μ l Internal Control RNA, 500 μ l Yellow Template Dilution Buffer, 250 μ l ROX Reference Dye, 1.9 μ l RNase-Free Water	208352
QuantiNova Probe PCR Kit (100)*	For 100 x 20 μ l reactions: 1 ml 2x QuantiNova Probe PCR Master Mix, 500 μ l QuantiNova Yellow Template Dilution Buffer, 250 μ l QN ROX Reference Dye, 1.9 ml Water	208252

*Larger kit sizes available.

The QuantiNova LNA Probe PCR Focus Panels are intended for molecular biology applications. These products are not intended for the diagnosis, prevention or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, LNA®, QuantiNova®, Sample to Insight® (QIAGEN Group); SYBR® (Life Technologies Corp.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

09/2019 © 2019 QIAGEN, all rights reserved.