

QuantiNova® LNA® Probe PCR Focus Panels (Rotor-Gene® Format)

Rat Fatty Acid Metabolism

Cat. no. 249955 UPRN-007ZR

For study focus gene expression analysis

Shipping and storage

QuantiNova LNA Probe PCR Focus Panels are shipped at room temperature. Immediately upon receipt, they should be stored protected from light at 2–8°C for short term storage or at –30°C to –15°C for long time storage. Under these conditions, all components are stable for at least 12 months.

Note: Open the package and store the products appropriately immediately upon receipt.

For optimal performance, QuantiNova LNA Probe PCR Focus Panels should be used together with the QuantiNova Reverse Transcription Kit for cDNA synthesis and the QuantiNova Probe PCR Kit (Mastermix) for PCR.

Panel layout (Rotor-Gene): QuantiNova LNA Probe PCR Focus Panel

The 96 real-time assays in the Rotor-Gene format are located in wells 1–96 of the Rotor-Disc® (plate A1–A12=Rotor-Disc 1–12, plate B1–B12=Rotor-Disc 13–24, etc.). To maintain data analysis compatibility, wells 97–100 do not contain real-time assays but will contain master mix to account for weight balance. Refer to the QuantiNova LNA Probe PCR Handbook at www.qiagen.com for further details.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Acaa1a	Acaa2	Acad10	Acad11	Acad9	Acadl	Acadm	Acads	Acadslb	Acadvl	Acat1	Acat2
B	Acat12	Acat2	Acat3	Acat7	Acat8	Acat9	Acox1	Acox2	Acox3	Acsbg1	Acsbg2	Acs1
C	Acs13	Acs14	Acs15	Acs16	Acsm2	Acsm3	Acsm4	Acsm5	Aldh2	Bdh1	Bdh2	Cpt1a
D	Cpt1b	Cpt1c	Cpt2	Crot	Crot	Decr1	Decr2	Echs1	Eci3	Ehhadh	Fabp1	Fabp2
E	Fabp3	Fabp4	Fabp5	Fabp6	Fabp7	Gcdh	Gk	Gk2	Gpd1	Gpd2	Hadh	Hmgcl
F	Hmgcs1	Hmgcs2	Lipe	Lpl	Mce	Mmut	Oxcl	Pecr	Ppa1	Prkaa1	Prkaa2	Prkab1
G	Prkab2	Prkaca	Prkacb	Kmt2d	Prkag2	Prkag3	Slc27a1	Slc27a2	Slc27a3	Slc27a4	Slc27a5	Slc27a6
H	Actb	B2m	Hprt1	Ldha	Rplp1	RGDC	QIC	QIC	QIC	PPC	PPC	PPC

Gene table: QuantiNova LNA Probe PCR Focus Panel

Position	Assay	Name	Symbol	Ensembl ID	Description
A01	UPFR1034645	ENSRNOT00000045049.5	Acaa1a	ENSRNOG0000032908	acetyl-CoA acyltransferase 1A Source RGD Symbol Acc 67379
A02	UPFR1053491	ENSRNOT00000085601.1	Acaa2	ENSRNOG0000013766	acetyl-CoA acyltransferase 2 Source RGD Symbol Acc 620482
A03	UPFR1013600	ENSRNOT00000073125.3	Acad10	ENSRNOG0000037815	acyl-CoA dehydrogenase family, member 10 Source RGD Symbol Acc 1310159
A04	UPFR1077165	ENSRNOT00000081281.1	Acad11	ENSRNOG0000010940	acyl-CoA dehydrogenase family, member 11 Source RGD Symbol Acc 1306270
A05	UPFR1029190	ENSRNOT00000089353.1	Acad9	ENSRNOG0000014178	acyl-CoA dehydrogenase family, member 9 Source RGD Symbol Acc 727973
A06	UPFR1080148	ENSRNOT00000017686.7	Acadl	ENSRNOG0000012966	acyl-CoA dehydrogenase, long chain Source RGD Symbol Acc 2011
A07	UPFR1090017	ENSRNOT00000013238.3	Acadm	ENSRNOG0000009845	acyl-CoA dehydrogenase medium chain Source RGD Symbol Acc 2012
A08	UPFR1113351	ENSRNOT00000001556.4	Acads	ENSRNOG0000001177	acyl-CoA dehydrogenase short chain Source RGD Symbol Acc 620514
A09	UPFR1033662	ENSRNOT00000027999.4	Acadsb	ENSRNOG0000020624	acyl-CoA dehydrogenase, short/branched chain Source RGD Symbol Acc 2013
A10	UPFR1056845	ENSRNOT00000024973.6	Acadvl	ENSRNOG0000018114	acyl-CoA dehydrogenase, very long chain Source RGD Symbol Acc 2014
A11	UPFR1024227	ENSRNOT00000010573.5	Acat1	ENSRNOG0000007862	acetyl-CoA acetyltransferase 1 Source RGD Symbol Acc 2016
A12	UPFR1013855	ENSRNOT00000077707.1	Acat2	ENSRNOG0000019189	acetyl-CoA acetyltransferase 2 Source RGD Symbol Acc 1359366
B01	UPFR1025316	ENSRNOT00000078919.1	Acot12	ENSRNOG0000061414	acyl-CoA thioesterase 12 Source RGD Symbol Acc 619752
B02	UPFR1101630	ENSRNOT00000013515.5	Acot2	ENSRNOG0000010134	acyl-CoA thioesterase 2 Source RGD Symbol Acc 621055
B03	UPFR1046042	ENSRNOT00000077555.1	Acot3	ENSRNOG0000053460	acyl-CoA thioesterase 3 Source RGD Symbol Acc 1564089
B04	UPFR1095676	ENSRNOT00000000079.8	Acot7	ENSRNOG0000010580	acyl-CoA thioesterase 7 Source RGD Symbol Acc 628856
B05	UPFR1033619	ENSRNOT00000020740.6	Acot8	ENSRNOG0000015187	acyl-CoA thioesterase 8 Source RGD Symbol Acc 70368
B06	UPFR1070636	ENSRNOT00000005033.5	Acot9	ENSRNOG0000003782	acyl-CoA thioesterase 9 Source RGD Symbol Acc 1359188
B07	UPFR1092066	ENSRNOT00000046754.5	Acox1	ENSRNOG0000008755	acyl-CoA oxidase 1 Source RGD Symbol Acc 619757
B08	UPFR1033324	ENSRNOT00000010260.4	Acox2	ENSRNOG0000007378	acyl-CoA oxidase 2 Source RGD Symbol Acc 628684
B09	UPFR1066235	ENSRNOT00000049798.3	Acox3	ENSRNOG0000008474	acyl-CoA oxidase 3, pristanoyl Source RGD Symbol Acc 69245
B10	UPFR1114754	ENSRNOT00000088208.1	Acsbg1	ENSRNOG0000011381	acyl-CoA synthetase bubblegum family member 1 Source RGD Symbol Acc 708557
B11	UPFR1022555	ENSRNOT00000074555.2	Acsbg2	ENSRNOG0000045947	acyl-CoA synthetase bubblegum family member 2 Source RGD Symbol Acc 1588580
B12	UPFR1104275	ENSRNOT00000089501.1	Acs1	ENSRNOG0000010633	acyl-CoA synthetase long-chain family member 1 Source RGD Symbol Acc 2015
C01	UPFR1023034	ENSRNOT00000020161.2	Acs13	ENSRNOG0000014718	acyl-CoA synthetase long-chain family member 3 Source RGD Symbol Acc 70552
C02	UPFR1082627	ENSRNOT00000026057.6	Acs14	ENSRNOG0000019180	acyl-CoA synthetase long-chain family member 4 Source RGD Symbol Acc 69401
C03	UPFR1034409	ENSRNOT00000022126.6	Acs15	ENSRNOG0000016265	acyl-CoA synthetase long-chain family member 5 Source RGD Symbol Acc 69402
C04	UPFR1061142	ENSRNOT00000058817.5	Acs16	ENSRNOG0000026745	acyl-CoA synthetase long-chain family member 6 Source RGD Symbol Acc 69403
C05	UPFR1030059	ENSRNOT00000084260.1	Acs2	ENSRNOG0000042084	acyl-CoA synthetase medium-chain family member 1 Source RGD Symbol Acc 1306813
C06	UPFR1076274	ENSRNOT00000020039.6	Acs3	ENSRNOG0000032246	acyl-CoA synthetase medium-chain family member 3 Source RGD Symbol Acc 62086
C07	UPFR1049002	ENSRNOT00000034610.4	Acs4	ENSRNOG0000014726	acyl-CoA synthetase medium-chain family member 4 Source RGD Symbol Acc 727928
C08	UPFR1108543	ENSRNOT00000046025.4	Acs5	ENSRNOG0000031211	acyl-CoA synthetase medium-chain family member 5 Source RGD Symbol Acc 1309578
C09	UPFR1023720	ENSRNOT00000001816.6	Aldh2	ENSRNOG0000001344	aldehyde dehydrogenase 2 family (mitochondrial) Source RGD Symbol Acc 69219
C10	UPFR1086005	ENSRNOT00000002366.5	Bdh1	ENSRNOG0000001736	3-hydroxybutyrate dehydrogenase 1 Source RGD Symbol Acc 620131
		ENSRNOT000000		ENSRNOG00	

Position	Assay	Name	Symbol	Ensembl ID	Description
C11	UPFR1064382	019507.5	Bdh2	000014490	3-hydroxybutyrate dehydrogenase 2 Source RGD Symbol Acc 1309898
C12	UPFR1101909	ENSRNOT00000 019652.3	Cpt1a	ENSRNOG00 000014254	carnitine palmitoyltransferase 1A Source RGD Symbol Acc 2396
D01	UPFR1021462	ENSRNOT00000 013985.3	Cpt1b	ENSRNOG00 000010438	carnitine palmitoyltransferase 1B Source RGD Symbol Acc 2397
D02	UPFR1052989	ENSRNOT00000 035908.4	Cpt1c	ENSRNOG00 000026163	carnitine palmitoyltransferase 1c Source RGD Symbol Acc 1305384
D03	UPFR1106957	ENSRNOT00000 016954.3	Cpt2	ENSRNOG00 000012443	carnitine palmitoyltransferase 2 Source RGD Symbol Acc 2398
D04	UPFR1065455	ENSRNOT00000 086717.1	Crat	ENSRNOG00 000018145	carnitine O-acetyltransferase Source RGD Symbol Acc 1303031
D05	UPFR1073310	ENSRNOT00000 090968.1	Crot	ENSRNOG00 000006779	carnitine O-octanoyltransferase Source RGD Symbol Acc 70908
D06	UPFR1028455	ENSRNOT00000 011330.6	Decr1	ENSRNOG00 000008236	2,4-dienoyl-CoA reductase 1 Source RGD Symbol Acc 70999
D07	UPFR1023701	ENSRNOT00000 077826.1	Decr2	ENSRNOG00 000050424	2,4-dienoyl-CoA reductase 2 Source RGD Symbol Acc 71002
D08	UPFR1017313	ENSRNOT00000 025446.6	Echs1	ENSRNOG00 000018522	enoyl-CoA hydratase, mitochondrial-like Source RGD Symbol Acc 6494098
D09	UPFR1076302	ENSRNOT00000 022022.6	Eci3	ENSRNOG00 000029549	enoyl-CoA delta isomerase 2 Source RGD Symbol Acc 1359427
D10	UPFR1042056	ENSRNOT00000 002410.4	Ehhadh	ENSRNOG00 000001770	enoyl-CoA hydratase and 3-hydroxyacyl CoA dehydrogenase Source RGD Symbol Acc 621441
D11	UPFR1068726	ENSRNOT00000 008840.5	Fabp1	ENSRNOG00 000006675	fatty acid binding protein 1 Source RGD Symbol Acc 2590
D12	UPFR1059794	ENSRNOT00000 029871.3	Fabp2	ENSRNOG00 000024947	fatty acid binding protein 2 Source RGD Symbol Acc 2591
E01	UPFR1083061	ENSRNOT00000 017325.5	Fabp3	ENSRNOG00 000012879	fatty acid binding protein 3 Source RGD Symbol Acc 69048
E02	UPFR1043524	ENSRNOT00000 014701.7	Fabp4	ENSRNOG00 000010805	fatty acid binding protein 4 Source RGD Symbol Acc 69309
E03	UPFR1101216	ENSRNOT00000 075493.2	Fabp5	ENSRNOG00 000049075	fatty acid binding protein 5 Source RGD Symbol Acc 70997
E04	UPFR1083830	ENSRNOT00000 083910.1	Fabp6	ENSRNOG00 000003902	fatty acid binding protein 6 Source RGD Symbol Acc 2592
E05	UPFR1075327	ENSRNOT00000 001079.5	Fabp7	ENSRNOG00 000000814	fatty acid binding protein 7 Source RGD Symbol Acc 69312
E06	UPFR1024909	ENSRNOT00000 004570.5	Gcdh	ENSRNOG00 000003307	glutaryl-CoA dehydrogenase Source RGD Symbol Acc 1308829
E07	UPFR1054947	ENSRNOT00000 064340.3	Gk	ENSRNOG00 000034116	glycerol kinase Source RGD Symbol Acc 70893
E08	UPFR1054058	ENSRNOT00000 044244.4	Gk2	ENSRNOG00 000029397	glycerol kinase 2 Source RGD Symbol Acc 1303263
E09	UPFR1095776	ENSRNOT00000 087662.1	Gpd1	ENSRNOG00 000056457	glycerol-3-phosphate dehydrogenase 1 Source RGD Symbol Acc 621381
E10	UPFR1065011	ENSRNOT00000 088589.1	Gpd2	ENSRNOG00 000033824	glycerol-3-phosphate dehydrogenase 2 Source RGD Symbol Acc 2726
E11	UPFR1051970	ENSRNOT00000 038649.6	Hadha	ENSRNOG00 000024629	hydroxyacyl-CoA dehydrogenase trifunctional multienzyme complex subunit alpha Source RGD Symbol Acc 620512
E12	UPFR1097191	ENSRNOT00000 012853.6	Hmgcl	ENSRNOG00 000009422	3-hydroxy-3-methylglutaryl-CoA lyase Source RGD Symbol Acc 620554
F01	UPFR1013564	ENSRNOT00000 023554.6	Hmgcs1	ENSRNOG00 000016552	3-hydroxy-3-methylglutaryl-CoA synthase 1 Source RGD Symbol Acc 70970
F02	UPFR1099961	ENSRNOT00000 026121.7	Hmgcs2	ENSRNOG00 000019120	3-hydroxy-3-methylglutaryl-CoA synthase 2 Source RGD Symbol Acc 2804
F03	UPFR1057468	ENSRNOT00000 027910.2	Lipe	ENSRNOG00 000020546	lipase E, hormone sensitive type Source RGD Symbol Acc 3010
F04	UPFR1016556	ENSRNOT00000 016543.3	Lpl	ENSRNOG00 000012181	lipoprotein lipase Source RGD Symbol Acc 3017
F05	UPFR1085181	ENSRNOT00000 021884.6	Mcee	ENSRNOG00 000016327	methylmalonyl CoA epimerase Source RGD Symbol Acc 1309966
F06	UPFR1076744	ENSRNOT00000 075279.2	Mmut	ENSRNOG00 000050843	methylmalonyl-CoA mutase Source RGD Symbol Acc 1587662
F07	UPFR1066906	ENSRNOT00000 065313.3	Oxct1	ENSRNOG00 000043094	3-oxoacid CoA transferase 1 Source RGD Symbol Acc 1584008
F08	UPFR1101536	ENSRNOT00000 089283.1	Pecr	ENSRNOG00 000055295	peroxisomal trans-2-enoyl-CoA reductase Source RGD Symbol Acc 70925
F09	UPFR1054159	ENSRNOT00000 000674.6	Ppa1	ENSRNOG00 000000557	pyrophosphatase (inorganic) 1 Source RGD Symbol Acc 1589773
F10	UPFR1073875	ENSRNOT00000 017626.5	Prkaa1	ENSRNOG00 000012799	protein kinase AMP-activated catalytic subunit alpha 1 Source RGD Symbol Acc 3387

Position	Assay	Name	Symbol	Ensembl ID	Description
F11	UPFR1026650	ENSRNOT00000010680.3	Prkaa2	ENSRNOG0000007706	protein kinase AMP-activated catalytic subunit alpha 2 Source RGD Symbol Acc 620893
F12	UPFR1066685	ENSRNOT00000001508.3	Prkab1	ENSRNOG0000001142	protein kinase AMP-activated non-catalytic subunit beta 1 Source RGD Symbol Acc 71057
G01	UPFR1070754	ENSRNOT00000086588.1	Prkab2	ENSRNOG0000018166	protein kinase AMP-activated non-catalytic subunit beta 2 Source RGD Symbol Acc 620905
G02	UPFR1048640	ENSRNOT00000041717.4	Prkaca	ENSRNOG0000005257	protein kinase cAMP-activated catalytic subunit alpha Source RGD Symbol Acc 3389
G03	UPFR1055817	ENSRNOT00000068739.1	Prkacb	ENSRNOG0000004978	protein kinase cAMP-activated catalytic subunit beta Source RGD Symbol Acc 1310574
G04	UPFR1057561	ENSRNOT00000077502.1	Kmt2d	ENSRNOG00000061499	protein kinase AMP-activated non-catalytic subunit gamma 1 Source RGD Symbol Acc 3388
G05	UPFR1041033	ENSRNOT00000067183.3	Prkag2	ENSRNOG0000009085	protein kinase AMP-activated non-catalytic subunit gamma 2 Source RGD Symbol Acc 727782
G06	UPFR1119563	ENSRNOT00000023365.5	Prkag3	ENSRNOG0000017248	protein kinase AMP-activated non-catalytic subunit gamma 3 Source RGD Symbol Acc 1308698
G07	UPFR1014919	ENSRNOT00000024659.4	Slc27a1	ENSRNOG0000018170	solute carrier family 27 member 1 Source RGD Symbol Acc 620927
G08	UPFR1017946	ENSRNOT00000067700.3	Slc27a2	ENSRNOG0000010128	solute carrier family 27 member 2 Source RGD Symbol Acc 71103
G09	UPFR1116141	ENSRNOT00000020815.5	Slc27a3	ENSRNOG0000015421	solute carrier family 27 member 3 Source RGD Symbol Acc 1310605
G10	UPFR1026149	ENSRNOT00000019500.6	Slc27a4	ENSRNOG0000014369	solute carrier family 27 member 4 Source RGD Symbol Acc 1307383
G11	UPFR1095460	ENSRNOT00000026575.4	Slc27a5	ENSRNOG0000019626	solute carrier family 27 member 5 Source RGD Symbol Acc 708535
G12	UPFR1020985	ENSRNOT00000077443.1	Slc27a6	ENSRNOG0000058193	solute carrier family 27 member 6 Source RGD Symbol Acc 1311977
H01	UPFR1132952	ENSRNOT00000080216.1	Actb	ENSRNOG00000034254	actin, beta Source RGD Symbol Acc 628837
H02	UPFR1132953	ENSRNOT00000023017.5	B2m	ENSRNOG0000017123	beta-2 microglobulin Source RGD Symbol Acc 2189
H03	UPFR1132959	ENSRNOT00000065935.3	Hprt1	ENSRNOG0000048561	hypoxanthine phosphoribosyltransferase 1 Source RGD Symbol Acc 2826
H04	UPFR1018740	ENSRNOT00000017468.2	Ldha	ENSRNOG0000013009	lactate dehydrogenase A Source RGD Symbol Acc 2996
H05	UPFR1132958	ENSRNOT00000018820.5	Rplp1	ENSRNOG0000013874	ribosomal protein lateral stalk subunit P1 Source RGD Symbol Acc 621774
H06	UPFR1126610	UPL_RGDC	RGDC	UPL_RGDC	Rat Genomic DNA Contamination
H07	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H08	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H09	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H10	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control
H11	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control
H12	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control

Related products

Product	Contents	Cat. no.
QuantiNova LNA Probe PCR QC Panel	These panels are designed to assess the quality of RNA samples before characterization using QuantiNova LNA Probe PCR Focus Panels; available in 96-well, 384-well, and Rotor-Disc 100 formats	249945
QuantiNova Reverse Transcription Kit (10)*	For 10 x 20 µl reactions: 20 µl 8x gDNA Removal Mix, 10 µl Reverse Transcription Enzyme, 40 µl Reverse Transcription Mix (containing RT primers), 20 µl Internal Control RNA, 1.9 ml RNase-Free Water	205410
QuantiNova Probe RT-PCR Kit (100)*	For 100 x 20 µl reactions: 1 ml QuantiNova Probe RT-PCR Master Mix, 20 µl QuantiNova Probe RT Mix, 20 µl Internal Control RNA, 500 µl Yellow Template Dilution Buffer, 250 µl ROX Reference Dye, 1.9 µl RNase-Free Water	208352
QuantiNova Probe PCR Kit (100)*	For 100 x 20 µl reactions: 1 ml 2x QuantiNova Probe PCR Master Mix, 500 µl QuantiNova Yellow Template Dilution Buffer, 250 µl QN ROX Reference Dye, 1.9 ml Water	208252

*Larger kit sizes available.

The QuantiNova LNA Probe PCR Focus Panels are intended for molecular biology applications. These products are not intended for the diagnosis, prevention or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, LNA®, QuantiNova®, Sample to Insight® (QIAGEN Group); SYBR® (Life Technologies Corp.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

09/2019 © 2019 QIAGEN, all rights reserved.