

QuantiNova® LNA® Probe PCR Focus Panels (Rotor-Gene® Format)

Rat Osteogenesis

Cat. no. 249955 UPRN-026ZR

For study focus gene expression analysis

Shipping and storage

QuantiNova LNA Probe PCR Focus Panels are shipped at room temperature. Immediately upon receipt, they should be stored protected from light at 2–8°C for short term storage or at –30°C to –15°C for long time storage. Under these conditions, all components are stable for at least 12 months.

Note: Open the package and store the products appropriately immediately upon receipt.

For optimal performance, QuantiNova LNA Probe PCR Focus Panels should be used together with the QuantiNova Reverse Transcription Kit for cDNA synthesis and the QuantiNova Probe PCR Kit (Mastermix) for PCR.

Panel layout (Rotor-Gene): QuantiNova LNA Probe PCR Focus Panel

The 96 real-time assays in the Rotor-Gene format are located in wells 1–96 of the Rotor-Disc® (plate A1–A12=Rotor-Disc 1–12, plate B1–B12=Rotor-Disc 13–24, etc.). To maintain data analysis compatibility, wells 97–100 do not contain real-time assays but will contain master mix to account for weight balance. Refer to the QuantiNova LNA Probe PCR Handbook at www.qiagen.com for further details.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Acvr1	AABR0703463 2.1	Alpl	Anxa5	Bglap	Bgn	Bmp1	Bmp2	Bmp3	Bmp4	Bmp5	Bmp6
B	Bmp7	Bmpr1a	Bmpr1b	Bmpr2	RGD1565355	Cdh11	Chrd	Col10a1	Col14a1	Col1a1	Col1a2	Col2a1
C	Col3a1	Col4a1	Col5a1	Col6a1	Comp	Csf1	Csf2	Csf3	Ctsk	Dlx5	Egf	Fgf1
D	Fgf2	Fgfr1	Fgfr2	Fth1	Fn1	Gdf10	Gli1	Icam1	Igf1	Igf1r	Ilh	Itga2
E	Itga3	Itgad	Itgav	Itgb1	Mmp10	Mmp2	Mmp8	Mmp9	Nfkb1	Nog	Pdgfra	Phex
F	Runx2	Serpinh1	Smad1	Smad2	Smad3	Smad4	Smad5	Sost	Sox9	Sp7	Spp1	Tgfb1
G	Tgfb2	Tgfb3	Tgfb1	Tgfb2	Tgfb3	LOC1036943 80	LOC1036899 74	Twist1	Vcam1	Vdr	Vegfa	Vegfb
H	Actb	B2m	Hprt1	Ldha	Rplp1	RGDC	QIC	QIC	QIC	PPC	PPC	PPC

Gene table: QuantiNova LNA Probe PCR Focus Panel

Position	Assay	Name	Symbol	Ensembl ID	Description
A01	UPFR1015339	ENSRNOT00000087071.1	Acvr1	ENSRNOG0000005033	activin A receptor type 1 Source RGD Symbol Acc 620200
A02	UPFR1058823	ENSRNOT00000078203.1	AABR07034632.1	ENSRNOG0000038370	
A03	UPFR1066331	ENSRNOT00000078851.1	Alpl	ENSRNOG0000013954	alkaline phosphatase, biomineralization associated Source RGD Symbol Acc 2100
A04	UPFR1079810	ENSRNOT00000019554.6	Anxa5	ENSRNOG0000014453	annexin A5 Source RGD Symbol Acc 2120
A05	UPFR1053867	ENSRNOT00000026530.3	Bglap	ENSRNOG0000019607	bone gamma-carboxyglutamate protein Source RGD Symbol Acc 2206
A06	UPFR1062701	ENSRNOT00000085201.1	Bgn	ENSRNOG0000055962	biglycan Source RGD Symbol Acc 2207
A07	UPFR1057170	ENSRNOT00000046054.4	Bmp1	ENSRNOG0000010890	bone morphogenetic protein 1 Source RGD Symbol Acc 620739
A08	UPFR1023820	ENSRNOT00000028904.4	Bmp2	ENSRNOG0000021276	bone morphogenetic protein 2 Source RGD Symbol Acc 2211
A09	UPFR1067060	ENSRNOT00000038872.3	Bmp3	ENSRNOG0000002381	bone morphogenetic protein 3 Source RGD Symbol Acc 2212
A10	UPFR1013450	ENSRNOT00000083268.1	Bmp4	ENSRNOG0000009694	bone morphogenetic protein 4 Source RGD Symbol Acc 2213
A11	UPFR1079381	ENSRNOT00000014846.5	Bmp5	ENSRNOG0000010917	bone morphogenetic protein 5 Source RGD Symbol Acc 1305979
A12	UPFR1048825	ENSRNOT00000018359.7	Bmp6	ENSRNOG0000013717	bone morphogenetic protein 6 Source RGD Symbol Acc 2214
B01	UPFR1105126	ENSRNOT00000084990.7	Bmp7	ENSRNOG0000053384	bone morphogenetic protein 7 Source RGD Symbol Acc 620743
B02	UPFR1025391	ENSRNOT00000079554.1	Bmpr1a	ENSRNOG0000052469	bone morphogenetic protein receptor type 1A Source RGD Symbol Acc 70989
B03	UPFR1039895	ENSRNOT00000021963.7	Bmpr1b	ENSRNOG0000016289	bone morphogenetic protein receptor type 1B Source RGD Symbol Acc 1595863
B04	UPFR1094942	ENSRNOT00000035238.4	Bmpr2	ENSRNOG0000022196	bone morphogenetic protein receptor type 2 Source RGD Symbol Acc 71082
B05	UPFR1013797	ENSRNOT00000008319.8	RGD1565355	ENSRNOG0000040108	similar to fatty acid translocase/CD36 Source RGD Symbol Acc 1565355
B06	UPFR1062734	ENSRNOT00000018535.5	Cdh11	ENSRNOG0000013481	cadherin 11 Source RGD Symbol Acc 619742
B07	UPFR1063759	ENSRNOT00000002394.7	Chrd	ENSRNOG0000001750	chordin Source RGD Symbol Acc 620181
B08	UPFR1057251	ENSRNOT00000081117.1	Col10a1	ENSRNOG0000051399	collagen type X alpha 1 chain Source RGD Symbol Acc 2371
B09	UPFR1100234	ENSRNOT00000067649.3	Col14a1	ENSRNOG0000026415	collagen type XIV alpha 1 chain Source RGD Symbol Acc 1305806
B10	UPFR1053469	ENSRNOT00000005311.6	Col1a1	ENSRNOG0000003897	collagen type I alpha 1 chain Source RGD Symbol Acc 61817
B11	UPFR1109083	ENSRNOT00000089292.1	Col1a2	ENSRNOG0000011292	collagen type I alpha 2 chain Source RGD Symbol Acc 621351
B12	UPFR1079135	ENSRNOT00000086062.1	Col2a1	ENSRNOG0000058560	collagen type II alpha 1 chain Source RGD Symbol Acc 2375
C01	UPFR1066020	ENSRNOT00000004956.4	Col3a1	ENSRNOG0000003357	collagen type III alpha 1 chain Source RGD Symbol Acc 71029
C02	UPFR1018897	ENSRNOT00000057386.2	Col4a1	ENSRNOG0000016281	collagen type IV alpha 1 chain Source RGD Symbol Acc 1307148
C03	UPFR1087101	ENSRNOT00000012334.7	Col5a1	ENSRNOG0000008749	collagen type V alpha 1 chain Source RGD Symbol Acc 70920
C04	UPFR1116990	ENSRNOT00000072951.2	Col6a1	ENSRNOG0000001249	collagen type VI alpha 1 chain Source RGD Symbol Acc 1565398
C05	UPFR1109491	ENSRNOT00000072536.2	Comp	ENSRNOG0000048472	cartilage oligomeric matrix protein Source RGD Symbol Acc 2378
C06	UPFR1022984	ENSRNOT00000081835.1	Csf1	ENSRNOG0000018659	colony stimulating factor 1 Source RGD Symbol Acc 621063
C07	UPFR1084630	ENSRNOT00000032333.3	Csf2	ENSRNOG0000026805	colony stimulating factor 2 Source RGD Symbol Acc 621065
C08	UPFR1037448	ENSRNOT00000011509.3	Csf3	ENSRNOG0000008525	colony stimulating factor 3 Source RGD Symbol Acc 2426
C09	UPFR1063624	ENSRNOT00000028730.6	Ctsk	ENSRNOG0000021155	cathepsin K Source RGD Symbol Acc 61810
C10	UPFR1057465	ENSRNOT00000014946.3	Dlx5	ENSRNOG0000010905	distal-less homeobox 5 Source RGD Symbol Acc 2506
		ENSRNOT000000		ENSRNOG00	

Position	Assay	Name	Symbol	Ensembl ID	Description
C11	UPFR1045222	087007.1	Egf	000053979	epidermal growth factor Source RGD Symbol Acc 2542
C12	UPFR1031298	ENSRNOT00000018577.3	Fgf1	ENSRNOG0000013867	fibroblast growth factor 1 Source RGD Symbol Acc 2605
D01	UPFR1033121	ENSRNOT00000023388.5	Fgf2	ENSRNOG0000017392	fibroblast growth factor 2 Source RGD Symbol Acc 2609
D02	UPFR1108643	ENSRNOT00000050980.5	Fgfr1	ENSRNOG0000016050	Fibroblast growth factor receptor 1 Source RGD Symbol Acc 620713
D03	UPFR1031158	ENSRNOT00000068183.2	Fgfr2	ENSRNOG0000016374	fibroblast growth factor receptor 2 Source RGD Symbol Acc 2611
D04	UPFR1070633	ENSRNOT00000001248.3	FIt1	ENSRNOG0000000940	FMS-related tyrosine kinase 1 Source RGD Symbol Acc 2621
D05	UPFR1089380	ENSRNOT00000057585.4	Fn1	ENSRNOG0000014288	fibronectin 1 Source RGD Symbol Acc 2624
D06	UPFR1064137	ENSRNOT00000081183.1	Gdf10	ENSRNOG0000051993	growth differentiation factor 10 Source RGD Symbol Acc 621512
D07	UPFR1075414	ENSRNOT00000009803.7	Gli1	ENSRNOG0000025120	GLI family zinc finger 1 Source RGD Symbol Acc 621673
D08	UPFR1031281	ENSRNOT00000028066.5	Icam1	ENSRNOG0000020679	intercellular adhesion molecule 1 Source RGD Symbol Acc 2857
D09	UPFR1050099	ENSRNOT00000038780.6	Igf1	ENSRNOG0000004517	insulin-like growth factor 1 Source RGD Symbol Acc 2868
D10	UPFR1018620	ENSRNOT00000019267.6	Igf1r	ENSRNOG0000014187	insulin-like growth factor 1 receptor Source RGD Symbol Acc 2869
D11	UPFR1064618	ENSRNOT00000024419.4	Ihh	ENSRNOG0000018059	Indian hedgehog signaling molecule Source RGD Symbol Acc 620021
D12	UPFR1068108	ENSRNOT00000077401.1	Itga2	ENSRNOG0000058111	integrin subunit alpha 2 Source RGD Symbol Acc 621632
E01	UPFR1093170	ENSRNOT00000005644.6	Itga3	ENSRNOG0000004276	integrin subunit alpha 3 Source RGD Symbol Acc 1310333
E02	UPFR1089159	ENSRNOT00000026748.4	Itgad	ENSRNOG0000019728	integrin subunit alpha M Source RGD Symbol Acc 2926
E03	UPFR1097241	ENSRNOT00000006961.7	Itgav	ENSRNOG0000004912	integrin subunit alpha V Source RGD Symbol Acc 1310613
E04	UPFR1091655	ENSRNOT00000014785.5	Itgb1	ENSRNOG0000010966	integrin subunit beta 1 Source RGD Symbol Acc 2927
E05	UPFR1089526	ENSRNOT00000013119.4	Mmp10	ENSRNOG0000032832	matrix metalloproteinase 10 Source RGD Symbol Acc 620192
E06	UPFR1080151	ENSRNOT00000022679.6	Mmp2	ENSRNOG0000016695	matrix metalloproteinase 2 Source RGD Symbol Acc 621316
E07	UPFR1084685	ENSRNOT00000013936.3	Mmp8	ENSRNOG0000009907	matrix metalloproteinase 8 Source RGD Symbol Acc 631408
E08	UPFR1032957	ENSRNOT00000023965.3	Mmp9	ENSRNOG0000017539	matrix metalloproteinase 9 Source RGD Symbol Acc 621320
E09	UPFR1116688	ENSRNOT00000036838.4	Nfkb1	ENSRNOG0000023258	nuclear factor kappa B subunit 1 Source RGD Symbol Acc 70498
E10	UPFR1059037	ENSRNOT00000032388.4	Nog	ENSRNOG0000023683	noggin Source RGD Symbol Acc 3183
E11	UPFR1035905	ENSRNOT00000042117.6	Pdgfra	ENSRNOG0000001312	platelet derived growth factor subunit A Source RGD Symbol Acc 3282
E12	UPFR1086497	ENSRNOT00000010970.7	Phex	ENSRNOG0000023440	phosphate regulating endopeptidase homolog, X-linked Source RGD Symbol Acc 3323
F01	UPFR1051266	ENSRNOT00000061014.4	Runx2	ENSRNOG0000020193	runt-related transcription factor 2 Source RGD Symbol Acc 2282
F02	UPFR1044270	ENSRNOT00000086091.1	Serpinh1	ENSRNOG0000016831	serpin family H member 1 Source RGD Symbol Acc 69302
F03	UPFR1071849	ENSRNOT00000025079.6	Smad1	ENSRNOG0000018483	SMAD family member 1 Source RGD Symbol Acc 3030
F04	UPFR1057367	ENSRNOT00000046847.4	Smad2	ENSRNOG0000018140	SMAD family member 2 Source RGD Symbol Acc 3031
F05	UPFR1033812	ENSRNOT00000039730.4	Smad3	ENSRNOG0000008620	SMAD family member 3 Source RGD Symbol Acc 3032
F06	UPFR1022605	ENSRNOT00000082484.1	Smad4	ENSRNOG0000051965	SMAD family member 4 Source RGD Symbol Acc 3033
F07	UPFR1115197	ENSRNOT00000016704.4	Smad5	ENSRNOG0000022870	SMAD family member 5 Source RGD Symbol Acc 620158
F08	UPFR1064691	ENSRNOT00000028238.2	Sost	ENSRNOG0000020805	sclerostin Source RGD Symbol Acc 69358
F09	UPFR1013060	ENSRNOT00000003511.6	Sox9	ENSRNOG0000002607	SRY box 9 Source RGD Symbol Acc 620474
F10	UPFR1014231	ENSRNOT00000018828.3	Sp7	ENSRNOG0000014082	Sp7 transcription factor Source RGD Symbol Acc 631377

Position	Assay	Name	Symbol	Ensembl ID	Description
F11	UPFR1100110	ENSRNOT00000067875.4	Spp1	ENSRNOG00000043451	secreted phosphoprotein 1 Source RGD Symbol Acc 3752
F12	UPFR1103881	ENSRNOT00000028051.4	Tgfb1	ENSRNOG00000020652	transforming growth factor, beta 1 Source RGD Symbol Acc 69051
G01	UPFR1046456	ENSRNOT00000087023.1	Tgfb2	ENSRNOG00000002418	transforming growth factor, beta 2 Source RGD Symbol Acc 70491
G02	UPFR1030751	ENSRNOT00000013516.5	Tgfb3	ENSRNOG0000009867	transforming growth factor, beta 3 Source RGD Symbol Acc 3851
G03	UPFR1035135	ENSRNOT00000081090.1	Tgfbr1	ENSRNOG0000007036	transforming growth factor, beta receptor 1 Source RGD Symbol Acc 3852
G04	UPFR1107633	ENSRNOT00000037883.4	Tgfbr2	ENSRNOG00000013265	transforming growth factor, beta receptor 2 Source RGD Symbol Acc 69651
G05	UPFR1058624	ENSRNOT00000002867.4	Tgfbr3	ENSRNOG00000002093	transforming growth factor beta receptor 3 Source RGD Symbol Acc 61821
G06	UPFR1098022	ENSRNOT00000079677.1	LOC103694380	ENSRNOG00000055156	tumor necrosis factor-like Source RGD Symbol Acc 9404643
G07	UPFR1105006	ENSRNOT00000049056.4	LOC103689974	ENSRNOG00000049547	TNF superfamily member 11 Source RGD Symbol Acc 620784
G08	UPFR1038411	ENSRNOT00000014763.6	Twist1	ENSRNOG00000011101	twist family bHLH transcription factor 1 Source RGD Symbol Acc 621455
G09	UPFR1101375	ENSRNOT00000019377.5	Vcam1	ENSRNOG00000014333	vascular cell adhesion molecule 1 Source RGD Symbol Acc 3952
G10	UPFR1048607	ENSRNOT00000082429.1	Vdr	ENSRNOG00000054420	vitamin D receptor Source RGD Symbol Acc 3959
G11	UPFR1117908	ENSRNOT00000026559.5	Vegfa	ENSRNOG00000019598	vascular endothelial growth factor A Source RGD Symbol Acc 619991
G12	UPFR1063692	ENSRNOT00000050891.3	Vegfb	ENSRNOG00000021156	vascular endothelial growth factor B Source RGD Symbol Acc 619799
H01	UPFR1132952	ENSRNOT00000080216.1	Actb	ENSRNOG00000034254	actin, beta Source RGD Symbol Acc 628837
H02	UPFR1132953	ENSRNOT00000023017.5	B2m	ENSRNOG00000017123	beta-2 microglobulin Source RGD Symbol Acc 2189
H03	UPFR1132959	ENSRNOT00000065935.3	Hprt1	ENSRNOG00000048561	hypoxanthine phosphoribosyltransferase 1 Source RGD Symbol Acc 2826
H04	UPFR1018740	ENSRNOT00000017468.2	Ldha	ENSRNOG00000013009	lactate dehydrogenase A Source RGD Symbol Acc 2996
H05	UPFR1132958	ENSRNOT00000018820.5	Rplp1	ENSRNOG00000013874	ribosomal protein lateral stalk subunit P1 Source RGD Symbol Acc 621774
H06	UPFR1126610	UPL_RGDC	RGDC	UPL_RGDC	Rat Genomic DNA Contamination
H07	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H08	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H09	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H10	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control
H11	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control
H12	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control

Related products

Product	Contents	Cat. no.
QuantiNova LNA Probe PCR QC Panel	These panels are designed to assess the quality of RNA samples before characterization using QuantiNova LNA Probe PCR Focus Panels; available in 96-well, 384-well, and Rotor-Disc 100 formats	249945
QuantiNova Reverse Transcription Kit (10)*	For 10 x 20 µl reactions: 20 µl 8x gDNA Removal Mix, 10 µl Reverse Transcription Enzyme, 40 µl Reverse Transcription Mix (containing RT primers), 20 µl Internal Control RNA, 1.9 ml RNase-Free Water	205410
QuantiNova Probe RT-PCR Kit (100)*	For 100 x 20 µl reactions: 1 ml QuantiNova Probe RT-PCR Master Mix, 20 µl QuantiNova Probe RT Mix, 20 µl Internal Control RNA, 500 µl Yellow Template Dilution Buffer, 250 µl ROX Reference Dye, 1.9 µl RNase-Free Water	208352
QuantiNova Probe PCR Kit (100)*	For 100 x 20 µl reactions: 1 ml 2x QuantiNova Probe PCR Master Mix, 500 µl QuantiNova Yellow Template Dilution Buffer, 250 µl QN ROX Reference Dye, 1.9 ml Water	208252

*Larger kit sizes available.

The QuantiNova LNA Probe PCR Focus Panels are intended for molecular biology applications. These products are not intended for the diagnosis, prevention or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, LNA®, QuantiNova®, Sample to Insight® (QIAGEN Group); SYBR® (Life Technologies Corp.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

09/2019 © 2019 QIAGEN, all rights reserved.