

QuantiNova® LNA® Probe PCR Focus Panels (Rotor-Gene® Format)

Rat EGF / PDGF Signaling Pathway

Cat. no. 249955 UPRN-040ZR

For study focus gene expression analysis

Shipping and storage

QuantiNova LNA Probe PCR Focus Panels are shipped at room temperature. Immediately upon receipt, they should be stored protected from light at 2–8°C for short term storage or at –30°C to –15°C for long time storage. Under these conditions, all components are stable for at least 12 months.

Note: Open the package and store the products appropriately immediately upon receipt.

For optimal performance, QuantiNova LNA Probe PCR Focus Panels should be used together with the QuantiNova Reverse Transcription Kit for cDNA synthesis and the QuantiNova Probe PCR Kit (Mastermix) for PCR.

Panel layout (Rotor-Gene): QuantiNova LNA Probe PCR Focus Panel

The 96 real-time assays in the Rotor-Gene format are located in wells 1–96 of the Rotor-Disc® (plate A1–A12=Rotor-Disc 1–12, plate B1–B12=Rotor-Disc 13–24, etc.). To maintain data analysis compatibility, wells 97–100 do not contain real-time assays but will contain master mix to account for weight balance. Refer to the QuantiNova LNA Probe PCR Handbook at www.qiagen.com for further details.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Actr2	Akt1	Akt2	Akt3	Araf	Aif1	Aif2	Bad	Bcar1	Bcl2	Braf	Casp3
B	Casp9	Cblb	Ccnd1	Chuk	Col1a1	Creb1	Csnk2a1	Csnk2b	Dusp1	Dusp6	Egf	Egfr
C	Egr1	Eif4e	Eps8	Faslg	Fn1	Fos	Foxo3	Gab1	Grb2	Gsk3a	Gsk3b	Hbegf
D	Hras	Ilkbb	Il2	Jak1	Jun	Kcnh8	Kras	Lta	Map2k1	Map2k4	Map2k7	Map3k2
E	Mapk1	Mapk10	Mapk3	Mapk8	Mapk9	Mknk1	Mmp7	Nck2	Nfatc3	Nfkb1	Nras	Nup62
F	Pdgfra	Pdgfrb	Pdgfra	Pdpk1	Pik3ca	Pik3r1	Pik3r2	Plat	Plcg1	LOC103694903	Prkca	Plen
G	Raf1	Rap1a	Rasa1	Rhoa	Rps6ka5	Rps6kb1	Shc1	Src	Stat1	Stat3	Stat5a	Tp53
H	Actb	B2m	Hprt1	Ldha	Rplp1	RGDC	QIC	QIC	QIC	PPC	PPC	PPC

Gene table: QuantiNova LNA Probe PCR Focus Panel

Position	Assay	Name	Symbol	Ensembl ID	Description
A01	UPFR1021386	ENSRNOT00000006607.5	Acr2	ENSRNOG0000004959	ARP2 actin related protein 2 homolog Source RGD Symbol Acc 1310826
A02	UPFR1083196	ENSRNOT00000031164.3	Akt1	ENSRNOG0000028629	AKT serine/threonine kinase 1 Source RGD Symbol Acc 2081
A03	UPFR1088137	ENSRNOT00000025303.3	Akt2	ENSRNOG0000018677	AKT serine/threonine kinase 2 Source RGD Symbol Acc 2082
A04	UPFR1105050	ENSRNOT00000085648.1	Akt3	ENSRNOG0000021497	AKT serine/threonine kinase 3 Source RGD Symbol Acc 62390
A05	UPFR1060431	ENSRNOT00000039911.3	Araf	ENSRNOG0000010838	A-Raf proto-oncogene, serine/threonine kinase Source RGD Symbol Acc 2148
A06	UPFR1117151	ENSRNOT00000079011.1	Atf1	ENSRNOG0000061088	activating transcription factor 1 Source RGD Symbol Acc 1307360
A07	UPFR1067157	ENSRNOT00000002174.7	Atf2	ENSRNOG0000001597	activating transcription factor 2 Source RGD Symbol Acc 621862
A08	UPFR1109796	ENSRNOT00000028712.5	Bad	ENSRNOG0000021147	BCL2-associated agonist of cell death Source RGD Symbol Acc 620103
A09	UPFR1037655	ENSRNOT00000058891.4	Bcar1	ENSRNOG0000019253	BCAR1, Cas family scaffold protein Source RGD Symbol Acc 2406
A10	UPFR1082998	ENSRNOT00000003768.2	Bcl2	ENSRNOG0000002791	BCL2, apoptosis regulator Source RGD Symbol Acc 2199
A11	UPFR1075695	ENSRNOT00000014604.7	Braf	ENSRNOG0000010957	B-Raf proto-oncogene, serine/threonine kinase Source RGD Symbol Acc 619908
A12	UPFR1091664	ENSRNOT00000014095.5	Casp3	ENSRNOG0000010475	caspase 3 Source RGD Symbol Acc 2275
B01	UPFR1054782	ENSRNOT00000085378.1	Casp9	ENSRNOG0000012944	caspase 9 Source RGD Symbol Acc 61867
B02	UPFR1107987	ENSRNOT00000077220.1	Cblb	ENSRNOG0000001982	Cbl proto-oncogene B Source RGD Symbol Acc 620535
B03	UPFR1042457	ENSRNOT00000088588.1	Ccnd1	ENSRNOG0000020918	cyclin D1 Source RGD Symbol Acc 68384
B04	UPFR1110725	ENSRNOT00000030782.4	Chuk	ENSRNOG0000022485	conserved helix-loop-helix ubiquitous kinase Source RGD Symbol Acc 1306661
B05	UPFR1053469	ENSRNOT00000005311.6	Col1a1	ENSRNOG0000003897	collagen type I alpha 1 chain Source RGD Symbol Acc 61817
B06	UPFR1042088	ENSRNOT00000018326.5	Creb1	ENSRNOG0000013412	cAMP responsive element binding protein 1 Source RGD Symbol Acc 620218
B07	UPFR1071105	ENSRNOT00000007558.5	Csnk2a1	ENSRNOG0000005276	casein kinase 2 alpha 1 Source RGD Symbol Acc 621663
B08	UPFR1014922	ENSRNOT00000076528.1	Csnk2b	ENSRNOG0000000847	casein kinase 2 beta Source RGD Symbol Acc 619978
B09	UPFR1053812	ENSRNOT00000005383.5	Dusp1	ENSRNOG0000003977	dual specificity phosphatase 1 Source RGD Symbol Acc 620897
B10	UPFR1071958	ENSRNOT00000037844.3	Dusp6	ENSRNOG0000023896	dual specificity phosphatase 6 Source RGD Symbol Acc 70978
B11	UPFR1045222	ENSRNOT00000087007.1	Egf	ENSRNOG0000053979	epidermal growth factor Source RGD Symbol Acc 2542
B12	UPFR1091984	ENSRNOT00000006087.2	Egfr	ENSRNOG0000004332	epidermal growth factor receptor Source RGD Symbol Acc 2543
C01	UPFR1034050	ENSRNOT00000026303.4	Egr1	ENSRNOG0000019422	early growth response 1 Source RGD Symbol Acc 2544
C02	UPFR1096247	ENSRNOT00000090526.2	Eif4e	ENSRNOG0000052343	eukaryotic translation initiation factor 4E Source RGD Symbol Acc 69647
C03	UPFR1046275	ENSRNOT00000009328.7	Eps8	ENSRNOG0000007047	epidermal growth factor receptor pathway substrate 8 Source RGD Symbol Acc 1310590
C04	UPFR1108356	ENSRNOT00000003998.2	Faslg	ENSRNOG0000002978	Fas ligand Source RGD Symbol Acc 3880
C05	UPFR1089380	ENSRNOT00000057585.4	Fn1	ENSRNOG0000014288	fibronectin 1 Source RGD Symbol Acc 2624
C06	UPFR1045359	ENSRNOT00000010712.3	Fos	ENSRNOG0000008015	Fos proto-oncogene, AP-1 transcription factor subunit Source RGD Symbol Acc 2626
C07	UPFR1095067	ENSRNOT00000000327.5	Foxo3	ENSRNOG0000000299	forkhead box O3 Source RGD Symbol Acc 1309196
C08	UPFR1032475	ENSRNOT00000083219.1	Gab1	ENSRNOG0000017879	GRB2-associated binding protein 1 Source RGD Symbol Acc 1311085
C09	UPFR1087543	ENSRNOT00000005347.5	Grb2	ENSRNOG0000003990	growth factor receptor bound protein 2 Source RGD Symbol Acc 619758
C10	UPFR1081475	ENSRNOT00000027677.5	Gsk3a	ENSRNOG0000020417	glycogen synthase kinase 3 alpha Source RGD Symbol Acc 620351
		ENSRNOT000000		ENSRNOG00	

Position	Assay	Name	Symbol	Ensembl ID	Description
C11	UPFR1109626	077612.1	Gsk3b	000002833	glycogen synthase kinase 3 beta Source RGD Symbol Acc 70982
C12	UPFR1043686	ENSRNOT00000025157.6	Hbegf	ENSRNOG0000018646	heparin-binding EGF-like growth factor Source RGD Symbol Acc 2526
D01	UPFR1076419	ENSRNOT00000022363.6	Hras	ENSRNOG0000016611	Hras proto-oncogene, GTPase Source RGD Symbol Acc 2827
D02	UPFR1104220	ENSRNOT00000025851.4	Ikbkb	ENSRNOG0000019073	inhibitor of nuclear factor kappa B kinase subunit beta Source RGD Symbol Acc 621375
D03	UPFR1070791	ENSRNOT00000023327.2	Il2	ENSRNOG0000017348	interleukin 2 Source RGD Symbol Acc 620047
D04	UPFR1050030	ENSRNOT00000064542.1	Jak1	ENSRNOG0000011157	Janus kinase 1 Source RGD Symbol Acc 69056
D05	UPFR1037970	ENSRNOT00000011731.3	Jun	ENSRNOG0000026293	Jun proto-oncogene, AP-1 transcription factor subunit Source RGD Symbol Acc 2943
D06	UPFR1013403	ENSRNOT00000087331.1	Kcnh8	ENSRNOG0000058626	potassium voltage-gated channel subfamily H member 8 Source RGD Symbol Acc 2549
D07	UPFR1088461	ENSRNOT00000012588.4	Kras	ENSRNOG0000009338	KRAS proto-oncogene, GTPase Source RGD Symbol Acc 2981
D08	UPFR1116910	ENSRNOT00000088936.1	Lta	ENSRNOG0000000838	lymphotoxin alpha Source RGD Symbol Acc 3020
D09	UPFR1053192	ENSRNOT00000013933.6	Map2k1	ENSRNOG0000010176	mitogen activated protein kinase kinase 1 Source RGD Symbol Acc 70495
D10	UPFR1096951	ENSRNOT00000005293.7	Map2k4	ENSRNOG0000003834	mitogen activated protein kinase kinase 4 Source RGD Symbol Acc 1307040
D11	UPFR1059686	ENSRNOT00000084460.1	Map2k7	ENSRNOG0000001047	mitogen activated protein kinase kinase 7 Source RGD Symbol Acc 1560043
D12	UPFR1053149	ENSRNOT00000060996.3	Map3k2	ENSRNOG0000014089	mitogen activated protein kinase kinase 2 Source RGD Symbol Acc 620967
E01	UPFR1055882	ENSRNOT00000002533.7	Mapk1	ENSRNOG0000001849	mitogen activated protein kinase 1 Source RGD Symbol Acc 70500
E02	UPFR1062108	ENSRNOT00000092395.1	Mapk10	ENSRNOG0000002079	mitogen activated protein kinase 10 Source RGD Symbol Acc 3663
E03	UPFR1049954	ENSRNOT00000087625.1	Mapk3	ENSRNOG0000053583	mitogen activated protein kinase 3 Source RGD Symbol Acc 3046
E04	UPFR1048766	ENSRNOT00000065216.3	Mapk8	ENSRNOG0000020155	mitogen-activated protein kinase 8 Source RGD Symbol Acc 621506
E05	UPFR1038865	ENSRNOT00000083688.1	Mapk9	ENSRNOG0000002823	mitogen-activated protein kinase 9 Source RGD Symbol Acc 628847
E06	UPFR1119759	ENSRNOT00000067655.2	Mknk1	ENSRNOG0000010381	MAP kinase-interacting serine/threonine kinase 1 Source RGD Symbol Acc 1559603
E07	UPFR1034539	ENSRNOT00000014041.5	Mmp7	ENSRNOG0000010507	matrix metalloproteinase 7 Source RGD Symbol Acc 3100
E08	UPFR1044689	ENSRNOT00000023029.7	Nck2	ENSRNOG0000017155	NCK adaptor protein 2 Source RGD Symbol Acc 1305965
E09	UPFR1047494	ENSRNOT00000087491.1	Nfatc3	ENSRNOG0000054264	nuclear factor of activated T-cells 3 Source RGD Symbol Acc 1308692
E10	UPFR1116688	ENSRNOT00000036838.4	Nfkb1	ENSRNOG0000023258	nuclear factor kappa B subunit 1 Source RGD Symbol Acc 70498
E11	UPFR1038224	ENSRNOT00000039572.3	Nras	ENSRNOG0000023079	NRAS proto-oncogene, GTPase Source RGD Symbol Acc 3205
E12	UPFR1110039	ENSRNOT00000074847.2	Nup62	ENSRNOG0000048733	nucleoporin 62 Source RGD Symbol Acc 619938
F01	UPFR1035905	ENSRNOT00000042117.6	Pdgfra	ENSRNOG0000001312	platelet derived growth factor subunit A Source RGD Symbol Acc 3282
F02	UPFR1088903	ENSRNOT00000023196.5	Pdgfb	ENSRNOG0000017197	platelet derived growth factor subunit B Source RGD Symbol Acc 3283
F03	UPFR1089489	ENSRNOT00000003077.5	Pdgfra	ENSRNOG0000002244	platelet derived growth factor receptor alpha Source RGD Symbol Acc 3284
F04	UPFR1093188	ENSRNOT00000067660.1	Pdpk1	ENSRNOG0000006136	3-phosphoinositide dependent protein kinase-1 Source RGD Symbol Acc 620307
F05	UPFR1035210	ENSRNOT00000083720.1	Pik3ca	ENSRNOG0000056371	phosphatidylinositol-4,5-bisphosphate 3-kinase, catalytic subunit alpha Source RGD Symbol Acc 620916
F06	UPFR1035231	ENSRNOT00000075057.1	Pik3r1	ENSRNOG0000018903	phosphoinositide-3-kinase regulatory subunit 1 Source RGD Symbol Acc 3329
F07	UPFR1126563	ENSRNOT00000026210.5	Pik3r2	ENSRNOG0000019228	phosphoinositide-3-kinase regulatory subunit 2 Source RGD Symbol Acc 68341
F08	UPFR1085291	ENSRNOT00000025763.6	Plat	ENSRNOG0000019018	plasminogen activator, tissue type Source RGD Symbol Acc 3342
F09	UPFR1035470	ENSRNOT00000078909.1	Plcg1	ENSRNOG0000051490	phospholipase C, gamma 1 Source RGD Symbol Acc 3347
F10	UPFR1111771	ENSRNOT00000007621.5	LOC103694903	ENSRNOG0000005389	serine/threonine-protein phosphatase 2A catalytic subunit alpha isoform Source RGD Symbol Acc 9257421

Position	Assay	Name	Symbol	Ensembl ID	Description
F11	UPFR1034906	ENSRNOT00000004699.8	Prkca	ENSRNOG0000003491	protein kinase C, alpha Source RGD Symbol Acc 3395
F12	UPFR1043890	ENSRNOT00000028143.3	Pten	ENSRNOG0000020723	phosphatase and tensin homolog Source RGD Symbol Acc 61995
G01	UPFR1087799	ENSRNOT00000013831.6	Raf1	ENSRNOG0000010153	Raf-1 proto-oncogene, serine/threonine kinase Source RGD Symbol Acc 3531
G02	UPFR1051806	ENSRNOT00000047836.4	Rap1a	ENSRNOG0000032463	RAP1A, member of RAS oncogene family Source RGD Symbol Acc 1359694
G03	UPFR1085068	ENSRNOT00000046251.2	Rasa1	ENSRNOG0000029185	RAS p21 protein activator 1 Source RGD Symbol Acc 3537
G04	UPFR1119707	ENSRNOT00000071664.2	Rhoa	ENSRNOG0000050519	ras homolog family member A Source RGD Symbol Acc 619921
G05	UPFR1013406	ENSRNOT00000091693.1	Rps6ka5	ENSRNOG0000004362	ribosomal protein S6 kinase A5 Source RGD Symbol Acc 1308336
G06	UPFR1107808	ENSRNOT00000005226.4	Rps6kb1	ENSRNOG0000003919	ribosomal protein S6 kinase B1 Source RGD Symbol Acc 620683
G07	UPFR1103707	ENSRNOT00000028038.6	Shc1	ENSRNOG0000020657	SHC adaptor protein 1 Source RGD Symbol Acc 620446
G08	UPFR1050158	ENSRNOT00000012739.4	Src	ENSRNOG0000009495	SRC proto-oncogene, non-receptor tyrosine kinase Source RGD Symbol Acc 620795
G09	UPFR1068760	ENSRNOT00000083514.1	Stat1	ENSRNOG0000014079	signal transducer and activator of transcription 4 Source RGD Symbol Acc 1305747
G10	UPFR1040072	ENSRNOT00000026760.4	Stat3	ENSRNOG0000019742	signal transducer and activator of transcription 3 Source RGD Symbol Acc 3772
G11	UPFR1042053	ENSRNOT00000026662.3	Stat5a	ENSRNOG0000019496	signal transducer and activator of transcription 5A Source RGD Symbol Acc 3773
G12	UPFR1095007	ENSRNOT00000085115.1	Tp53	ENSRNOG0000010756	tumor protein p53 Source RGD Symbol Acc 3889
H01	UPFR1132952	ENSRNOT00000080216.1	Actb	ENSRNOG0000034254	actin, beta Source RGD Symbol Acc 628837
H02	UPFR1132953	ENSRNOT00000023017.5	B2m	ENSRNOG0000017123	beta-2 microglobulin Source RGD Symbol Acc 2189
H03	UPFR1132959	ENSRNOT00000065935.3	Hprt1	ENSRNOG0000048561	hypoxanthine phosphoribosyltransferase 1 Source RGD Symbol Acc 2826
H04	UPFR1018740	ENSRNOT00000017468.2	Ldha	ENSRNOG0000013009	lactate dehydrogenase A Source RGD Symbol Acc 2996
H05	UPFR1132958	ENSRNOT00000018820.5	Rplp1	ENSRNOG0000013874	ribosomal protein lateral stalk subunit P1 Source RGD Symbol Acc 621774
H06	UPFR1126610	UPL_RGDC	RGDC	UPL_RGDC	Rat Genomic DNA Contamination
H07	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H08	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H09	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H10	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control
H11	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control
H12	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control

Related products

Product	Contents	Cat. no.
QuantiNova LNA Probe PCR QC Panel	These panels are designed to assess the quality of RNA samples before characterization using QuantiNova LNA Probe PCR Focus Panels; available in 96-well, 384-well, and Rotor-Disc 100 formats	249945
QuantiNova Reverse Transcription Kit (10)*	For 10 x 20 µl reactions: 20 µl 8x gDNA Removal Mix, 10 µl Reverse Transcription Enzyme, 40 µl Reverse Transcription Mix (containing RT primers), 20 µl Internal Control RNA, 1.9 ml RNase-Free Water	205410
QuantiNova Probe RT-PCR Kit (100)*	For 100 x 20 µl reactions: 1 ml QuantiNova Probe RT-PCR Master Mix, 20 µl QuantiNova Probe RT Mix, 20 µl Internal Control RNA, 500 µl Yellow Template Dilution Buffer, 250 µl ROX Reference Dye, 1.9 µl RNase-Free Water	208352
QuantiNova Probe PCR Kit (100)*	For 100 x 20 µl reactions: 1 ml 2x QuantiNova Probe PCR Master Mix, 500 µl QuantiNova Yellow Template Dilution Buffer, 250 µl QN ROX Reference Dye, 1.9 ml Water	208252

*Larger kit sizes available.

The QuantiNova LNA Probe PCR Focus Panels are intended for molecular biology applications. These products are not intended for the diagnosis, prevention or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, LNA®, QuantiNova®, Sample to Insight® (QIAGEN Group); SYBR® (Life Technologies Corp.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

09/2019 © 2019 QIAGEN, all rights reserved.