

QuantiNova® LNA® Probe PCR Focus Panels (Rotor-Gene® Format)

Rat Parkinson's Disease

Cat. no. 249955 UPRN-124ZR

For study focus gene expression analysis

Shipping and storage

QuantiNova LNA Probe PCR Focus Panels are shipped at room temperature. Immediately upon receipt, they should be stored protected from light at 2–8°C for short term storage or at –30°C to –15°C for long time storage. Under these conditions, all components are stable for at least 12 months.

Note: Open the package and store the products appropriately immediately upon receipt.

For optimal performance, QuantiNova LNA Probe PCR Focus Panels should be used together with the QuantiNova Reverse Transcription Kit for cDNA synthesis and the QuantiNova Probe PCR Kit (Mastermix) for PCR.

Panel layout (Rotor-Gene): QuantiNova LNA Probe PCR Focus Panel

The 96 real-time assays in the Rotor-Gene format are located in wells 1–96 of the Rotor-Disc® (plate A1–A12=Rotor-Disc 1–12, plate B1–B12=Rotor-Disc 13–24, etc.). To maintain data analysis compatibility, wells 97–100 do not contain real-time assays but will contain master mix to account for weight balance. Refer to the QuantiNova LNA Probe PCR Handbook at www.qiagen.com for further details.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Aldh1a1	Apc	App	Atp2b2	Atxn2	Atxn3	Basp1	Bdnf	Cadps	Casp1	Casp3	Casp7
B	Casp8	Casp9	Cd44	Cdc27	Cdc42	Cdh8	Chgb	Cirbp	Cul2	Cxcl1	Ddc	Dlk1
C	Drd2	Egln1	Fbxo9	Fgf13	Fpx1	Ftn1	Gabbg2	Gpr37	Gria3	Hmgcs1	Hspa4	Htr2a
D	Kcnj6	Mapk9	Mapt	Nefl	Nfasc	Nr4a2	Nrxn3	Nsf	Nsg1	Ntrk2	Opa1	Pan2
E	AABR0700151 9.1	Park7	Pink1	Ppid	Ppp2r1b	Prdx2	Psen2	Plen	Rgs4	Rtn1	S100b	Skp1
F	Slc18a2	Slc25a4	Slc6a3	Slt1	Sncg	Spen	Srsf7	Stub1	AC126641.1	Syng3	Syt1	Syt11
G	Tcf7l2	Th	Tpbg	Uba1	Ubb	Ube2i	Ube2k	Ube2l3	Uchl1	Vamp1	Vdac3	Ywhaz
H	Actb	B2m	Hprt1	Ldha	Rplp1	RGDC	QIC	QIC	QIC	PPC	PPC	PPC

Gene table: QuantiNova LNA Probe PCR Focus Panel

Position	Assay	Name	Symbol	Ensembl ID	Description
A01	UPFR1084336	ENSRNOT00000024000.6	Aldh1a1	ENSRNOG0000017619	aldehyde dehydrogenase 1 family, member A1 Source RGD Symbol Acc 2087
A02	UPFR1090312	ENSRNOT00000090264.1	Apc	ENSRNOG0000020423	APC, WNT signaling pathway regulator Source RGD Symbol Acc 2123
A03	UPFR1021296	ENSRNOT00000044081.5	App	ENSRNOG0000006997	amyloid beta precursor protein Source RGD Symbol Acc 2139
A04	UPFR1055524	ENSRNOT00000067931.3	Atp2b2	ENSRNOG0000030269	ATPase plasma membrane Ca2+ transporting 2 Source RGD Symbol Acc 2176
A05	UPFR1025098	ENSRNOT00000001691.8	Atxn2	ENSRNOG0000001256	ataxin 2 Source RGD Symbol Acc 1306637
A06	UPFR1071648	ENSRNOT00000007505.2	Atxn3	ENSRNOG0000005470	ataxin 3 Source RGD Symbol Acc 621567
A07	UPFR1054974	ENSRNOT00000073915.2	Basp1	ENSRNOG0000046313	brain abundant, membrane attached signal protein 1 Source RGD Symbol Acc 621491
A08	UPFR1014257	ENSRNOT00000078543.1	Bdnf	ENSRNOG0000047466	brain-derived neurotrophic factor Source RGD Symbol Acc 2202
A09	UPFR1104126	ENSRNOT00000012153.6	Cadps	ENSRNOG0000008570	calcium dependent secretion activator Source RGD Symbol Acc 708573
A10	UPFR1052168	ENSRNOT00000009993.6	Casp1	ENSRNOG0000007372	caspase 1 Source RGD Symbol Acc 2274
A11	UPFR1091664	ENSRNOT00000014095.5	Casp3	ENSRNOG0000010475	caspase 3 Source RGD Symbol Acc 2275
A12	UPFR1051332	ENSRNOT000000080511.1	Casp7	ENSRNOG00000056216	caspase 7 Source RGD Symbol Acc 620944
B01	UPFR1030636	ENSRNOT00000016613.4	Casp8	ENSRNOG0000012331	caspase 8 Source RGD Symbol Acc 620945
B02	UPFR1054782	ENSRNOT000000085378.1	Casp9	ENSRNOG0000012944	caspase 9 Source RGD Symbol Acc 61867
B03	UPFR1040363	ENSRNOT00000009073.7	Cd44	ENSRNOG0000006094	CD44 molecule (Indian blood group) Source RGD Symbol Acc 2307
B04	UPFR1087432	ENSRNOT00000007963.5	Cdc27	ENSRNOG0000005904	cell division cycle 27 Source RGD Symbol Acc 1304921
B05	UPFR1103160	ENSRNOT00000018118.5	Cdc42	ENSRNOG0000013536	cell division cycle 42 Source RGD Symbol Acc 71043
B06	UPFR1070588	ENSRNOT000000084126.1	Cdh8	ENSRNOG00000056643	cadherin 8 Source RGD Symbol Acc 69286
B07	UPFR1114232	ENSRNOT00000028892.4	Chgb	ENSRNOG00000021269	chromogranin B Source RGD Symbol Acc 2339
B08	UPFR1097525	ENSRNOT00000021648.6	Cirbp	ENSRNOG0000015999	cold inducible RNA binding protein Source RGD Symbol Acc 620756
B09	UPFR1119212	ENSRNOT00000079607.1	Cul2	ENSRNOG0000015292	cullin 2 Source RGD Symbol Acc 1310644
B10	UPFR1078963	ENSRNOT00000020472.6	Cxxc1	ENSRNOG0000014614	CXXC finger protein 1 Source RGD Symbol Acc 1310755
B11	UPFR1086801	ENSRNOT000000057051.4	Ddc	ENSRNOG0000004327	dopa decarboxylase Source RGD Symbol Acc 2494
B12	UPFR1103659	ENSRNOT00000006339.6	Dlk1	ENSRNOG0000019584	delta like non-canonical Notch ligand 1 Source RGD Symbol Acc 619931
C01	UPFR1040837	ENSRNOT00000045944.5	Drd2	ENSRNOG0000008428	dopamine receptor D2 Source RGD Symbol Acc 2520
C02	UPFR1097483	ENSRNOT00000026767.7	Egln1	ENSRNOG0000019773	egl-9 family hypoxia-inducible factor 1 Source RGD Symbol Acc 631375
C03	UPFR1100246	ENSRNOT00000010961.5	Fbxo9	ENSRNOG0000008214	f-box protein 9 Source RGD Symbol Acc 1310374
C04	UPFR1116340	ENSRNOT000000091330.1	Fgf13	ENSRNOG0000042753	fibroblast growth factor 13 Source RGD Symbol Acc 620164
C05	UPFR1113935	ENSRNOT00000007077.5	Fjx1	ENSRNOG0000005347	four-jointed box kinase 1 Source RGD Symbol Acc 1306113
C06	UPFR1089380	ENSRNOT000000057585.4	Fn1	ENSRNOG0000014288	fibronectin 1 Source RGD Symbol Acc 2624
C07	UPFR1038575	ENSRNOT000000111573.4	Gabbr2	ENSRNOG0000008431	gamma-aminobutyric acid type B receptor subunit 2 Source RGD Symbol Acc 619864
C08	UPFR1100665	ENSRNOT00000003593.6	Gpr37	ENSRNOG0000002524	G protein-coupled receptor 37 Source RGD Symbol Acc 619848
C09	UPFR1047643	ENSRNOT00000010367.6	Gria3	ENSRNOG0000007682	glutamate ionotropic receptor AMPA type subunit 3 Source RGD Symbol Acc 70958
C10	UPFR1013564	ENSRNOT00000023554.6	Hmgcs1	ENSRNOG0000016552	3-hydroxy-3-methylglutaryl-CoA synthase 1 Source RGD Symbol Acc 70970
		ENSRNOT000000		ENSRNOG00	

Position	Assay	Name	Symbol	Ensembl ID	Description
C11	UPFR1042317	023628.3	Hspa4	000016596	heat shock protein family A member 4 Source RGD Symbol Acc 628878
C12	UPFR1042475	ENSRNOT00000013408.6	Hlr2a	ENSRNOG0000010063	5-hydroxytryptamine receptor 2A Source RGD Symbol Acc 61800
D01	UPFR1047292	ENSRNOT00000079955.1	Kcnj6	ENSRNOG0000001658	potassium voltage-gated channel subfamily J member 6 Source RGD Symbol Acc 2959
D02	UPFR1038865	ENSRNOT00000083688.1	Mapk9	ENSRNOG0000002823	mitogen-activated protein kinase 9 Source RGD Symbol Acc 628847
D03	UPFR1041157	ENSRNOT00000045127.6	Mapt	ENSRNOG0000005133	microtubule-associated protein tau Source RGD Symbol Acc 69329
D04	UPFR1018074	ENSRNOT00000018599.3	Nefl	ENSRNOG0000013658	neurofilament light Source RGD Symbol Acc 621458
D05	UPFR1061749	ENSRNOT00000051512.5	Nfasc	ENSRNOG0000030515	neurofascin Source RGD Symbol Acc 620911
D06	UPFR1088324	ENSRNOT00000078001.1	Nr4a2	ENSRNOG0000005600	nuclear receptor subfamily 4, group A, member 2 Source RGD Symbol Acc 3202
D07	UPFR1114223	ENSRNOT00000031205.5	Nrxn3	ENSRNOG0000047574	neurexin 3 Source RGD Symbol Acc 620212
D08	UPFR1036287	ENSRNOT00000006361.6	Nsf	ENSRNOG0000003905	N-ethylmaleimide sensitive factor, vesicle fusing ATPase Source RGD Symbol Acc 621594
D09	UPFR1063359	ENSRNOT00000008030.7	Nsg1	ENSRNOG0000005700	neuronal vesicle trafficking associated 1 Source RGD Symbol Acc 2222
D10	UPFR1075665	ENSRNOT00000042145.2	Ntrk2	ENSRNOG0000018839	neurotrophic receptor tyrosine kinase 2 Source RGD Symbol Acc 3213
D11	UPFR1114549	ENSRNOT00000002343.7	Opa1	ENSRNOG0000001717	OPA1, mitochondrial dynamin like GTPase Source RGD Symbol Acc 708423
D12	UPFR1075051	ENSRNOT00000050792.4	Pan2	ENSRNOG0000032441	poly(A) specific ribonuclease subunit PAN2 Source RGD Symbol Acc 1303301
E01	UPFR1080845	ENSRNOT00000086127.1	AABR07001519.1	ENSRNOG0000055547	
E02	UPFR1045325	ENSRNOT00000091206.1	Park7	ENSRNOG0000018289	Parkinsonism associated deglycase Source RGD Symbol Acc 621808
E03	UPFR1021822	ENSRNOT00000083474.1	Pink1	ENSRNOG0000015385	PTEN induced putative kinase 1 Source RGD Symbol Acc 1305769
E04	UPFR1030191	ENSRNOT00000037890.4	Ppid	ENSRNOG0000027408	peptidylprolyl isomerase D Source RGD Symbol Acc 1303174
E05	UPFR1072545	ENSRNOT00000064290.2	Ppp2r1b	ENSRNOG0000010922	protein phosphatase 2 scaffold subunit A beta Source RGD Symbol Acc 1304764
E06	UPFR1094945	ENSRNOT00000004799.5	Prdx2	ENSRNOG0000003520	peroxiredoxin 2 Source RGD Symbol Acc 3838
E07	UPFR1089192	ENSRNOT00000040203.4	Psen2	ENSRNOG0000002879	presenilin 2 Source RGD Symbol Acc 621060
E08	UPFR1043890	ENSRNOT00000028143.3	Pten	ENSRNOG0000020723	phosphatase and tensin homolog Source RGD Symbol Acc 61995
E09	UPFR1118025	ENSRNOT00000003774.2	Rgs4	ENSRNOG0000002773	regulator of G-protein signaling 4 Source RGD Symbol Acc 3567
E10	UPFR1030105	ENSRNOT00000006373.4	Rtn1	ENSRNOG0000004794	reticulin 1 Source RGD Symbol Acc 620986
E11	UPFR1048763	ENSRNOT00000001743.2	S100b	ENSRNOG0000001295	S100 calcium binding protein B Source RGD Symbol Acc 3615
E12	UPFR1117229	ENSRNOT00000007676.7	Skp1	ENSRNOG0000005828	S-phase kinase-associated protein 1 Source RGD Symbol Acc 1359648
F01	UPFR1012902	ENSRNOT00000011983.2	Slc18a2	ENSRNOG0000008890	solute carrier family 18 member A2 Source RGD Symbol Acc 3694
F02	UPFR1053922	ENSRNOT00000014704.4	Slc25a4	ENSRNOG0000010830	solute carrier family 25 member 4 Source RGD Symbol Acc 620352
F03	UPFR1108894	ENSRNOT00000040291.1	Slc6a3	ENSRNOG0000017302	solute carrier family 6 member 3 Source RGD Symbol Acc 3715
F04	UPFR1064093	ENSRNOT00000035415.6	Slit1	ENSRNOG0000026065	slit guidance ligand 1 Source RGD Symbol Acc 69307
F05	UPFR1068205	ENSRNOT00000039247.4	Snca	ENSRNOG0000008656	synuclein alpha Source RGD Symbol Acc 3729
F06	UPFR1103041	ENSRNOT00000067517.3	Spen	ENSRNOG0000033556	spen family transcriptional repressor Source RGD Symbol Acc 1589867
F07	UPFR1062562	ENSRNOT00000038448.5	Srsf7	ENSRNOG0000027360	serine and arginine rich splicing factor 7 Source RGD Symbol Acc 1307425
F08	UPFR1091827	ENSRNOT00000026921.6	Stub1	ENSRNOG0000019798	STIP1 homology and U-box containing protein 1 Source RGD Symbol Acc 1306167
F09	UPFR1110731	ENSRNOT00000015084.6	AC126641.1	ENSRNOG0000011160	synaptic vesicle glycoprotein 2b Source NCBI gene Acc 117556
F10	UPFR1118921	ENSRNOT00000017833.5	Syngr3	ENSRNOG0000012620	synaptogyrin 3 Source RGD Symbol Acc 1311801

Position	Assay	Name	Symbol	Ensembl ID	Description
F11	UPFR1065705	ENSRNOT00000088950.1	Syt1	ENSRNOG0000006426	synaptotagmin 1 Source RGD Symbol Acc 3803
F12	UPFR1077061	ENSRNOT00000027474.3	Syt11	ENSRNOG0000020279	synaptotagmin 11 Source RGD Symbol Acc 62042
G01	UPFR1096861	ENSRNOT00000019681.6	Tcf7l2	ENSRNOG0000049232	transcription factor 7 like 2 Source RGD Symbol Acc 1583621
G02	UPFR1024503	ENSRNOT00000027682.5	Th	ENSRNOG0000020410	tyrosine hydroxylase Source RGD Symbol Acc 3853
G03	UPFR1041949	ENSRNOT00000089491.1	Tpbp	ENSRNOG0000010694	trophoblast glycoprotein Source RGD Symbol Acc 621453
G04	UPFR1120795	ENSRNOT00000031115.4	Uba1	ENSRNOG0000019164	ubiquitin-like modifier activating enzyme 1 Source RGD Symbol Acc 1359327
G05	UPFR1051928	ENSRNOT00000066885.3	Ubb	ENSRNOG0000042271	ubiquitin B Source RGD Symbol Acc 621562
G06	UPFR1022099	ENSRNOT00000024406.6	Ube2i	ENSRNOG0000017907	ubiquitin-conjugating enzyme E2I Source RGD Symbol Acc 3926
G07	UPFR1077947	ENSRNOT00000030485.5	Ube2k	ENSRNOG0000027088	ubiquitin-conjugating enzyme E2K Source RGD Symbol Acc 1311626
G08	UPFR1017507	ENSRNOT00000046159.4	Ube2l3	ENSRNOG000001862	ubiquitin-conjugating enzyme E2L 3 Source RGD Symbol Acc 1308361
G09	UPFR1110967	ENSRNOT00000003248.5	Uchl1	ENSRNOG0000002343	ubiquitin C-terminal hydrolase L1 Source RGD Symbol Acc 3928
G10	UPFR1029388	ENSRNOT00000025081.5	Vamp1	ENSRNOG0000019219	vesicle-associated membrane protein 1 Source RGD Symbol Acc 3948
G11	UPFR1052803	ENSRNOT00000026197.6	Vdac3	ENSRNOG0000019277	voltage-dependent anion channel 3 Source RGD Symbol Acc 621577
G12	UPFR1099982	ENSRNOT00000080676.1	Ywhaz	ENSRNOG0000008195	tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, zeta Source RGD Symbol Acc 3980
H01	UPFR1132952	ENSRNOT00000080216.1	Actb	ENSRNOG0000034254	actin, beta Source RGD Symbol Acc 628837
H02	UPFR1132953	ENSRNOT00000023017.5	B2m	ENSRNOG0000017123	beta-2 microglobulin Source RGD Symbol Acc 2189
H03	UPFR1132959	ENSRNOT00000065935.3	Hprt1	ENSRNOG0000048561	hypoxanthine phosphoribosyltransferase 1 Source RGD Symbol Acc 2826
H04	UPFR1018740	ENSRNOT00000017468.2	Ldha	ENSRNOG0000013009	lactate dehydrogenase A Source RGD Symbol Acc 2996
H05	UPFR1132958	ENSRNOT00000018820.5	Rplp1	ENSRNOG0000013874	ribosomal protein lateral stalk subunit P1 Source RGD Symbol Acc 621774
H06	UPFR1126610	UPL_RGDC	RGDC	UPL_RGDC	Rat Genomic DNA Contamination
H07	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H08	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H09	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H10	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control
H11	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control
H12	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control

Related products

Product	Contents	Cat. no.
QuantiNova LNA Probe PCR QC Panel	These panels are designed to assess the quality of RNA samples before characterization using QuantiNova LNA Probe PCR Focus Panels; available in 96-well, 384-well, and Rotor-Disc 100 formats	249945
QuantiNova Reverse Transcription Kit (10)*	For 10 x 20 μ l reactions: 20 μ l 8x gDNA Removal Mix, 10 μ l Reverse Transcription Enzyme, 40 μ l Reverse Transcription Mix (containing RT primers), 20 μ l Internal Control RNA, 1.9 ml RNase-Free Water	205410
QuantiNova Probe RT-PCR Kit (100)*	For 100 x 20 μ l reactions: 1 ml QuantiNova Probe RT-PCR Master Mix, 20 μ l QuantiNova Probe RT Mix, 20 μ l Internal Control RNA, 500 μ l Yellow Template Dilution Buffer, 250 μ l ROX Reference Dye, 1.9 μ l RNase-Free Water	208352
QuantiNova Probe PCR Kit (100)*	For 100 x 20 μ l reactions: 1 ml 2x QuantiNova Probe PCR Master Mix, 500 μ l QuantiNova Yellow Template Dilution Buffer, 250 μ l QN ROX Reference Dye, 1.9 ml Water	208252

*Larger kit sizes available.

The QuantiNova LNA Probe PCR Focus Panels are intended for molecular biology applications. These products are not intended for the diagnosis, prevention or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, LNA®, QuantiNova®, Sample to Insight® (QIAGEN Group); SYBR® (Life Technologies Corp.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

09/2019 © 2019 QIAGEN, all rights reserved.