

QuantiNova® LNA® Probe PCR Focus Panels (Rotor-Gene® Format)

Rat Cytokines & Chemokines

Cat. no. 249955 UPRN-150ZR

For study focus gene expression analysis

Shipping and storage

QuantiNova LNA Probe PCR Focus Panels are shipped at room temperature. Immediately upon receipt, they should be stored protected from light at 2–8°C for short term storage or at –30°C to –15°C for long time storage. Under these conditions, all components are stable for at least 12 months.

Note: Open the package and store the products appropriately immediately upon receipt.

For optimal performance, QuantiNova LNA Probe PCR Focus Panels should be used together with the QuantiNova Reverse Transcription Kit for cDNA synthesis and the QuantiNova Probe PCR Kit (Mastermix) for PCR.

Panel layout (Rotor-Gene): QuantiNova LNA Probe PCR Focus Panel

The 96 real-time assays in the Rotor-Gene format are located in wells 1–96 of the Rotor-Disc® (plate A1–A12=Rotor-Disc 1–12, plate B1–B12=Rotor-Disc 13–24, etc.). To maintain data analysis compatibility, wells 97–100 do not contain real-time assays but will contain master mix to account for weight balance. Refer to the QuantiNova LNA Probe PCR Handbook at www.qiagen.com for further details.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Adipoq	Bmp2	Bmp4	Bmp6	Bmp7	C5	Cd1	Cd11	Cd12	Cd17	Cd19	Cd2
B	Cd20	Cd21	Cd22	Cd24	Cd3	Cd4	Cd5	Cd7	Cd40lg	Cd70	Cntf	Csf1
C	Csf2	Csf3	Cf1	Cx3cl1	Cxcl1	Cxcl10	Cxcl11	Cxcl12	Cxcl13	Cxcl16	Cxcl3	Cxcl9
D	Faslg	Gdf15	Gpi	Ifna2	Ifng	Il10	Il11	Il12a	Il12b	Il13	Il15	Il16
E	Il17a	Il17f	Il18	Il1a	Il1b	Il1m	Il2	Il21	Il22	Il23a	Il24	Il27
F	Il3	Il4	Il5	Il6	Il7	Il9	Lif	Lta	LOC1036943 81	LOC1036948 77	Mtn	Nodal
G	Osm	Pf4	Pbp	Spp1	Tgfb2	Thpo1	LOC1036943 80	Tnfrsf11b	Tnfrsf10	LOC1036899 74	Vegfa	Xcl1
H	Actb	B2m	Hprt1	Ldha	Rplp1	RGDC	QIC	QIC	QIC	PPC	PPC	PPC

Gene table: QuantiNova LNA Probe PCR Focus Panel

Position	Assay	Name	Symbol	Ensembl ID	Description
A01	UPFR1094195	ENSRNOT00000089988.1	Adipoq	ENSRNOG0000001821	adiponectin, C1Q and collagen domain containing Source RGD Symbol Acc 628748
A02	UPFR1023820	ENSRNOT00000028904.4	Bmp2	ENSRNOG00000021276	bone morphogenetic protein 2 Source RGD Symbol Acc 2211
A03	UPFR1013450	ENSRNOT00000083268.1	Bmp4	ENSRNOG0000009694	bone morphogenetic protein 4 Source RGD Symbol Acc 2213
A04	UPFR1048825	ENSRNOT00000018359.7	Bmp6	ENSRNOG0000013717	bone morphogenetic protein 6 Source RGD Symbol Acc 2214
A05	UPFR1105126	ENSRNOT00000084990.1	Bmp7	ENSRNOG0000053384	bone morphogenetic protein 7 Source RGD Symbol Acc 620743
A06	UPFR1034130	ENSRNOT00000025534.7	C5	ENSRNOG0000018899	complement C5 Source RGD Symbol Acc 2237
A07	UPFR1032143	ENSRNOT00000031360.3	Ccl1	ENSRNOG00000021851	C-C motif chemokine ligand 1 Source RGD Symbol Acc 1586495
A08	UPFR1034165	ENSRNOT00000009756.5	Ccl11	ENSRNOG0000007335	C-C motif chemokine ligand 11 Source RGD Symbol Acc 3644
A09	UPFR1080637	ENSRNOT00000031626.5	Ccl12	ENSRNOG0000029768	chemokine (C-C motif) ligand 12 Source RGD Symbol Acc 1309255
A10	UPFR1085781	ENSRNOT00000021842.2	Ccl17	ENSRNOG0000016278	C-C motif chemokine ligand 17 Source RGD Symbol Acc 619924
A11	UPFR1074099	ENSRNOT00000020982.5	Ccl19	ENSRNOG0000015668	C-C motif chemokine ligand 19 Source RGD Symbol Acc 1310336
A12	UPFR1014050	ENSRNOT00000009448.3	Ccl2	ENSRNOG0000007159	C-C motif chemokine ligand 2 Source RGD Symbol Acc 3645
B01	UPFR1113129	ENSRNOT00000021730.3	Ccl20	ENSRNOG0000015992	C-C motif chemokine ligand 20 Source RGD Symbol Acc 3646
B02	UPFR1022549	ENSRNOT00000023951.6	Ccl21	ENSRNOG00000034290	C-C motif chemokine ligand 21 Source RGD Symbol Acc 1311996
B03	UPFR1070569	ENSRNOT00000022167.3	Ccl22	ENSRNOG0000016535	C-C motif chemokine ligand 22 Source RGD Symbol Acc 619933
B04	UPFR1058995	ENSRNOT00000045157.3	Ccl24	ENSRNOG00000031162	C-C motif chemokine ligand 24 Source RGD Symbol Acc 1310245
B05	UPFR1057376	ENSRNOT00000015139.3	Ccl3	ENSRNOG0000011205	C-C motif chemokine ligand 3 Source RGD Symbol Acc 3647
B06	UPFR1110783	ENSRNOT00000015199.3	Ccl4	ENSRNOG0000011406	C-C motif chemokine ligand 4 Source RGD Symbol Acc 620441
B07	UPFR1019085	ENSRNOT00000014865.6	Ccl5	ENSRNOG0000010906	C-C motif chemokine ligand 5 Source RGD Symbol Acc 69069
B08	UPFR1083626	ENSRNOT00000000256.4	Ccl7	ENSRNOG0000000239	C-C motif chemokine ligand 7 Source RGD Symbol Acc 1359152
B09	UPFR1038741	ENSRNOT00000001162.4	Cd40lg	ENSRNOG0000000871	CD40 ligand Source RGD Symbol Acc 708418
B10	UPFR1090606	ENSRNOT00000070806.1	Cd70	ENSRNOG0000051015	Cd70 molecule Source RGD Symbol Acc 1588573
B11	UPFR1111031	ENSRNOT00000016690.2	Cnlf	ENSRNOG0000012460	ciliary neurotrophic factor Source RGD Symbol Acc 2370
B12	UPFR1022984	ENSRNOT00000081835.1	Csf1	ENSRNOG0000018659	colony stimulating factor 1 Source RGD Symbol Acc 621063
C01	UPFR1084630	ENSRNOT00000032333.3	Csf2	ENSRNOG0000026805	colony stimulating factor 2 Source RGD Symbol Acc 621065
C02	UPFR1037448	ENSRNOT00000011509.3	Csf3	ENSRNOG0000008525	colony stimulating factor 3 Source RGD Symbol Acc 2426
C03	UPFR1081638	ENSRNOT00000025646.2	Ctf1	ENSRNOG0000018962	cardiotrophin 1 Source RGD Symbol Acc 2442
C04	UPFR1070853	ENSRNOT00000022128.5	Cx3cl1	ENSRNOG0000016326	C-X3-C motif chemokine ligand 1 Source RGD Symbol Acc 620458
C05	UPFR1056296	ENSRNOT00000003778.3	Cxcl1	ENSRNOG0000002802	C-X-C motif chemokine ligand 1 Source RGD Symbol Acc 619869
C06	UPFR1049597	ENSRNOT00000003075.7	Cxcl10	ENSRNOG0000022256	C-X-C motif chemokine ligand 10 Source RGD Symbol Acc 620209
C07	UPFR1033909	ENSRNOT00000031667.4	Cxcl11	ENSRNOG0000022298	C-X-C motif chemokine ligand 11 Source RGD Symbol Acc 727827
C08	UPFR1044005	ENSRNOT00000066670.3	Cxcl12	ENSRNOG0000013589	C-X-C motif chemokine ligand 12 Source RGD Symbol Acc 3651
C09	UPFR1043327	ENSRNOT00000039383.4	Cxcl13	ENSRNOG0000024899	C-X-C motif chemokine ligand 13 Source RGD Symbol Acc 1564256
C10	UPFR1076114	ENSRNOT00000032926.4	Cxcl16	ENSRNOG0000026647	C-X-C motif chemokine ligand 16 Source RGD Symbol Acc 1562928
		ENSRNOT000000		ENSRNOG00	

Position	Assay	Name	Symbol	Ensembl ID	Description
C11	UPFR1076775	034090.3	Cxcl3	000028043	chemokine (C-X-C motif) ligand 3 Source RGD Symbol Acc 621812
C12	UPFR1089821	ENSRNOT00000003082.7	Cxcl9	ENSRNOG000002242	C-X-C motif chemokine ligand 9 Source RGD Symbol Acc 628798
D01	UPFR1108356	ENSRNOT00000003998.2	Faslg	ENSRNOG0000002978	Fas ligand Source RGD Symbol Acc 3880
D02	UPFR1099357	ENSRNOT00000026652.3	Gdf15	ENSRNOG0000019661	growth differentiation factor 15 Source RGD Symbol Acc 2674
D03	UPFR1020266	ENSRNOT00000032613.4	Gpi	ENSRNOG0000023150	glucose-6-phosphate isomerase Source RGD Symbol Acc 2727
D04	UPFR1086334	ENSRNOT00000045433.4	Ifna2	ENSRNOG0000058365	interferon alpha 2 Source RGD Symbol Acc 1311425
D05	UPFR1084134	ENSRNOT00000009919.2	Ifng	ENSRNOG0000007468	interferon gamma Source RGD Symbol Acc 2866
D06	UPFR1113235	ENSRNOT00000006246.5	Il10	ENSRNOG0000004647	interleukin 10 Source RGD Symbol Acc 2886
D07	UPFR1065504	ENSRNOT00000023489.5	Il11	ENSRNOG0000017386	interleukin 11 Source RGD Symbol Acc 621475
D08	UPFR1065603	ENSRNOT00000012831.3	Il12a	ENSRNOG0000009468	interleukin 12A Source RGD Symbol Acc 620562
D09	UPFR1060283	ENSRNOT00000005898.6	Il12b	ENSRNOG0000004380	interleukin 12B Source RGD Symbol Acc 628704
D10	UPFR1038252	ENSRNOT00000010121.2	Il13	ENSRNOG0000007652	interleukin 13 Source RGD Symbol Acc 68949
D11	UPFR1102147	ENSRNOT00000004590.6	Il15	ENSRNOG0000003439	interleukin 15 Source RGD Symbol Acc 2887
D12	UPFR1122227	ENSRNOT00000016289.4	Il16	ENSRNOG0000011680	interleukin 16 Source RGD Symbol Acc 1583872
E01	UPFR1062668	ENSRNOT00000016664.4	Il17a	ENSRNOG0000012467	interleukin 17A Source RGD Symbol Acc 2888
E02	UPFR1042305	ENSRNOT00000016717.5	Il17f	ENSRNOG0000012509	interleukin 17F Source RGD Symbol Acc 1311398
E03	UPFR1053861	ENSRNOT00000013093.6	Il18	ENSRNOG0000009848	interleukin 18 Source RGD Symbol Acc 2889
E04	UPFR1066959	ENSRNOT00000080818.1	Il1a	ENSRNOG0000004575	interleukin 1 alpha Source RGD Symbol Acc 2890
E05	UPFR1121351	ENSRNOT00000006308.4	Il1b	ENSRNOG0000004649	interleukin 1 beta Source RGD Symbol Acc 2891
E06	UPFR1086251	ENSRNOT00000007949.6	Il1rn	ENSRNOG0000005871	interleukin 1 receptor antagonist Source RGD Symbol Acc 621159
E07	UPFR1070791	ENSRNOT00000023327.2	Il2	ENSRNOG0000017348	interleukin 2 Source RGD Symbol Acc 620047
E08	UPFR1051958	ENSRNOT00000023348.4	Il21	ENSRNOG0000017376	interleukin 21 Source RGD Symbol Acc 1307384
E09	UPFR1058955	ENSRNOT00000009776.4	Il22	ENSRNOG0000007365	interleukin 22 Source RGD Symbol Acc 1561292
E10	UPFR1075411	ENSRNOT00000004363.3	Il23a	ENSRNOG0000003254	interleukin 23 subunit alpha Source RGD Symbol Acc 620873
E11	UPFR1023656	ENSRNOT00000006106.3	Il24	ENSRNOG0000004470	interleukin 24 Source RGD Symbol Acc 621484
E12	UPFR1104994	ENSRNOT00000025950.5	Il27	ENSRNOG0000019170	interleukin 27 Source RGD Symbol Acc 1561420
F01	UPFR1073307	ENSRNOT00000031762.1	Il3	ENSRNOG0000026786	interleukin 3 Source RGD Symbol Acc 2897
F02	UPFR1099283	ENSRNOT00000010029.3	Il4	ENSRNOG0000007624	interleukin 4 Source RGD Symbol Acc 2898
F03	UPFR1086524	ENSRNOT00000010729.4	Il5	ENSRNOG0000008111	interleukin 5 Source RGD Symbol Acc 2900
F04	UPFR1098910	ENSRNOT00000013732.6	Il6	ENSRNOG0000010278	interleukin 6 Source RGD Symbol Acc 2901
F05	UPFR1101991	ENSRNOT00000016091.7	Il7	ENSRNOG0000011973	interleukin 7 Source RGD Symbol Acc 2904
F06	UPFR1021636	ENSRNOT00000016187.2	Il9	ENSRNOG0000012131	interleukin 9 Source RGD Symbol Acc 620049
F07	UPFR1017798	ENSRNOT00000009313.4	Lif	ENSRNOG0000007002	LIF, interleukin 6 family cytokine Source RGD Symbol Acc 620865
F08	UPFR1116910	ENSRNOT00000088936.1	Lta	ENSRNOG0000000838	lymphotoxin alpha Source RGD Symbol Acc 3020
F09	UPFR1050435	ENSRNOT00000089411.1	LOC103694381	ENSRNOG0000058566	lymphotoxin-beta Source RGD Symbol Acc 9083628
F10	UPFR1126542	ENSRNOT00000008608.4	LOC103694877	ENSRNOG0000006589	macrophage migration inhibitory factor Source RGD Symbol Acc 621163

Position	Assay	Name	Symbol	Ensembl ID	Description
F11	UPFR1114916	ENSRNOT00000038093.3	Mstn	ENSRNOG00000021294	myostatin Source RGD Symbol Acc 3115
F12	UPFR1057301	ENSRNOT00000000672.4	Nodal	ENSRNOG0000000556	nodal growth differentiation factor Source RGD Symbol Acc 1305994
G01	UPFR1058346	ENSRNOT00000035591.3	Osm	ENSRNOG00000024390	oncostatin M Source RGD Symbol Acc 1585012
G02	UPFR1099077	ENSRNOT00000033406.1	Pf4	ENSRNOG00000028015	platelet factor 4 Source RGD Symbol Acc 3305
G03	UPFR1100499	ENSRNOT00000003794.3	Pbp	ENSRNOG0000002829	pro-platelet basic protein Source RGD Symbol Acc 70548
G04	UPFR1100110	ENSRNOT00000067875.4	Spp1	ENSRNOG00000043451	secreted phosphoprotein 1 Source RGD Symbol Acc 3752
G05	UPFR1046456	ENSRNOT00000087023.1	Tgfb2	ENSRNOG0000002418	transforming growth factor, beta 2 Source RGD Symbol Acc 70491
G06	UPFR1100728	ENSRNOT00000073743.2	Thpo1	ENSRNOG00000050172	thrombopoietin like 1 Source RGD Symbol Acc 11455920
G07	UPFR1098022	ENSRNOT00000079677.1	LOC103694380	ENSRNOG00000055156	tumor necrosis factor-like Source RGD Symbol Acc 9404643
G08	UPFR1015258	ENSRNOT00000011344.5	Tnfrsf11b	ENSRNOG0000008336	TNF receptor superfamily member 11B Source RGD Symbol Acc 619802
G09	UPFR1107113	ENSRNOT00000017758.5	Tnfsf10	ENSRNOG00000013269	TNF superfamily member 10 Source RGD Symbol Acc 628734
G10	UPFR1105006	ENSRNOT00000049056.4	LOC103689974	ENSRNOG00000049547	TNF superfamily member 11 Source RGD Symbol Acc 620784
G11	UPFR1117908	ENSRNOT00000026559.5	Vegfa	ENSRNOG00000019598	vascular endothelial growth factor A Source RGD Symbol Acc 619991
G12	UPFR1029900	ENSRNOT00000003976.3	Xcl1	ENSRNOG0000002964	X-C motif chemokine ligand 1 Source RGD Symbol Acc 620452
H01	UPFR1132952	ENSRNOT00000080216.1	Actb	ENSRNOG00000034254	actin, beta Source RGD Symbol Acc 628837
H02	UPFR1132953	ENSRNOT00000023017.5	B2m	ENSRNOG00000017123	beta-2 microglobulin Source RGD Symbol Acc 2189
H03	UPFR1132959	ENSRNOT00000065935.3	Hprt1	ENSRNOG00000048561	hypoxanthine phosphoribosyltransferase 1 Source RGD Symbol Acc 2826
H04	UPFR1018740	ENSRNOT00000017468.2	Ldha	ENSRNOG00000013009	lactate dehydrogenase A Source RGD Symbol Acc 2996
H05	UPFR1132958	ENSRNOT00000018820.5	Rplp1	ENSRNOG00000013874	ribosomal protein lateral stalk subunit P1 Source RGD Symbol Acc 621774
H06	UPFR1126610	UPL_RGDC	RGDC	UPL_RGDC	Rat Genomic DNA Contamination
H07	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H08	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H09	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H10	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control
H11	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control
H12	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control

Related products

Product	Contents	Cat. no.
QuantiNova LNA Probe PCR QC Panel	These panels are designed to assess the quality of RNA samples before characterization using QuantiNova LNA Probe PCR Focus Panels; available in 96-well, 384-well, and Rotor-Disc 100 formats	249945
QuantiNova Reverse Transcription Kit (10)*	For 10 x 20 μ l reactions: 20 μ l 8x gDNA Removal Mix, 10 μ l Reverse Transcription Enzyme, 40 μ l Reverse Transcription Mix (containing RT primers), 20 μ l Internal Control RNA, 1.9 ml RNase-Free Water	205410
QuantiNova Probe RT-PCR Kit (100)*	For 100 x 20 μ l reactions: 1 ml QuantiNova Probe RT-PCR Master Mix, 20 μ l QuantiNova Probe RT Mix, 20 μ l Internal Control RNA, 500 μ l Yellow Template Dilution Buffer, 250 μ l ROX Reference Dye, 1.9 μ l RNase-Free Water	208352
QuantiNova Probe PCR Kit (100)*	For 100 x 20 μ l reactions: 1 ml 2x QuantiNova Probe PCR Master Mix, 500 μ l QuantiNova Yellow Template Dilution Buffer, 250 μ l QN ROX Reference Dye, 1.9 ml Water	208252

*Larger kit sizes available.

The QuantiNova LNA Probe PCR Focus Panels are intended for molecular biology applications. These products are not intended for the diagnosis, prevention or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, LNA®, QuantiNova®, Sample to Insight® (QIAGEN Group); SYBR® (Life Technologies Corp.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

09/2019 © 2019 QIAGEN, all rights reserved.