

QuantiNova® LNA® Probe PCR Focus Panels (96-Well Format and 384-Well [4 x 96] Format)

Rat WNT Signaling Targets

Cat. no. 249955 UPRN-243ZA

For study focus gene expression analysis

Shipping and storage

QuantiNova LNA Probe PCR Focus Panels are shipped at room temperature. Immediately upon receipt, they should be stored protected from light at 2–8°C for short term storage or at –30°C to –15°C for long time storage. Under these conditions, all components are stable for at least 12 months.

Note: Open the package and store the products appropriately immediately upon receipt.

For optimal performance, QuantiNova LNA Probe PCR Focus Panels should be used together with the QuantiNova Reverse Transcription Kit for cDNA synthesis and the QuantiNova Probe PCR Kit (Mastermix) for PCR.

Panel layout (96-well): QuantiNova LNA Probe PCR Focus Panel

For the 384-well (4 × 96) PCR panels, genes are present in a staggered format. Refer to the QuantiNova LNA Probe PCR Handbook at www.qiagen.com for further details.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Abcb1a	Ahr	Angptl4	Anthr1	Axin2	Bglap	Birc5	Bmp4	Bltrc	Cacna2d3	Ccnd1	Ccnd2
B	Cd44	Cdh1	Cdkn2a	Cdon	Cebpd	Ccn2	Cubn	Dab2	Dkk1	Dlk1	Dpp10	Efnb1
C	Egfr	Egr1	Enpp2	Ets2	Fgf20	Fgf4	Fgf7	Fgf9	Fn1	Fosl1	Fstl	Gdf5
D	Gdnf	Gja1	Id2	Igf1	Igf2	Il6	Irs1	Jag1	Klf5	Lef1	Lrp1	Met
E	Mmp2	Mmp7	Mmp9	Myc	Nanog	Nrcam	Nrp1	Ntrk2	Pdgfra	Pitx2	Plaur	Pou5f1
F	Ppp3	Ppard	Ptch1	Ptgs2	Runx2	Sfrp2	Six1	Smo	Sox2	Sox9	Tbst	Tcf7l1
G	Tcf4	Tcf7	Tcf7l2	Tgfb3	Tle1	Twist1	Vegfa	Ccn4	Ccn5	Wnt3a	Wnt5a	Wnt9a
H	Actb	B2m	Hprt1	Ldha	Rplp1	RGDC	QIC	QIC	QIC	PPC	PPC	PPC

Gene table: QuantiNova LNA Probe PCR Focus Panel

Position	Assay	Name	Symbol	Ensembl ID	Description
A01	UPFR1045141	ENSRNOT00000037712.5	Abcb1a	ENSRNOG0000008012	ATP binding cassette subfamily B member 4 Source RGD Symbol Acc 620248
A02	UPFR1105470	ENSRNOT00000006618.7	Ahr	ENSRNOG0000004342	aryl hydrocarbon receptor Source RGD Symbol Acc 2074
A03	UPFR1032714	ENSRNOT00000010031.6	Angptl4	ENSRNOG0000007545	angiopoietin-like 4 Source RGD Symbol Acc 735058
A04	UPFR1041604	ENSRNOT00000011675.4	Antxr1	ENSRNOG0000008678	ANTXR cell adhesion molecule 1 Source RGD Symbol Acc 1307144
A05	UPFR1102592	ENSRNOT00000088599.1	Axin2	ENSRNOG0000055010	axin 2 Source RGD Symbol Acc 69259
A06	UPFR1053867	ENSRNOT00000026530.3	Bglap	ENSRNOG0000019607	bone gamma-carboxylglutamate protein Source RGD Symbol Acc 2206
A07	UPFR1053500	ENSRNOT00000072387.1	Birc5	ENSRNOG0000050819	baculoviral IAP repeat-containing 5 Source RGD Symbol Acc 70499
A08	UPFR1013450	ENSRNOT00000083268.1	Bmp4	ENSRNOG0000009694	bone morphogenetic protein 4 Source RGD Symbol Acc 2213
A09	UPFR1066655	ENSRNOT00000022426.7	Birc	ENSRNOG0000016280	beta-transducin repeat containing E3 ubiquitin protein ligase Source RGD Symbol Acc 1359721
A10	UPFR1023487	ENSRNOT00000048043.3	Cacna2d3	ENSRNOG0000031287	calcium voltage-gated channel auxiliary subunit alpha2delta 3 Source RGD Symbol Acc 631361
A11	UPFR1042457	ENSRNOT00000088588.1	Ccnd1	ENSRNOG0000020918	cyclin D1 Source RGD Symbol Acc 68384
A12	UPFR1088052	ENSRNOT00000086440.1	Ccnd2	ENSRNOG0000057710	cyclin D2 Source RGD Symbol Acc 621083
B01	UPFR1040363	ENSRNOT00000009073.7	Cd44	ENSRNOG0000006094	CD44 molecule (Indian blood group) Source RGD Symbol Acc 2307
B02	UPFR1037315	ENSRNOT00000027346.2	Cdh1	ENSRNOG0000020151	cadherin 1 Source RGD Symbol Acc 69279
B03	UPFR1104432	ENSRNOT00000079251.1	Cdkn2a	ENSRNOG0000059837	cyclin-dependent kinase inhibitor 2A Source RGD Symbol Acc 2323
B04	UPFR1021608	ENSRNOT00000016248.4	Cdon	ENSRNOG0000011789	cell adhesion associated, oncogene regulated Source RGD Symbol Acc 708433
B05	UPFR1073105	ENSRNOT00000074586.2	Cebpd	ENSRNOG0000050869	CCAAT/enhancer binding protein delta Source RGD Symbol Acc 2328
B06	UPFR1047569	ENSRNOT00000089196.1	Ccn2	ENSRNOG0000015036	cellular communication network factor 2 Source RGD Symbol Acc 621392
B07	UPFR1021724	ENSRNOT00000090064.1	Cubn	ENSRNOG0000029047	cubilin Source RGD Symbol Acc 68355
B08	UPFR1110844	ENSRNOT00000050655.5	Dab2	ENSRNOG0000028930	DAB2, clathrin adaptor protein Source RGD Symbol Acc 621007
B09	UPFR1075162	ENSRNOT00000015771.7	Dkk1	ENSRNOG0000011692	dickkopf WNT signaling pathway inhibitor 1 Source RGD Symbol Acc 1307313
B10	UPFR1103659	ENSRNOT00000006339.6	Dlk1	ENSRNOG0000019584	delta like non-canonical Notch ligand 1 Source RGD Symbol Acc 619931
B11	UPFR1115875	ENSRNOT00000003527.5	Dpp10	ENSRNOG0000002595	dipeptidylpeptidase 10 Source RGD Symbol Acc 1306427
B12	UPFR1117680	ENSRNOT00000009635.6	Efnb1	ENSRNOG0000006877	ephrin B1 Source RGD Symbol Acc 2540
C01	UPFR1091984	ENSRNOT00000006087.2	Egfr	ENSRNOG0000004332	epidermal growth factor receptor Source RGD Symbol Acc 2543
C02	UPFR1034050	ENSRNOT00000026303.4	Egr1	ENSRNOG0000019422	early growth response 1 Source RGD Symbol Acc 2544
C03	UPFR1030168	ENSRNOT00000051139.5	Enpp2	ENSRNOG0000004089	ectonucleotide pyrophosphatase/phosphodiesterase 2 Source RGD Symbol Acc 69298
C04	UPFR1103791	ENSRNOT00000002247.6	Ets2	ENSRNOG0000001647	ETS proto-oncogene 2, transcription factor Source RGD Symbol Acc 1584977
C05	UPFR1083724	ENSRNOT00000000121.2	Fgf20	ENSRNOG0000000109	fibroblast growth factor 20 Source RGD Symbol Acc 71068
C06	UPFR1084715	ENSRNOT00000028355.5	Fgf4	ENSRNOG0000020890	fibroblast growth factor 4 Source RGD Symbol Acc 620127
C07	UPFR1021146	ENSRNOT00000012700.4	Fgf7	ENSRNOG0000009425	fibroblast growth factor 7 Source RGD Symbol Acc 61805
C08	UPFR1020537	ENSRNOT00000015367.5	Fgf9	ENSRNOG0000011471	fibroblast growth factor 9 Source RGD Symbol Acc 2610
C09	UPFR1089380	ENSRNOT00000057585.4	Fn1	ENSRNOG0000014288	fibronectin 1 Source RGD Symbol Acc 2624
C10	UPFR1040543	ENSRNOT00000027891.2	Fosl1	ENSRNOG0000020552	FOS like 1, AP-1 transcription factor subunit Source RGD Symbol Acc 2627
		ENSRNOT000000		ENSRNOG00	

Position	Assay	Name	Symbol	Ensembl ID	Description
C11	UPFR1120134	015680.4	Fst	000011631	folliculin Source RGD Symbol Acc 2633
C12	UPFR1064525	ENSRNOT00000073736.3	Gdf5	ENSRNOG0000050123	growth differentiation factor 5 Source RGD Symbol Acc 620102
D01	UPFR1111724	ENSRNOT00000041288.2	Gdnf	ENSRNOG0000012819	glial cell derived neurotrophic factor Source RGD Symbol Acc 2677
D02	UPFR1122020	ENSRNOT00000001054.4	Gja1	ENSRNOG0000000805	gap junction protein, alpha 1 Source RGD Symbol Acc 2690
D03	UPFR1117393	ENSRNOT00000009491.6	Id2	ENSRNOG0000007237	inhibitor of DNA binding 2 Source RGD Symbol Acc 2859
D04	UPFR1050099	ENSRNOT00000038780.6	Igf1	ENSRNOG0000004517	insulin-like growth factor 1 Source RGD Symbol Acc 2868
D05	UPFR1085899	ENSRNOT00000089171.1	Igf2	ENSRNOG0000020369	insulin-like growth factor 2 Source RGD Symbol Acc 2870
D06	UPFR1098910	ENSRNOT00000013732.6	Il6	ENSRNOG0000010278	interleukin 6 Source RGD Symbol Acc 2901
D07	UPFR1099858	ENSRNOT00000019579.5	Irs1	ENSRNOG0000014597	insulin receptor substrate 1 Source RGD Symbol Acc 2922
D08	UPFR1082062	ENSRNOT00000010638.7	Jag1	ENSRNOG0000007443	jagged 1 Source RGD Symbol Acc 2937
D09	UPFR1080710	ENSRNOT00000057843.2	Klf5	ENSRNOG0000008785	Kruppel-like factor 5 Source RGD Symbol Acc 621446
D10	UPFR1076567	ENSRNOT00000013694.3	Lef1	ENSRNOG0000010121	lymphoid enhancer binding factor 1 Source RGD Symbol Acc 620241
D11	UPFR1046102	ENSRNOT00000031005.5	Lrp1	ENSRNOG0000025053	LDL receptor related protein 1 Source RGD Symbol Acc 1307535
D12	UPFR1071342	ENSRNOT00000086902.1	Met	ENSRNOG0000052745	MET proto-oncogene, receptor tyrosine kinase Source RGD Symbol Acc 3082
E01	UPFR1080151	ENSRNOT00000022679.6	Mmp2	ENSRNOG0000016695	matrix metalloproteinase 2 Source RGD Symbol Acc 621316
E02	UPFR1034539	ENSRNOT00000014041.5	Mmp7	ENSRNOG0000010507	matrix metalloproteinase 7 Source RGD Symbol Acc 3100
E03	UPFR1032957	ENSRNOT00000023965.3	Mmp9	ENSRNOG0000017539	matrix metalloproteinase 9 Source RGD Symbol Acc 621320
E04	UPFR1101439	ENSRNOT00000006188.5	Myc	ENSRNOG0000004500	MYC proto-oncogene, bHLH transcription factor Source RGD Symbol Acc 3130
E05	UPFR1091249	ENSRNOT00000060937.3	Nanog	ENSRNOG0000008368	Nanog homeobox Source RGD Symbol Acc 1303178
E06	UPFR1059699	ENSRNOT00000083088.2	Nrcam	ENSRNOG0000004067	neuronal cell adhesion molecule Source RGD Symbol Acc 3209
E07	UPFR1025300	ENSRNOT00000080465.1	Nrp1	ENSRNOG0000010744	neuropilin 1 Source RGD Symbol Acc 621588
E08	UPFR1075665	ENSRNOT00000042145.2	Ntrk2	ENSRNOG0000018839	neurotrophic receptor tyrosine kinase 2 Source RGD Symbol Acc 3213
E09	UPFR1089489	ENSRNOT00000003077.5	Pdgfra	ENSRNOG0000002244	platelet derived growth factor receptor alpha Source RGD Symbol Acc 3284
E10	UPFR1089140	ENSRNOT00000051009.4	Pitx2	ENSRNOG0000010681	paired-like homeodomain 2 Source RGD Symbol Acc 3331
E11	UPFR1041940	ENSRNOT00000075167.3	Plaur	ENSRNOG0000037931	plasminogen activator, urokinase receptor Source RGD Symbol Acc 620597
E12	UPFR1099582	ENSRNOT00000072288.2	Pou5f1	ENSRNOG0000046487	POU class 5 homeobox 1 Source RGD Symbol Acc 1359491
F01	UPFR1081880	ENSRNOT00000011237.5	Plpp3	ENSRNOG0000008116	phospholipid phosphatase 3 Source RGD Symbol Acc 620454
F02	UPFR1082868	ENSRNOT00000042539.5	Ppard	ENSRNOG0000000503	peroxisome proliferator-activated receptor delta Source RGD Symbol Acc 3370
F03	UPFR1054272	ENSRNOT00000086130.1	Ptch1	ENSRNOG0000019354	patched 1 Source RGD Symbol Acc 621425
F04	UPFR1035449	ENSRNOT00000003567.4	Ptgs2	ENSRNOG0000002525	prostaglandin-endoperoxide synthase 2 Source RGD Symbol Acc 620349
F05	UPFR1051266	ENSRNOT00000061014.4	Runx2	ENSRNOG0000020193	runt-related transcription factor 2 Source RGD Symbol Acc 2282
F06	UPFR1086653	ENSRNOT00000012658.6	Sfrp2	ENSRNOG0000009465	secreted frizzled-related protein 2 Source RGD Symbol Acc 735163
F07	UPFR1053538	ENSRNOT00000034338.4	Six1	ENSRNOG0000022777	SIX homeobox 1 Source RGD Symbol Acc 620906
F08	UPFR1032074	ENSRNOT00000080223.1	Smo	ENSRNOG0000008332	smoothed, frizzled class receptor Source RGD Symbol Acc 3726
F09	UPFR1042145	ENSRNOT00000016236.3	Sox2	ENSRNOG0000012199	SRY box 2 Source RGD Symbol Acc 1565646
F10	UPFR1013060	ENSRNOT00000003511.6	Sox9	ENSRNOG0000002607	SRY box 9 Source RGD Symbol Acc 620474

Position	Assay	Name	Symbol	Ensembl ID	Description
F11	UPFR1021899	ENSRNOT00000016471.3	Tbxt	ENSRNOG0000012229	T-box transcription factor T Source RGD Symbol Acc 1310141
F12	UPFR1026727	ENSRNOT00000020005.7	Tcf7l1	ENSRNOG0000014753	transcription factor 7 like 1 Source RGD Symbol Acc 1311671
G01	UPFR1028213	ENSRNOT00000081797.1	Tcf4	ENSRNOG0000012405	transcription factor 4 Source RGD Symbol Acc 69271
G02	UPFR1041848	ENSRNOT00000008022.7	Tcf7	ENSRNOG0000005872	transcription factor 7 Source RGD Symbol Acc 1305894
G03	UPFR1096861	ENSRNOT00000019681.6	Tcf7l2	ENSRNOG0000049232	transcription factor 7 like 2 Source RGD Symbol Acc 1583621
G04	UPFR1030751	ENSRNOT00000013516.5	Tgfb3	ENSRNOG0000009867	transforming growth factor, beta 3 Source RGD Symbol Acc 3851
G05	UPFR1113015	ENSRNOT00000058906.3	Tle1	ENSRNOG0000005882	TLE family member 1, transcriptional corepressor Source RGD Symbol Acc 1309237
G06	UPFR1038411	ENSRNOT00000014763.6	Twist1	ENSRNOG0000011101	twist family bHLH transcription factor 1 Source RGD Symbol Acc 621455
G07	UPFR1117908	ENSRNOT00000026559.5	Vegfa	ENSRNOG0000019598	vascular endothelial growth factor A Source RGD Symbol Acc 619991
G08	UPFR1109359	ENSRNOT00000009673.7	Ccn4	ENSRNOG0000007078	cellular communication network factor 4 Source RGD Symbol Acc 69431
G09	UPFR1121952	ENSRNOT00000014346.6	Ccn5	ENSRNOG0000010666	cellular communication network factor 5 Source RGD Symbol Acc 621867
G10	UPFR1043063	ENSRNOT00000064505.3	Wnt3a	ENSRNOG0000003039	Wnt family member 3A Source RGD Symbol Acc 1308057
G11	UPFR1100946	ENSRNOT00000021164.3	Wnt5a	ENSRNOG0000015618	Wnt family member 5A Source RGD Symbol Acc 69250
G12	UPFR1073122	ENSRNOT00000058327.3	Wnt9a	ENSRNOG0000003066	Wnt family member 9A Source RGD Symbol Acc 1305018
H01	UPFR1132952	ENSRNOT00000080216.1	Actb	ENSRNOG00000034254	actin, beta Source RGD Symbol Acc 628837
H02	UPFR1132953	ENSRNOT00000023017.5	B2m	ENSRNOG0000017123	beta-2 microglobulin Source RGD Symbol Acc 2189
H03	UPFR1132959	ENSRNOT00000065935.3	Hprt1	ENSRNOG0000048561	hypoxanthine phosphoribosyltransferase 1 Source RGD Symbol Acc 2826
H04	UPFR1018740	ENSRNOT00000017468.2	Ldha	ENSRNOG0000013009	lactate dehydrogenase A Source RGD Symbol Acc 2996
H05	UPFR1132958	ENSRNOT00000018820.5	Rplp1	ENSRNOG0000013874	ribosomal protein lateral stalk subunit P1 Source RGD Symbol Acc 621774
H06	UPFR1126610	UPL_RGDC	RGDC	UPL_RGDC	Rat Genomic DNA Contamination
H07	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H08	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H09	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H10	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control
H11	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control
H12	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control

Related products

Product	Contents	Cat. no.
QuantiNova LNA Probe PCR QC Panel	These panels are designed to assess the quality of RNA samples before characterization using QuantiNova LNA Probe PCR Focus Panels; available in 96-well, 384-well, and Rotor-Disc 100 formats	249945
QuantiNova Reverse Transcription Kit (10)*	For 10 x 20 μ l reactions: 20 μ l 8x gDNA Removal Mix, 10 μ l Reverse Transcription Enzyme, 40 μ l Reverse Transcription Mix (containing RT primers), 20 μ l Internal Control RNA, 1.9 ml RNase-Free Water	205410
QuantiNova Probe RT-PCR Kit (100)*	For 100 x 20 μ l reactions: 1 ml QuantiNova Probe RT-PCR Master Mix, 20 μ l QuantiNova Probe RT Mix, 20 μ l Internal Control RNA, 500 μ l Yellow Template Dilution Buffer, 250 μ l ROX Reference Dye, 1.9 μ l RNase-Free Water	208352
QuantiNova Probe PCR Kit (100)*	For 100 x 20 μ l reactions: 1 ml 2x QuantiNova Probe PCR Master Mix, 500 μ l QuantiNova Yellow Template Dilution Buffer, 250 μ l QN ROX Reference Dye, 1.9 ml Water	208252

*Larger kit sizes available.

The QuantiNova LNA Probe PCR Focus Panels are intended for molecular biology applications. These products are not intended for the diagnosis, prevention or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, LNA®, QuantiNova®, Sample to Insight® (QIAGEN Group); SYBR® (Life Technologies Corp.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

09/2019 © 2019 QIAGEN, all rights reserved.