

miRCURY LNA™ miRNA Focus PCR Panels

Human Brain Cancer Product Data Sheet

Cat. no. 339325 YAHS-208Z

For mature miRNA expression profiling using real-time PCR

Format	Suitable real-time cyclers	Plate	Cat. no.
A	Applied Biosystems® models 5700, 7000, 7300, 7500, 7700, 7900HT, ViiA™ 7 (96-well block); Bio-Rad® models iCycler®, iQ™5, MyiQ™, MyiQ2; Bio-Rad/MJ Research Chromo4™; Eppendorf® MasterCycler® ep realplex models 2, 2s, 4, 4s; Stratagene® models Mx3005P®, Mx3000P®; Takara: TP-800	96-well	YAHS-208ZA
C	Applied Biosystems models 7500 (Fast block), 7900HT (Fast block), StepOnePlus™, ViiA 7 (Fast block)	96-well	YAHS-208ZC
D	Bio-Rad CFX96™; Bio-Rad/MJ Research models DNA Engine Opticon®, DNA Engine Opticon 2; Stratagene Mx4000®	96-well	YAHS-208ZD
E	Applied Biosystems® models 7900HT (384-well block), ViiA™ 7 (384-well block); Bio-Rad CFX384™	384-well	YAHS-208ZE
F	Roche® LightCycler® 480 (96-well block)	96-well	YAHS-208ZF
G	Roche® LightCycler® 480 (384-well block)	384-well	YAHS-208ZG

Description

The Human Brain Cancer miRCURY LNA™ miRNA Focus PCR Panel profiles the expression of 84 miRNAs known to alter their expression during nervous system-related carcinogenesis. This array provides cancer researchers with a convenient way to quickly analyze the miRNAs most relevant to brain and central nervous system tumorigenesis. Brain tumors, which cause uncontrolled cellular growth in the cranium or spinal column, involve any of the cell types found in the brain, such as neurons, glial cells, astrocytes, or cells of the meninges. By attenuating gene expression, miRNAs regulate many biological functions, and some, such as those in the miR-17-92 cluster, have known oncogenic properties. Expression studies of different brain tumor types have identified dysregulated miRNA; however, few are functionally linked to pathogenic mechanisms. Research into the functional consequences of miRNA dysregulation during neurocarcinogenesis may yield insights into the molecular mechanisms behind the pathogenesis of brain tumors. The array contains miRNAs with known functional consequences as well as miRNAs with a significant dysregulation but unknown function during oncogenesis. The profiling results from this array may yield insights into the molecular mechanisms behind the pathogenesis of multiple nervous system cancers. A set of controls present on this array enables data analysis using the $\Delta\Delta CT$ method of relative quantification, assessment of reverse transcription performance, and assessment of PCR performance. Using SYBR Green-based real-time PCR, the expression of a focused panel of miRNAs related to brain cancer can be easily and reliably analyzed with this miRCURY LNA™ miRNA Focus PCR Panel.

For further details, consult the *miRCURY LNA™ miRNA Focus PCR Panels Handbook*.

Array Layout

	1	2	3	4	5	6	7	8	9	10	11	12
A	hsa-lef-7b-5p	hsa-miR-101-3p	hsa-miR-106b-5p	hsa-miR-107	hsa-miR-10b-5p	hsa-miR-124-3p	hsa-miR-125a-5p	hsa-miR-125b-5p	hsa-miR-127-5p	hsa-miR-128-3p	hsa-miR-129-5p	hsa-miR-130a-3p
B	hsa-miR-130b-3p	hsa-miR-132-3p	hsa-miR-133a-3p	hsa-miR-133b	hsa-miR-137-3p	hsa-miR-138-5p	hsa-miR-141-3p	hsa-miR-144-3p	hsa-miR-146a-5p	hsa-miR-146b-5p	hsa-miR-148a-3p	hsa-miR-149-5p
C	hsa-miR-150-5p	hsa-miR-153-3p	hsa-miR-15a-5p	hsa-miR-15b-5p	hsa-miR-16-5p	hsa-miR-17-5p	hsa-miR-17-3p	hsa-miR-181a-5p	hsa-miR-181b-5p	hsa-miR-182-5p	hsa-miR-183-5p	hsa-miR-184
D	hsa-miR-185-5p	hsa-miR-187-3p	hsa-miR-18a-5p	hsa-miR-190a-5p	hsa-miR-191-5p	hsa-miR-19a-3p	hsa-miR-19b-3p	hsa-miR-200a-3p	hsa-miR-203a-3p	hsa-miR-20a-5p	hsa-miR-21-5p	hsa-miR-210-3p
E	hsa-miR-216a-5p	hsa-miR-217-5p	hsa-miR-218-5p	hsa-miR-221-3p	hsa-miR-222-3p	hsa-miR-222-5p	hsa-miR-23a-3p	hsa-miR-23b-3p	hsa-miR-24-3p	hsa-miR-25-3p	hsa-miR-26a-5p	hsa-miR-27a-3p
F	hsa-miR-296-5p	hsa-miR-29a-3p	hsa-miR-29b-3p	hsa-miR-29c-3p	hsa-miR-30b-5p	hsa-miR-30c-5p	hsa-miR-31-5p	hsa-miR-320a-3p	hsa-miR-323a-5p	hsa-miR-324-5p	hsa-miR-326	hsa-miR-328-3p
G	hsa-miR-331-5p	hsa-miR-335-5p	hsa-miR-34a-5p	hsa-miR-425-5p	hsa-miR-451a	hsa-miR-486-5p	hsa-miR-7-5p	hsa-miR-9-5p	hsa-miR-9-3p	hsa-miR-92a-3p	hsa-miR-93-5p	hsa-miR-96-5p
H	cel-miR-39-3p	cel-miR-39-3p	SNORD44 (hsa)	SNORD388 (hsa)	SNORD49A (hsa)	U6 snRNA (v2)	UniSp2	UniSp4	UniSp5	UniSp6	UniSp3	UniSp3

miRNA Table

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
A01	hsa-let-7b-5p	MIMAT0000063	YP00204750	E01	hsa-miR-216a-5p	MIMAT0000273	YP00204167
A02	hsa-miR-101-3p	MIMAT0000099	YP00204786	E02	hsa-miR-217-5p	MIMAT0000274	YP00204010
A03	hsa-miR-106b-5p	MIMAT0000680	YP00205884	E03	hsa-miR-218-5p	MIMAT0000275	YP00206034
A04	hsa-miR-107	MIMAT0000104	YP00204468	E04	hsa-miR-221-3p	MIMAT0000278	YP00204532
A05	hsa-miR-10b-5p	MIMAT0000254	YP00205637	E05	hsa-miR-222-3p	MIMAT0000279	YP00204551
A06	hsa-miR-124-3p	MIMAT0000422	YP00206026	E06	hsa-miR-222-5p	MIMAT0004569	YP00204314
A07	hsa-miR-125a-5p	MIMAT0000443	YP00204339	E07	hsa-miR-23a-3p	MIMAT0000078	YP00204772
A08	hsa-miR-125b-5p	MIMAT0000423	YP00205713	E08	hsa-miR-23b-3p	MIMAT0000418	YP02119314
A09	hsa-miR-127-5p	MIMAT0004604	YP00204161	E09	hsa-miR-24-3p	MIMAT0000080	YP00204260
A10	hsa-miR-128-3p	MIMAT0000424	YP00205995	E10	hsa-miR-25-3p	MIMAT0000081	YP00204361
A11	hsa-miR-129-5p	MIMAT0000242	YP00204534	E11	hsa-miR-26a-5p	MIMAT0000082	YP00206023
A12	hsa-miR-130a-3p	MIMAT0000425	YP00204658	E12	hsa-miR-27a-3p	MIMAT0000084	YP00206038
B01	hsa-miR-130b-3p	MIMAT0000691	YP00204317	F01	hsa-miR-296-5p	MIMAT0000690	YP00204436
B02	hsa-miR-132-3p	MIMAT0000426	YP00206035	F02	hsa-miR-29a-3p	MIMAT0000086	YP00204698
B03	hsa-miR-133a-3p	MIMAT0000427	YP00204788	F03	hsa-miR-29b-3p	MIMAT0000100	YP00204679
B04	hsa-miR-133b	MIMAT0000770	YP00206058	F04	hsa-miR-29c-3p	MIMAT0000681	YP00204729
B05	hsa-miR-137-3p	MIMAT0000429	YP00206062	F05	hsa-miR-30b-5p	MIMAT0000420	YP00204765
B06	hsa-miR-138-5p	MIMAT0000430	YP00206078	F06	hsa-miR-30c-5p	MIMAT0000244	YP00204783
B07	hsa-miR-141-3p	MIMAT0000432	YP00204504	F07	hsa-miR-31-5p	MIMAT0000089	YP00204236
B08	hsa-miR-144-3p	MIMAT0000436	YP00204754	F08	hsa-miR-320a-3p	MIMAT0000510	YP00206042
B09	hsa-miR-146a-5p	MIMAT0000449	YP00204688	F09	hsa-miR-323a-5p	MIMAT0004696	YP00206080
B10	hsa-miR-146b-5p	MIMAT0002809	YP02119310	F10	hsa-miR-324-5p	MIMAT0000761	YP02119473
B11	hsa-miR-148a-3p	MIMAT0000243	YP00205867	F11	hsa-miR-326	MIMAT0000756	YP00204512
B12	hsa-miR-149-5p	MIMAT0000450	YP00204321	F12	hsa-miR-328-3p	MIMAT0000752	YP00204364
C01	hsa-miR-150-5p	MIMAT0000451	YP00204660	G01	hsa-miR-331-5p	MIMAT0004700	YP00204423
C02	hsa-miR-153-3p	MIMAT0000439	YP00204338	G02	hsa-miR-335-5p	MIMAT0000765	YP02119293
C03	hsa-miR-15a-5p	MIMAT0000068	YP00204066	G03	hsa-miR-34a-5p	MIMAT0000255	YP00204486
C04	hsa-miR-15b-5p	MIMAT0000417	YP00204243	G04	hsa-miR-425-5p	MIMAT0003393	YP00204337
C05	hsa-miR-16-5p	MIMAT0000069	YP00205702	G05	hsa-miR-451a	MIMAT0001631	YP02119305
C06	hsa-miR-17-5p	MIMAT0000070	YP02119304	G06	hsa-miR-486-5p	MIMAT0002177	YP00204001
C07	hsa-miR-17-3p	MIMAT0000071	YP00206008	G07	hsa-miR-7-5p	MIMAT0000252	YP02119317
C08	hsa-miR-181a-5p	MIMAT0000256	YP00206081	G08	hsa-miR-9-5p	MIMAT0000441	YP00204513
C09	hsa-miR-181b-5p	MIMAT0000257	YP00204530	G09	hsa-miR-9-3p	MIMAT0000442	YP00204620
C10	hsa-miR-182-5p	MIMAT0000259	YP00206070	G10	hsa-miR-92a-3p	MIMAT0000092	YP00204258
C11	hsa-miR-183-5p	MIMAT0000261	YP00206030	G11	hsa-miR-93-5p	MIMAT0000093	YP00204715
C12	hsa-miR-184	MIMAT0000454	YP00204601	G12	hsa-miR-96-5p	MIMAT0000095	YP00204417
D01	hsa-miR-185-5p	MIMAT0000455	YP00206037	H01	cel-miR-39-3p	MIMAT0000010	YP00203952
D02	hsa-miR-187-3p	MIMAT0000262	YP00204018	H02	cel-miR-39-3p	MIMAT0000010	YP00203952
D03	hsa-miR-18a-5p	MIMAT0000072	YP00204207	H03	SNORD44 (hsa)	N/A	YP00203902
D04	hsa-miR-190a-5p	MIMAT0000458	YP00204763	H04	SNORD38B (hsa)	N/A	YP00203901
D05	hsa-miR-191-5p	MIMAT0000440	YP00204306	H05	SNORD49A (hsa)	N/A	YP00203904
D06	hsa-miR-19a-3p	MIMAT0000073	YP00205862	H06	U6 snRNA (v2)	N/A	YP02119464
D07	hsa-miR-19b-3p	MIMAT0000074	YP00204450	H07	UniSp2	N/A	YP00203950
D08	hsa-miR-200a-3p	MIMAT0000682	YP00204707	H08	UniSp4	N/A	YP00203953
D09	hsa-miR-203a-3p	MIMAT0000264	YP00205914	H09	UniSp5	N/A	YP00203955
D10	hsa-miR-20a-5p	MIMAT0000075	YP00204292	H10	UniSp6	N/A	YP00203954
D11	hsa-miR-21-5p	MIMAT0000076	YP00204230	H11	UniSP3	N/A	YP02119288
D12	hsa-miR-210-3p	MIMAT0000267	YP00204333	H12	UniSP3	N/A	YP02119288

Functional Groupings

Anaplastic Astrocytomas: hsa-miR-107, hsa-miR-10b-5p, hsa-miR-124-3p, hsa-miR-125b-5p, hsa-miR-127-5p, hsa-miR-129-5p, hsa-miR-137-3p, hsa-miR-138-5p, hsa-miR-187-3p, hsa-miR-203a-3p, hsa-miR-218-5p, hsa-miR-296-5p, hsa-miR-31-5p, hsa-miR-320a-3p, hsa-miR-7-5p.

Glioblastomas: hsa-miR-101-3p, hsa-miR-107, hsa-miR-10b-5p, hsa-miR-124-3p, hsa-miR-128-3p, hsa-miR-129-5p, hsa-miR-130a-3p, hsa-miR-132-3p, hsa-miR-133a-3p, hsa-miR-133b, hsa-miR-137-3p, hsa-miR-138-5p, hsa-miR-149-5p, hsa-miR-153-3p, hsa-miR-16-5p, hsa-miR-181a-5p, hsa-miR-181b-5p, hsa-miR-185-5p, hsa-miR-187-3p, hsa-miR-191-5p, hsa-miR-203a-3p, hsa-miR-21-5p, hsa-miR-210-3p, hsa-miR-218-5p, hsa-miR-221-3p, hsa-miR-222-3p, hsa-miR-222-5p, hsa-miR-23a-3p, hsa-miR-23b-3p, hsa-miR-25-3p, hsa-miR-29b-3p, hsa-miR-31-5p, hsa-miR-323a-5p, hsa-miR-328-3p, hsa-miR-34a-5p, hsa-miR-425-5p, hsa-miR-451a, hsa-miR-486-5p, hsa-miR-7-5p.

Gliomas: hsa-miR-106b-5p, hsa-miR-10b-5p, hsa-miR-128-3p, hsa-miR-130b-3p, hsa-miR-144-3p, hsa-miR-146b-5p, hsa-miR-148a-3p, hsa-miR-15a-5p, hsa-miR-15b-5p, hsa-miR-16-5p, hsa-miR-17-5p, hsa-miR-181a-5p, hsa-miR-181b-5p, hsa-miR-182-5p, hsa-miR-183-5p, hsa-miR-184, hsa-miR-19a-3p, hsa-miR-200a-3p, hsa-miR-20a-5p, hsa-miR-21-5p, hsa-miR-210-3p, hsa-miR-221-3p, hsa-miR-222-3p, hsa-miR-23a-3p, hsa-miR-25-3p, hsa-miR-26a-5p, hsa-miR-27a-3p, hsa-miR-296-5p, hsa-miR-328-3p, hsa-miR-34a-5p, hsa-miR-451a, hsa-miR-486-5p, hsa-miR-9-3p, hsa-miR-9-5p, hsa-miR-93-5p.

Medulloblastomas: hsa-miR-106b-5p, hsa-miR-124-3p, hsa-miR-125a-5p, hsa-miR-125b-5p, hsa-miR-127-5p, hsa-miR-128-3p, hsa-miR-129-5p, hsa-miR-130a-3p, hsa-miR-130b-3p, hsa-miR-133b, hsa-miR-137-3p, hsa-miR-138-5p, hsa-miR-148a-3p, hsa-miR-150-5p, hsa-miR-153-3p, hsa-miR-17-3p, hsa-miR-17-5p, hsa-miR-187-3p, hsa-miR-18a-5p, hsa-miR-190a-5p, hsa-miR-191-5p, hsa-miR-19a-3p, hsa-miR-19b-3p, hsa-miR-203a-3p, hsa-miR-20a-5p, hsa-miR-210-3p, hsa-miR-216a-5p, hsa-miR-217-5p, hsa-miR-218-5p, hsa-miR-25-3p, hsa-miR-29a-3p, hsa-miR-29b-3p, hsa-miR-30b-5p, hsa-miR-30c-5p, hsa-miR-320a-3p, hsa-miR-323a-5p, hsa-miR-324-5p, hsa-miR-326, hsa-miR-328-3p, hsa-miR-331-5p, hsa-miR-335-5p, hsa-miR-34a-5p, hsa-miR-425-5p, hsa-miR-7-5p, hsa-miR-9-3p, hsa-miR-92a-3p, hsa-miR-93-5p, hsa-miR-96-5p.

Meningiomas: hsa-let-7b-5p, hsa-miR-106b-5p, hsa-miR-125a-5p, hsa-miR-153-3p, hsa-miR-15a-5p, hsa-miR-17-5p, hsa-miR-181a-5p, hsa-miR-185-5p, hsa-miR-19b-3p, hsa-miR-200a-3p, hsa-miR-217-5p, hsa-miR-23b-3p, hsa-miR-26a-5p, hsa-miR-27a-3p, hsa-miR-29a-3p, hsa-miR-335-5p, hsa-miR-451a.

Neuroblastomas: hsa-let-7b-5p, hsa-miR-107, hsa-miR-125a-5p, hsa-miR-125b-5p, hsa-miR-128-3p, hsa-miR-129-5p, hsa-miR-137-3p, hsa-miR-141-3p, hsa-miR-146a-5p, hsa-miR-149-5p, hsa-miR-150-5p, hsa-miR-17-3p, hsa-miR-17-5p, hsa-miR-181a-5p, hsa-miR-181b-5p, hsa-miR-184, hsa-miR-187-3p, hsa-miR-190a-5p, hsa-miR-200a-3p, hsa-miR-216a-5p, hsa-miR-221-3p, hsa-miR-23b-3p, hsa-miR-24-3p, hsa-miR-29a-3p, hsa-miR-29b-3p, hsa-miR-29c-3p, hsa-miR-30b-5p, hsa-miR-30c-5p, hsa-miR-320a-3p, hsa-miR-323a-5p, hsa-miR-324-5p, hsa-miR-326, hsa-miR-331-5p, hsa-miR-335-5p, hsa-miR-34a-5p, hsa-miR-9-5p, hsa-miR-92a-3p, hsa-miR-93-5p.

Retinoblastomas: hsa-miR-320a-3p.

Ordering Information

Product	Contents	Cat. no.
miRCURY LNA miRNA Focus PCR Panels	miRCURY LNA miRNA PCR Panels for application-based miRNome profiling, available in 96-well or 384-well format; for SYBR® Green-based detection	339325
miRCURY LNA miRNA miRNome PCR Panels	miRCURY LNA miRNA PCR Panels for PCR-based miRNome profiling, available in 384-well format; for SYBR® Green-based detection	339322
miRCURY LNA miRNA QC PCR Panel	miRCURY LNA miRNA PCR Panel of quality control assays, available in 96-well or 384-well format; for SYBR® Green-based detection	339331
miRCURY LNA miRNA Custom PCR Panels	8 identical, ready-to-use 96- or 384-well plates; each well contains primers sufficient for one 10 µl reaction; for SYBR® Green-based detection	339330
miRCURY LNA Custom PCR Panel Additional Plate	Additional miRCURY LNA Custom PCR Panel plates; set of 4 plates; only available in addition to the base plates ordered through the core product (cat. no. 339330)	339332
miRCURY LNA miRNA PCR Assays	Contains forward and reverse primers for 200 SYBR® Green-based, real-time qPCR reactions, 166 EvaGreen-based digital PCR reactions for Nanoplate 8.5k or 50 EvaGreen-based digital PCR reactions for Nanoplate 26k	339306
miRCURY LNA miRNA Probe PCR Assays	Complete premixed assays containing LNA-enhanced target-specific forward primer and probe. For 200 reactions.	339350
miRCURY LNA miRNA Custom Probe PCR Assays	Custom-designed, target-specific forward primer and probe for any user-defined miRNA target. Complete premixed assay for 200 reactions.	339351

Related Products

Product	Contents	Cat. no.
miRCURY LNA RT Kit	For 8–64 cDNA synthesis reactions: 5x RT SYBR Green Reaction Buffer, 5x RT Probe Reaction Buffer, 10x RT Enzyme Mix, UniSp6, RNA Spike-in template, RNase-Free Water	339340
RNA Spike-In Kit, For RT	Contains the UniSp2, UniSp4, and UniSp5 RNA Spike-in Template Mix and the cel-miR-39-3p RNA Spike-in Template	339390
miRCURY LNA SYBR® Green PCR Kits (200)	For 200 reactions: 2X miRCURY SYBR Green Master Mix, RNase-Free Water	339345
miRCURY LNA SYBR® Green PCR Kits (600)	For 600 reactions: 2X miRCURY SYBR Green Master Mix, RNase-Free Water	339346
miRCURY LNA SYBR® Green PCR Kits (4000)	For 4000 reactions: 2X miRCURY SYBR Green Master Mix, RNase-Free Water	339347
miRCURY LNA Probe PCR Kit (200)	For 200 reactions: 2X QuantiNova Probe Master Mix, 10X miRCURY Probe Univ. Primer, Rox Reference Dye, RNase-Free Water	339371
miRCURY LNA Probe PCR Kit (800)	For 800 reactions: 2X QuantiNova Probe Master Mix, 10X miRCURY Probe Univ. Primer, Rox Reference Dye, RNase-Free Water	339372
miRCURY LNA Probe PCR Kit (4000)	For 4000 reactions: 2X QuantiNova Probe Master Mix, 10X miRCURY Probe Univ. Primer, Rox Reference Dye, RNase-Free Water	339373
miRCURY LNA miRNA PCR Starter Kit	Two miRCURY LNA PCR Assays of your choice, spike-in control Assay (UniSp6), one candidate endogenous control assay (miR-103-3p) and all reagents for 20 reverse transcription reactions and 100 PCR amplifications; for SYBR® Green-based qPCR detection	339320

Important

The miRCURY LNA™ miRNA PCR Panels are Ready-to-Use and designed for optimal performance with the miRCURY LNA RT Kit and the miRCURY LNA SYBR® Green PCR Kit. The performance of the primer sets will be affected when used in combination with less than optimal reagents. Use the miRCURY LNA miRNA PCR Panel Handbook for experiment setup. RNA work requires specific handling and precautions should be taken to prevent RNase contamination and degradation of the RNA sample and reagents.

Shipping and storage

The Ready-to-Use miRCURY LNA miRNA PCR Panels are shipped at room temperature and can be stored at 4°C for at least 6 months. For long term storage, it is recommended to place the panels at -20°C. Under these conditions, the LNA PCR primers are stable for at least 6 months after receipt.

Intended use

miRCURY LNA miRNA PCR Assays are intended for molecular biology applications. This product is not intended for the diagnosis, prevention or treatment of a disease.

Safety information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). These are available online in convenient and compact PDF format at www.qiagen.com/safety where you can find, view and print the SDS for each QIAGEN kit and kit component.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, Sample to Insight®, GeneGlobe®, LNA™ (QIAGEN Group), SYBR® (Life Technologies Corporation). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law. 1107212 05/2017 HB-2376-001 © 2017 QIAGEN, all rights reserved.
