

miRCURY LNA miRNA miRNome PCR Panels

Human Panel 1-7, V6 Product Data Sheet

Cat. no. 339322 YAHS-370Z

For mature miRNA expression profiling using real-time PCR

Format	Suitable real-time cyclers	Plate	Cat. no.
E	Applied Biosystems® models 7900HT (384-well block), ViiA™ 7 (384-well block); Bio-Rad CFX384™	384-well	YAHS-370ZE
G	Roche® LightCycler® 480 (384-well block)	384-well	YAHS-370ZG

Description

The miRCURY LNA miRNA miRNome PCR Panels profile the expression of hundreds of the most abundantly expressed and best characterized miRNA sequences in the a given species-specific miRNome as annotated in miRBase (www.mirbase.org). Although they are well characterized, each of these miRNA sequences can regulate one to several messenger RNA transcripts, and conversely one mRNA can be regulated by one to several miRNA sequences. Therefore, the complex role that any given known miRNA sequence plays has yet to be completely defined. Use of the miRCURY LNA miRNA miRNome PCR Panels maximizes the likelihood of discovering miRNA sequences whose expression patterns correlate with the biological phenotypes under study. A set of controls included on each plate enables data analysis using the $\Delta\Delta CT$ method of relative quantification, assessment of reverse transcription performance, and assessment of PCR performance. The Human miRNome miScript miRNA PCR Array enables easy and reliable, SYBR Green-based real-time PCR analysis of the expression of 1008 miRNA sequences and discovery of those most important to the research area of interest.

For further details, consult the *miRCURY LNA™ miRCURY LNA™ miRNA Focus PCR Panels Handbook*.

miRNA Table: Plate 1

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
A01	hsa-miR-7-5p	MIMAT0000252	YP02119317	I01	hsa-miR-19b-3p	MIMAT0000074	YP00204450
A02	hsa-miR-217-5p	MIMAT0000274	YP00204010	I02	hsa-miR-526b-5p	MIMAT0002835	YP00204259
A03	hsa-miR-337-5p	MIMAT0004695	YP00204391	I03	hsa-miR-215-5p	MIMAT0000272	YP00204598
A04	hsa-miR-328-3p	MIMAT0000752	YP00204364	I04	hsa-miR-30b-5p	MIMAT0000420	YP00204765
A05	hsa-miR-374b-3p	MIMAT0004956	YP00205917	I05	hsa-miR-184	MIMAT0000454	YP00204601
A06	hsa-miR-143-3p	MIMAT0000435	YP00205992	I06	hsa-miR-422a	MIMAT0001339	YP00204164
A07	hsa-miR-623	MIMAT0003292	YP00204019	I07	hsa-miR-199a-3p	MIMAT0000232	YP00204536
A08	hsa-miR-520c-3p	MIMAT0002846	YP00204497	I08	hsa-miR-335-5p	MIMAT0000765	YP02119293
A09	hsa-miR-557	MIMAT0003221	YP00204346	I09	hsa-miR-519a-3p	MIMAT0002869	YP00205919
A10	hsa-miR-218-5p	MIMAT0000275	YP00206034	I10	hsa-miR-21-5p	MIMAT0000076	YP00204230
A11	hsa-miR-136-5p	MIMAT0000448	YP00204779	I11	hsa-miR-129-2-3p	MIMAT0004605	YP00206067
A12	hsa-miR-127-5p	MIMAT0004604	YP00204161	I12	hsa-miR-26b-5p	MIMAT0000083	YP00204172
A13	hsa-miR-140-5p	MIMAT0000431	YP00204540	I13	hsa-miR-214-3p	MIMAT0000271	YP00204510
A14	hsa-miR-31-3p	MIMAT0004504	YP00204079	I14	hsa-miR-32-5p	MIMAT0000090	YP00204792
A15	hsa-miR-20b-3p	MIMAT0004752	YP00204569	I15	hsa-miR-324-3p	MIMAT0000762	YP00204303
A16	hsa-miR-325	MIMAT0000771	YP00204464	I16	hsa-miR-488-3p	MIMAT0004763	YP00204469
A17	hsa-miR-509-3-5p	MIMAT0004975	YP00204503	I17	hsa-miR-371a-5p	MIMAT0004687	YP00204493
A18	hsa-miR-210-3p	MIMAT0000267	YP00204333	I18	hsa-miR-455-5p	MIMAT0003150	YP00204363
A19	hsa-miR-199b-5p	MIMAT0000263	YP00204152	I19	hsa-miR-891a-5p	MIMAT0004902	YP00204220
A20	hsa-miR-194-5p	MIMAT0000460	YP00204080	I20	hsa-miR-549a-3p	MIMAT0003333	YP00205975
A21	hsa-let-7g-5p	MIMAT0000414	YP00204565	I21	hsa-miR-205-5p	MIMAT0000266	YP00204487
A22	hsa-miR-203a-3p	MIMAT0000264	YP00205914	I22	hsa-miR-518b	MIMAT0002844	YP00204405
A23	hsa-miR-181a-3p	MIMAT0000270	YP00204110	I23	hsa-miR-19a-3p	MIMAT0000073	YP00205862
A24	hsa-miR-137-3p	MIMAT0000429	YP00206062	I24	hsa-miR-150-5p	MIMAT0000451	YP00204660
B01	hsa-miR-551b-3p	MIMAT0003233	YP00204067	J01	hsa-miR-15a-5p	MIMAT0000068	YP00204066
B02	BlankSpot	N/A	YP02119291	J02	hsa-let-7d-3p	MIMAT0004484	YP00205627
B03	hsa-miR-524-3p	MIMAT0002850	YP00204030	J03	hsa-miR-608	MIMAT0003276	YP00204286
B04	UniSp2	N/A	YP00203950	J04	hsa-miR-671-5p	MIMAT0003880	YP00205649
B05	hsa-miR-486-5p	MIMAT0002177	YP00204001	J05	hsa-miR-497-5p	MIMAT0002820	YP00204354
B06	UniSp4	N/A	YP00203953	J06	hsa-miR-877-5p	MIMAT0004949	YP00205626
B07	UniSp3	N/A	YP02119288	J07	hsa-miR-187-5p	MIMAT0004561	YP00205920
B08	UniSp5	N/A	YP00203955	J08	hsa-miR-10b-5p	MIMAT0000254	YP00205637
B09	hsa-miR-329-3p	MIMAT0001629	YP00206052	J09	hsa-let-7i-5p	MIMAT0000415	YP00204394
B10	UniSp6	N/A	YP00203954	J10	hsa-miR-202-5p	MIMAT0002810	YP00205616
B11	hsa-miR-487b-3p	MIMAT0003180	YP00204489	J11	hsa-miR-652-3p	MIMAT0003322	YP00204387
B12	cel-miR-39-3p	MIMAT0000010	YP00203952	J12	hsa-miR-126-5p	MIMAT0000444	YP00206010
B13	hsa-miR-138-5p	MIMAT0000430	YP00206078	J13	hsa-miR-30e-3p	MIMAT0000693	YP00204410
B14	hsa-miR-191-5p	MIMAT0000440	YP00204306	J14	hsa-miR-181c-5p	MIMAT0000258	YP00204683
B15	hsa-miR-378a-3p	MIMAT0000732	YP00205946	J15	hsa-miR-9-3p	MIMAT0000442	YP00204620
B16	hsa-miR-103a-3p	MIMAT0000101	YP00204063	J16	hsa-miR-548c-3p	MIMAT0003285	YP00204697
B17	hsa-miR-890	MIMAT0004912	YP00205859	J17	hsa-miR-152-3p	MIMAT0000438	YP00204294
B18	hsa-miR-423-5p	MIMAT0004748	YP00205624	J18	hsa-miR-93-5p	MIMAT0000093	YP00204715
B19	hsa-miR-221-3p	MIMAT0000278	YP00204532	J19	hsa-miR-365a-3p	MIMAT0000710	YP00204622
B20	SNORD38B (hsa)	N/A	YP00203901	J20	hsa-miR-29c-3p	MIMAT0000681	YP00204729
B21	hsa-miR-301b-3p	MIMAT0004958	YP00204390	J21	hsa-miR-372-3p	MIMAT0000724	YP00204137
B22	SNORD49A (hsa)	N/A	YP00203904	J22	hsa-miR-133a-3p	MIMAT0000427	YP00204788
B23	hsa-miR-550a-5p	MIMAT0004800	YP00204638	J23	hsa-miR-124-3p	MIMAT0000422	YP00206026
B24	U6 snRNA (v2)	N/A	YP02119464	J24	hsa-miR-190a-5p	MIMAT0000458	YP00204763
C01	hsa-miR-532-5p	MIMAT0002888	YP00204221	K01	hsa-miR-302a-3p	MIMAT0000684	YP00206059

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
C02	hsa-miR-99a-5p	MIMAT0000097	YP00204521	K02	hsa-miR-595	MIMAT0003263	YP00204070
C03	hsa-miR-16-5p	MIMAT0000069	YP00205702	K03	hsa-miR-602	MIMAT0003270	YP00204316
C04	hsa-miR-98-5p	MIMAT0000096	YP00204640	K04	hsa-miR-223-3p	MIMAT0000280	YP00205986
C05	hsa-miR-185-5p	MIMAT0000455	YP00206037	K05	hsa-miR-627-5p	MIMAT0003296	YP00205979
C06	hsa-miR-25-3p	MIMAT0000081	YP00204361	K06	hsa-miR-34b-3p	MIMAT0004676	YP00204005
C07	UniSP3	N/A	YP02119288	K07	hsa-miR-410-3p	MIMAT0002171	YP00204042
C08	UniSP3	N/A	YP02119288	K08	hsa-miR-17-5p	MIMAT0000070	YP02119304
C09	hsa-miR-765	MIMAT0003945	YP00204071	K09	hsa-miR-376a-3p	MIMAT0000729	YP00204508
C10	hsa-miR-24-3p	MIMAT0000080	YP00204260	K10	hsa-miR-514a-3p	MIMAT0002883	YP00205931
C11	hsa-miR-369-5p	MIMAT0001621	YP00206014	K11	hsa-miR-512-5p	MIMAT0002822	YP00204134
C12	hsa-miR-425-5p	MIMAT0003393	YP00204337	K12	hsa-miR-449a	MIMAT0001541	YP00204481
C13	hsa-miR-590-5p	MIMAT0003258	YP00204222	K13	hsa-miR-498-5p	MIMAT0002824	YP00204150
C14	hsa-miR-760	MIMAT0004957	YP00204549	K14	hsa-miR-148b-3p	MIMAT0000759	YP00204047
C15	hsa-miR-574-3p	MIMAT0003239	YP00206011	K15	hsa-miR-127-3p	MIMAT0000446	YP00204048
C16	hsa-miR-130b-3p	MIMAT0000691	YP00204317	K16	hsa-miR-598-3p	MIMAT0003266	YP00204320
C17	hsa-miR-30c-5p	MIMAT0000244	YP00204783	K17	hsa-miR-96-5p	MIMAT0000095	YP00204417
C18	hsa-miR-133b	MIMAT0000770	YP00206058	K18	hsa-let-7d-5p	MIMAT0000065	YP00204124
C19	hsa-miR-524-5p	MIMAT0002849	YP00204135	K19	hsa-miR-135b-5p	MIMAT0000758	YP00204130
C20	hsa-miR-23a-3p	MIMAT0000078	YP00204772	K20	hsa-miR-495-3p	MIMAT0002817	YP00206015
C21	hsa-miR-193b-3p	MIMAT0002819	YP00204226	K21	hsa-miR-299-5p	MIMAT0002890	YP00204544
C22	hsa-miR-501-5p	MIMAT0002872	YP00204648	K22	hsa-miR-34c-3p	MIMAT0004677	YP00204373
C23	hsa-miR-518c-5p	MIMAT0002847	YP00204280	K23	hsa-miR-596	MIMAT0003264	YP00204002
C24	hsa-miR-130a-3p	MIMAT0000425	YP00204658	K24	hsa-miR-744-5p	MIMAT0004945	YP00204663
D01	hsa-miR-933	MIMAT0004976	YP00205996	L01	hsa-miR-145-5p	MIMAT0000437	YP00204483
D02	hsa-miR-379-5p	MIMAT0000733	YP00205658	L02	hsa-miR-622	MIMAT0003291	YP00205615
D03	hsa-miR-452-5p	MIMAT0001635	YP00204301	L03	hsa-miR-516a-5p	MIMAT0004770	YP00204471
D04	hsa-miR-589-5p	MIMAT0004799	YP00205675	L04	hsa-let-7a-5p	MIMAT0000062	YP00205727
D05	hsa-miR-141-3p	MIMAT0000432	YP00204504	L05	hsa-miR-96-3p	MIMAT0004510	YP00205899
D06	hsa-miR-342-3p	MIMAT0000753	YP00205625	L06	hsa-miR-185-3p	MIMAT0004611	YP00205710
D07	hsa-miR-668-3p	MIMAT0003881	YP00204206	L07	hsa-miR-615-3p	MIMAT0003283	YP00204453
D08	hsa-miR-934	MIMAT0004977	YP02119292	L08	hsa-miR-128-3p	MIMAT0000424	YP00205995
D09	hsa-miR-101-3p	MIMAT0000099	YP00204786	L09	hsa-miR-766-3p	MIMAT0003888	YP00204499
D10	hsa-miR-539-5p	MIMAT0003163	YP00205656	L10	hsa-miR-206	MIMAT0000462	YP00206073
D11	hsa-miR-331-3p	MIMAT0000760	YP00206046	L11	hsa-miR-298	MIMAT0004901	YP00204115
D12	hsa-miR-499a-5p	MIMAT0002870	YP00205935	L12	hsa-miR-193a-5p	MIMAT0004614	YP00204665
D13	hsa-miR-196a-5p	MIMAT0000226	YP00204386	L13	hsa-miR-449b-5p	MIMAT0003327	YP00204751
D14	hsa-miR-888-5p	MIMAT0004916	YP00206024	L14	hsa-miR-520d-5p	MIMAT0002855	YP00204684
D15	hsa-miR-330-3p	MIMAT0000751	YP00204017	L15	hsa-miR-192-5p	MIMAT0000222	YP00204099
D16	hsa-miR-570-3p	MIMAT0003235	YP00204694	L16	hsa-miR-29a-3p	MIMAT0000086	YP00204698
D17	hsa-miR-518c-3p	MIMAT0002848	YP00204554	L17	hsa-miR-18a-3p	MIMAT0002891	YP00204523
D18	hsa-miR-200a-3p	MIMAT0000682	YP00204707	L18	hsa-miR-383-5p	MIMAT0000738	YP00205904
D19	hsa-miR-188-5p	MIMAT0000457	YP00204239	L19	hsa-miR-9-5p	MIMAT0000441	YP00204513
D20	hsa-miR-26a-5p	MIMAT0000082	YP00206023	L20	hsa-miR-142-5p	MIMAT0000433	YP00204722
D21	hsa-miR-99b-5p	MIMAT0000689	YP00205983	L21	hsa-miR-363-5p	MIMAT0003385	YP00204173
D22	hsa-miR-431-5p	MIMAT0001625	YP00204737	L22	hsa-miR-147b-3p	MIMAT0004928	YP00204368
D23	hsa-miR-23b-3p	MIMAT0000418	YP02119314	L23	hsa-miR-197-3p	MIMAT0000227	YP00204380
D24	hsa-miR-367-3p	MIMAT0000719	YP00204784	L24	hsa-miR-597-5p	MIMAT0003265	YP00204769
E01	hsa-miR-505-3p	MIMAT0002876	YP00204214	M01	hsa-miR-326	MIMAT0000756	YP00204512
E02	hsa-miR-18a-5p	MIMAT0000072	YP00204207	M02	hsa-miR-15b-5p	MIMAT0000417	YP00204243
E03	hsa-miR-92a-3p	MIMAT0000092	YP00204258	M03	hsa-miR-105-5p	MIMAT0000102	YP00204389

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
E04	hsa-miR-500a-5p	MIMAT0004773	YP00204794	M04	hsa-miR-196b-5p	MIMAT0001080	YP00204555
E05	hsa-miR-887-3p	MIMAT0004951	YP00204472	M05	hsa-miR-296-5p	MIMAT0000690	YP00204436
E06	hsa-miR-491-3p	MIMAT0004765	YP00205977	M06	hsa-miR-20b-5p	MIMAT0001413	YP00204755
E07	hsa-miR-423-3p	MIMAT0001340	YP00204488	M07	hsa-miR-147a	MIMAT0000251	YP00204398
E08	hsa-miR-126-3p	MIMAT0000445	YP00204227	M08	hsa-miR-198	MIMAT0000228	YP00204341
E09	hsa-miR-421	MIMAT0003339	YP00204603	M09	hsa-miR-375-3p	MIMAT0000728	YP00204362
E10	hsa-miR-376b-3p	MIMAT0002172	YP00204218	M10	hsa-miR-517a-3p	MIMAT0002852	YP00206019
E11	hsa-miR-302c-3p	MIMAT0000717	YP00204403	M11	hsa-miR-361-3p	MIMAT0004682	YP00204008
E12	hsa-miR-625-3p	MIMAT0004808	YP00204647	M12	hsa-miR-21-3p	MIMAT0004494	YP00204302
E13	hsa-miR-339-5p	MIMAT0000764	YP00206007	M13	hsa-miR-373-3p	MIMAT0000726	YP00204604
E14	hsa-miR-873-5p	MIMAT0004953	YP00204175	M14	hsa-miR-518f-3p	MIMAT0002842	YP00204140
E15	hsa-miR-323a-3p	MIMAT0000755	YP00204278	M15	hsa-miR-222-3p	MIMAT0000279	YP00204551
E16	hsa-miR-181d-5p	MIMAT0002821	YP00204789	M16	hsa-miR-617	MIMAT0003286	YP00204564
E17	hsa-miR-125a-5p	MIMAT0000443	YP00204339	M17	hsa-miR-154-5p	MIMAT0000452	YP00204518
E18	hsa-miR-129-5p	MIMAT0000242	YP00204534	M18	hsa-miR-708-5p	MIMAT0004926	YP00204490
E19	hsa-miR-492	MIMAT0002812	YP00204053	M19	hsa-let-7b-5p	MIMAT0000063	YP00204750
E20	hsa-miR-20a-5p	MIMAT0000075	YP00204292	M20	hsa-miR-95-3p	MIMAT0000094	YP00204288
E21	hsa-miR-374b-5p	MIMAT0004955	YP00204608	M21	hsa-miR-517c-3p	MIMAT0002866	YP00204177
E22	hsa-miR-302d-3p	MIMAT0000718	YP00204311	M22	hsa-miR-151a-5p	MIMAT0004697	YP00204007
E23	hsa-miR-346	MIMAT0000773	YP00206009	M23	hsa-miR-502-5p	MIMAT0002873	YP00204449
E24	hsa-miR-151a-3p	MIMAT0000757	YP00204576	M24	hsa-miR-345-5p	MIMAT0000772	YP00206006
F01	hsa-miR-493-3p	MIMAT0003161	YP00204557	N01	hsa-miR-509-3p	MIMAT0002881	YP00204458
F02	hsa-miR-122-5p	MIMAT0000421	YP00205664	N02	hsa-miR-134-5p	MIMAT0000447	YP00205989
F03	hsa-miR-99a-3p	MIMAT0004511	YP00204520	N03	hsa-miR-382-5p	MIMAT0000737	YP00204169
F04	hsa-miR-361-5p	MIMAT0000703	YP00206054	N04	hsa-miR-490-3p	MIMAT0002806	YP00205999
F05	hsa-miR-202-3p	MIMAT0002811	YP00205990	N05	hsa-miR-200c-3p	MIMAT0000617	YP00204482
F06	hsa-miR-125b-5p	MIMAT0000423	YP00205713	N06	hsa-miR-30a-5p	MIMAT0000087	YP00205695
F07	hsa-miR-503-5p	MIMAT0002874	YP00204334	N07	hsa-miR-181b-5p	MIMAT0000257	YP00204530
F08	hsa-miR-204-5p	MIMAT0000265	YP00206072	N08	hsa-miR-33a-5p	MIMAT0000091	YP00205690
F09	hsa-miR-30d-5p	MIMAT0000245	YP00206047	N09	hsa-miR-195-5p	MIMAT0000461	YP00205869
F10	hsa-miR-301a-3p	MIMAT0000688	YP00205601	N10	hsa-miR-874-3p	MIMAT0004911	YP00204761
F11	hsa-miR-362-5p	MIMAT0000705	YP00204618	N11	hsa-miR-135a-5p	MIMAT0000428	YP00204762
F12	hsa-miR-30b-3p	MIMAT0004589	YP00206003	N12	hsa-miR-26a-2-3p	MIMAT0004681	YP00204676
F13	hsa-miR-654-5p	MIMAT0003330	YP00204439	N13	hsa-miR-146b-5p	MIMAT0002809	YP02119310
F14	hsa-miR-545-3p	MIMAT0003165	YP00206087	N14	hsa-miR-412-3p	MIMAT0002170	YP00204460
F15	hsa-miR-29b-2-5p	MIMAT0004515	YP00204208	N15	hsa-miR-1-3p	MIMAT0000416	YP00204344
F16	hsa-miR-491-5p	MIMAT0002807	YP00204695	N16	hsa-miR-299-3p	MIMAT0000687	YP00204702
F17	hsa-miR-92b-3p	MIMAT0003218	YP00204384	N17	hsa-miR-142-3p	MIMAT0000434	YP00204291
F18	hsa-miR-665	MIMAT0004952	YP00204710	N18	hsa-miR-338-3p	MIMAT0000763	YP00204719
F19	hsa-miR-506-3p	MIMAT0002878	YP00204539	N19	hsa-miR-584-5p	MIMAT0003249	YP00204568
F20	hsa-miR-363-3p	MIMAT0000707	YP00204726	N20	hsa-miR-377-3p	MIMAT0000730	YP00204733
F21	hsa-miR-132-3p	MIMAT0000426	YP00206035	N21	hsa-miR-216a-5p	MIMAT0000273	YP00204167
F22	hsa-miR-651-5p	MIMAT0003321	YP00205864	N22	hsa-miR-424-5p	MIMAT0001341	YP00204736
F23	hsa-miR-628-3p	MIMAT0003297	YP00206057	N23	hsa-miR-921	MIMAT0004971	YP00204281
F24	hsa-miR-432-5p	MIMAT0002814	YP00204776	N24	hsa-miR-513a-5p	MIMAT0002877	YP00205900
G01	hsa-miR-154-3p	MIMAT0000453	YP00204096	O01	hsa-miR-140-3p	MIMAT0004597	YP00204304
G02	hsa-miR-27a-3p	MIMAT0000084	YP00206038	O02	hsa-miR-181a-5p	MIMAT0000256	YP00206081
G03	hsa-miR-376c-3p	MIMAT0000720	YP00204442	O03	hsa-miR-10a-5p	MIMAT0000253	YP00204778
G04	hsa-miR-940	MIMAT0004983	YP00204094	O04	hsa-miR-106a-5p	MIMAT0000103	YP00204563
G05	hsa-miR-22-5p	MIMAT0004495	YP00204255	O05	hsa-miR-182-5p	MIMAT0000259	YP00206070

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
G06	hsa-miR-224-5p	MIMAT0000281	YP02119313	O06	hsa-miR-370-3p	MIMAT0000722	YP00204011
G07	hsa-miR-885-5p	MIMAT0004947	YP00204473	O07	hsa-miR-576-5p	MIMAT0003241	YP00206064
G08	hsa-miR-320a-3p	MIMAT0000510	YP00206042	O08	hsa-miR-425-3p	MIMAT0001343	YP00204038
G09	hsa-miR-18b-5p	MIMAT0001412	YP00204084	O09	hsa-miR-450a-5p	MIMAT0001545	YP00206085
G10	hsa-miR-187-3p	MIMAT0000262	YP00204018	O10	hsa-miR-411-5p	MIMAT0003329	YP00204531
G11	hsa-miR-516b-5p	MIMAT0002859	YP00204558	O11	hsa-miR-216b-5p	MIMAT0004959	YP00204289
G12	hsa-miR-302c-5p	MIMAT0000716	YP00205940	O12	hsa-miR-106b-5p	MIMAT0000680	YP00205884
G13	hsa-miR-548b-3p	MIMAT0003254	YP00204485	O13	hsa-miR-22-3p	MIMAT0000077	YP00204606
G14	hsa-miR-186-5p	MIMAT0000456	YP00206053	O14	hsa-miR-510-5p	MIMAT0002882	YP00204349
G15	hsa-miR-199a-5p	MIMAT0000231	YP00204494	O15	hsa-miR-212-3p	MIMAT0000269	YP00204170
G16	hsa-miR-155-5p	MIMAT0000646	YP02119311	O16	hsa-miR-525-5p	MIMAT0002838	YP00204355
G17	hsa-miR-107	MIMAT0000104	YP00204468	O17	hsa-miR-542-5p	MIMAT0003340	YP00204198
G18	hsa-miR-302b-3p	MIMAT0000715	YP00204773	O18	hsa-miR-576-3p	MIMAT0004796	YP00204195
G19	hsa-miR-662	MIMAT0003325	YP00204180	O19	hsa-miR-583	MIMAT0003248	YP00204021
G20	hsa-miR-519d-3p	MIMAT0002853	YP00204062	O20	hsa-miR-483-3p	MIMAT0002173	YP00204012
G21	hsa-miR-485-3p	MIMAT0002176	YP00206055	O21	hsa-miR-582-5p	MIMAT0003247	YP00204254
G22	hsa-miR-200b-3p	MIMAT0000318	YP00206071	O22	hsa-miR-183-5p	MIMAT0000261	YP00206030
G23	hsa-miR-337-3p	MIMAT0000754	YP00205938	O23	hsa-miR-33b-5p	MIMAT0003301	YP00205860
G24	hsa-miR-494-3p	MIMAT0002816	YP00204579	O24	hsa-miR-193a-3p	MIMAT0000459	YP00204591
H01	hsa-miR-371a-3p	MIMAT0000723	YP00204299	P01	hsa-miR-153-3p	MIMAT0000439	YP00204338
H02	hsa-miR-637	MIMAT0003307	YP00205679	P02	hsa-let-7e-5p	MIMAT0000066	YP00205711
H03	hsa-miR-144-3p	MIMAT0000436	YP00204754	P03	hsa-miR-409-3p	MIMAT0001639	YP00204358
H04	hsa-miR-161-3p	MIMAT0004489	YP00206012	P04	hsa-miR-100-5p	MIMAT0000098	YP00205689
H05	hsa-miR-631	MIMAT0003300	YP00204312	P05	hsa-miR-629-5p	MIMAT0004810	YP00204370
H06	hsa-miR-34c-5p	MIMAT0000686	YP00205659	P06	hsa-miR-484	MIMAT0002174	YP00205636
H07	hsa-miR-211-5p	MIMAT0000268	YP00204009	P07	hsa-miR-429	MIMAT0001536	YP00205901
H08	hsa-miR-454-3p	MIMAT0003885	YP00205663	P08	hsa-miR-30c-2-3p	MIMAT0004550	YP00205632
H09	hsa-let-7f-5p	MIMAT0000067	YP00204359	P09	hsa-miR-518a-3p	MIMAT0002863	YP00206069
H10	hsa-miR-30e-5p	MIMAT0000692	YP00204714	P10	hsa-miR-340-5p	MIMAT0004692	YP00206068
H11	hsa-miR-34a-5p	MIMAT0000255	YP00204486	P11	hsa-miR-508-3p	MIMAT0002880	YP00204480
H12	hsa-miR-663a	MIMAT0003326	YP00204284	P12	hsa-miR-381-3p	MIMAT0000736	YP00205887
H13	hsa-miR-518e-3p	MIMAT0002861	YP00204051	P13	hsa-miR-148a-3p	MIMAT0000243	YP00205867
H14	hsa-miR-29b-3p	MIMAT0000100	YP00204679	P14	hsa-miR-146a-5p	MIMAT0000449	YP00204688
H15	hsa-miR-658	MIMAT0003336	YP00204522	P15	hsa-miR-139-5p	MIMAT0000250	YP00205874
H16	hsa-miR-572	MIMAT0003237	YP00204696	P16	hsa-miR-373-5p	MIMAT0000725	YP00204703
H17	hsa-miR-802	MIMAT0004185	YP00205980	P17	hsa-miR-149-5p	MIMAT0000450	YP00204321
H18	hsa-miR-521	MIMAT0002854	YP00205890	P18	hsa-miR-642a-5p	MIMAT0003312	YP00204793
H19	hsa-miR-433-3p	MIMAT0001627	YP00204036	P19	hsa-miR-31-5p	MIMAT0000089	YP00204236
H20	hsa-miR-660-5p	MIMAT0003338	YP00205911	P20	hsa-miR-451a	MIMAT0001631	YP02119305
H21	hsa-let-7c-5p	MIMAT0000064	YP00204767	P21	hsa-miR-620	MIMAT0003289	YP00204637
H22	hsa-miR-28-5p	MIMAT0000085	YP00204322	P22	hsa-miR-27b-3p	MIMAT0000419	YP00205915
H23	hsa-miR-324-5p	MIMAT0000761	YP02119473	P23	hsa-miR-523-3p	MIMAT0002840	YP00204452
H24	hsa-miR-219a-5p	MIMAT0000276	YP00204780	P24	hsa-miR-374a-5p	MIMAT0000727	YP00204758

miRNA Table: Plate 2

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
A01	hsa-miR-92a-1-5p	MIMAT0004507	YP00204560	I01	hsa-miR-519c-3p	MIMAT0002832	YP00204016
A02	hsa-miR-219a-1-	MIMAT0004567	YP00205880	I02	hsa-miR-554	MIMAT0003217	YP00204556
A03	hsa-miR-1913	MIMAT0007888	YP00204673	I03	hsa-miR-938	MIMAT0004981	YP00205912
A04	hsa-miR-1245a	MIMAT0005897	YP00204229	I04	hsa-miR-1243	MIMAT0005894	YP00204211
A05	hsa-miR-522-3p	MIMAT0002868	YP00205913	I05	hsa-miR-708-3p	MIMAT0004927	YP00204681
A06	hsa-miR-571	MIMAT0003236	YP00206065	I06	hsa-miR-1185-5p	MIMAT0005798	YP00205926
A07	hsa-miR-323a-5p	MIMAT0004696	YP00206080	I07	hsa-miR-512-3p	MIMAT0002823	YP00204382
A08	hsa-miR-592	MIMAT0003260	YP00204375	I08	hsa-miR-587	MIMAT0003253	YP00205876
A09	hsa-miR-487a-3p	MIMAT0002178	YP00204381	I09	hsa-miR-603	MIMAT0003271	YP00204112
A10	hsa-miR-1249-3p	MIMAT0005901	YP00204122	I10	hsa-miR-1184	MIMAT0005829	YP00204330
A11	hsa-miR-25-5p	MIMAT0004498	YP00204031	I11	hsa-miR-20a-3p	MIMAT0004493	YP00204052
A12	hsa-miR-922	MIMAT0004972	YP00204478	I12	hsa-miR-588	MIMAT0003255	YP00205885
A13	hsa-miR-124-5p	MIMAT0004591	YP00204266	I13	hsa-miR-455-3p	MIMAT0004784	YP00204035
A14	hsa-miR-1264	MIMAT0005791	YP00204056	I14	hsa-miR-582-3p	MIMAT0004797	YP00204072
A15	hsa-miR-504-5p	MIMAT0002875	YP00204396	I15	hsa-miR-409-5p	MIMAT0001638	YP00204014
A16	hsa-miR-138-1-3p	MIMAT0004607	YP00205881	I16	hsa-miR-452-3p	MIMAT0001636	YP00204201
A17	hsa-miR-502-3p	MIMAT0004775	YP00204043	I17	hsa-miR-19b-1-5p	MIMAT0004491	YP00204310
A18	hsa-miR-490-5p	MIMAT0004764	YP00206077	I18	hsa-miR-610	MIMAT0003278	YP00204492
A19	hsa-miR-567	MIMAT0003231	YP00204165	I19	hsa-miR-511-5p	MIMAT0002808	YP00204046
A20	hsa-miR-18b-3p	MIMAT0004751	YP00204089	I20	hsa-miR-200c-5p	MIMAT0004657	YP02119294
A21	hsa-miR-125a-3p	MIMAT0004602	YP00204446	I21	hsa-let-7a-3p	MIMAT0004481	YP00206084
A22	hsa-miR-653-5p	MIMAT0003328	YP00204669	I22	hsa-miR-135a-3p	MIMAT0004595	YP00204022
A23	hsa-miR-891b	MIMAT0004913	YP00204109	I23	hsa-miR-520a-5p	MIMAT0002833	YP00204340
A24	hsa-miR-144-5p	MIMAT0004600	YP00204670	I24	hsa-miR-1468-5p	MIMAT0006789	YP00204577
B01	hsa-miR-1538	MIMAT0007400	YP00204615	J01	hsa-miR-628-5p	MIMAT0004809	YP00205893
B02	hsa-miR-384	MIMAT0001075	YP00205879	J02	hsa-miR-552-3p	MIMAT0003215	YP00206032
B03	hsa-miR-196b-3p	MIMAT0009201	YP00206018	J03	hsa-miR-145-3p	MIMAT0004601	YP00204192
B04	hsa-miR-649	MIMAT0003319	YP00204146	J04	BlankSpot	N/A	YP02119291
B05	hsa-miR-143-5p	MIMAT0004599	YP00204570	J05	hsa-miR-378a-5p	MIMAT0000731	YP00204347
B06	hsa-miR-1207-5p	MIMAT0005871	YP00204693	J06	hsa-miR-7-1-3p	MIMAT0004553	YP00205888
B07	hsa-miR-943	MIMAT0004986	YP00204545	J07	hsa-miR-181c-3p	MIMAT0004559	YP00204213
B08	hsa-miR-675-3p	MIMAT0006790	YP00204732	J08	hsa-miR-195-3p	MIMAT0004615	YP00204128
B09	hsa-miR-200b-5p	MIMAT0004571	YP00204144	J09	hsa-miR-578	MIMAT0003243	YP00204269
B10	hsa-miR-519e-5p	MIMAT0002828	YP00204029	J10	hsa-miR-505-5p	MIMAT0004776	YP00205657
B11	hsa-miR-942-5p	MIMAT0004985	YP00204440	J11	hsa-miR-875-3p	MIMAT0004923	YP00204231
B12	hsa-miR-450b-3p	MIMAT0004910	YP00205683	J12	hsa-miR-450b-5p	MIMAT0004909	YP00205607
B13	hsa-miR-553	MIMAT0003216	YP00204612	J13	hsa-miR-876-5p	MIMAT0004924	YP00204527
B14	hsa-miR-605-5p	MIMAT0003273	YP00205594	J14	hsa-miR-362-3p	MIMAT0004683	YP00205612
B15	hsa-miR-24-2-5p	MIMAT0004497	YP00204187	J15	hsa-miR-624-5p	MIMAT0003293	YP00204625
B16	hsa-miR-23a-5p	MIMAT0004496	YP00205631	J16	hsa-miR-27a-5p	MIMAT0004501	YP00206021
B17	hsa-miR-27b-5p	MIMAT0004588	YP00204219	J17	hsa-miR-744-3p	MIMAT0004946	YP00206082
B18	hsa-miR-759	MIMAT0010497	YP00206000	J18	hsa-miR-139-3p	MIMAT0004552	YP00205661
B19	hsa-miR-770-5p	MIMAT0003948	YP00204197	J19	hsa-miR-138-2-3p	MIMAT0004596	YP00204477
B20	hsa-miR-585-3p	MIMAT0003250	YP00205701	J20	hsa-miR-655-3p	MIMAT0003331	YP00205645
B21	hsa-miR-376a-5p	MIMAT0003386	YP00205988	J21	hsa-miR-99b-3p	MIMAT0004678	YP00204064
B22	hsa-miR-507	MIMAT0002879	YP00205662	J22	hsa-miR-581	MIMAT0003246	YP00205696
B23	hsa-miR-520b-3p	MIMAT0002843	YP00205939	J23	hsa-miR-191-3p	MIMAT0001618	YP00204196
B24	hsa-miR-302f	MIMAT0005932	YP00205634	J24	hsa-miR-32-3p	MIMAT0004505	YP00205648
C01	hsa-miR-28-3p	MIMAT0004502	YP00204119	K01	hsa-miR-1204	MIMAT0005868	YP00204123

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
C02	hsa-miR-875-5p	MIMAT0004922	YP00205978	K02	hsa-miR-548j-5p	MIMAT0005875	YP00205861
C03	hsa-miR-219a-2-	MIMAT0004675	YP00204674	K03	hsa-miR-555	MIMAT0003219	YP00204377
C04	hsa-miR-1183	MIMAT0005828	YP00204176	K04	hsa-miR-1224-3p	MIMAT0005459	YP00204045
C05	hsa-miR-758-3p	MIMAT0003879	YP00204353	K05	hsa-miR-1539	MIMAT0007401	YP00204682
C06	hsa-miR-1244	MIMAT0005896	YP00204245	K06	hsa-miR-663b	MIMAT0005867	YP00205985
C07	hsa-miR-4315	MIMAT0016866	YP02101265	K07	UniSP3	N/A	YP02119288
C08	hsa-miR-1256	MIMAT0005907	YP00204059	K08	hsa-miR-1248	MIMAT0005900	YP00204253
C09	hsa-miR-516a-3p	MIMAT0006778	YP00204305	K09	hsa-miR-889-3p	MIMAT0004921	YP00204113
C10	hsa-miR-548c-5p	MIMAT0004806	YP00205882	K10	hsa-miR-1227-3p	MIMAT0005580	YP00204200
C11	hsa-miR-496	MIMAT0002818	YP00204441	K11	hsa-miR-548h-5p	MIMAT0005928	YP00204777
C12	hsa-miR-876-3p	MIMAT0004925	YP00205927	K12	hsa-miR-1255b-5p	MIMAT0005945	YP00205878
C13	hsa-miR-532-3p	MIMAT0004780	YP00204003	K13	hsa-miR-330-5p	MIMAT0004693	YP00204372
C14	hsa-miR-654-3p	MIMAT0004814	YP00204054	K14	hsa-miR-1238-3p	MIMAT0005593	YP00204004
C15	hsa-miR-659-3p	MIMAT0003337	YP00204199	K15	hsa-miR-188-3p	MIMAT0004613	YP02119296
C16	hsa-miR-135b-3p	MIMAT0004698	YP00204427	K16	hsa-miR-589-3p	MIMAT0003256	YP00206048
C17	hsa-miR-641	MIMAT0003311	YP00204411	K17	hsa-miR-125b-2-3p	MIMAT0004603	YP00204295
C18	hsa-miR-2113	MIMAT0009206	YP00204686	K18	hsa-miR-16-2-3p	MIMAT0004518	YP00204309
C19	hsa-miR-644a	MIMAT0003314	YP02118540	K19	hsa-miR-515-5p	MIMAT0002826	YP00204431
C20	hsa-miR-661	MIMAT0003324	YP00204657	K20	hsa-miR-340-3p	MIMAT0000750	YP00204250
C21	hsa-miR-892a	MIMAT0004907	YP00204092	K21	hsa-miR-513a-3p	MIMAT0004777	YP00204597
C22	hsa-miR-10b-3p	MIMAT0004556	YP00205909	K22	hsa-miR-34a-3p	MIMAT0004557	YP00206061
C23	hsa-miR-122-3p	MIMAT0004590	YP00204437	K23	hsa-miR-342-5p	MIMAT0004694	YP00204516
C24	hsa-miR-100-3p	MIMAT0004512	YP00204535	K24	hsa-miR-639	MIMAT0003309	YP02119298
D01	hsa-miR-769-3p	MIMAT0003887	YP00204104	L01	hsa-let-7i-3p	MIMAT0004585	YP00204247
D02	hsa-miR-300	MIMAT0004903	YP00204756	L02	hsa-miR-543	MIMAT0004954	YP00204447
D03	hsa-miR-518e-5p	MIMAT0005450	YP00206040	L03	hsa-miR-645	MIMAT0003315	YP00206076
D04	hsa-miR-489-3p	MIMAT0002805	YP00204395	L04	hsa-miR-548d-5p	MIMAT0004812	YP00205907
D05	hsa-miR-937-3p	MIMAT0004980	YP00204060	L05	hsa-miR-33a-3p	MIMAT0004506	YP00204086
D06	hsa-miR-381-5p	MIMAT0022862	YP00206075	L06	hsa-miR-664a-3p	MIMAT0005949	YP00204725
D07	hsa-miR-640	MIMAT0003310	YP00204533	L07	UniSP3	N/A	YP02119288
D08	hsa-miR-148b-5p	MIMAT0004699	YP00204735	L08	UniSP3	N/A	YP02119288
D09	hsa-miR-29c-5p	MIMAT0004673	YP00204132	L09	hsa-miR-379-3p	MIMAT0004690	YP00204345
D10	hsa-miR-499a-3p	MIMAT0004772	YP00204438	L10	hsa-miR-556-3p	MIMAT0004793	YP00205717
D11	hsa-let-7f-1-3p	MIMAT0004486	YP00204323	L11	hsa-miR-614	MIMAT0003282	YP00204273
D12	hsa-miR-382-3p	MIMAT0022697	YP00206066	L12	hsa-miR-616-5p	MIMAT0003284	YP00205611
D13	hsa-miR-609	MIMAT0003277	YP00204613	L13	hsa-miR-93-3p	MIMAT0004509	YP00204470
D14	hsa-miR-10a-3p	MIMAT0004555	YP00205688	L14	hsa-miR-1972	MIMAT0009447	YP00205605
D15	hsa-miR-106a-3p	MIMAT0004517	YP00204443	L15	hsa-miR-616-3p	MIMAT0004805	YP00204116
D16	hsa-let-7e-3p	MIMAT0004485	YP00205301	L16	hsa-miR-369-3p	MIMAT0000721	YP00206028
D17	hsa-miR-580-3p	MIMAT0003245	YP00204505	L17	hsa-miR-2110	MIMAT0010133	YP00204328
D18	hsa-miR-761	MIMAT0010364	YP00205475	L18	hsa-miR-548a-3p	MIMAT0003251	YP00205650
D19	hsa-miR-643	MIMAT0003313	YP00204275	L19	hsa-miR-634	MIMAT0003304	YP00204228
D20	hsa-miR-618	MIMAT0003287	YP00205674	L20	hsa-miR-320c	MIMAT0005793	YP00205706
D21	hsa-miR-221-5p	MIMAT0004568	YP00204032	L21	hsa-miR-636	MIMAT0003306	YP00204298
D22	hsa-miR-513b-5p	MIMAT0005788	YP00205620	L22	hsa-miR-606	MIMAT0003274	YP00205629
D23	hsa-miR-411-3p	MIMAT0004813	YP00204434	L23	hsa-miR-208b-3p	MIMAT0004960	YP00204636
D24	hsa-miR-19a-5p	MIMAT0004490	YP00205604	L24	hsa-miR-367-5p	MIMAT0004686	YP00205652
E01	hsa-miR-338-5p	MIMAT0004701	YP00204114	M01	hsa-miR-520d-3p	MIMAT0002856	YP00204040
E02	hsa-miR-1914-3p	MIMAT0007890	YP00204644	M02	hsa-miR-1265	MIMAT0005918	YP00204082
E03	hsa-miR-323b-5p	MIMAT0001630	YP00204796	M03	hsa-miR-1203	MIMAT0005866	YP00204414

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
E04	hsa-miR-548i	MIMAT0005935	YP00204103	M04	hsa-miR-548k	MIMAT0005882	YP00205906
E05	hsa-miR-541-3p	MIMAT0004920	YP00204276	M05	hsa-miR-548a-5p	MIMAT0004803	YP00204416
E06	hsa-miR-1272	MIMAT0005925	YP00205863	M06	hsa-miR-1253	MIMAT0005904	YP00204525
E07	hsa-miR-1205	MIMAT0005869	YP00204143	M07	hsa-miR-615-5p	MIMAT0004804	YP00204445
E08	hsa-miR-544a	MIMAT0003164	YP00204646	M08	hsa-miR-607	MIMAT0003275	YP00204649
E09	hsa-miR-431-3p	MIMAT0004757	YP00204567	M09	hsa-miR-1208	MIMAT0005873	YP00205937
E10	hsa-miR-621	MIMAT0003290	YP00204307	M10	hsa-miR-302e	MIMAT0005931	YP00205858
E11	hsa-miR-556-5p	MIMAT0003220	YP00204240	M11	hsa-miR-1206	MIMAT0005870	YP00204147
E12	hsa-miR-1267	MIMAT0005921	YP00204543	M12	hsa-miR-1270	MIMAT0005924	YP00204025
E13	hsa-miR-141-5p	MIMAT0004598	YP00206088	M13	hsa-miR-525-3p	MIMAT0002839	YP00204111
E14	hsa-miR-1269a	MIMAT0005923	YP00204455	M14	hsa-miR-1200	MIMAT0005863	YP00204044
E15	hsa-miR-501-3p	MIMAT0004774	YP00204178	M15	hsa-miR-1911-5p	MIMAT0007885	YP00204595
E16	hsa-miR-15b-3p	MIMAT0004586	YP00205898	M16	hsa-miR-33b-3p	MIMAT0004811	YP00204462
E17	hsa-miR-146b-3p	MIMAT0004766	YP02119309	M17	hsa-miR-223-5p	MIMAT0004570	YP00204529
E18	hsa-miR-222-5p	MIMAT0004569	YP00204314	M18	hsa-miR-34b-5p	MIMAT0000685	YP00204424
E19	hsa-miR-601	MIMAT0003269	YP00204426	M19	hsa-miR-888-3p	MIMAT0004917	YP00205934
E20	hsa-miR-924	MIMAT0004974	YP02119297	M20	hsa-miR-424-3p	MIMAT0004749	YP00205918
E21	hsa-miR-29a-5p	MIMAT0004503	YP00204430	M21	hsa-miR-339-3p	MIMAT0004702	YP00204160
E22	hsa-let-7a-2-3p	MIMAT0010195	YP00206045	M22	hsa-miR-380-5p	MIMAT0000734	YP00205903
E23	hsa-miR-520f-3p	MIMAT0002830	YP00204283	M23	hsa-miR-647	MIMAT0003317	YP00204027
E24	hsa-miR-101-5p	MIMAT0004513	YP00204379	M24	hsa-miR-518f-5p	MIMAT0002841	YP00204671
F01	hsa-miR-520a-3p	MIMAT0002834	YP00204074	N01	hsa-miR-92b-5p	MIMAT0004792	YP00204415
F02	hsa-miR-548m	MIMAT0005917	YP00205929	N02	hsa-miR-551a	MIMAT0003214	YP00204331
F03	hsa-miR-517-5p	MIMAT0002851	YP00204336	N03	hsa-miR-146a-3p	MIMAT0004608	YP00204401
F04	hsa-miR-448	MIMAT0001532	YP00204643	N04	hsa-miR-218-1-3p	MIMAT0004565	YP00204691
F05	hsa-miR-1296-5p	MIMAT0005794	YP00206020	N05	hsa-miR-593-5p	MIMAT0003261	YP00204550
F06	hsa-miR-1537-3p	MIMAT0007399	YP00205994	N06	hsa-miR-561-3p	MIMAT0003225	YP00205895
F07	hsa-miR-920	MIMAT0004970	YP00204248	N07	hsa-miR-767-3p	MIMAT0003883	YP02119295
F08	hsa-miR-1247-5p	MIMAT0005899	YP00204770	N08	hsa-miR-526b-3p	MIMAT0002836	YP00205984
F09	hsa-miR-19b-2-5p	MIMAT0004492	YP00205897	N09	hsa-miR-24-1-5p	MIMAT0000079	YP00204357
F10	hsa-miR-558	MIMAT0003222	YP00205670	N10	hsa-let-7b-3p	MIMAT0004482	YP00205653
F11	hsa-miR-106b-3p	MIMAT0004672	YP00204020	N11	hsa-miR-193b-5p	MIMAT0004767	YP00204065
F12	hsa-miR-1258	MIMAT0005909	YP00205709	N12	hsa-miR-335-3p	MIMAT0004703	YP00205613
F13	hsa-miR-619-3p	MIMAT0003288	YP00204120	N13	hsa-miR-541-5p	MIMAT0004919	YP00204118
F14	hsa-miR-208a-3p	MIMAT0000241	YP00205619	N14	hsa-miR-30c-1-3p	MIMAT0004674	YP00205707
F15	hsa-miR-17-3p	MIMAT0000071	YP00206008	N15	hsa-miR-629-3p	MIMAT0003298	YP00204210
F16	hsa-miR-136-3p	MIMAT0004606	YP00205503	N16	hsa-miR-377-5p	MIMAT0004689	YP00205428
F17	hsa-miR-877-3p	MIMAT0004950	YP00204263	N17	hsa-miR-630	MIMAT0003299	YP00204392
F18	hsa-miR-935	MIMAT0004978	YP00205579	N18	hsa-miR-548d-3p	MIMAT0003323	YP00205608
F19	hsa-miR-224-3p	MIMAT0009198	YP00204629	N19	hsa-miR-885-3p	MIMAT0004948	YP00204136
F20	hsa-miR-624-3p	MIMAT0004807	YP00205684	N20	hsa-miR-320d	MIMAT0006764	YP00205667
F21	hsa-miR-767-5p	MIMAT0003882	YP00204238	N21	hsa-miR-2053	MIMAT0009978	YP00204605
F22	hsa-miR-559	MIMAT0003223	YP00205883	N22	hsa-miR-675-5p	MIMAT0004284	YP00205687
F23	hsa-miR-449b-3p	MIMAT0009203	YP00204631	N23	hsa-miR-1252-5p	MIMAT0005944	YP00204701
F24	hsa-miR-205-3p	MIMAT0009197	YP00205602	N24	hsa-miR-548e-3p	MIMAT0005874	YP00206005
G01	hsa-miR-604	MIMAT0003272	YP00204075	O01	hsa-miR-1914-5p	MIMAT0007889	YP00204672
G02	hsa-miR-130b-5p	MIMAT0004680	YP00204456	O02	hsa-miR-513c-5p	MIMAT0005789	YP00205928
G03	hsa-miR-149-3p	MIMAT0004609	YP00204093	O03	hsa-miR-331-5p	MIMAT0004700	YP00204423
G04	hsa-miR-1271-5p	MIMAT0005796	YP00204351	O04	hsa-miR-1182	MIMAT0005827	YP00204466
G05	hsa-miR-520h	MIMAT0002867	YP00206050	O05	hsa-miR-611	MIMAT0003279	YP00204041

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
G06	hsa-miR-769-5p	MIMAT0003886	YP00204270	O06	hsa-miR-1181	MIMAT0005826	YP00204050
G07	hsa-miR-612	MIMAT0003280	YP00204068	O07	hsa-miR-638	MIMAT0003308	YP00204234
G08	hsa-miR-1237-3p	MIMAT0005592	YP00204069	O08	hsa-miR-515-3p	MIMAT0002827	YP00204651
G09	hsa-miR-1908-5p	MIMAT0007881	YP00204587	O09	hsa-miR-650	MIMAT0003320	YP00204233
G10	hsa-miR-1260a	MIMAT0005911	YP00205892	O10	hsa-miR-1178-3p	MIMAT0005823	YP00205870
G11	hsa-miR-182-3p	MIMAT0000260	YP00204098	O11	hsa-miR-600	MIMAT0003268	YP00204076
G12	hsa-miR-365b-5p	MIMAT0022833	YP00204654	O12	hsa-miR-599	MIMAT0003267	YP00205908
G13	hsa-miR-508-5p	MIMAT0004778	YP00204077	O13	hsa-miR-520g-3p	MIMAT0002858	YP00204055
G14	hsa-miR-671-3p	MIMAT0004819	YP00204024	O14	hsa-miR-564	MIMAT0003228	YP00205925
G15	hsa-miR-941	MIMAT0004984	YP00204574	O15	hsa-miR-132-5p	MIMAT0004594	YP00204552
G16	hsa-miR-23b-5p	MIMAT0004587	YP00205873	O16	hsa-miR-577	MIMAT0003242	YP00205998
G17	hsa-miR-591	MIMAT0003259	YP00205872	O17	hsa-miR-1911-3p	MIMAT0007886	YP00204599
G18	hsa-miR-26b-3p	MIMAT0004500	YP02119315	O18	hsa-let-7f-2-3p	MIMAT0004487	YP00204095
G19	hsa-miR-519b-3p	MIMAT0002837	YP00204332	O19	hsa-miR-155-3p	MIMAT0004658	YP00204000
G20	hsa-miR-30d-3p	MIMAT0004551	YP00204023	O20	hsa-miR-105-3p	MIMAT0004516	YP00204193
G21	hsa-miR-518d-5p	MIMAT0005456	YP00204461	O21	hsa-miR-486-3p	MIMAT0004762	YP00204107
G22	hsa-miR-212-5p	MIMAT0022695	YP00205401	O22	hsa-miR-320b	MIMAT0005792	YP02119299
G23	hsa-miR-520e-3p	MIMAT0002825	YP00204013	O23	hsa-miR-296-3p	MIMAT0004679	YP00204393
G24	hsa-miR-646	MIMAT0003316	YP00204546	O24	hsa-miR-7-2-3p	MIMAT0004554	YP00204581
H01	hsa-miR-519e-3p	MIMAT0002829	YP00204491	P01	hsa-miR-550a-3p	MIMAT0003257	YP00204795
H02	hsa-miR-626	MIMAT0003295	YP00204712	P02	hsa-miR-380-3p	MIMAT0000735	YP00204586
H03	hsa-miR-26a-1-3p	MIMAT0004499	YP00204184	P03	hsa-miR-593-3p	MIMAT0004802	YP00204083
H04	hsa-miR-190b-5p	MIMAT0004929	YP02119312	P04	hsa-miR-1912-3p	MIMAT0007887	YP00205866
H05	hsa-miR-1471	MIMAT0007349	YP00204623	P05	hsa-miR-493-5p	MIMAT0002813	YP00204166
H06	hsa-miR-548l	MIMAT0005889	YP00205891	P06	hsa-miR-432-3p	MIMAT0002815	YP00205894
H07	hsa-miR-586	MIMAT0003252	YP00204205	P07	hsa-miR-454-5p	MIMAT0003884	YP00204279
H08	hsa-miR-103b	MIMAT0007402	YP00204740	P08	hsa-miR-936	MIMAT0004979	YP00204185
H09	hsa-miR-488-5p	MIMAT0002804	YP00205981	P09	hsa-miR-30a-3p	MIMAT0000088	YP00204457
H10	hsa-miR-129-1-3p	MIMAT0004548	YP00206074	P10	hsa-let-7g-3p	MIMAT0004584	YP00205865
H11	hsa-miR-192-3p	MIMAT0004543	YP00204272	P11	hsa-miR-214-5p	MIMAT0004564	YP00204575
H12	hsa-miR-632	MIMAT0003302	YP00205600	P12	hsa-miR-183-3p	MIMAT0004560	YP00205673
H13	hsa-miR-181a-2-	MIMAT0004558	YP00204142	P13	hsa-miR-1179	MIMAT0005824	YP00204624
H14	hsa-miR-1909-3p	MIMAT0007883	YP00205697	P14	hsa-miR-562	MIMAT0003226	YP00205941
H15	hsa-miR-573	MIMAT0003238	YP00204216	P15	hsa-miR-579-3p	MIMAT0003244	YP00204627
H16	hsa-miR-302d-5p	MIMAT0004685	YP00204033	P16	hsa-miR-590-3p	MIMAT0004801	YP00205448
H17	hsa-miR-194-3p	MIMAT0004671	YP00204204	P17	hsa-miR-130a-5p	MIMAT0004593	YP02119308
H18	hsa-miR-302b-5p	MIMAT0000714	YP00205676	P18	hsa-miR-563	MIMAT0003227	YP00205618
H19	hsa-miR-551b-5p	MIMAT0004794	YP00204630	P19	hsa-miR-200a-5p	MIMAT0001620	YP00206063
H20	hsa-miR-635	MIMAT0003305	YP00205936	P20	hsa-miR-483-5p	MIMAT0004761	YP00205693
H21	hsa-miR-518d-3p	MIMAT0002864	YP00204628	P21	hsa-miR-15a-3p	MIMAT0004488	YP00204435
H22	hsa-miR-569	MIMAT0003234	YP00205666	P22	hsa-miR-944	MIMAT0004987	YP00205700
H23	hsa-miR-125b-1-	MIMAT0004592	YP00204400	P23	hsa-miR-92a-2-5p	MIMAT0004508	YP00204277
H24	hsa-miR-218-2-3p	MIMAT0004566	YP00205639	P24	hsa-miR-548n	MIMAT0005916	YP00205705

miRNA Table: Plate 3

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
A01	hsa-miR-1305	MIMAT0005893	YP02118778	I01	hsa-miR-1294	MIMAT0005884	YP02109036
A02	hsa-miR-1306-3p	MIMAT0005950	YP02119520	I02	hsa-miR-4289	MIMAT0016920	YP02113941
A03	hsa-miR-4452	MIMAT0018974	YP02105054	I03	hsa-miR-1295a	MIMAT0005885	YP02113131
A04	hsa-miR-3173-3p	MIMAT0015048	YP02114341	I04	hsa-miR-3687	MIMAT0018115	YP02114671
A05	hsa-miR-320e	MIMAT0015072	YP02103235	I05	hsa-miR-1261	MIMAT0005913	YP00205716
A06	hsa-miR-5706	MIMAT0022500	YP02113815	I06	hsa-miR-1297	MIMAT0005886	YP02100221
A07	hsa-miR-3613-3p	MIMAT0017991	YP02100646	I07	hsa-miR-4305	MIMAT0016857	YP02112386
A08	hsa-miR-3154	MIMAT0015028	YP02114160	I08	hsa-miR-3692-3p	MIMAT0018122	YP02113419
A09	hsa-miR-3913-5p	MIMAT0018187	YP02103624	I09	hsa-miR-1298-5p	MIMAT0005800	YP02108491
A10	hsa-miR-3670	MIMAT0018093	YP02104621	I10	hsa-miR-2052	MIMAT0009977	YP02116015
A11	hsa-miR-2355-3p	MIMAT0017950	YP02118968	I11	hsa-miR-103a-2-5p	MIMAT0009196	YP00204609
A12	hsa-miR-4659a-5p	MIMAT0019726	YP02114501	I12	hsa-miR-3130-3p	MIMAT0014994	YP02108951
A13	hsa-miR-29b-1-5p	MIMAT0004514	YP00205623	I13	hsa-miR-23c	MIMAT0018000	YP02105008
A14	hsa-miR-3667-5p	MIMAT0018089	YP02108319	I14	hsa-miR-1197	MIMAT0005955	YP02113584
A15	hsa-miR-4276	MIMAT0016904	YP02106443	I15	hsa-miR-3119	MIMAT0014981	YP02104128
A16	hsa-miR-302a-5p	MIMAT0000683	YP00204699	I16	hsa-miR-3135b	MIMAT0018985	YP02114177
A17	hsa-miR-3612	MIMAT0017989	YP02111812	I17	hsa-miR-3200-3p	MIMAT0015085	YP02112132
A18	hsa-miR-4254	MIMAT0016884	YP02112721	I18	hsa-miR-1231	MIMAT0005586	YP02119498
A19	hsa-miR-4260	MIMAT0016881	YP02101761	I19	hsa-miR-378c	MIMAT0016847	YP02109675
A20	hsa-miR-4324	MIMAT0016876	YP02108405	I20	hsa-miR-1233-3p	MIMAT0005588	YP02116905
A21	hsa-miR-3938	MIMAT0018353	YP02116726	I21	hsa-miR-3168	MIMAT0015043	YP02106500
A22	hsa-miR-642b-3p	MIMAT0018444	YP02103759	I22	hsa-miR-1236-3p	MIMAT0005591	YP02111648
A23	hsa-miR-548v	MIMAT0015020	YP02107607	I23	hsa-miR-4518	MIMAT0019055	YP02106110
A24	hsa-miR-4262	MIMAT0016894	YP02118051	I24	hsa-miR-3939	MIMAT0018355	YP02115536
B01	hsa-miR-3118	MIMAT0014980	YP02106587	J01	hsa-miR-2115-3p	MIMAT0011159	YP02103136
B02	hsa-miR-3919	MIMAT0018193	YP02109375	J02	hsa-miR-3117-3p	MIMAT0014979	YP02108014
B03	hsa-miR-1307-3p	MIMAT0005951	YP02103132	J03	hsa-miR-4477b	MIMAT0019005	YP02105940
B04	hsa-miR-3137	MIMAT0015005	YP02109715	J04	hsa-miR-1250-5p	MIMAT0005902	YP00205651
B05	hsa-miR-3605-3p	MIMAT0017982	YP02113576	J05	hsa-miR-3622a-5p	MIMAT0018003	YP02108403
B06	hsa-miR-3666	MIMAT0018088	YP02101822	J06	hsa-miR-3198	MIMAT0015083	YP02119447
B07	hsa-miR-3149	MIMAT0015022	YP02113078	J07	hsa-miR-1251-5p	MIMAT0005903	YP00205715
B08	hsa-miR-323b-3p	MIMAT0015050	YP02100923	J08	hsa-miR-3201	MIMAT0015086	YP02119150
B09	hsa-miR-4255	MIMAT0016885	YP02110254	J09	hsa-miR-3922-3p	MIMAT0018197	YP02117005
B10	hsa-miR-3912-3p	MIMAT0018186	YP02102954	J10	hsa-miR-3605-5p	MIMAT0017981	YP02103942
B11	hsa-miR-3925-5p	MIMAT0018200	YP02106663	J11	hsa-miR-148a-5p	MIMAT0004549	YP00204188
B12	hsa-miR-4285	MIMAT0016913	YP02104979	J12	hsa-miR-3909	MIMAT0018183	YP02103937
B13	hsa-miR-4252	MIMAT0016886	YP02105111	J13	hsa-miR-4310	MIMAT0016862	YP02119555
B14	hsa-miR-1193	MIMAT0015049	YP02115297	J14	hsa-miR-3914	MIMAT0018188	YP02105750
B15	hsa-miR-3944-3p	MIMAT0018360	YP02119421	J15	hsa-miR-3678-3p	MIMAT0018103	YP02115318
B16	hsa-miR-509-5p	MIMAT0004779	YP02103795	J16	hsa-miR-4319	MIMAT0016870	YP02111871
B17	hsa-miR-4732-3p	MIMAT0019856	YP02106069	J17	hsa-miR-3679-5p	MIMAT0018104	YP02109537
B18	hsa-miR-4536-3p	MIMAT0020959	YP02102426	J18	hsa-miR-4300	MIMAT0016853	YP02117892
B19	hsa-miR-1260b	MIMAT0015041	YP02106579	J19	hsa-miR-4265	MIMAT0016891	YP02118209
B20	hsa-miR-3667-3p	MIMAT0018090	YP02112332	J20	hsa-miR-3150b-3p	MIMAT0018194	YP02109258
B21	hsa-miR-4312	MIMAT0016864	YP02107126	J21	hsa-miR-3661	MIMAT0018082	YP02103023
B22	hsa-miR-4326	MIMAT0016888	YP02103754	J22	hsa-miR-2116-3p	MIMAT0011161	YP02104585
B23	hsa-miR-4304	MIMAT0016854	YP02103864	J23	hsa-miR-1293	MIMAT0005883	YP02110424
B24	hsa-miR-3918	MIMAT0018192	YP02119250	J24	hsa-miR-4279	MIMAT0016909	YP02106388
C01	hsa-miR-3674	MIMAT0018097	YP02109756	K01	hsa-miR-4303	MIMAT0016856	YP02101703

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
C02	hsa-miR-3174	MIMAT0015051	YP02117538	K02	hsa-miR-3126-3p	MIMAT0015377	YP02115717
C03	hsa-miR-3671	MIMAT0018094	YP02109323	K03	hsa-miR-3176	MIMAT0015053	YP02113244
C04	hsa-miR-466	MIMAT0015002	YP02113885	K04	hsa-miR-676-5p	MIMAT0018203	YP02114028
C05	hsa-miR-3175	MIMAT0015052	YP02112727	K05	hsa-miR-374a-3p	MIMAT0004688	YP00205692
C06	hsa-miR-4294	MIMAT0016849	YP02111533	K06	hsa-miR-3167	MIMAT0015042	YP02101060
C07	hsa-miR-3182	MIMAT0015062	YP02115918	K07	hsa-miR-497-3p	MIMAT0004768	YP00205660
C08	hsa-miR-548aa	MIMAT0018447	YP02106297	K08	hsa-miR-3166	MIMAT0015040	YP02104301
C09	hsa-miR-3678-5p	MIMAT0018102	YP02113776	K09	hsa-miR-3927-3p	MIMAT0018202	YP02106115
C10	hsa-miR-297	MIMAT0004450	YP02119527	K10	hsa-miR-3126-5p	MIMAT0014989	YP02110514
C11	hsa-miR-3124-5p	MIMAT0014986	YP02119537	K11	hsa-miR-3184-5p	MIMAT0015064	YP02114720
C12	hsa-miR-3162-5p	MIMAT0015036	YP02107673	K12	hsa-miR-545-5p	MIMAT0004785	YP00205682
C13	hsa-miR-3684	MIMAT0018112	YP02109780	K13	hsa-miR-3065-5p	MIMAT0015066	YP02106617
C14	hsa-miR-3659	MIMAT0018080	YP02116948	K14	hsa-miR-550b-3p	MIMAT0018445	YP02103887
C15	hsa-miR-4325	MIMAT0016887	YP02102566	K15	hsa-miR-4732-3p	MIMAT0019856	YP02106069
C16	hsa-miR-2355-5p	MIMAT0016895	YP02100607	K16	hsa-miR-3617-5p	MIMAT0017997	YP02116991
C17	hsa-miR-3198	MIMAT0015083	YP02119447	K17	hsa-miR-3649	MIMAT0018069	YP02109330
C18	hsa-miR-676-3p	MIMAT0018204	YP02112616	K18	hsa-miR-4799-5p	MIMAT0019976	YP02112015
C19	hsa-miR-3199	MIMAT0015084	YP02116364	K19	hsa-miR-3135b	MIMAT0018985	YP02114177
C20	hsa-miR-1321	MIMAT0005952	YP02118243	K20	hsa-miR-3650	MIMAT0018070	YP02102887
C21	hsa-miR-670-5p	MIMAT0010357	YP02109934	K21	hsa-miR-4290	MIMAT0016921	YP02113894
C22	hsa-miR-3130-5p	MIMAT0014995	YP02101081	K22	hsa-miR-3161	MIMAT0015035	YP02117476
C23	hsa-miR-4293	MIMAT0016848	YP02119330	K23	hsa-miR-3128	MIMAT0014991	YP02119445
C24	hsa-miR-3921	MIMAT0018196	YP02109256	K24	hsa-miR-4267	MIMAT0016893	YP02108101
D01	hsa-miR-3682-3p	MIMAT0018110	YP02111637	L01	hsa-miR-4299	MIMAT0016851	YP02103113
D02	hsa-miR-548s	MIMAT0014987	YP02117763	L02	hsa-miR-3908	MIMAT0018182	YP02102183
D03	hsa-miR-633	MIMAT0003303	YP02106059	L03	hsa-miR-3681-5p	MIMAT0018108	YP02119120
D04	hsa-miR-3907	MIMAT0018179	YP02118102	L04	hsa-miR-3139	MIMAT0015007	YP02104758
D05	hsa-miR-485-5p	MIMAT0002175	YP02112548	L05	hsa-miR-656-3p	MIMAT0003332	YP02119466
D06	hsa-miR-3943	MIMAT0018359	YP02109349	L06	hsa-miR-3664-5p	MIMAT0018086	YP02105587
D07	hsa-miR-3675-5p	MIMAT0018098	YP02104697	L07	hsa-miR-664a-5p	MIMAT0005948	YP02115612
D08	hsa-miR-378b	MIMAT0014999	YP02119100	L08	hsa-miR-550a-3p	MIMAT0003257	YP00204795
D09	hsa-miR-3685	MIMAT0018113	YP02106245	L09	hsa-miR-3121-3p	MIMAT0014983	YP02119245
D10	hsa-miR-3606-5p	MIMAT0017983	YP02113117	L10	hsa-miR-1825	MIMAT0006765	YP02119443
D11	hsa-miR-1299	MIMAT0005887	YP02114891	L11	hsa-miR-3135a	MIMAT0015001	YP02107752
D12	hsa-miR-3120-3p	MIMAT0014982	YP02118957	L12	hsa-miR-514b-5p	MIMAT0015087	YP02115289
D13	hsa-miR-3194-5p	MIMAT0015078	YP02101080	L13	hsa-miR-1255a	MIMAT0005906	YP00205597
D14	hsa-miR-3164	MIMAT0015038	YP02113491	L14	hsa-miR-3652	MIMAT0018072	YP02117517
D15	hsa-miR-4261	MIMAT0016890	YP02115003	L15	hsa-miR-3714	MIMAT0018165	YP02112326
D16	hsa-miR-1301-3p	MIMAT0005797	YP02111482	L16	hsa-miR-1257	MIMAT0005908	YP02113278
D17	hsa-miR-3681-3p	MIMAT0018109	YP02112938	L17	hsa-miR-3186-5p	MIMAT0015067	YP02101723
D18	hsa-miR-575	MIMAT0003240	YP02114484	L18	hsa-miR-4306	MIMAT0016858	YP02105455
D19	hsa-miR-3713	MIMAT0018164	YP02117308	L19	hsa-miR-3125	MIMAT0014988	YP02118320
D20	hsa-miR-3616-3p	MIMAT0017996	YP02114542	L20	hsa-miR-1262	MIMAT0005914	YP00205672
D21	hsa-miR-3616-5p	MIMAT0017995	YP02112077	L21	hsa-miR-4251	MIMAT0016883	YP02100490
D22	hsa-miR-4288	MIMAT0016918	YP02102818	L22	hsa-miR-1263	MIMAT0005915	YP00204174
D23	hsa-miR-1302	MIMAT0005890	YP02115302	L23	hsa-miR-449c-5p	MIMAT0010251	YP02119451
D24	hsa-miR-4273	MIMAT0016903	YP02105868	L24	hsa-miR-613	MIMAT0003281	YP02107726
E01	hsa-miR-1322	MIMAT0005953	YP02119522	M01	hsa-miR-4258	MIMAT0016879	YP02119162
E02	hsa-miR-4291	MIMAT0016922	YP02115206	M02	hsa-miR-1266-5p	MIMAT0005920	YP00206044
E03	hsa-miR-1278	MIMAT0005936	YP02109979	M03	hsa-miR-4309	MIMAT0016859	YP02105234

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
E04	hsa-miR-764	MIMAT0010367	YP02110982	M04	hsa-miR-3177-3p	MIMAT0015054	YP02103017
E05	hsa-miR-4799-5p	MIMAT0019976	YP02112015	M05	hsa-miR-3668	MIMAT0018091	YP02112764
E06	hsa-miR-3158-5p	MIMAT0019211	YP02116711	M06	hsa-miR-4323	MIMAT0016875	YP02110325
E07	hsa-miR-3945	MIMAT0018361	YP02115442	M07	hsa-miR-3660	MIMAT0018081	YP02106800
E08	hsa-miR-1281	MIMAT0005939	YP02119361	M08	hsa-miR-2278	MIMAT0011778	YP02105280
E09	hsa-miR-2117	MIMAT0011162	YP02103840	M09	hsa-miR-3156-5p	MIMAT0015030	YP02106011
E10	hsa-miR-1282	MIMAT0005940	YP02104324	M10	hsa-miR-1276	MIMAT0005930	YP02107183
E11	hsa-miR-4283	MIMAT0016914	YP02101408	M11	hsa-miR-4275	MIMAT0016905	YP02101219
E12	hsa-miR-1283	MIMAT0005799	YP02104470	M12	hsa-miR-1277-3p	MIMAT0005933	YP02114217
E13	hsa-miR-3186-3p	MIMAT0015068	YP02108664	M13	hsa-miR-3074-3p	MIMAT0015027	YP02116315
E14	hsa-miR-542-3p	MIMAT0003389	YP00205444	M14	hsa-miR-544b	MIMAT0015004	YP02109575
E15	hsa-miR-1284	MIMAT0005941	YP02107484	M15	hsa-miR-3920	MIMAT0018195	YP02101091
E16	hsa-miR-3662	MIMAT0018083	YP02118684	M16	hsa-miR-3193	MIMAT0015077	YP02100644
E17	hsa-miR-1323	MIMAT0005795	YP02103741	M17	hsa-miR-4296	MIMAT0016845	YP02100857
E18	hsa-miR-4311	MIMAT0016863	YP02108881	M18	hsa-miR-4263	MIMAT0016898	YP02119396
E19	hsa-miR-1285-3p	MIMAT0005876	YP02113268	M19	hsa-miR-3159	MIMAT0015033	YP02101997
E20	hsa-miR-3163	MIMAT0015037	YP02107780	M20	hsa-miR-2115-5p	MIMAT0011158	YP02114843
E21	hsa-miR-1286	MIMAT0005877	YP02111048	M21	hsa-miR-4766-3p	MIMAT0019918	YP02104335
E22	hsa-miR-3935	MIMAT0018350	YP02108832	M22	hsa-miR-3170	MIMAT0015045	YP02103985
E23	hsa-miR-1287-5p	MIMAT0005878	YP02107699	M23	hsa-miR-548u	MIMAT0015013	YP02104219
E24	hsa-miR-2276-3p	MIMAT0011775	YP02115594	M24	hsa-miR-3190-3p	MIMAT0022839	YP02105534
F01	hsa-miR-1288-3p	MIMAT0005942	YP02101879	N01	hsa-miR-3171	MIMAT0015046	YP02105869
F02	hsa-miR-514b-3p	MIMAT0015088	YP02108297	N02	hsa-miR-3690	MIMAT0018119	YP02115943
F03	hsa-miR-1303	MIMAT0005891	YP02114522	N03	hsa-miR-4330	MIMAT0016924	YP02105730
F04	hsa-miR-4307	MIMAT0016860	YP02110717	N04	hsa-miR-4766-3p	MIMAT0019918	YP02104335
F05	hsa-miR-1289	MIMAT0005879	YP02109506	N05	hsa-miR-3150a-3p	MIMAT0015023	YP02113040
F06	hsa-miR-1827	MIMAT0006767	YP02112497	N06	hsa-miR-3153	MIMAT0015026	YP02113354
F07	hsa-miR-648	MIMAT0003318	YP02105325	N07	hsa-miR-3134	MIMAT0015000	YP02115497
F08	hsa-miR-1228-5p	MIMAT0005582	YP02112262	N08	hsa-miR-4278	MIMAT0016910	YP02103415
F09	hsa-miR-4328	MIMAT0016926	YP02117568	N09	hsa-miR-3115	MIMAT0014977	YP02104625
F10	hsa-miR-1290	MIMAT0005880	YP02118634	N10	hsa-miR-3151-5p	MIMAT0015024	YP02113164
F11	hsa-miR-3200-5p	MIMAT0017392	YP02109903	N11	hsa-miR-3157-5p	MIMAT0015031	YP02114115
F12	hsa-miR-3122	MIMAT0014984	YP02112390	N12	hsa-miR-4287	MIMAT0016917	YP02111033
F13	hsa-miR-3654	MIMAT0018074	YP02116828	N13	hsa-miR-3148	MIMAT0015021	YP02119446
F14	hsa-miR-5187-3p	MIMAT0021118	YP02115122	N14	hsa-miR-4317	MIMAT0016872	YP02117292
F15	hsa-miR-1180-3p	MIMAT0005825	YP02102171	N15	hsa-miR-3133	MIMAT0014998	YP02105647
F16	hsa-miR-4282	MIMAT0016912	YP02100348	N16	hsa-miR-3136-5p	MIMAT0015003	YP02107031
F17	hsa-miR-374c-5p	MIMAT0018443	YP02114350	N17	hsa-miR-3614-3p	MIMAT0017993	YP02116362
F18	hsa-miR-3146	MIMAT0015018	YP02108476	N18	hsa-miR-3131	MIMAT0014996	YP02105160
F19	hsa-miR-4280	MIMAT0016911	YP02118556	N19	hsa-miR-3155a	MIMAT0015029	YP02100663
F20	hsa-miR-548w	MIMAT0015060	YP02107802	N20	hsa-miR-4297	MIMAT0016846	YP02112513
F21	hsa-miR-3147	MIMAT0015019	YP02106354	N21	hsa-miR-548ah-3p	MIMAT0020957	YP02116102
F22	hsa-miR-3143	MIMAT0015012	YP02113902	N22	hsa-miR-3689e	MIMAT0019009	YP02113857
F23	hsa-miR-1207-3p	MIMAT0005872	YP00205610	N23	hsa-miR-5706	MIMAT0022500	YP02113815
F24	hsa-miR-4268	MIMAT0016896	YP02104995	N24	hsa-miR-3129-5p	MIMAT0014992	YP02114960
G01	hsa-miR-3680-5p	MIMAT0018106	YP02115428	O01	hsa-miR-3191-3p	MIMAT0015075	YP02114802
G02	hsa-miR-1224-5p	MIMAT0005458	YP02108797	O02	hsa-miR-3916	MIMAT0018190	YP02108269
G03	hsa-miR-1226-3p	MIMAT0005577	YP02102736	O03	hsa-miR-2114-5p	MIMAT0011156	YP02109221
G04	hsa-miR-4277	MIMAT0016908	YP02116118	O04	hsa-miR-4298	MIMAT0016852	YP02119449
G05	hsa-miR-3917	MIMAT0018191	YP02106642	O05	hsa-miR-3924	MIMAT0018199	YP02104754

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
G06	hsa-miR-3185	MIMAT0015065	YP02107693	O06	hsa-miR-3655	MIMAT0018075	YP02118256
G07	hsa-miR-500a-3p	MIMAT0002871	YP00205982	O07	hsa-miR-3152-3p	MIMAT0015025	YP02118670
G08	hsa-miR-186-3p	MIMAT0004612	YP02110599	O08	hsa-miR-3144-5p	MIMAT0015014	YP02116220
G09	hsa-miR-3926	MIMAT0018201	YP02110637	O09	hsa-miR-548q	MIMAT0011163	YP02103735
G10	hsa-miR-892b	MIMAT0004918	YP02110497	O10	hsa-miR-3180-3p	MIMAT0015058	YP02104649
G11	hsa-miR-3140-3p	MIMAT0015008	YP02119541	O11	hsa-miR-3145-3p	MIMAT0015016	YP02113486
G12	hsa-miR-3611	MIMAT0017988	YP02106378	O12	hsa-miR-3657	MIMAT0018077	YP02109393
G13	hsa-miR-3142	MIMAT0015011	YP02110330	O13	hsa-miR-3941	MIMAT0018357	YP02108607
G14	hsa-miR-1909-5p	MIMAT0007882	YP00205703	O14	hsa-miR-4264	MIMAT0016899	YP02118385
G15	hsa-miR-4714-3p	MIMAT0019823	YP02102316	O15	hsa-miR-3923	MIMAT0018198	YP02111856
G16	hsa-miR-4302	MIMAT0016855	YP02108009	O16	hsa-miR-548f-5p	MIMAT0015009	YP02116139
G17	hsa-miR-2114-3p	MIMAT0011157	YP02116512	O17	hsa-miR-4271	MIMAT0016901	YP02103404
G18	hsa-miR-4286	MIMAT0016916	YP02112356	O18	hsa-miR-4318	MIMAT0016869	YP02112808
G19	hsa-miR-3618	MIMAT0017998	YP02118233	O19	hsa-miR-3116	MIMAT0014978	YP02119047
G20	hsa-miR-3683	MIMAT0018111	YP02106322	O20	hsa-miR-3065-3p	MIMAT0015378	YP02114101
G21	hsa-miR-548a-3p	MIMAT0005919	YP00204752	O21	hsa-miR-196a-3p	MIMAT0004562	YP02117597
G22	hsa-miR-4314	MIMAT0016868	YP02113557	O22	hsa-miR-3141	MIMAT0015010	YP02107812
G23	hsa-miR-3132	MIMAT0014997	YP02113218	O23	hsa-miR-3179	MIMAT0015056	YP02113829
G24	hsa-miR-1910-5p	MIMAT0007884	YP00206033	O24	hsa-miR-1976	MIMAT0009451	YP02112898
H01	hsa-miR-1324	MIMAT0005956	YP02104152	P01	hsa-miR-4274	MIMAT0016906	YP02109237
H02	hsa-miR-3679-3p	MIMAT0018105	YP02118543	P02	hsa-miR-3165	MIMAT0015039	YP02118817
H03	hsa-miR-3658	MIMAT0018078	YP02119545	P03	hsa-miR-4308	MIMAT0016861	YP02119335
H04	hsa-miR-3915	MIMAT0018189	YP02110570	P04	hsa-miR-3160-3p	MIMAT0015034	YP02114929
H05	hsa-miR-3688-3p	MIMAT0018116	YP02107387	P05	hsa-miR-4316	MIMAT0016867	YP02116593
H06	hsa-miR-4253	MIMAT0016882	YP02119435	P06	hsa-miR-3144-3p	MIMAT0015015	YP02109767
H07	hsa-miR-4301	MIMAT0016850	YP02116818	P07	hsa-miR-942-3p	MIMAT0026734	YP02118338
H08	hsa-miR-3158-3p	MIMAT0015032	YP02107763	P08	hsa-miR-95-5p	MIMAT0026473	YP02101365
H09	hsa-miR-4999-5p	MIMAT0021017	YP02111952	P09	hsa-miR-3910	MIMAT0018184	YP02108262
H10	hsa-miR-3202	MIMAT0015089	YP02116603	P10	hsa-miR-1304-5p	MIMAT0005892	YP02110866
H11	hsa-miR-3127-5p	MIMAT0014990	YP02104545	P11	hsa-miR-301b-3p	MIMAT0004958	YP00204390
H12	hsa-miR-4470	MIMAT0018997	YP02118874	P12	BlankSpot	N/A	YP02119291
H13	hsa-miR-449c-3p	MIMAT0013771	YP02100963	P13	SNORD44 (hsa)	N/A	YP00203902
H14	hsa-miR-1291	MIMAT0005881	YP02107164	P14	SNORD38B (hsa)	N/A	YP00203901
H15	hsa-miR-3942-5p	MIMAT0018358	YP02101651	P15	SNORD49A (hsa)	N/A	YP00203904
H16	hsa-miR-3691-5p	MIMAT0018120	YP02113779	P16	U6 snRNA (v2)	N/A	YP02119464
H17	hsa-miR-657	MIMAT0003335	YP02108736	P17	cel-miR-39-3p	MIMAT0000010	YP00203952
H18	hsa-miR-3934-5p	MIMAT0018349	YP02109655	P18	UniSp2	N/A	YP00203950
H19	hsa-miR-4257	MIMAT0016878	YP02105566	P19	UniSp4	N/A	YP00203953
H20	hsa-miR-3672	MIMAT0018095	YP02107559	P20	UniSp5	N/A	YP00203955
H21	hsa-miR-3651	MIMAT0018071	YP02113836	P21	UniSp6	N/A	YP00203954
H22	hsa-miR-2116-5p	MIMAT0011160	YP02110412	P22	UniSP3	N/A	YP02119288
H23	hsa-miR-2909	MIMAT0013863	YP02108840	P23	UniSP3	N/A	YP02119288
H24	hsa-miR-3622b-3p	MIMAT0018006	YP02118315	P24	UniSP3	N/A	YP02119288

miRNA Table: Plate 4

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
A01	hsa-miR-5192	MIMAT0021123	YP02105152	I01	hsa-miR-4259	MIMAT0016880	YP02108171
A02	hsa-miR-181b-3p	MIMAT0022692	YP02106895	I02	hsa-miR-3619-3p	MIMAT0019219	YP02106037
A03	hsa-miR-4429	MIMAT0018944	YP02110346	I03	hsa-miR-3688-5p	MIMAT0019223	YP02115698
A04	hsa-miR-548ah-	MIMAT0018972	YP02108508	I04	hsa-miR-4437	MIMAT0018953	YP02105246
A05	hsa-miR-4786-5p	MIMAT0019954	YP02115269	I05	hsa-miR-4726-5p	MIMAT0019845	YP02113284
A06	hsa-miR-4740-3p	MIMAT0019870	YP02107508	I06	hsa-miR-3664-3p	MIMAT0019220	YP02109280
A07	hsa-miR-4790-3p	MIMAT0019962	YP02100995	I07	hsa-miR-4639-3p	MIMAT0019698	YP02116422
A08	hsa-miR-4649-5p	MIMAT0019711	YP02116842	I08	hsa-miR-5682	MIMAT0022470	YP02102240
A09	hsa-miR-4499	MIMAT0019035	YP02119336	I09	hsa-miR-4538	MIMAT0019081	YP02104874
A10	hsa-miR-4436b-5p	MIMAT0019940	YP02112379	I10	hsa-miR-4703-3p	MIMAT0019802	YP02114058
A11	hsa-miR-2861	MIMAT0013802	YP02119565	I11	hsa-miR-371b-5p	MIMAT0019892	YP02108253
A12	hsa-miR-4751	MIMAT0019888	YP02114926	I12	hsa-miR-4771	MIMAT0019925	YP02101627
A13	hsa-miR-4458	MIMAT0018980	YP02117234	I13	hsa-miR-450a-2-3p	MIMAT0031074	YP02107563
A14	hsa-miR-4744	MIMAT0019875	YP02102581	I14	hsa-miR-4714-5p	MIMAT0019822	YP02103832
A15	hsa-miR-4641	MIMAT0019701	YP02102652	I15	hsa-miR-1304-3p	MIMAT0022720	YP02106758
A16	hsa-miR-4460	MIMAT0018982	YP02118742	I16	hsa-miR-5698	MIMAT0022491	YP02109293
A17	hsa-miR-4487	MIMAT0019021	YP02103660	I17	hsa-miR-203b-5p	MIMAT0019813	YP02114497
A18	hsa-miR-5579-3p	MIMAT0022270	YP02119457	I18	hsa-miR-4501	MIMAT0019037	YP02111910
A19	hsa-miR-4482-5p	MIMAT0019016	YP02110162	I19	hsa-miR-561-5p	MIMAT0022706	YP02117324
A20	hsa-miR-4788	MIMAT0019958	YP02108226	I20	hsa-miR-5007-5p	MIMAT0021035	YP02112000
A21	hsa-miR-3977	MIMAT0019362	YP02117768	I21	hsa-miR-4467	MIMAT0018994	YP02116307
A22	hsa-miR-4531	MIMAT0019070	YP02112923	I22	hsa-miR-3074-5p	MIMAT0019208	YP02110958
A23	hsa-miR-5011-5p	MIMAT0021045	YP02103375	I23	hsa-miR-4266	MIMAT0016892	YP02119331
A24	hsa-miR-4535	MIMAT0019075	YP02117966	I24	hsa-miR-4676-3p	MIMAT0019759	YP02102548
B01	hsa-miR-4480	MIMAT0019014	YP02100758	J01	hsa-miR-1247-3p	MIMAT0022721	YP02103125
B02	hsa-miR-4733-5p	MIMAT0019857	YP02103988	J02	hsa-miR-4521	MIMAT0019058	YP02109260
B03	hsa-miR-568	MIMAT0003232	YP02119300	J03	hsa-miR-203b-3p	MIMAT0019814	YP02103829
B04	hsa-miR-548ao-	MIMAT0021029	YP02119456	J04	hsa-miR-4666b	MIMAT0022485	YP02119223
B05	hsa-miR-4661-5p	MIMAT0019729	YP02115180	J05	hsa-miR-4749-5p	MIMAT0019885	YP02119410
B06	hsa-miR-4735-5p	MIMAT0019860	YP02113099	J06	hsa-miR-4534	MIMAT0019073	YP02107694
B07	hsa-miR-642b-5p	MIMAT0022736	YP02109732	J07	hsa-miR-4477a	MIMAT0019004	YP02118912
B08	hsa-miR-4720-5p	MIMAT0019833	YP02115260	J08	hsa-miR-451b	MIMAT0019840	YP02103713
B09	hsa-miR-539-3p	MIMAT0022705	YP02101713	J09	hsa-miR-4670-3p	MIMAT0019751	YP02110005
B10	hsa-miR-4767	MIMAT0019919	YP02119597	J10	hsa-miR-4527	MIMAT0019066	YP02105459
B11	hsa-miR-4649-3p	MIMAT0019712	YP02112040	J11	hsa-miR-4430	MIMAT0018945	YP02114986
B12	hsa-miR-4431	MIMAT0018947	YP02110679	J12	hsa-miR-548aj-3p	MIMAT0018990	YP02100008
B13	hsa-miR-5047	MIMAT0020541	YP02114300	J13	hsa-miR-5091	MIMAT0021083	YP02110257
B14	hsa-miR-4659b-5p	MIMAT0019733	YP02107053	J14	hsa-miR-3646	MIMAT0018065	YP02116669
B15	hsa-miR-5588-5p	MIMAT0022295	YP02101379	J15	hsa-miR-4483	MIMAT0019017	YP02102993
B16	hsa-miR-506-5p	MIMAT0022701	YP02119376	J16	hsa-miR-4443	MIMAT0018961	YP02104824
B17	hsa-miR-4769-5p	MIMAT0019922	YP02105772	J17	hsa-miR-3689f	MIMAT0019010	YP02106748
B18	hsa-miR-4770	MIMAT0019924	YP02113433	J18	hsa-miR-4723-5p	MIMAT0019838	YP02101281
B19	hsa-miR-5587-3p	MIMAT0022290	YP02107712	J19	hsa-miR-548g-5p	MIMAT0022722	YP02114059
B20	hsa-miR-3177-5p	MIMAT0019215	YP02106664	J20	hsa-miR-4683	MIMAT0019768	YP02118119
B21	hsa-miR-4732-5p	MIMAT0019855	YP02113664	J21	hsa-miR-4786-3p	MIMAT0019955	YP02101935
B22	hsa-miR-4693-3p	MIMAT0019785	YP02109295	J22	hsa-miR-2682-3p	MIMAT0013518	YP02112755
B23	hsa-miR-4800-5p	MIMAT0019978	YP02115647	J23	hsa-miR-4804-3p	MIMAT0019985	YP02111044
B24	hsa-miR-4781-5p	MIMAT0019942	YP02100671	J24	hsa-miR-1343-3p	MIMAT0019776	YP02110797
C01	hsa-miR-4659b-3p	MIMAT0019734	YP02100566	K01	hsa-miR-345-3p	MIMAT0022698	YP02101775

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
C02	hsa-miR-4448	MIMAT0018967	YP02108427	K02	hsa-miR-378g	MIMAT0018937	YP02110329
C03	hsa-miR-4471	MIMAT0018998	YP02100524	K03	hsa-miR-4502	MIMAT0019038	YP02119075
C04	hsa-miR-211-3p	MIMAT0022694	YP02119691	K04	hsa-miR-4678	MIMAT0019762	YP02108423
C05	hsa-miR-4426	MIMAT0018941	YP02112311	K05	hsa-miR-4491	MIMAT0019026	YP02104904
C06	hsa-miR-4436a	MIMAT0018952	YP02103033	K06	hsa-miR-4795-3p	MIMAT0019969	YP02102711
C07	hsa-miR-3976	MIMAT0019361	YP02110827	K07	hsa-miR-548ae-3p	MIMAT0018954	YP02104542
C08	hsa-miR-4777-5p	MIMAT0019934	YP02111756	K08	hsa-miR-1238-5p	MIMAT0022947	YP02107444
C09	hsa-miR-4690-5p	MIMAT0019779	YP02100659	K09	hsa-miR-5571-3p	MIMAT0022258	YP02111694
C10	hsa-miR-1185-1-	MIMAT0022838	YP02103158	K10	hsa-miR-5680	MIMAT0022468	YP02108554
C11	hsa-miR-1245b-5p	MIMAT0019950	YP02102232	K11	hsa-miR-5580-3p	MIMAT0022274	YP02102108
C12	hsa-miR-3940-3p	MIMAT0018356	YP02119252	K12	hsa-miR-5591-5p	MIMAT0022301	YP02106649
C13	hsa-miR-548as-5p	MIMAT0022267	YP02106262	K13	hsa-miR-4802-3p	MIMAT0019982	YP02114328
C14	hsa-miR-3621	MIMAT0018002	YP02107851	K14	hsa-miR-548ai	MIMAT0019024	YP02104648
C15	hsa-miR-5582-5p	MIMAT0022279	YP02113310	K15	hsa-miR-518e-5p	MIMAT0005450	YP00206040
C16	hsa-miR-4666a-5p	MIMAT0019741	YP02103330	K16	hsa-miR-518d-5p	MIMAT0005456	YP00204461
C17	hsa-miR-4717-5p	MIMAT0019829	YP02109154	K17	hsa-miR-5683	MIMAT0022472	YP02104302
C18	hsa-miR-3196	MIMAT0015080	YP02119217	K18	hsa-miR-3189-3p	MIMAT0015071	YP02113106
C19	hsa-miR-1268b	MIMAT0018925	YP02100365	K19	hsa-miR-4478	MIMAT0019006	YP02108950
C20	hsa-miR-4709-3p	MIMAT0019812	YP02115820	K20	hsa-miR-4449	MIMAT0018968	YP02101057
C21	hsa-miR-4464	MIMAT0018988	YP02116618	K21	hsa-miR-5694	MIMAT0022487	YP02101037
C22	hsa-miR-4737	MIMAT0019863	YP02115554	K22	hsa-miR-4520-3p	MIMAT0019057	YP02115335
C23	hsa-miR-664b-5p	MIMAT0022271	YP02109261	K23	hsa-miR-4759	MIMAT0019905	YP02116285
C24	hsa-miR-4679	MIMAT0019763	YP02115588	K24	hsa-miR-4506	MIMAT0019042	YP02116642
D01	hsa-miR-4638-3p	MIMAT0019696	YP02119408	L01	hsa-miR-5087	MIMAT0021079	YP02103586
D02	hsa-miR-4765	MIMAT0019916	YP02102134	L02	hsa-miR-3152-5p	MIMAT0019207	YP02106814
D03	hsa-miR-5696	MIMAT0022489	YP02103209	L03	hsa-miR-5702	MIMAT0022495	YP02111721
D04	hsa-miR-4522	MIMAT0019060	YP02110473	L04	hsa-miR-4717-3p	MIMAT0019830	YP02102971
D05	hsa-miR-4731-3p	MIMAT0019854	YP02115342	L05	hsa-miR-513c-3p	MIMAT0022728	YP02100149
D06	hsa-miR-5004-5p	MIMAT0021027	YP02103151	L06	hsa-miR-4671-3p	MIMAT0019753	YP02108547
D07	hsa-miR-4764-5p	MIMAT0019914	YP02118604	L07	hsa-miR-548aq-5p	MIMAT0022263	YP02113472
D08	hsa-miR-3183	MIMAT0015063	YP02109689	L08	hsa-miR-4329	MIMAT0016923	YP02101238
D09	hsa-miR-4436b-3p	MIMAT0019941	YP02117489	L09	hsa-miR-4511	MIMAT0019048	YP02117909
D10	hsa-miR-4639-5p	MIMAT0019697	YP02110154	L10	hsa-miR-3689a-3p	MIMAT0018118	YP02115246
D11	hsa-miR-5707	MIMAT0022501	YP02112846	L11	hsa-miR-4539	MIMAT0019082	YP02111967
D12	hsa-miR-4780	MIMAT0019939	YP02101575	L12	hsa-miR-4774-3p	MIMAT0019930	YP02118334
D13	hsa-miR-4701-5p	MIMAT0019798	YP02108244	L13	hsa-miR-5093	MIMAT0021085	YP02107366
D14	hsa-miR-2681-5p	MIMAT0013515	YP02108377	L14	hsa-miR-216a-3p	MIMAT0022844	YP02116567
D15	hsa-miR-5590-3p	MIMAT0022300	YP02115723	L15	hsa-miR-4474-5p	MIMAT0019234	YP02106856
D16	hsa-miR-5196-3p	MIMAT0021129	YP02100723	L16	hsa-miR-378f	MIMAT0018932	YP02105315
D17	hsa-miR-5701	MIMAT0022494	YP02106962	L17	hsa-miR-4781-3p	MIMAT0019943	YP02109946
D18	hsa-miR-4524b-5p	MIMAT0022255	YP02116497	L18	hsa-miR-4524b-3p	MIMAT0022256	YP02111618
D19	hsa-miR-4640-3p	MIMAT0019700	YP02107962	L19	hsa-miR-4668-3p	MIMAT0019746	YP02115463
D20	hsa-miR-548an	MIMAT0019079	YP02108431	L20	hsa-miR-5191	MIMAT0021122	YP02108188
D21	hsa-miR-4800-3p	MIMAT0019979	YP02111766	L21	hsa-miR-5004-3p	MIMAT0021028	YP02113266
D22	hsa-miR-219b-3p	MIMAT0019748	YP02102754	L22	hsa-miR-4650-5p	MIMAT0019713	YP02115568
D23	hsa-miR-4481	MIMAT0019015	YP02119450	L23	hsa-miR-5580-5p	MIMAT0022273	YP02108180
D24	hsa-miR-4636	MIMAT0019693	YP02112644	L24	hsa-miR-4735-3p	MIMAT0019861	YP02103941
E01	hsa-miR-5581-3p	MIMAT0022276	YP02117351	M01	hsa-miR-4445-5p	MIMAT0018963	YP02115139
E02	hsa-miR-5006-3p	MIMAT0021034	YP02107641	M02	hsa-miR-3160-5p	MIMAT0019212	YP02110391
E03	hsa-miR-1228-3p	MIMAT0005583	YP02119681	M03	hsa-miR-5690	MIMAT0022482	YP02110964

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
E04	hsa-miR-4653-5p	MIMAT0019718	YP02104820	M04	hsa-miR-374c-3p	MIMAT0022735	YP02105481
E05	hsa-miR-4711-5p	MIMAT0019816	YP02103794	M05	hsa-miR-5009-3p	MIMAT0021042	YP02109300
E06	hsa-miR-1587	MIMAT0019077	YP02119501	M06	hsa-miR-4440	MIMAT0018958	YP02112148
E07	hsa-miR-3680-3p	MIMAT0018107	YP02102098	M07	hsa-miR-3064-5p	MIMAT0019864	YP02107200
E08	hsa-miR-4668-5p	MIMAT0019745	YP02114533	M08	hsa-miR-4789-5p	MIMAT0019959	YP02118087
E09	hsa-miR-4667-5p	MIMAT0019743	YP02111936	M09	hsa-miR-4453	MIMAT0018975	YP02119375
E10	hsa-miR-4754	MIMAT0019894	YP02107017	M10	hsa-miR-625-5p	MIMAT0003294	YP02119372
E11	hsa-miR-3682-5p	MIMAT0019222	YP02103591	M11	hsa-miR-4423-3p	MIMAT0018936	YP02111002
E12	hsa-miR-1295b-3p	MIMAT0022294	YP02111901	M12	hsa-miR-4256	MIMAT0016877	YP02100638
E13	hsa-miR-4648	MIMAT0019710	YP02103519	M13	hsa-miR-4743-5p	MIMAT0019874	YP02115797
E14	hsa-miR-4447	MIMAT0018966	YP02119381	M14	hsa-miR-4646-3p	MIMAT0019708	YP02111753
E15	hsa-miR-3123	MIMAT0014985	YP02119567	M15	hsa-miR-4644	MIMAT0019704	YP02115350
E16	hsa-miR-4482-3p	MIMAT0020958	YP02115702	M16	hsa-miR-5100	MIMAT0022259	YP02112520
E17	hsa-miR-4716-5p	MIMAT0019826	YP02119441	M17	hsa-miR-4708-5p	MIMAT0019809	YP02113054
E18	hsa-miR-4456	MIMAT0018978	YP02109793	M18	hsa-miR-4646-5p	MIMAT0019707	YP02118828
E19	hsa-miR-4442	MIMAT0018960	YP02115257	M19	hsa-miR-518a-5p	MIMAT0005457	YP00205718
E20	hsa-miR-5193	MIMAT0021124	YP02103663	M20	hsa-miR-4686	MIMAT0019773	YP02106034
E21	hsa-miR-4656	MIMAT0019723	YP02102459	M21	hsa-miR-499b-3p	MIMAT0019898	YP02109070
E22	hsa-miR-1915-5p	MIMAT0007891	YP02119687	M22	hsa-miR-4795-5p	MIMAT0019968	YP02119413
E23	hsa-miR-4748	MIMAT0019884	YP02114288	M23	hsa-miR-4424	MIMAT0018939	YP02107719
E24	hsa-miR-4725-3p	MIMAT0019844	YP02106604	M24	hsa-miR-4652-5p	MIMAT0019716	YP02105494
F01	hsa-miR-4736	MIMAT0019862	YP02119341	N01	hsa-miR-2681-3p	MIMAT0013516	YP02114773
F02	hsa-miR-4523	MIMAT0019061	YP02102394	N02	hsa-miR-4760-3p	MIMAT0019907	YP02101287
F03	hsa-miR-5571-5p	MIMAT0022257	YP02109692	N03	hsa-miR-548a-3p	MIMAT0022266	YP02100148
F04	hsa-miR-4645-3p	MIMAT0019706	YP02111813	N04	hsa-miR-4799-3p	MIMAT0019977	YP02101030
F05	hsa-miR-4692	MIMAT0019783	YP02104557	N05	hsa-miR-4724-5p	MIMAT0019841	YP02105974
F06	hsa-miR-4475	MIMAT0019002	YP02101426	N06	hsa-miR-3529-5p	MIMAT0019828	YP02110253
F07	hsa-miR-301a-5p	MIMAT0022696	YP02109957	N07	hsa-miR-5704	MIMAT0022498	YP02110143
F08	hsa-miR-873-3p	MIMAT0022717	YP02109020	N08	hsa-miR-4533	MIMAT0019072	YP02105683
F09	hsa-miR-4745-3p	MIMAT0019879	YP02109728	N09	hsa-miR-4321	MIMAT0016874	YP02117073
F10	hsa-miR-5089-5p	MIMAT0021081	YP02104012	N10	hsa-miR-5194	MIMAT0021125	YP02108038
F11	hsa-miR-2467-5p	MIMAT0019952	YP02101409	N11	hsa-miR-4528	MIMAT0019067	YP02104107
F12	hsa-miR-3117-5p	MIMAT0019197	YP02103204	N12	hsa-miR-5000-5p	MIMAT0021019	YP02102468
F13	hsa-miR-4664-3p	MIMAT0019738	YP02119369	N13	hsa-miR-5700	MIMAT0022493	YP02118073
F14	hsa-miR-4721	MIMAT0019835	YP02113239	N14	hsa-miR-4729	MIMAT0019851	YP02101851
F15	hsa-miR-4762-5p	MIMAT0019910	YP02118723	N15	hsa-miR-4766-5p	MIMAT0019917	YP02108430
F16	hsa-miR-4675	MIMAT0019757	YP02119422	N16	hsa-miR-5196-5p	MIMAT0021128	YP02102885
F17	hsa-miR-4680-5p	MIMAT0019764	YP02101985	N17	hsa-miR-5687	MIMAT0022478	YP02104733
F18	hsa-miR-4755-5p	MIMAT0019895	YP02103815	N18	hsa-miR-4782-5p	MIMAT0019944	YP02113456
F19	hsa-miR-3173-5p	MIMAT0019214	YP02119133	N19	hsa-miR-4422	MIMAT0018935	YP02111931
F20	hsa-miR-4689	MIMAT0019778	YP02119407	N20	hsa-miR-3686	MIMAT0018114	YP02111718
F21	hsa-miR-3121-5p	MIMAT0019199	YP02115858	N21	hsa-miR-3192-5p	MIMAT0015076	YP02118077
F22	hsa-miR-4652-3p	MIMAT0019717	YP02106951	N22	hsa-miR-4777-3p	MIMAT0019935	YP02107689
F23	hsa-miR-4803	MIMAT0019983	YP02108716	N23	hsa-miR-4710	MIMAT0019815	YP02113415
F24	hsa-miR-1185-2-	MIMAT0022713	YP02111736	N24	hsa-miR-4446-5p	MIMAT0019233	YP02107725
G01	hsa-miR-122b-5p	MIMAT0019876	YP02118784	O01	hsa-miR-4804-5p	MIMAT0019984	YP02101634
G02	hsa-miR-5587-5p	MIMAT0022289	YP02110826	O02	hsa-miR-4749-3p	MIMAT0019886	YP02103226
G03	hsa-miR-548a	MIMAT0019013	YP02100504	O03	hsa-miR-4715-3p	MIMAT0019825	YP02106824
G04	hsa-miR-2467-3p	MIMAT0019953	YP02116006	O04	hsa-miR-4441	MIMAT0018959	YP02112472
G05	hsa-miR-4761-3p	MIMAT0019909	YP02116589	O05	hsa-miR-3677-5p	MIMAT0019221	YP02118473

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
G06	hsa-miR-5693	MIMAT0022486	YP02108370	O06	hsa-miR-4752	MIMAT0019889	YP02111471
G07	hsa-miR-5092	MIMAT0021084	YP02114487	O07	hsa-miR-3622a-3p	MIMAT0018004	YP02114267
G08	hsa-miR-4500	MIMAT0019036	YP02119475	O08	hsa-miR-4774-5p	MIMAT0019929	YP02105955
G09	hsa-miR-5681b	MIMAT0022480	YP02106375	O09	hsa-miR-4665-3p	MIMAT0019740	YP02119591
G10	hsa-miR-4707-3p	MIMAT0019808	YP02114088	O10	hsa-miR-4704-3p	MIMAT0019804	YP02110762
G11	hsa-miR-4719	MIMAT0019832	YP02100839	O11	hsa-miR-4433a-3p	MIMAT0018949	YP02116443
G12	hsa-miR-4661-3p	MIMAT0019730	YP02105803	O12	hsa-miR-4747-5p	MIMAT0019882	YP02110146
G13	hsa-miR-652-5p	MIMAT0022709	YP02102689	O13	hsa-miR-4700-5p	MIMAT0019796	YP02113820
G14	hsa-miR-5006-5p	MIMAT0021033	YP02103769	O14	hsa-miR-5590-5p	MIMAT0022299	YP02103565
G15	hsa-miR-4693-5p	MIMAT0019784	YP02102077	O15	hsa-miR-4667-3p	MIMAT0019744	YP02113337
G16	hsa-miR-3181	MIMAT0015061	YP02119571	O16	hsa-miR-4537	MIMAT0019080	YP02116832
G17	hsa-miR-4715-5p	MIMAT0019824	YP02118306	O17	hsa-miR-4451	MIMAT0018973	YP02102301
G18	hsa-miR-5697	MIMAT0022490	YP02107133	O18	hsa-miR-4517	MIMAT0019054	YP02107279
G19	hsa-miR-4722-3p	MIMAT0019837	YP02117952	O19	hsa-miR-1469	MIMAT0007347	YP02119684
G20	hsa-miR-514a-5p	MIMAT0022702	YP02101601	O20	hsa-miR-4520-5p	MIMAT0019235	YP02116120
G21	hsa-miR-4772-3p	MIMAT0019927	YP02102250	O21	hsa-miR-3610	MIMAT0017987	YP02119507
G22	hsa-miR-1295b-5p	MIMAT0022293	YP02110835	O22	hsa-miR-5008-5p	MIMAT0021039	YP02109719
G23	hsa-miR-4750-5p	MIMAT0019887	YP02101986	O23	hsa-miR-4681	MIMAT0019766	YP02101916
G24	hsa-miR-4455	MIMAT0018977	YP02105370	O24	hsa-miR-4695-3p	MIMAT0019789	YP02115642
H01	hsa-miR-4783-3p	MIMAT0019947	YP02111616	P01	hsa-miR-4746-3p	MIMAT0019881	YP02119631
H02	hsa-miR-4536-5p	MIMAT0019078	YP02114967	P02	hsa-miR-3619-5p	MIMAT0017999	YP02114593
H03	hsa-miR-5094	MIMAT0021086	YP02119090	P03	hsa-miR-4474-3p	MIMAT0019001	YP02111455
H04	hsa-miR-4655-5p	MIMAT0019721	YP02107460	P04	hsa-miR-4479	MIMAT0019011	YP02107943
H05	hsa-miR-718	MIMAT0012735	YP02119603	P05	hsa-miR-5586-5p	MIMAT0022287	YP02103469
H06	hsa-miR-4514	MIMAT0019051	YP02100591	P06	hsa-miR-4733-3p	MIMAT0019858	YP02114943
H07	hsa-miR-4633-5p	MIMAT0019689	YP02110587	P07	hsa-miR-548a1-3p	MIMAT0022278	YP02116407
H08	hsa-miR-378e	MIMAT0018927	YP02103282	P08	hsa-miR-3150a-5p	MIMAT0019206	YP02107066
H09	hsa-miR-4691-5p	MIMAT0019781	YP02117647	P09	hsa-miR-3138	MIMAT0015006	YP02103329
H10	hsa-miR-4676-5p	MIMAT0019758	YP02111160	P10	hsa-miR-4712-5p	MIMAT0019818	YP02102710
H11	hsa-miR-548aw	MIMAT0022471	YP02106950	P11	hsa-miR-4671-5p	MIMAT0019752	YP02110416
H12	hsa-miR-4510	MIMAT0019047	YP02103091	P12	BlankSpot	N/A	YP02119291
H13	hsa-miR-766-5p	MIMAT0022714	YP02113546	P13	SNORD44 (hsa)	N/A	YP00203902
H14	hsa-miR-3124-3p	MIMAT0019200	YP02104478	P14	SNORD38B (hsa)	N/A	YP00203901
H15	hsa-miR-4802-5p	MIMAT0019981	YP02118807	P15	SNORD49A (hsa)	N/A	YP00203904
H16	hsa-miR-4697-3p	MIMAT0019792	YP02115344	P16	U6 snRNA (v2)	N/A	YP02119464
H17	hsa-miR-4476	MIMAT0019003	YP02103882	P17	cel-miR-39-3p	MIMAT0000010	YP00203952
H18	hsa-miR-3677-3p	MIMAT0018101	YP02119574	P18	UniSp2	N/A	YP00203950
H19	hsa-miR-5708	MIMAT0022502	YP02107170	P19	UniSp4	N/A	YP00203953
H20	hsa-miR-4727-3p	MIMAT0019848	YP02108445	P20	UniSp5	N/A	YP00203955
H21	hsa-miR-3191-5p	MIMAT0022732	YP02106558	P21	UniSp6	N/A	YP00203954
H22	hsa-miR-4740-5p	MIMAT0019869	YP02101587	P22	UniSP3	N/A	YP02119288
H23	hsa-miR-4540	MIMAT0019083	YP02114695	P23	UniSP3	N/A	YP02119288
H24	hsa-miR-4421	MIMAT0018934	YP02117627	P24	UniSP3	N/A	YP02119288

miRNA Table: Plate 5

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
A01	hsa-miR-4757-3p	MIMAT0019902	YP02118705	I01	hsa-miR-204-3p	MIMAT0022693	YP02113689
A02	hsa-miR-5703	MIMAT0022496	YP02113499	I02	hsa-miR-4450	MIMAT0018971	YP02107801
A03	hsa-miR-5692b	MIMAT0022497	YP02102673	I03	hsa-miR-4504	MIMAT0019040	YP02100624
A04	hsa-miR-4773	MIMAT0019928	YP02113507	I04	hsa-miR-5189-5p	MIMAT0021120	YP02101590
A05	hsa-miR-4677-3p	MIMAT0019761	YP02114408	I05	hsa-miR-548ac	MIMAT0018938	YP02118887
A06	hsa-miR-4432	MIMAT0018948	YP02106170	I06	hsa-miR-548az-3p	MIMAT0025457	YP02119349
A07	hsa-miR-4718	MIMAT0019831	YP02110067	I07	hsa-miR-548ab	MIMAT0018928	YP02109941
A08	hsa-miR-1285-5p	MIMAT0022719	YP02118966	I08	hsa-miR-3615	MIMAT0017994	YP02106599
A09	hsa-miR-4728-3p	MIMAT0019850	YP02118404	I09	hsa-miR-4789-3p	MIMAT0019960	YP02110273
A10	hsa-miR-1306-5p	MIMAT0022726	YP02105660	I10	hsa-miR-4677-5p	MIMAT0019760	YP02110460
A11	hsa-miR-499b-5p	MIMAT0019897	YP02116465	I11	hsa-miR-4797-3p	MIMAT0019973	YP02116034
A12	hsa-miR-3940-5p	MIMAT0019229	YP02114897	I12	hsa-miR-548az-5p	MIMAT0025456	YP02117220
A13	hsa-miR-5681a	MIMAT0022469	YP02117877	I13	hsa-miR-6505-5p	MIMAT0025466	YP02109966
A14	hsa-miR-4682	MIMAT0019767	YP02105579	I14	hsa-miR-3922-5p	MIMAT0019227	YP02113921
A15	hsa-miR-4746-5p	MIMAT0019880	YP02113549	I15	hsa-miR-5739	MIMAT0023116	YP02118381
A16	hsa-miR-5583-5p	MIMAT0022281	YP02119415	I16	hsa-miR-4670-5p	MIMAT0019750	YP02110349
A17	hsa-miR-3606-3p	MIMAT0022965	YP02108816	I17	hsa-miR-4647	MIMAT0019709	YP02103494
A18	hsa-miR-3180	MIMAT0018178	YP02100211	I18	hsa-miR-4281	MIMAT0016907	YP02119448
A19	hsa-miR-3911	MIMAT0018185	YP02104224	I19	hsa-miR-1255b-2-	MIMAT0022725	YP02104394
A20	hsa-miR-4505	MIMAT0019041	YP02100715	I20	hsa-miR-4663	MIMAT0019735	YP02115216
A21	hsa-miR-6505-3p	MIMAT0025467	YP02104161	I21	hsa-miR-4662a-5p	MIMAT0019731	YP02115325
A22	hsa-miR-3691-3p	MIMAT0019224	YP02101388	I22	hsa-miR-4654	MIMAT0019720	YP02113601
A23	hsa-miR-4727-5p	MIMAT0019847	YP02110633	I23	hsa-miR-548ao-3p	MIMAT0021030	YP02119455
A24	hsa-miR-4444	MIMAT0018962	YP02101842	I24	hsa-miR-4713-5p	MIMAT0019820	YP02111735
B01	hsa-miR-3120-5p	MIMAT0019198	YP02108916	J01	hsa-miR-1225-3p	MIMAT0005573	YP02114067
B02	hsa-miR-5689	MIMAT0022481	YP02110675	J02	hsa-miR-4690-3p	MIMAT0019780	YP02108184
B03	hsa-miR-4722-5p	MIMAT0019836	YP02103297	J03	hsa-miR-5589-3p	MIMAT0022298	YP02103230
B04	hsa-miR-4446-3p	MIMAT0018965	YP02102613	J04	hsa-miR-3689c	MIMAT0019007	YP02113684
B05	hsa-miR-3978	MIMAT0019363	YP02117756	J05	hsa-miR-4796-5p	MIMAT0019970	YP02101503
B06	hsa-miR-4712-3p	MIMAT0019819	YP02102929	J06	hsa-miR-2682-5p	MIMAT0013517	YP02101607
B07	hsa-miR-4776-3p	MIMAT0019933	YP02103595	J07	hsa-miR-4764-3p	MIMAT0019915	YP02104048
B08	hsa-miR-4423-5p	MIMAT0019232	YP02116196	J08	hsa-miR-4445-3p	MIMAT0018964	YP02115453
B09	hsa-miR-4633-3p	MIMAT0019690	YP02102111	J09	hsa-miR-4711-3p	MIMAT0019817	YP02108632
B10	hsa-miR-3190-5p	MIMAT0015073	YP02104632	J10	hsa-miR-1245b-3p	MIMAT0019951	YP02111711
B11	hsa-miR-5090	MIMAT0021082	YP02109448	J11	hsa-miR-5787	MIMAT0023252	YP02119640
B12	hsa-miR-5584-5p	MIMAT0022283	YP02105488	J12	hsa-miR-4696	MIMAT0019790	YP02116605
B13	hsa-miR-5195-5p	MIMAT0021126	YP02107913	J13	hsa-miR-5581-5p	MIMAT0022275	YP02116140
B14	hsa-miR-5583-3p	MIMAT0022282	YP02103724	J14	hsa-miR-4465	MIMAT0018992	YP02116480
B15	hsa-miR-5088-5p	MIMAT0021080	YP02102004	J15	hsa-miR-4796-3p	MIMAT0019971	YP02118178
B16	hsa-miR-3617-3p	MIMAT0022966	YP02111742	J16	hsa-miR-6068	MIMAT0023693	YP02105247
B17	hsa-miR-376a-2-	MIMAT0022928	YP02103763	J17	hsa-miR-5007-3p	MIMAT0021036	YP02105202
B18	hsa-miR-4776-5p	MIMAT0019932	YP02118311	J18	hsa-miR-5588-3p	MIMAT0022296	YP02115355
B19	hsa-miR-3180-5p	MIMAT0015057	YP02113327	J19	hsa-miR-3150b-5p	MIMAT0019226	YP02108329
B20	hsa-miR-4520-2-	MIMAT0020300	YP02118548	J20	hsa-miR-3157-3p	MIMAT0019210	YP02105632
B21	hsa-miR-4779	MIMAT0019938	YP02109279	J21	hsa-miR-4643	MIMAT0019703	YP02113643
B22	hsa-miR-4516	MIMAT0019053	YP02112882	J22	hsa-miR-4761-5p	MIMAT0019908	YP02108787
B23	hsa-miR-4525	MIMAT0019064	YP02111531	J23	hsa-miR-4758-3p	MIMAT0019904	YP02118014
B24	hsa-miR-4530	MIMAT0019069	YP02105012	J24	hsa-miR-4720-3p	MIMAT0019834	YP02108093
C01	hsa-miR-5688	MIMAT0022479	YP02115742	K01	hsa-miR-548ax	MIMAT0022474	YP02118201

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
C02	hsa-miR-5008-3p	MIMAT0021040	YP02112563	K02	hsa-miR-4704-5p	MIMAT0019803	YP02115935
C03	hsa-miR-3929	MIMAT0018206	YP02116898	K03	hsa-miR-3972	MIMAT0019357	YP02105947
C04	hsa-miR-4742-3p	MIMAT0019873	YP02102892	K04	hsa-miR-5011-3p	MIMAT0021046	YP02113764
C05	hsa-miR-4684-5p	MIMAT0019769	YP02107585	K05	hsa-miR-6506-3p	MIMAT0025469	YP02109545
C06	hsa-miR-4791	MIMAT0019963	YP02104047	K06	hsa-miR-4723-3p	MIMAT0019839	YP02105274
C07	hsa-miR-5699-3p	MIMAT0022492	YP02119458	K07	hsa-miR-4658	MIMAT0019725	YP02119620
C08	hsa-miR-4775	MIMAT0019931	YP02113269	K08	hsa-miR-4790-5p	MIMAT0019961	YP02107814
C09	hsa-miR-3136-3p	MIMAT0019203	YP02103774	K09	hsa-miR-4496	MIMAT0019031	YP02108046
C10	hsa-miR-4701-3p	MIMAT0019799	YP02119454	K10	hsa-miR-4768-5p	MIMAT0019920	YP02110525
C11	hsa-miR-3194-3p	MIMAT0019218	YP02119067	K11	hsa-miR-5002-3p	MIMAT0021024	YP02109038
C12	hsa-miR-3184-3p	MIMAT0022731	YP02115490	K12	hsa-miR-5188	MIMAT0021119	YP02119346
C13	hsa-miR-548at-5p	MIMAT0022277	YP02108691	K13	hsa-miR-548ag	MIMAT0018969	YP02102770
C14	hsa-miR-4439	MIMAT0018957	YP02113184	K14	hsa-miR-4320	MIMAT0016871	YP02107197
C15	hsa-miR-3975	MIMAT0019360	YP02107044	K15	hsa-miR-4769-3p	MIMAT0019923	YP02104105
C16	hsa-miR-3925-3p	MIMAT0019228	YP02105559	K16	hsa-miR-6506-5p	MIMAT0025468	YP02100686
C17	hsa-miR-4738-3p	MIMAT0019867	YP02118789	K17	hsa-miR-4680-3p	MIMAT0019765	YP02101178
C18	hsa-miR-376b-5p	MIMAT0022923	YP02109956	K18	hsa-miR-4772-5p	MIMAT0019926	YP02115719
C19	hsa-miR-550a-3-	MIMAT0020925	YP02100573	K19	hsa-miR-584-3p	MIMAT0022708	YP02109832
C20	hsa-miR-3973	MIMAT0019358	YP02105898	K20	hsa-miR-6069	MIMAT0023694	YP02119643
C21	hsa-miR-4454	MIMAT0018976	YP02114119	K21	hsa-miR-3944-5p	MIMAT0019231	YP02110628
C22	hsa-miR-4468	MIMAT0018995	YP02119414	K22	hsa-miR-6070	MIMAT0023695	YP02116476
C23	hsa-miR-4708-3p	MIMAT0019810	YP02113933	K23	hsa-miR-4638-5p	MIMAT0019695	YP02117139
C24	hsa-miR-3913-3p	MIMAT0019225	YP02113656	K24	hsa-miR-4742-5p	MIMAT0019872	YP02103400
D01	hsa-miR-1269b	MIMAT0019059	YP02105174	L01	hsa-miR-3663-3p	MIMAT0018085	YP02119418
D02	hsa-miR-5187-5p	MIMAT0021117	YP02116675	L02	hsa-miR-4659a-3p	MIMAT0019727	YP02100136
D03	hsa-miR-4753-5p	MIMAT0019890	YP02101410	L03	hsa-miR-4699-5p	MIMAT0019794	YP02108533
D04	hsa-miR-376c-5p	MIMAT0022861	YP02105282	L04	hsa-miR-5009-5p	MIMAT0021041	YP02107327
D05	hsa-miR-3156-3p	MIMAT0019209	YP02106491	L05	hsa-miR-5591-3p	MIMAT0022302	YP02117651
D06	hsa-miR-4691-3p	MIMAT0019782	YP02102608	L06	hsa-miR-1229-3p	MIMAT0005584	YP02119607
D07	hsa-miR-5003-5p	MIMAT0021025	YP02110201	L07	hsa-miR-4524a-3p	MIMAT0019063	YP02109562
D08	hsa-miR-4513	MIMAT0019050	YP02110923	L08	hsa-miR-4798-3p	MIMAT0019975	YP02105074
D09	hsa-miR-378j	MIMAT0024612	YP02117697	L09	hsa-miR-3974	MIMAT0019359	YP02103868
D10	hsa-miR-3187-5p	MIMAT0019216	YP02110153	L10	hsa-miR-3127-3p	MIMAT0019201	YP02100657
D11	hsa-miR-4798-5p	MIMAT0019974	YP02104836	L11	hsa-miR-3178	MIMAT0015055	YP02100320
D12	hsa-miR-3622b-5p	MIMAT0018005	YP02113236	L12	hsa-miR-5584-3p	MIMAT0022284	YP02102452
D13	hsa-miR-762	MIMAT0010313	YP02119423	L13	hsa-miR-151b	MIMAT0010214	YP02110935
D14	hsa-miR-4756-5p	MIMAT0019899	YP02105920	L14	hsa-miR-3928-3p	MIMAT0018205	YP02111343
D15	hsa-miR-4529-5p	MIMAT0019236	YP02102902	L15	hsa-miR-4512	MIMAT0019049	YP02113646
D16	hsa-miR-4457	MIMAT0018979	YP02103909	L16	hsa-miR-4495	MIMAT0019030	YP02114014
D17	hsa-miR-5010-5p	MIMAT0021043	YP02106494	L17	hsa-miR-5586-3p	MIMAT0022288	YP02109152
D18	hsa-miR-548as-3p	MIMAT0022268	YP02108615	L18	hsa-miR-378i	MIMAT0019074	YP02119362
D19	hsa-miR-4420	MIMAT0018933	YP02117276	L19	hsa-miR-378h	MIMAT0018984	YP02101913
D20	hsa-miR-1275	MIMAT0005929	YP02101914	L20	hsa-miR-548ad-3p	MIMAT0018946	YP02111820
D21	hsa-miR-5190	MIMAT0021121	YP02119377	L21	hsa-miR-4724-3p	MIMAT0019842	YP02116086
D22	hsa-miR-4489	MIMAT0019023	YP02101469	L22	hsa-miR-4741	MIMAT0019871	YP02106249
D23	hsa-miR-4662b	MIMAT0019736	YP02113149	L23	hsa-miR-4662a-3p	MIMAT0019732	YP02108595
D24	hsa-miR-4725-5p	MIMAT0019843	YP02114387	L24	hsa-miR-4637	MIMAT0019694	YP02101574
E01	hsa-miR-5000-3p	MIMAT0021020	YP02106843	M01	hsa-miR-4660	MIMAT0019728	YP02106968
E02	hsa-miR-4739	MIMAT0019868	YP02110996	M02	hsa-miR-4493	MIMAT0019028	YP02110202
E03	hsa-miR-4433a-5p	MIMAT0020956	YP02114971	M03	hsa-miR-4734	MIMAT0019859	YP02108428

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
E04	hsa-miR-3927-5p	MIMAT0022970	YP02116343	M04	hsa-miR-6507-5p	MIMAT0025470	YP02101463
E05	hsa-miR-4473	MIMAT0019000	YP02106035	M05	hsa-miR-4698	MIMAT0019793	YP02110425
E06	hsa-miR-3609	MIMAT0017986	YP02107756	M06	hsa-miR-4657	MIMAT0019724	YP02109302
E07	hsa-miR-4529-3p	MIMAT0019068	YP02118987	M07	hsa-miR-6071	MIMAT0023696	YP02100434
E08	hsa-miR-4503	MIMAT0019039	YP02109937	M08	hsa-miR-4632-3p	MIMAT0019688	YP02116820
E09	hsa-miR-4509	MIMAT0019046	YP02102368	M09	hsa-miR-4664-5p	MIMAT0019737	YP02104688
E10	hsa-miR-3936	MIMAT0018351	YP02104742	M10	hsa-miR-5003-3p	MIMAT0021026	YP02103453
E11	hsa-miR-4435	MIMAT0018951	YP02115241	M11	hsa-miR-365a-5p	MIMAT0009199	YP02107446
E12	hsa-miR-4726-3p	MIMAT0019846	YP02103341	M12	hsa-miR-5579-5p	MIMAT0022269	YP02103427
E13	hsa-miR-4642	MIMAT0019702	YP02107627	M13	hsa-miR-6072	MIMAT0023697	YP02105896
E14	hsa-miR-3529-3p	MIMAT0022741	YP02106413	M14	hsa-miR-4760-5p	MIMAT0019906	YP02102642
E15	hsa-miR-4700-3p	MIMAT0019797	YP02113909	M15	hsa-miR-5585-5p	MIMAT0022285	YP02116681
E16	hsa-miR-4705	MIMAT0019805	YP02102778	M16	hsa-miR-5589-5p	MIMAT0022297	YP02117154
E17	hsa-miR-4490	MIMAT0019025	YP02104753	M17	hsa-miR-4762-3p	MIMAT0019911	YP02101891
E18	hsa-miR-3189-5p	MIMAT0019217	YP02112919	M18	hsa-miR-550b-2-5p	MIMAT0022737	YP02105365
E19	hsa-miR-5002-5p	MIMAT0021023	YP02115089	M19	hsa-miR-4694-3p	MIMAT0019787	YP02108694
E20	hsa-miR-4653-3p	MIMAT0019719	YP02116403	M20	hsa-miR-4666a-3p	MIMAT0019742	YP02102707
E21	hsa-miR-4685-5p	MIMAT0019771	YP02113728	M21	hsa-miR-5691	MIMAT0022483	YP02102410
E22	hsa-miR-3934-3p	MIMAT0022975	YP02108518	M22	hsa-miR-5197-5p	MIMAT0021130	YP02106539
E23	hsa-miR-548av-5p	MIMAT0022303	YP02114904	M23	hsa-miR-2392	MIMAT0019043	YP02104616
E24	hsa-miR-4703-5p	MIMAT0019801	YP02114283	M24	hsa-miR-548z	MIMAT0018446	YP02104422
F01	hsa-miR-4747-3p	MIMAT0019883	YP02105244	N01	hsa-miR-4684-3p	MIMAT0019770	YP02108035
F02	hsa-miR-4728-5p	MIMAT0019849	YP02105760	N02	hsa-miR-4688	MIMAT0019777	YP02106632
F03	hsa-miR-3145-5p	MIMAT0019205	YP02117024	N03	hsa-miR-4695-5p	MIMAT0019788	YP02106722
F04	hsa-miR-4706	MIMAT0019806	YP02102599	N04	hsa-miR-3188	MIMAT0015070	YP02106155
F05	hsa-miR-660-3p	MIMAT0022711	YP02112277	N05	hsa-miR-4507	MIMAT0019044	YP02119452
F06	hsa-miR-4743-3p	MIMAT0022978	YP02112978	N06	hsa-miR-4697-5p	MIMAT0019791	YP02110543
F07	hsa-miR-4494	MIMAT0019029	YP02103631	N07	hsa-miR-1178-5p	MIMAT0022940	YP02112008
F08	hsa-miR-4750-3p	MIMAT0022979	YP02109761	N08	hsa-miR-1229-5p	MIMAT0022942	YP02107908
F09	hsa-miR-495-5p	MIMAT0022924	YP02116565	N09	hsa-miR-1233-5p	MIMAT0022943	YP02105468
F10	hsa-miR-4797-5p	MIMAT0019972	YP02102015	N10	hsa-miR-1236-5p	MIMAT0022945	YP02109228
F11	hsa-miR-664b-3p	MIMAT0022272	YP02111078	N11	hsa-miR-1237-5p	MIMAT0022946	YP02119348
F12	hsa-miR-503-3p	MIMAT0022925	YP02109611	N12	hsa-miR-6073	MIMAT0023698	YP02114330
F13	hsa-miR-4999-3p	MIMAT0021018	YP02116831	N13	hsa-miR-6074	MIMAT0023699	YP02106450
F14	hsa-miR-3692-5p	MIMAT0018121	YP02111748	N14	hsa-miR-6075	MIMAT0023700	YP02118362
F15	hsa-miR-5585-3p	MIMAT0022286	YP02115478	N15	hsa-miR-6076	MIMAT0023701	YP02111947
F16	hsa-miR-548au-	MIMAT0022292	YP02117096	N16	hsa-miR-6077	MIMAT0023702	YP02103758
F17	hsa-miR-4694-5p	MIMAT0019786	YP02101085	N17	hsa-miR-6078	MIMAT0023703	YP02104937
F18	hsa-miR-3129-3p	MIMAT0019202	YP02106263	N18	hsa-miR-6079	MIMAT0023704	YP02105595
F19	hsa-miR-4756-3p	MIMAT0019900	YP02109527	N19	hsa-miR-6080	MIMAT0023705	YP02117171
F20	hsa-miR-5685	MIMAT0022475	YP02101480	N20	hsa-miR-6508-3p	MIMAT0025473	YP02112150
F21	hsa-miR-3620-3p	MIMAT0018001	YP02109143	N21	hsa-miR-6509-3p	MIMAT0025475	YP02118144
F22	hsa-miR-4713-3p	MIMAT0019821	YP02117930	N22	hsa-miR-6086	MIMAT0023711	YP02107036
F23	hsa-miR-4778-3p	MIMAT0019937	YP02113533	N23	hsa-miR-6124	MIMAT0024597	YP02117732
F24	hsa-miR-4685-3p	MIMAT0019772	YP02115818	N24	hsa-miR-6126	MIMAT0024599	YP02119651
G01	hsa-miR-122b-3p	MIMAT0019877	YP02119380	O01	hsa-miR-6127	MIMAT0024610	YP02114505
G02	hsa-miR-4673	MIMAT0019755	YP02102822	O02	hsa-miR-6128	MIMAT0024611	YP02119353
G03	hsa-miR-4768-3p	MIMAT0019921	YP02110163	O03	hsa-miR-6129	MIMAT0024613	YP02113680
G04	hsa-miR-3064-3p	MIMAT0019865	YP02110730	O04	hsa-miR-6130	MIMAT0024614	YP02109513
G05	hsa-miR-5089-3p	MIMAT0022984	YP02117169	O05	hsa-miR-6131	MIMAT0024615	YP02102492

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
G06	hsa-miR-548a-5p	MIMAT0022265	YP02106657	O06	hsa-miR-6132	MIMAT0024616	YP02111067
G07	hsa-miR-548a-5p	MIMAT0021037	YP02100835	O07	hsa-miR-6133	MIMAT0024617	YP02111821
G08	hsa-miR-4418	MIMAT0018930	YP02117035	O08	hsa-miR-6134	MIMAT0024618	YP02112380
G09	hsa-miR-5692a	MIMAT0022484	YP02119363	O09	hsa-miR-6165	MIMAT0024782	YP02117564
G10	hsa-miR-5692c	MIMAT0022476	YP02106294	O10	hsa-miR-6499-3p	MIMAT0025451	YP02104793
G11	hsa-miR-2277-5p	MIMAT0017352	YP02104427	O11	hsa-miR-6499-5p	MIMAT0025450	YP02117976
G12	hsa-miR-4699-3p	MIMAT0019795	YP02119478	O12	hsa-miR-6500-3p	MIMAT0025455	YP02101176
G13	hsa-miR-4434	MIMAT0018950	YP02115659	O13	hsa-miR-6500-5p	MIMAT0025454	YP02108457
G14	hsa-miR-4755-3p	MIMAT0019896	YP02105715	O14	hsa-miR-6501-3p	MIMAT0025459	YP02108442
G15	hsa-miR-4645-5p	MIMAT0019705	YP02108522	O15	hsa-miR-6501-5p	MIMAT0025458	YP02108075
G16	hsa-miR-548a-5p	MIMAT0025453	YP02110459	O16	hsa-miR-6502-3p	MIMAT0025461	YP02116807
G17	hsa-miR-5705	MIMAT0022499	YP02113091	O17	hsa-miR-6503-3p	MIMAT0025463	YP02119048
G18	hsa-miR-4428	MIMAT0018943	YP02119279	O18	hsa-miR-6503-5p	MIMAT0025462	YP02106551
G19	hsa-miR-5010-3p	MIMAT0021044	YP02109984	O19	hsa-miR-6504-3p	MIMAT0025465	YP02119354
G20	hsa-miR-3614-5p	MIMAT0017992	YP02108616	O20	hsa-miR-6504-5p	MIMAT0025464	YP02116070
G21	hsa-miR-4438	MIMAT0018956	YP02119615	O21	hsa-miR-6509-5p	MIMAT0025474	YP02102580
G22	hsa-miR-4793-5p	MIMAT0019965	YP02111009	O22	hsa-miR-6510-5p	MIMAT0025476	YP02101889
G23	hsa-miR-2277-3p	MIMAT0011777	YP02110810	O23	hsa-miR-6511b-5p	MIMAT0025847	YP02108808
G24	hsa-miR-5001-3p	MIMAT0021022	YP02101303	O24	hsa-miR-6512-5p	MIMAT0025480	YP02114608
H01	hsa-miR-4793-3p	MIMAT0019966	YP02113347	P01	hsa-miR-6513-3p	MIMAT0025483	YP02115098
H02	hsa-miR-4709-5p	MIMAT0019811	YP02105739	P02	hsa-miR-6513-5p	MIMAT0025482	YP02117963
H03	hsa-miR-4738-5p	MIMAT0019866	YP02112441	P03	hsa-miR-6514-3p	MIMAT0025485	YP02117257
H04	hsa-miR-4753-3p	MIMAT0019891	YP02108145	P04	hsa-miR-6514-5p	MIMAT0025484	YP02109317
H05	hsa-miR-4650-3p	MIMAT0019714	YP02116929	P05	hsa-miR-6515-3p	MIMAT0025487	YP02118952
H06	hsa-miR-4783-5p	MIMAT0019946	YP02108366	P06	hsa-miR-6515-5p	MIMAT0025486	YP02114457
H07	hsa-miR-371b-3p	MIMAT0019893	YP02108243	P07	hsa-miR-6715a-3p	MIMAT0025841	YP02109588
H08	hsa-miR-642a-3p	MIMAT0020924	YP02112940	P08	hsa-miR-6715b-5p	MIMAT0025842	YP02102216
H09	hsa-miR-3140-5p	MIMAT0019204	YP02105747	P09	hsa-miR-6716-3p	MIMAT0025845	YP02110812
H10	hsa-miR-5582-3p	MIMAT0022280	YP02111017	P10	hsa-miR-6716-5p	MIMAT0025844	YP02115284
H11	hsa-miR-548a-5p	MIMAT0018989	YP02115184	P11	hsa-miR-6717-5p	MIMAT0025846	YP02102075
H12	hsa-miR-3663-5p	MIMAT0018084	YP02109334	P12	BlankSpot	N/A	YP02119291
H13	hsa-miR-4519	MIMAT0019056	YP02113093	P13	SNORD44 (hsa)	N/A	YP00203902
H14	hsa-miR-3197	MIMAT0015082	YP02118682	P14	SNORD38B (hsa)	N/A	YP00203901
H15	hsa-miR-4425	MIMAT0018940	YP02116995	P15	SNORD49A (hsa)	N/A	YP00203904
H16	hsa-miR-4794	MIMAT0019967	YP02114578	P16	U6 snRNA (v2)	N/A	YP02119464
H17	hsa-miR-4763-5p	MIMAT0019912	YP02110500	P17	cel-miR-39-3p	MIMAT0000010	YP00203952
H18	hsa-miR-4801	MIMAT0019980	YP02101187	P18	UniSp2	N/A	YP00203950
H19	hsa-miR-4757-5p	MIMAT0019901	YP02117442	P19	UniSp4	N/A	YP00203953
H20	hsa-miR-548a-5p	MIMAT0025452	YP02100340	P20	UniSp5	N/A	YP00203955
H21	hsa-miR-659-5p	MIMAT0022710	YP02108749	P21	UniSp6	N/A	YP00203954
H22	hsa-miR-5695	MIMAT0022488	YP02115028	P22	UniSP3	N/A	YP02119288
H23	hsa-miR-378d	MIMAT0018926	YP02100497	P23	UniSP3	N/A	YP02119288
H24	hsa-miR-1271-3p	MIMAT0022712	YP02112839	P24	UniSP3	N/A	YP02119288

miRNA Table: Plate 6

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
A01	hsa-miR-6718-5p	MIMAT0025849	YP02113259	I01	hsa-miR-6796-5p	MIMAT0027492	YP02112426
A02	hsa-miR-6719-3p	MIMAT0025850	YP02111364	I02	hsa-miR-6797-5p	MIMAT0027494	YP02118253
A03	hsa-miR-6720-3p	MIMAT0025851	YP02114592	I03	hsa-miR-6798-5p	MIMAT0027496	YP02105728
A04	hsa-miR-6721-5p	MIMAT0025852	YP02110957	I04	hsa-miR-6799-3p	MIMAT0027499	YP02101545
A05	hsa-miR-6722-3p	MIMAT0025854	YP02105759	I05	hsa-miR-6799-5p	MIMAT0027498	YP02100691
A06	hsa-miR-6722-5p	MIMAT0025853	YP02113668	I06	hsa-miR-6800-3p	MIMAT0027501	YP02107834
A07	hsa-miR-758-5p	MIMAT0022929	YP02116137	I07	hsa-miR-6800-5p	MIMAT0027500	YP02112192
A08	hsa-miR-892c-3p	MIMAT0025858	YP02117971	I08	hsa-miR-6801-3p	MIMAT0027503	YP02101519
A09	hsa-miR-892c-5p	MIMAT0025857	YP02114947	I09	hsa-miR-6801-5p	MIMAT0027502	YP02115784
A10	hsa-miR-937-5p	MIMAT0022938	YP02117509	I10	hsa-miR-6802-3p	MIMAT0027505	YP02112736
A11	hsa-miR-98-3p	MIMAT0022842	YP02106353	I11	hsa-miR-6802-5p	MIMAT0027504	YP02119388
A12	hsa-let-7c-3p	MIMAT0026472	YP00205202	I12	hsa-miR-6803-3p	MIMAT0027507	YP02115732
A13	hsa-miR-1180-5p	MIMAT0026735	YP02100476	I13	hsa-miR-6804-3p	MIMAT0027509	YP02108252
A14	hsa-miR-1199-3p	MIMAT0031120	YP02102766	I14	hsa-miR-6804-5p	MIMAT0027508	YP02119084
A15	hsa-miR-1250-3p	MIMAT0026740	YP02110614	I15	hsa-miR-6805-5p	MIMAT0027510	YP02112777
A16	hsa-miR-1251-3p	MIMAT0026741	YP02119283	I16	hsa-miR-6806-3p	MIMAT0027513	YP02110431
A17	hsa-miR-1252-3p	MIMAT0026744	YP02114253	I17	hsa-miR-6806-5p	MIMAT0027512	YP02115030
A18	hsa-miR-1266-3p	MIMAT0026742	YP02109017	I18	hsa-miR-6807-3p	MIMAT0027515	YP02105778
A19	hsa-miR-1273h-3p	MIMAT0030416	YP02111365	I19	hsa-miR-6807-5p	MIMAT0027514	YP02119459
A20	hsa-miR-128-1-5p	MIMAT0026477	YP02117210	I20	hsa-miR-6808-3p	MIMAT0027517	YP02118534
A21	hsa-miR-128-2-5p	MIMAT0031095	YP02110351	I21	hsa-miR-6808-5p	MIMAT0027516	YP02106511
A22	hsa-miR-1288-5p	MIMAT0026743	YP02110095	I22	hsa-miR-6809-3p	MIMAT0027519	YP02111411
A23	hsa-miR-1296-3p	MIMAT0026637	YP02102262	I23	hsa-miR-6809-5p	MIMAT0027518	YP02101251
A24	hsa-miR-1298-3p	MIMAT0026641	YP02108006	I24	hsa-miR-6810-5p	MIMAT0027520	YP02101678
B01	hsa-miR-1301-5p	MIMAT0026639	YP02118512	J01	hsa-miR-6811-3p	MIMAT0027523	YP02111254
B02	hsa-miR-133a-5p	MIMAT0026478	YP00205501	J02	hsa-miR-6811-5p	MIMAT0027522	YP02114988
B03	hsa-miR-1343-5p	MIMAT0027038	YP02100503	J03	hsa-miR-6812-5p	MIMAT0027524	YP02110112
B04	hsa-miR-134-3p	MIMAT0026481	YP02114563	J04	hsa-miR-6813-5p	MIMAT0027526	YP02113477
B05	hsa-miR-1468-3p	MIMAT0026638	YP02110580	J05	hsa-miR-6814-5p	MIMAT0027528	YP02105238
B06	hsa-miR-152-5p	MIMAT0026479	YP02108962	J06	hsa-miR-6815-3p	MIMAT0027531	YP02102659
B07	hsa-miR-1537-5p	MIMAT0026765	YP02116182	J07	hsa-miR-6815-5p	MIMAT0027530	YP02106133
B08	hsa-miR-181d-3p	MIMAT0026608	YP02105350	J08	hsa-miR-6816-3p	MIMAT0027533	YP02109792
B09	hsa-miR-190a-3p	MIMAT0026482	YP02100791	J09	hsa-miR-6816-5p	MIMAT0027532	YP02118525
B10	hsa-miR-1910-3p	MIMAT0026917	YP02110313	J10	hsa-miR-6817-5p	MIMAT0027534	YP02119229
B11	hsa-miR-208a-5p	MIMAT0026474	YP00205398	J11	hsa-miR-6818-3p	MIMAT0027537	YP02116125
B12	hsa-miR-208b-5p	MIMAT0026722	YP02110745	J12	hsa-miR-6818-5p	MIMAT0027536	YP02111698
B13	hsa-miR-210-5p	MIMAT0026475	YP02104321	J13	hsa-miR-6819-5p	MIMAT0027538	YP02110218
B14	hsa-miR-215-3p	MIMAT0026476	YP02114252	J14	hsa-miR-6820-3p	MIMAT0027541	YP02101831
B15	hsa-miR-216b-3p	MIMAT0026721	YP02117492	J15	hsa-miR-6820-5p	MIMAT0027540	YP02116855
B16	hsa-miR-2276-5p	MIMAT0026921	YP02108484	J16	hsa-miR-6821-3p	MIMAT0027543	YP02115191
B17	hsa-miR-3151-3p	MIMAT0027026	YP02118703	J17	hsa-miR-6821-5p	MIMAT0027542	YP02108669
B18	hsa-miR-3192-3p	MIMAT0027027	YP02107026	J18	hsa-miR-6822-3p	MIMAT0027545	YP02107987
B19	hsa-miR-328-5p	MIMAT0026486	YP02113205	J19	hsa-miR-6822-5p	MIMAT0027544	YP02112115
B20	hsa-miR-329-5p	MIMAT0026555	YP02107092	J20	hsa-miR-6823-3p	MIMAT0027547	YP02118078
B21	hsa-miR-370-5p	MIMAT0026483	YP02107400	J21	hsa-miR-6823-5p	MIMAT0027546	YP02114286
B22	hsa-miR-372-5p	MIMAT0026484	YP02115668	J22	hsa-miR-6824-3p	MIMAT0027549	YP02111157
B23	hsa-miR-383-3p	MIMAT0026485	YP02108841	J23	hsa-miR-6824-5p	MIMAT0027548	YP02102777
B24	hsa-miR-3912-5p	MIMAT0027036	YP02119128	J24	hsa-miR-6825-3p	MIMAT0027551	YP02115710
C01	hsa-miR-3928-5p	MIMAT0027037	YP02116032	K01	hsa-miR-6825-5p	MIMAT0027550	YP02104934

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
C02	hsa-miR-410-5p	MIMAT0026558	YP02111083	K02	hsa-miR-6826-3p	MIMAT0027553	YP02109244
C03	hsa-miR-412-5p	MIMAT0026557	YP02115902	K03	hsa-miR-6826-5p	MIMAT0027552	YP02102567
C04	hsa-miR-433-5p	MIMAT0026554	YP00205430	K04	hsa-miR-6827-3p	MIMAT0027555	YP02118217
C05	hsa-miR-4433b-3p	MIMAT0030414	YP02115942	K05	hsa-miR-6827-5p	MIMAT0027554	YP02114959
C06	hsa-miR-4433b-5p	MIMAT0030413	YP02118899	K06	hsa-miR-6828-3p	MIMAT0027557	YP02117525
C07	hsa-miR-487a-5p	MIMAT0026559	YP02112873	K07	hsa-miR-6828-5p	MIMAT0027556	YP02101120
C08	hsa-miR-487b-5p	MIMAT0026614	YP02115384	K08	hsa-miR-6829-3p	MIMAT0027559	YP02109974
C09	hsa-miR-489-5p	MIMAT0026605	YP02107863	K09	hsa-miR-6829-5p	MIMAT0027558	YP02107920
C10	hsa-miR-494-5p	MIMAT0026607	YP02111151	K10	hsa-miR-6830-3p	MIMAT0027561	YP02108622
C11	hsa-miR-500b-3p	MIMAT0027032	YP02105743	K11	hsa-miR-6830-5p	MIMAT0027560	YP02116880
C12	hsa-miR-6889-3p	MIMAT0027679	YP02104407	K12	hsa-miR-6831-3p	MIMAT0027563	YP02109795
C13	hsa-miR-511-3p	MIMAT0026606	YP02108582	K13	hsa-miR-6831-5p	MIMAT0027562	YP02109034
C14	hsa-miR-513b-3p	MIMAT0026749	YP02102666	K14	hsa-miR-6832-3p	MIMAT0027565	YP02116095
C15	hsa-miR-5189-3p	MIMAT0027088	YP02106506	K15	hsa-miR-6832-5p	MIMAT0027564	YP02112352
C16	hsa-miR-520g-5p	MIMAT0026611	YP02118483	K16	hsa-miR-6833-3p	MIMAT0027567	YP02119433
C17	hsa-miR-548ba	MIMAT0031175	YP02113326	K17	hsa-miR-6833-5p	MIMAT0027566	YP02114308
C18	hsa-miR-548e-5p	MIMAT0026736	YP02117828	K18	hsa-miR-6834-3p	MIMAT0027569	YP02107263
C19	hsa-miR-548f-5p	MIMAT0026739	YP02112172	K19	hsa-miR-6834-5p	MIMAT0027568	YP02100744
C20	hsa-miR-548j-3p	MIMAT0026737	YP02119385	K20	hsa-miR-6835-3p	MIMAT0027571	YP02116069
C21	hsa-miR-5699-5p	MIMAT0027103	YP02113295	K21	hsa-miR-6835-5p	MIMAT0027570	YP02119658
C22	hsa-miR-579-5p	MIMAT0026616	YP02117190	K22	hsa-miR-6836-3p	MIMAT0027575	YP02119660
C23	hsa-miR-580-5p	MIMAT0026617	YP02104524	K23	hsa-miR-6836-5p	MIMAT0027574	YP02105775
C24	hsa-miR-585-5p	MIMAT0026618	YP02114778	K24	hsa-miR-6837-3p	MIMAT0027577	YP02111870
D01	hsa-miR-598-5p	MIMAT0026620	YP02113869	L01	hsa-miR-6837-5p	MIMAT0027576	YP02118732
D02	hsa-miR-605-3p	MIMAT0026621	YP02105161	L02	hsa-miR-6838-3p	MIMAT0027579	YP02108844
D03	hsa-miR-619-5p	MIMAT0026622	YP02118098	L03	hsa-miR-6838-5p	MIMAT0027578	YP02109467
D04	hsa-miR-651-3p	MIMAT0026624	YP02112173	L04	hsa-miR-6839-3p	MIMAT0027581	YP02100700
D05	hsa-miR-6516-3p	MIMAT0030418	YP02101857	L05	hsa-miR-6839-5p	MIMAT0027580	YP02118949
D06	hsa-miR-6516-5p	MIMAT0030417	YP02117973	L06	hsa-miR-6840-3p	MIMAT0027583	YP02101589
D07	hsa-miR-653-3p	MIMAT0026625	YP02102836	L07	hsa-miR-6840-5p	MIMAT0027582	YP02110643
D08	hsa-miR-655-5p	MIMAT0026626	YP02104256	L08	hsa-miR-6841-3p	MIMAT0027585	YP02118210
D09	hsa-miR-656-5p	MIMAT0026627	YP02101477	L09	hsa-miR-6841-5p	MIMAT0027584	YP02107704
D10	hsa-miR-668-5p	MIMAT0026636	YP02108164	L10	hsa-miR-6842-3p	MIMAT0027587	YP02115310
D11	hsa-miR-670-3p	MIMAT0026640	YP02118565	L11	hsa-miR-6842-5p	MIMAT0027586	YP02118059
D12	hsa-miR-6889-5p	MIMAT0027678	YP02102772	L12	hsa-miR-6843-3p	MIMAT0027588	YP02105311
D13	hsa-miR-6726-3p	MIMAT0027354	YP02105180	L13	hsa-miR-6890-3p	MIMAT0027681	YP02102771
D14	hsa-miR-6726-5p	MIMAT0027353	YP02117526	L14	hsa-miR-6845-5p	MIMAT0027590	YP02110607
D15	hsa-miR-6728-3p	MIMAT0027358	YP02114170	L15	hsa-miR-6846-3p	MIMAT0027593	YP02107207
D16	hsa-miR-6730-5p	MIMAT0027361	YP02108676	L16	hsa-miR-6846-5p	MIMAT0027592	YP02106013
D17	hsa-miR-6731-3p	MIMAT0027364	YP02119123	L17	hsa-miR-6847-3p	MIMAT0027595	YP02103437
D18	hsa-miR-6732-3p	MIMAT0027366	YP02114748	L18	hsa-miR-6847-5p	MIMAT0027594	YP02113204
D19	hsa-miR-6732-5p	MIMAT0027365	YP02119428	L19	hsa-miR-6848-3p	MIMAT0027597	YP02101654
D20	hsa-miR-6733-3p	MIMAT0027368	YP02115361	L20	hsa-miR-6848-5p	MIMAT0027596	YP02107929
D21	hsa-miR-6733-5p	MIMAT0027367	YP02110920	L21	hsa-miR-6849-3p	MIMAT0027599	YP02109972
D22	hsa-miR-6734-3p	MIMAT0027370	YP02108135	L22	hsa-miR-6849-5p	MIMAT0027598	YP02116737
D23	hsa-miR-6734-5p	MIMAT0027369	YP02117038	L23	hsa-miR-6850-3p	MIMAT0027601	YP02113302
D24	hsa-miR-6735-5p	MIMAT0027371	YP02109568	L24	hsa-miR-6851-3p	MIMAT0027603	YP02117632
E01	hsa-miR-6736-3p	MIMAT0027374	YP02117652	M01	hsa-miR-6851-5p	MIMAT0027602	YP02109092
E02	hsa-miR-6736-5p	MIMAT0027373	YP02107530	M02	hsa-miR-6852-3p	MIMAT0027605	YP02116751
E03	hsa-miR-6737-5p	MIMAT0027375	YP02103653	M03	hsa-miR-6852-5p	MIMAT0027604	YP02104799

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
E04	hsa-miR-6738-3p	MIMAT0027378	YP02101750	M04	hsa-miR-6853-3p	MIMAT0027607	YP02118391
E05	hsa-miR-6738-5p	MIMAT0027377	YP02117886	M05	hsa-miR-6853-5p	MIMAT0027606	YP02105449
E06	hsa-miR-6739-3p	MIMAT0027380	YP02102392	M06	hsa-miR-6854-3p	MIMAT0027609	YP02108291
E07	hsa-miR-6739-5p	MIMAT0027379	YP02104637	M07	hsa-miR-6854-5p	MIMAT0027608	YP02115043
E08	hsa-miR-6740-3p	MIMAT0027382	YP02115616	M08	hsa-miR-6855-3p	MIMAT0027611	YP02112157
E09	hsa-miR-6740-5p	MIMAT0027381	YP02111295	M09	hsa-miR-6855-5p	MIMAT0027610	YP02116164
E10	hsa-miR-6741-3p	MIMAT0027384	YP02119233	M10	hsa-miR-6856-3p	MIMAT0027613	YP02113125
E11	hsa-miR-6741-5p	MIMAT0027383	YP02108759	M11	hsa-miR-6857-3p	MIMAT0027615	YP02116660
E12	hsa-miR-6742-3p	MIMAT0027386	YP02110516	M12	hsa-miR-6857-5p	MIMAT0027614	YP02111562
E13	hsa-miR-6742-5p	MIMAT0027385	YP02119213	M13	hsa-miR-6858-5p	MIMAT0027616	YP02112076
E14	hsa-miR-6743-5p	MIMAT0027387	YP02117389	M14	hsa-miR-6859-5p	MIMAT0027618	YP02103332
E15	hsa-miR-6744-3p	MIMAT0027390	YP02114950	M15	hsa-miR-6860	MIMAT0027622	YP02105127
E16	hsa-miR-6744-5p	MIMAT0027389	YP02104363	M16	hsa-miR-6861-5p	MIMAT0027623	YP02119257
E17	hsa-miR-6746-3p	MIMAT0027393	YP02116628	M17	hsa-miR-6862-3p	MIMAT0027626	YP02116509
E18	hsa-miR-6747-3p	MIMAT0027395	YP02113949	M18	hsa-miR-6862-5p	MIMAT0027625	YP02116492
E19	hsa-miR-6748-3p	MIMAT0027397	YP02106726	M19	hsa-miR-6863	MIMAT0027627	YP02112191
E20	hsa-miR-6748-5p	MIMAT0027396	YP02103103	M20	hsa-miR-6864-3p	MIMAT0027629	YP02105403
E21	hsa-miR-6749-5p	MIMAT0027398	YP02108544	M21	hsa-miR-6864-5p	MIMAT0027628	YP02106668
E22	hsa-miR-6750-5p	MIMAT0027400	YP02111905	M22	hsa-miR-6865-3p	MIMAT0027631	YP02115783
E23	hsa-miR-6751-3p	MIMAT0027403	YP02106215	M23	hsa-miR-6865-5p	MIMAT0027630	YP02118138
E24	hsa-miR-6751-5p	MIMAT0027402	YP02118769	M24	hsa-miR-6866-3p	MIMAT0027633	YP02118145
F01	hsa-miR-6752-5p	MIMAT0027404	YP02110467	N01	hsa-miR-6866-5p	MIMAT0027632	YP02111715
F02	hsa-miR-6753-3p	MIMAT0027407	YP02118515	N02	hsa-miR-6867-3p	MIMAT0027635	YP02111051
F03	hsa-miR-6753-5p	MIMAT0027406	YP02105614	N03	hsa-miR-6867-5p	MIMAT0027634	YP02116720
F04	hsa-miR-6754-5p	MIMAT0027408	YP02116950	N04	hsa-miR-6868-3p	MIMAT0027637	YP02102933
F05	hsa-miR-6755-3p	MIMAT0027411	YP02117866	N05	hsa-miR-6868-5p	MIMAT0027636	YP02111215
F06	hsa-miR-6755-5p	MIMAT0027410	YP02109049	N06	hsa-miR-6869-3p	MIMAT0027639	YP02103497
F07	hsa-miR-6756-5p	MIMAT0027412	YP02119429	N07	hsa-miR-6870-3p	MIMAT0027641	YP02117175
F08	hsa-miR-6757-3p	MIMAT0027415	YP02107562	N08	hsa-miR-6871-3p	MIMAT0027643	YP02103093
F09	hsa-miR-6757-5p	MIMAT0027414	YP02109195	N09	hsa-miR-6871-5p	MIMAT0027642	YP02113232
F10	hsa-miR-6758-3p	MIMAT0027417	YP02109503	N10	hsa-miR-6872-5p	MIMAT0027644	YP02101577
F11	hsa-miR-6759-3p	MIMAT0027419	YP02115722	N11	hsa-miR-6873-3p	MIMAT0027647	YP02102681
F12	hsa-miR-6759-5p	MIMAT0027418	YP02107923	N12	hsa-miR-6873-5p	MIMAT0027646	YP02103510
F13	hsa-miR-6761-3p	MIMAT0027423	YP02106554	N13	hsa-miR-6874-3p	MIMAT0027649	YP02117460
F14	hsa-miR-6762-3p	MIMAT0027425	YP02106317	N14	hsa-miR-6874-5p	MIMAT0027648	YP02109241
F15	hsa-miR-6762-5p	MIMAT0027424	YP02105421	N15	hsa-miR-6875-3p	MIMAT0027651	YP02108527
F16	hsa-miR-6763-5p	MIMAT0027426	YP02118635	N16	hsa-miR-6875-5p	MIMAT0027650	YP02116574
F17	hsa-miR-6764-3p	MIMAT0027429	YP02117884	N17	hsa-miR-6876-3p	MIMAT0027653	YP02115467
F18	hsa-miR-6764-5p	MIMAT0027428	YP02109894	N18	hsa-miR-6876-5p	MIMAT0027652	YP02110823
F19	hsa-miR-6765-5p	MIMAT0027430	YP02114389	N19	hsa-miR-6877-5p	MIMAT0027654	YP02115155
F20	hsa-miR-6766-5p	MIMAT0027432	YP02102746	N20	hsa-miR-6878-3p	MIMAT0027657	YP02110113
F21	hsa-miR-6767-3p	MIMAT0027435	YP02109520	N21	hsa-miR-6878-5p	MIMAT0027656	YP02113438
F22	hsa-miR-6767-5p	MIMAT0027434	YP02118572	N22	hsa-miR-6879-5p	MIMAT0027658	YP02114612
F23	hsa-miR-6768-3p	MIMAT0027437	YP02102014	N23	hsa-miR-6880-5p	MIMAT0027660	YP02116477
F24	hsa-miR-6768-5p	MIMAT0027436	YP02110713	N24	hsa-miR-6881-3p	MIMAT0027663	YP02101083
G01	hsa-miR-6769a-5p	MIMAT0027438	YP02106693	O01	hsa-miR-6881-5p	MIMAT0027662	YP02109019
G02	hsa-miR-6769b-3p	MIMAT0027621	YP02101670	O02	hsa-miR-6882-3p	MIMAT0027665	YP02108189
G03	hsa-miR-6769b-5p	MIMAT0027620	YP02113037	O03	hsa-miR-6882-5p	MIMAT0027664	YP02106235
G04	hsa-miR-6770-3p	MIMAT0027441	YP02111845	O04	hsa-miR-6883-3p	MIMAT0027667	YP02117395
G05	hsa-miR-6770-5p	MIMAT0027440	YP02103217	O05	hsa-miR-6883-5p	MIMAT0027666	YP02118150

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
G06	hsa-miR-6771-3p	MIMAT0027443	YP02119355	O06	hsa-miR-6884-3p	MIMAT0027669	YP02107360
G07	hsa-miR-6771-5p	MIMAT0027442	YP02105212	O07	hsa-miR-6884-5p	MIMAT0027668	YP02118648
G08	hsa-miR-6772-3p	MIMAT0027445	YP02103982	O08	hsa-miR-6885-3p	MIMAT0027671	YP02101345
G09	hsa-miR-6772-5p	MIMAT0027444	YP02118541	O09	hsa-miR-6885-5p	MIMAT0027670	YP02119482
G10	hsa-miR-6773-3p	MIMAT0027447	YP02106901	O10	hsa-miR-6886-5p	MIMAT0027672	YP02105593
G11	hsa-miR-6773-5p	MIMAT0027446	YP02116984	O11	hsa-miR-6887-3p	MIMAT0027675	YP02116052
G12	hsa-miR-6774-3p	MIMAT0027449	YP02106422	O12	hsa-miR-6887-5p	MIMAT0027674	YP02113868
G13	hsa-miR-6774-5p	MIMAT0027448	YP02107421	O13	hsa-miR-6888-3p	MIMAT0027677	YP02111557
G14	hsa-miR-6775-3p	MIMAT0027451	YP02118761	O14	hsa-miR-6888-5p	MIMAT0027676	YP02104718
G15	hsa-miR-6775-5p	MIMAT0027450	YP02119417	O15	hsa-miR-6890-5p	MIMAT0027680	YP02110989
G16	hsa-miR-6776-3p	MIMAT0027453	YP02114654	O16	hsa-miR-6891-3p	MIMAT0027683	YP02106363
G17	hsa-miR-6776-5p	MIMAT0027452	YP02116812	O17	hsa-miR-6891-5p	MIMAT0027682	YP02102735
G18	hsa-miR-6777-5p	MIMAT0027454	YP02101311	O18	hsa-miR-6892-3p	MIMAT0027685	YP02114726
G19	hsa-miR-6778-3p	MIMAT0027457	YP02105003	O19	hsa-miR-6892-5p	MIMAT0027684	YP02113175
G20	hsa-miR-6778-5p	MIMAT0027456	YP02101666	O20	hsa-miR-6893-3p	MIMAT0027687	YP02119432
G21	hsa-miR-6779-3p	MIMAT0027459	YP02114622	O21	hsa-miR-6893-5p	MIMAT0027686	YP02116619
G22	hsa-miR-6779-5p	MIMAT0027458	YP02117843	O22	hsa-miR-7107-3p	MIMAT0028112	YP02101788
G23	hsa-miR-6780a-3p	MIMAT0027461	YP02117483	O23	hsa-miR-7108-3p	MIMAT0028114	YP02101501
G24	hsa-miR-6780a-5p	MIMAT0027460	YP02117364	O24	hsa-miR-7109-5p	MIMAT0028115	YP02111712
H01	hsa-miR-6780b-3p	MIMAT0027573	YP02100688	P01	hsa-miR-7110-3p	MIMAT0028118	YP02114244
H02	hsa-miR-6781-5p	MIMAT0027462	YP02116752	P02	hsa-miR-7111-3p	MIMAT0028120	YP02110682
H03	hsa-miR-6782-3p	MIMAT0027465	YP02104981	P03	hsa-miR-7112-3p	MIMAT0028122	YP02106857
H04	hsa-miR-6783-3p	MIMAT0027467	YP02111036	P04	hsa-miR-7112-5p	MIMAT0028121	YP02109316
H05	hsa-miR-6783-5p	MIMAT0027466	YP02113087	P05	hsa-miR-7113-3p	MIMAT0028124	YP02106063
H06	hsa-miR-6785-3p	MIMAT0027471	YP02118992	P06	hsa-miR-7113-5p	MIMAT0028123	YP02107502
H07	hsa-miR-6785-5p	MIMAT0027470	YP02103050	P07	hsa-miR-7114-3p	MIMAT0028126	YP02109524
H08	hsa-miR-6786-3p	MIMAT0027473	YP02115225	P08	hsa-miR-7114-5p	MIMAT0028125	YP02115348
H09	hsa-miR-6786-5p	MIMAT0027472	YP02107023	P09	hsa-miR-7150	MIMAT0028211	YP02109509
H10	hsa-miR-6787-3p	MIMAT0027475	YP02104871	P10	hsa-miR-7151-3p	MIMAT0028213	YP02108098
H11	hsa-miR-6787-5p	MIMAT0027474	YP02119387	P11	hsa-miR-7151-5p	MIMAT0028212	YP02113505
H12	hsa-miR-6788-3p	MIMAT0027477	YP02106685	P12	BlankSpot	N/A	YP02119291
H13	hsa-miR-6788-5p	MIMAT0027476	YP02107621	P13	SNORD44 (hsa)	N/A	YP00203902
H14	hsa-miR-6789-5p	MIMAT0027478	YP02100084	P14	SNORD38B (hsa)	N/A	YP00203901
H15	hsa-miR-6790-3p	MIMAT0027481	YP02105201	P15	SNORD49A (hsa)	N/A	YP00203904
H16	hsa-miR-6790-5p	MIMAT0027480	YP02105097	P16	U6 snRNA (v2)	N/A	YP02119464
H17	hsa-miR-6791-3p	MIMAT0027483	YP02113917	P17	cel-miR-39-3p	MIMAT0000010	YP00203952
H18	hsa-miR-6792-3p	MIMAT0027485	YP02115808	P18	UniSp2	N/A	YP00203950
H19	hsa-miR-6792-5p	MIMAT0027484	YP02117927	P19	UniSp4	N/A	YP00203953
H20	hsa-miR-6793-5p	MIMAT0027486	YP02114368	P20	UniSp5	N/A	YP00203955
H21	hsa-miR-6794-3p	MIMAT0027489	YP02102949	P21	UniSp6	N/A	YP00203954
H22	hsa-miR-6794-5p	MIMAT0027488	YP02116385	P22	UniSP3	N/A	YP02119288
H23	hsa-miR-6795-3p	MIMAT0027491	YP02119119	P23	UniSP3	N/A	YP02119288
H24	hsa-miR-6795-5p	MIMAT0027490	YP02103188	P24	UniSP3	N/A	YP02119288

miRNA Table: Plate 7

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
A01	hsa-miR-7152-3p	MIMAT0028215	YP02119389	I01	BlankSpot	N/A	YP02119291
A02	hsa-miR-7153-3p	MIMAT0028217	YP02113293	I02	BlankSpot	N/A	YP02119291
A03	hsa-miR-7153-5p	MIMAT0028216	YP02115850	I03	BlankSpot	N/A	YP02119291
A04	hsa-miR-7154-3p	MIMAT0028219	YP02113308	I04	BlankSpot	N/A	YP02119291
A05	hsa-miR-7154-5p	MIMAT0028218	YP02104853	I05	BlankSpot	N/A	YP02119291
A06	hsa-miR-7155-3p	MIMAT0028221	YP02107402	I06	BlankSpot	N/A	YP02119291
A07	hsa-miR-7155-5p	MIMAT0028220	YP02106993	I07	BlankSpot	N/A	YP02119291
A08	hsa-miR-7156-3p	MIMAT0028223	YP02111003	I08	BlankSpot	N/A	YP02119291
A09	hsa-miR-7156-5p	MIMAT0028222	YP02108860	I09	BlankSpot	N/A	YP02119291
A10	hsa-miR-7157-3p	MIMAT0028225	YP02109750	I10	BlankSpot	N/A	YP02119291
A11	hsa-miR-7157-5p	MIMAT0028224	YP02101903	I11	BlankSpot	N/A	YP02119291
A12	hsa-miR-7158-3p	MIMAT0028227	YP02100572	I12	BlankSpot	N/A	YP02119291
A13	hsa-miR-7158-5p	MIMAT0028226	YP02112447	I13	BlankSpot	N/A	YP02119291
A14	hsa-miR-7159-5p	MIMAT0028228	YP02107640	I14	BlankSpot	N/A	YP02119291
A15	hsa-miR-7160-3p	MIMAT0028231	YP02102653	I15	BlankSpot	N/A	YP02119291
A16	hsa-miR-7160-5p	MIMAT0028230	YP02104288	I16	BlankSpot	N/A	YP02119291
A17	hsa-miR-7161-3p	MIMAT0028233	YP02102404	I17	BlankSpot	N/A	YP02119291
A18	hsa-miR-7161-5p	MIMAT0028232	YP02114009	I18	BlankSpot	N/A	YP02119291
A19	hsa-miR-7162-3p	MIMAT0028235	YP02102435	I19	BlankSpot	N/A	YP02119291
A20	hsa-miR-7162-5p	MIMAT0028234	YP02114557	I20	BlankSpot	N/A	YP02119291
A21	hsa-miR-7515	MIMAT0029310	YP02111847	I21	BlankSpot	N/A	YP02119291
A22	hsa-miR-7702	MIMAT0030017	YP02111608	I22	BlankSpot	N/A	YP02119291
A23	hsa-miR-7703	MIMAT0030018	YP02107095	I23	BlankSpot	N/A	YP02119291
A24	hsa-miR-7705	MIMAT0030020	YP02108293	I24	BlankSpot	N/A	YP02119291
B01	hsa-miR-7706	MIMAT0030021	YP02104933	J01	BlankSpot	N/A	YP02119291
B02	hsa-miR-7843-3p	MIMAT0030412	YP02101312	J02	BlankSpot	N/A	YP02119291
B03	hsa-miR-7843-5p	MIMAT0030411	YP02108752	J03	BlankSpot	N/A	YP02119291
B04	hsa-miR-7844-5p	MIMAT0030419	YP02111467	J04	BlankSpot	N/A	YP02119291
B05	hsa-miR-7845-5p	MIMAT0030420	YP02113971	J05	BlankSpot	N/A	YP02119291
B06	hsa-miR-7846-3p	MIMAT0030421	YP02112396	J06	BlankSpot	N/A	YP02119291
B07	hsa-miR-7847-3p	MIMAT0030422	YP02105945	J07	BlankSpot	N/A	YP02119291
B08	hsa-miR-7848-3p	MIMAT0030423	YP02103827	J08	BlankSpot	N/A	YP02119291
B09	hsa-miR-7849-3p	MIMAT0030424	YP02108167	J09	BlankSpot	N/A	YP02119291
B10	hsa-miR-7850-5p	MIMAT0030425	YP02104739	J10	BlankSpot	N/A	YP02119291
B11	hsa-miR-7851-3p	MIMAT0030426	YP02107258	J11	BlankSpot	N/A	YP02119291
B12	hsa-miR-7853-5p	MIMAT0030428	YP02116597	J12	BlankSpot	N/A	YP02119291
B13	hsa-miR-7854-3p	MIMAT0030429	YP02117617	J13	BlankSpot	N/A	YP02119291
B14	hsa-miR-7856-5p	MIMAT0030431	YP02113465	J14	BlankSpot	N/A	YP02119291
B15	hsa-miR-7973	MIMAT0031176	YP02114193	J15	BlankSpot	N/A	YP02119291
B16	hsa-miR-7974	MIMAT0031177	YP02102804	J16	BlankSpot	N/A	YP02119291
B17	hsa-miR-7975	MIMAT0031178	YP02107687	J17	BlankSpot	N/A	YP02119291
B18	hsa-miR-7976	MIMAT0031179	YP02100942	J18	BlankSpot	N/A	YP02119291
B19	hsa-miR-7977	MIMAT0031180	YP02101949	J19	BlankSpot	N/A	YP02119291
B20	hsa-miR-7978	MIMAT0031181	YP02100721	J20	BlankSpot	N/A	YP02119291
B21	hsa-miR-8052	MIMAT0030979	YP02104941	J21	BlankSpot	N/A	YP02119291
B22	hsa-miR-8053	MIMAT0030980	YP02109702	J22	BlankSpot	N/A	YP02119291
B23	hsa-miR-8054	MIMAT0030981	YP02118531	J23	BlankSpot	N/A	YP02119291
B24	hsa-miR-8055	MIMAT0030982	YP02114781	J24	BlankSpot	N/A	YP02119291
C01	hsa-miR-8057	MIMAT0030984	YP02107662	K01	BlankSpot	N/A	YP02119291

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
C02	hsa-miR-8058	MIMAT0030985	YP02109831	K02	BlankSpot	N/A	YP02119291
C03	hsa-miR-8059	MIMAT0030986	YP02101572	K03	BlankSpot	N/A	YP02119291
C04	hsa-miR-8060	MIMAT0030987	YP02114499	K04	BlankSpot	N/A	YP02119291
C05	hsa-miR-8061	MIMAT0030988	YP02100921	K05	BlankSpot	N/A	YP02119291
C06	hsa-miR-8062	MIMAT0030989	YP02113053	K06	BlankSpot	N/A	YP02119291
C07	hsa-miR-8063	MIMAT0030990	YP02107624	K07	BlankSpot	N/A	YP02119291
C08	hsa-miR-8064	MIMAT0030991	YP02105507	K08	BlankSpot	N/A	YP02119291
C09	hsa-miR-8066	MIMAT0030993	YP02109835	K09	BlankSpot	N/A	YP02119291
C10	hsa-miR-8067	MIMAT0030994	YP02108809	K10	BlankSpot	N/A	YP02119291
C11	hsa-miR-8068	MIMAT0030995	YP02113074	K11	BlankSpot	N/A	YP02119291
C12	hsa-miR-8070	MIMAT0030997	YP02105844	K12	BlankSpot	N/A	YP02119291
C13	hsa-miR-8071	MIMAT0030998	YP02106533	K13	BlankSpot	N/A	YP02119291
C14	hsa-miR-8072	MIMAT0030999	YP02116569	K14	BlankSpot	N/A	YP02119291
C15	hsa-miR-8073	MIMAT0031000	YP02116426	K15	BlankSpot	N/A	YP02119291
C16	hsa-miR-8074	MIMAT0031001	YP02102467	K16	BlankSpot	N/A	YP02119291
C17	hsa-miR-8075	MIMAT0031002	YP02106929	K17	BlankSpot	N/A	YP02119291
C18	hsa-miR-8076	MIMAT0031003	YP02108200	K18	BlankSpot	N/A	YP02119291
C19	hsa-miR-8077	MIMAT0031004	YP02110898	K19	BlankSpot	N/A	YP02119291
C20	hsa-miR-8078	MIMAT0031005	YP02106446	K20	BlankSpot	N/A	YP02119291
C21	hsa-miR-8079	MIMAT0031006	YP02106474	K21	BlankSpot	N/A	YP02119291
C22	hsa-miR-8080	MIMAT0031007	YP02118960	K22	BlankSpot	N/A	YP02119291
C23	hsa-miR-8081	MIMAT0031008	YP02113595	K23	BlankSpot	N/A	YP02119291
C24	hsa-miR-8082	MIMAT0031009	YP02103264	K24	BlankSpot	N/A	YP02119291
D01	hsa-miR-8083	MIMAT0031010	YP02116867	L01	BlankSpot	N/A	YP02119291
D02	hsa-miR-8085	MIMAT0031012	YP02113598	L02	BlankSpot	N/A	YP02119291
D03	hsa-miR-8086	MIMAT0031013	YP02114444	L03	BlankSpot	N/A	YP02119291
D04	hsa-miR-8087	MIMAT0031014	YP02111259	L04	BlankSpot	N/A	YP02119291
D05	hsa-miR-8088	MIMAT0031015	YP02118508	L05	BlankSpot	N/A	YP02119291
D06	hsa-miR-8089	MIMAT0031016	YP02118075	L06	BlankSpot	N/A	YP02119291
D07	hsa-miR-874-5p	MIMAT0026718	YP02105779	L07	BlankSpot	N/A	YP02119291
D08	hsa-miR-887-5p	MIMAT0026720	YP02118645	L08	BlankSpot	N/A	YP02119291
D09	hsa-miR-889-5p	MIMAT0026719	YP02119061	L09	BlankSpot	N/A	YP02119291
D10	hsa-miR-891a-3p	MIMAT0026717	YP02118478	L10	BlankSpot	N/A	YP02119291
D11	BlankSpot	N/A	YP02119291	L11	BlankSpot	N/A	YP02119291
D12	BlankSpot	N/A	YP02119291	L12	BlankSpot	N/A	YP02119291
D13	BlankSpot	N/A	YP02119291	L13	BlankSpot	N/A	YP02119291
D14	BlankSpot	N/A	YP02119291	L14	BlankSpot	N/A	YP02119291
D15	BlankSpot	N/A	YP02119291	L15	BlankSpot	N/A	YP02119291
D16	BlankSpot	N/A	YP02119291	L16	BlankSpot	N/A	YP02119291
D17	BlankSpot	N/A	YP02119291	L17	BlankSpot	N/A	YP02119291
D18	BlankSpot	N/A	YP02119291	L18	BlankSpot	N/A	YP02119291
D19	BlankSpot	N/A	YP02119291	L19	BlankSpot	N/A	YP02119291
D20	BlankSpot	N/A	YP02119291	L20	BlankSpot	N/A	YP02119291
D21	BlankSpot	N/A	YP02119291	L21	BlankSpot	N/A	YP02119291
D22	BlankSpot	N/A	YP02119291	L22	BlankSpot	N/A	YP02119291
D23	BlankSpot	N/A	YP02119291	L23	BlankSpot	N/A	YP02119291
D24	BlankSpot	N/A	YP02119291	L24	BlankSpot	N/A	YP02119291
E01	BlankSpot	N/A	YP02119291	M01	BlankSpot	N/A	YP02119291
E02	BlankSpot	N/A	YP02119291	M02	BlankSpot	N/A	YP02119291
E03	BlankSpot	N/A	YP02119291	M03	BlankSpot	N/A	YP02119291

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
E04	BlankSpot	N/A	YP02119291	M04	BlankSpot	N/A	YP02119291
E05	BlankSpot	N/A	YP02119291	M05	BlankSpot	N/A	YP02119291
E06	BlankSpot	N/A	YP02119291	M06	BlankSpot	N/A	YP02119291
E07	BlankSpot	N/A	YP02119291	M07	BlankSpot	N/A	YP02119291
E08	BlankSpot	N/A	YP02119291	M08	BlankSpot	N/A	YP02119291
E09	BlankSpot	N/A	YP02119291	M09	BlankSpot	N/A	YP02119291
E10	BlankSpot	N/A	YP02119291	M10	BlankSpot	N/A	YP02119291
E11	BlankSpot	N/A	YP02119291	M11	BlankSpot	N/A	YP02119291
E12	BlankSpot	N/A	YP02119291	M12	BlankSpot	N/A	YP02119291
E13	BlankSpot	N/A	YP02119291	M13	BlankSpot	N/A	YP02119291
E14	BlankSpot	N/A	YP02119291	M14	BlankSpot	N/A	YP02119291
E15	BlankSpot	N/A	YP02119291	M15	BlankSpot	N/A	YP02119291
E16	BlankSpot	N/A	YP02119291	M16	BlankSpot	N/A	YP02119291
E17	BlankSpot	N/A	YP02119291	M17	BlankSpot	N/A	YP02119291
E18	BlankSpot	N/A	YP02119291	M18	BlankSpot	N/A	YP02119291
E19	BlankSpot	N/A	YP02119291	M19	BlankSpot	N/A	YP02119291
E20	BlankSpot	N/A	YP02119291	M20	BlankSpot	N/A	YP02119291
E21	BlankSpot	N/A	YP02119291	M21	BlankSpot	N/A	YP02119291
E22	BlankSpot	N/A	YP02119291	M22	BlankSpot	N/A	YP02119291
E23	BlankSpot	N/A	YP02119291	M23	BlankSpot	N/A	YP02119291
E24	BlankSpot	N/A	YP02119291	M24	BlankSpot	N/A	YP02119291
F01	BlankSpot	N/A	YP02119291	N01	BlankSpot	N/A	YP02119291
F02	BlankSpot	N/A	YP02119291	N02	BlankSpot	N/A	YP02119291
F03	BlankSpot	N/A	YP02119291	N03	BlankSpot	N/A	YP02119291
F04	BlankSpot	N/A	YP02119291	N04	BlankSpot	N/A	YP02119291
F05	BlankSpot	N/A	YP02119291	N05	BlankSpot	N/A	YP02119291
F06	BlankSpot	N/A	YP02119291	N06	BlankSpot	N/A	YP02119291
F07	BlankSpot	N/A	YP02119291	N07	BlankSpot	N/A	YP02119291
F08	BlankSpot	N/A	YP02119291	N08	BlankSpot	N/A	YP02119291
F09	BlankSpot	N/A	YP02119291	N09	BlankSpot	N/A	YP02119291
F10	BlankSpot	N/A	YP02119291	N10	BlankSpot	N/A	YP02119291
F11	BlankSpot	N/A	YP02119291	N11	BlankSpot	N/A	YP02119291
F12	BlankSpot	N/A	YP02119291	N12	BlankSpot	N/A	YP02119291
F13	BlankSpot	N/A	YP02119291	N13	BlankSpot	N/A	YP02119291
F14	BlankSpot	N/A	YP02119291	N14	BlankSpot	N/A	YP02119291
F15	BlankSpot	N/A	YP02119291	N15	BlankSpot	N/A	YP02119291
F16	BlankSpot	N/A	YP02119291	N16	BlankSpot	N/A	YP02119291
F17	BlankSpot	N/A	YP02119291	N17	BlankSpot	N/A	YP02119291
F18	BlankSpot	N/A	YP02119291	N18	BlankSpot	N/A	YP02119291
F19	BlankSpot	N/A	YP02119291	N19	BlankSpot	N/A	YP02119291
F20	BlankSpot	N/A	YP02119291	N20	BlankSpot	N/A	YP02119291
F21	BlankSpot	N/A	YP02119291	N21	BlankSpot	N/A	YP02119291
F22	BlankSpot	N/A	YP02119291	N22	BlankSpot	N/A	YP02119291
F23	BlankSpot	N/A	YP02119291	N23	BlankSpot	N/A	YP02119291
F24	BlankSpot	N/A	YP02119291	N24	BlankSpot	N/A	YP02119291
G01	BlankSpot	N/A	YP02119291	O01	BlankSpot	N/A	YP02119291
G02	BlankSpot	N/A	YP02119291	O02	BlankSpot	N/A	YP02119291
G03	BlankSpot	N/A	YP02119291	O03	BlankSpot	N/A	YP02119291
G04	BlankSpot	N/A	YP02119291	O04	BlankSpot	N/A	YP02119291
G05	BlankSpot	N/A	YP02119291	O05	BlankSpot	N/A	YP02119291

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
G06	BlankSpot	N/A	YP02119291	O06	BlankSpot	N/A	YP02119291
G07	BlankSpot	N/A	YP02119291	O07	BlankSpot	N/A	YP02119291
G08	BlankSpot	N/A	YP02119291	O08	BlankSpot	N/A	YP02119291
G09	BlankSpot	N/A	YP02119291	O09	BlankSpot	N/A	YP02119291
G10	BlankSpot	N/A	YP02119291	O10	BlankSpot	N/A	YP02119291
G11	BlankSpot	N/A	YP02119291	O11	BlankSpot	N/A	YP02119291
G12	BlankSpot	N/A	YP02119291	O12	BlankSpot	N/A	YP02119291
G13	BlankSpot	N/A	YP02119291	O13	BlankSpot	N/A	YP02119291
G14	BlankSpot	N/A	YP02119291	O14	BlankSpot	N/A	YP02119291
G15	BlankSpot	N/A	YP02119291	O15	BlankSpot	N/A	YP02119291
G16	BlankSpot	N/A	YP02119291	O16	BlankSpot	N/A	YP02119291
G17	BlankSpot	N/A	YP02119291	O17	BlankSpot	N/A	YP02119291
G18	BlankSpot	N/A	YP02119291	O18	BlankSpot	N/A	YP02119291
G19	BlankSpot	N/A	YP02119291	O19	BlankSpot	N/A	YP02119291
G20	BlankSpot	N/A	YP02119291	O20	BlankSpot	N/A	YP02119291
G21	BlankSpot	N/A	YP02119291	O21	BlankSpot	N/A	YP02119291
G22	BlankSpot	N/A	YP02119291	O22	BlankSpot	N/A	YP02119291
G23	BlankSpot	N/A	YP02119291	O23	BlankSpot	N/A	YP02119291
G24	BlankSpot	N/A	YP02119291	O24	BlankSpot	N/A	YP02119291
H01	BlankSpot	N/A	YP02119291	P01	BlankSpot	N/A	YP02119291
H02	BlankSpot	N/A	YP02119291	P02	BlankSpot	N/A	YP02119291
H03	BlankSpot	N/A	YP02119291	P03	BlankSpot	N/A	YP02119291
H04	BlankSpot	N/A	YP02119291	P04	BlankSpot	N/A	YP02119291
H05	BlankSpot	N/A	YP02119291	P05	BlankSpot	N/A	YP02119291
H06	BlankSpot	N/A	YP02119291	P06	BlankSpot	N/A	YP02119291
H07	BlankSpot	N/A	YP02119291	P07	BlankSpot	N/A	YP02119291
H08	BlankSpot	N/A	YP02119291	P08	BlankSpot	N/A	YP02119291
H09	BlankSpot	N/A	YP02119291	P09	BlankSpot	N/A	YP02119291
H10	BlankSpot	N/A	YP02119291	P10	BlankSpot	N/A	YP02119291
H11	BlankSpot	N/A	YP02119291	P11	BlankSpot	N/A	YP02119291
H12	BlankSpot	N/A	YP02119291	P12	BlankSpot	N/A	YP02119291
H13	BlankSpot	N/A	YP02119291	P13	SNORD44 (hsa)	N/A	YP00203902
H14	BlankSpot	N/A	YP02119291	P14	SNORD38B (hsa)	N/A	YP00203901
H15	BlankSpot	N/A	YP02119291	P15	SNORD49A (hsa)	N/A	YP00203904
H16	BlankSpot	N/A	YP02119291	P16	U6 snRNA (v2)	N/A	YP02119464
H17	BlankSpot	N/A	YP02119291	P17	cel-miR-39-3p	MIMAT0000010	YP00203952
H18	BlankSpot	N/A	YP02119291	P18	UniSp2	N/A	YP00203950
H19	BlankSpot	N/A	YP02119291	P19	UniSp4	N/A	YP00203953
H20	BlankSpot	N/A	YP02119291	P20	UniSp5	N/A	YP00203955
H21	BlankSpot	N/A	YP02119291	P21	UniSp6	N/A	YP00203954
H22	BlankSpot	N/A	YP02119291	P22	UniSP3	N/A	YP02119288
H23	BlankSpot	N/A	YP02119291	P23	UniSP3	N/A	YP02119288
H24	BlankSpot	N/A	YP02119291	P24	UniSP3	N/A	YP02119288

Ordering Information

Product	Contents	Cat. no.
miRCURY LNA miRNA Focus PCR Panels	miRCURY LNA miRNA PCR Panels for application-based miRNome profiling, available in 96-well or 384-well format; for SYBR® Green-based detection	339325
miRCURY LNA miRNA miRNome PCR Panels	miRCURY LNA miRNA PCR Panels for PCR-based miRNome profiling, available in 384-well format; for SYBR® Green-based detection	339322
miRCURY LNA miRNA QC PCR Panel	miRCURY LNA miRNA PCR Panel of quality control assays, available in 96-well or 384-well format; for SYBR® Green-based detection	339331
miRCURY LNA miRNA Custom PCR Panels	8 identical, ready-to-use 96- or 384-well plates; each well contains primers sufficient for one 10 µl reaction; for SYBR® Green-based detection	339330
miRCURY LNA Custom PCR Panel Additional Plate	Additional miRCURY LNA Custom PCR Panel plates; set of 4 plates; only available in addition to the base plates ordered through the core product (cat. no. 339330)	339332
miRCURY LNA miRNA PCR Assays	Contains forward and reverse primers for 200 SYBR® Green-based, real-time qPCR reactions, 166 EvaGreen-based digital PCR reactions for Nanoplate 8.5k or 50 EvaGreen-based digital PCR reactions for Nanoplate 26k	339306
miRCURY LNA miRNA Probe PCR Assays	Complete premixed assays containing LNA-enhanced target-specific forward primer and probe. For 200 reactions.	339350
miRCURY LNA miRNA Custom Probe PCR Assays	Custom-designed, target-specific forward primer and probe for any user-defined miRNA target. Complete premixed assay for 200 reactions.	339351

Related Products

Product	Contents	Cat. no.
miRCURY LNA RT Kit	For 8–64 cDNA synthesis reactions: 5x RT SYBR Green Reaction Buffer, 5x RT Probe Reaction Buffer, 10x RT Enzyme Mix, UniSp6, RNA Spike-in template, RNase-Free Water	339340
RNA Spike-In Kit, For RT	Contains the UniSp2, UniSp4, and UniSp5 RNA Spike-in Template Mix and the cel-miR-39-3p RNA Spike-in Template	339390
miRCURY LNA SYBR® Green PCR Kits (200)	For 200 reactions: 2X miRCURY SYBR Green Master Mix, RNase-Free Water	339345
miRCURY LNA SYBR® Green PCR Kits (600)	For 600 reactions: 2X miRCURY SYBR Green Master Mix, RNase-Free Water	339346
miRCURY LNA SYBR® Green PCR Kits (4000)	For 4000 reactions: 2X miRCURY SYBR Green Master Mix, RNase-Free Water	339347
miRCURY LNA Probe PCR Kit (200)	For 200 reactions: 2X QuantiNova Probe Master Mix, 10X miRCURY Probe Univ. Primer, Rox Reference Dye, RNase-Free Water	339371
miRCURY LNA Probe PCR Kit (800)	For 800 reactions: 2X QuantiNova Probe Master Mix, 10X miRCURY Probe Univ. Primer, Rox Reference Dye, RNase-Free Water	339372
miRCURY LNA Probe PCR Kit (4000)	For 4000 reactions: 2X QuantiNova Probe Master Mix, 10X miRCURY Probe Univ. Primer, Rox Reference Dye, RNase-Free Water	339373
miRCURY LNA miRNA PCR Starter Kit	Two miRCURY LNA PCR Assays of your choice, spike-in control Assay (UniSp6), one candidate endogenous control assay (miR-103-3p) and all reagents for 20 reverse transcription reactions and 100 PCR amplifications; for SYBR® Green-based qPCR detection	339320

Important

The miRCURY LNA™ miRNA PCR Panels are Ready-to-Use and designed for optimal performance with the miRCURY LNA RT Kit and the miRCURY LNA SYBR® Green PCR Kit. The performance of the primer sets will be affected when used in combination with less than optimal reagents. Use the miRCURY LNA miRNA PCR Panel Handbook for experiment setup. RNA work requires specific handling and precautions should be taken to prevent RNase contamination and degradation of the RNA sample and reagents.

Shipping and storage

The Ready-to-Use miRCURY LNA miRNA PCR Panels are shipped at room temperature and can be stored at 4°C for at least 6 months. For long term storage, it is recommended to place the panels at -20°C. Under these conditions, the LNA PCR primers are stable for at least 6 months after receipt.

Intended use

miRCURY LNA miRNA PCR Assays are intended for molecular biology applications. This product is not intended for the diagnosis, prevention or treatment of a disease.

Safety information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). These are available online in convenient and compact PDF format at www.qiagen.com/safety where you can find, view and print the SDS for each QIAGEN kit and kit component.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, Sample to Insight®, GeneGlobe®, LNA™ (QIAGEN Group), SYBR® (Life Technologies Corporation). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law. 1107212 05/2017 HB-2376-001 © 2017 QIAGEN, all rights reserved.
