

miRCURY LNA™ miRNA Focus PCR Panels

Rat Inflammatory Response & Autoimmunity V2

Product Data Sheet

Cat. no. 339325 YARN-205Y

For mature miRNA expression profiling using real-time PCR

Format	Suitable real-time cyclers	Plate	Cat. no.
A	Applied Biosystems® models 5700, 7000, 7300, 7500, 7700, 7900HT, ViiA™ 7 (96-well block); Bio-Rad® models iCycler®, iQ™5, MyiQ™, MyiQ2; Bio-Rad/MJ Research Chromo4™; Eppendorf® MasterCycler® ep realplex models 2, 2s, 4, 4s; Stratagene® models Mx3005P®, Mx3000P®; Takara: TP-800	96-well	YARN-205YA
C	Applied Biosystems models 7500 (Fast block), 7900HT (Fast block), StepOnePlus™, ViiA 7 (Fast block)	96-well	YARN-205YC
D	Bio-Rad CFX96™; Bio-Rad/MJ Research models DNA Engine Opticon®, DNA Engine Opticon 2; Stratagene Mx4000®	96-well	YARN-205YD
E	Applied Biosystems® models 7900HT (384-well block), ViiA™ 7 (384-well block); Bio-Rad CFX384™	384-well	YARN-205YE
F	Roche® LightCycler® 480 (96-well block)	96-well	YARN-205YF
G	Roche® LightCycler® 480 (384-well block)	384-well	YARN-205YG

Description

The Rat Inflammatory Response & Autoimmunity V2 miRCURY LNA™ miRNA Focus PCR Panel profiles the expression of 84 miRNAs predicted to regulate the expression of proinflammatory or antiinflammatory genes. This array provides researchers with a convenient way to quickly analyze the miRNAs most likely to be relevant to inflammatory and autoimmune disorders. These miRNAs have been carefully selected based on those predicted by bioinformatic algorithms and databases to regulate genes known to be relevant to inflammation. The profiling results from this array can serve as a useful molecular marker for the study of normal and pathological inflammation. The results may also provide a deeper understanding of the molecular mechanisms regulating inflammatory gene expression by miRNA. A set of controls present on this array enables data analysis using the $\Delta\Delta CT$ method of relative quantification, assessment of reverse transcription performance, and assessment of PCR performance. Using SYBR Green-based real-time PCR, the expression of a focused panel of miRNAs predicted to regulate inflammatory and autoimmune response genes can be easily and reliably analyzed with this miRCURY LNA™ miRNA Focus PCR Panel.

For further details, consult the *miRCURY LNA™ miRNA Focus PCR Panels Handbook*.

Array Layout

	1	2	3	4	5	6	7	8	9	10	11	12
A	rno-let-7a-5p	rno-let-7b-5p	rno-let-7c-5p	rno-let-7d-5p	rno-let-7e-5p	rno-let-7f-5p	rno-let-7i-5p	rno-miR-101a-3p	rno-miR-101b-3p	rno-miR-106b-5p	rno-miR-125a-5p	rno-miR-125b-5p
B	rno-miR-128-3p	rno-miR-136-5p	rno-miR-140-5p	rno-miR-141-3p	rno-miR-142-3p	rno-miR-144-3p	rno-miR-145-5p	rno-miR-148b-3p	rno-miR-152-3p	rno-miR-15b-5p	rno-miR-16-5p	rno-miR-17-5p
C	rno-miR-181a-5p	rno-miR-181b-5p	rno-miR-181c-5p	rno-miR-181d-5p	rno-miR-182	rno-miR-183-5p	rno-miR-186-5p	rno-miR-195-5p	rno-miR-19a-3p	rno-miR-19b-3p	rno-miR-200a-3p	rno-miR-200c-3p
D	rno-miR-203a-3p	rno-miR-205	rno-miR-20a-5p	rno-miR-20b-5p	rno-miR-21-5p	rno-miR-221-3p	rno-miR-222-3p	rno-miR-23a-3p	rno-miR-23b-3p	rno-miR-26a-5p	rno-miR-26b-5p	rno-miR-27a-3p
E	rno-miR-27b-3p	rno-miR-291a-3p	rno-miR-29a-3p	rno-miR-29b-3p	rno-miR-29c-3p	rno-miR-30a-5p	rno-miR-30b-5p	rno-miR-30c-5p	rno-miR-30d-5p	rno-miR-30e-5p	rno-miR-320-3p	rno-miR-322-5p
F	rno-miR-323-3p	rno-miR-325-3p	rno-miR-327	rno-miR-34a-5p	rno-miR-34c-5p	rno-miR-351-5p	rno-miR-369-3p	rno-miR-374-5p	rno-miR-381-3p	rno-miR-384-5p	rno-miR-410-3p	rno-miR-429
G	rno-miR-448-3p	rno-miR-449a-5p	rno-miR-495	rno-miR-497-5p	rno-miR-539-5p	rno-miR-664-3p	rno-miR-673-5p	rno-miR-743b-3p	rno-miR-878	rno-miR-9a-5p	rno-miR-93-5p	rno-miR-98-5p
H	cel-miR-39-3p	cel-miR-39-3p	U6 snRNA (V2)	5S rRNA	RNU5G	RNU1A1	UniSp2	UniSp4	UniSp5	UniSp6	UniSp3	UniSp3

miRNA Table

Well	miRNA ID	Accession #	Assay Catalog #	Well	miRNA ID	Accession #	Assay Catalog #
A01	rno-let-7a-5p	MIMAT0000062	YP00205727	E01	rno-miR-27b-3p	MIMAT0000419	YP00205915
A02	rno-let-7b-5p	MIMAT0000063	YP00204750	E02	rno-miR-291a-3p	MIMAT0000895	YP00206041
A03	rno-let-7c-5p	MIMAT0000064	YP00204767	E03	rno-miR-29a-3p	MIMAT0000086	YP00204698
A04	rno-let-7d-5p	MIMAT0000065	YP00204124	E04	rno-miR-29b-3p	MIMAT0000100	YP00204679
A05	rno-let-7e-5p	MIMAT0000066	YP00205711	E05	rno-miR-29c-3p	MIMAT0000681	YP00204729
A06	rno-let-7f-5p	MIMAT0000067	YP00204359	E06	rno-miR-30a-5p	MIMAT0000087	YP00205695
A07	rno-let-7i-5p	MIMAT0000415	YP00204394	E07	rno-miR-30b-5p	MIMAT0000420	YP00204765
A08	rno-miR-101a-3p	MIMAT0000099	YP00204786	E08	rno-miR-30c-5p	MIMAT0000244	YP00204783
A09	rno-miR-101b-3p	MIMAT0000615	YP02127460	E09	rno-miR-30d-5p	MIMAT0000245	YP00206047
A10	rno-miR-106b-5p	MIMAT0000680	YP00205884	E10	rno-miR-30e-5p	MIMAT0000692	YP00204714
A11	rno-miR-125a-5p	MIMAT0000443	YP00204339	E11	rno-miR-320-3p	MIMAT0000510	YP00206042
A12	rno-miR-125b-5p	MIMAT0000423	YP00205713	E12	rno-miR-322-5p	MIMAT0000548	YP00205182
B01	rno-miR-128-3p	MIMAT0000424	YP00205995	F01	rno-miR-323-3p	MIMAT0000755	YP00204278
B02	rno-miR-136-5p	MIMAT0000448	YP00204779	F02	rno-miR-325-3p	MIMAT0000460	YP00205088
B03	rno-miR-140-5p	MIMAT0000431	YP00204540	F03	rno-miR-327	MIMAT0000561	YP00205524
B04	rno-miR-141-3p	MIMAT0000432	YP00204504	F04	rno-miR-34a-5p	MIMAT0000255	YP00204486
B05	rno-miR-142-3p	MIMAT0000434	YP00204291	F05	rno-miR-34c-5p	MIMAT0000686	YP00205659
B06	rno-miR-144-3p	MIMAT0000436	YP00204754	F06	rno-miR-351-5p	MIMAT0000608	YP00205534
B07	rno-miR-145-5p	MIMAT0000437	YP00204483	F07	rno-miR-369-3p	MIMAT0000721	YP00206028
B08	rno-miR-148b-3p	MIMAT0000759	YP00204047	F08	rno-miR-374-5p	MIMAT00004955	YP00204608
B09	rno-miR-152-3p	MIMAT0000438	YP00204294	F09	rno-miR-381-3p	MIMAT0003199	YP00205540
B10	rno-miR-15b-5p	MIMAT0000417	YP00204243	F10	rno-miR-384-5p	MIMAT00004745	YP00205189
B11	rno-miR-16-5p	MIMAT0000069	YP00205702	F11	rno-miR-410-3p	MIMAT0002171	YP00204042
B12	rno-miR-17-5p	MIMAT0000070	YP02119304	F12	rno-miR-429	MIMAT0001537	YP00205068
C01	rno-miR-181a-5p	MIMAT0000256	YP00206081	G01	rno-miR-448-3p	MIMAT0001534	YP02102683
C02	rno-miR-181b-5p	MIMAT0000257	YP00204530	G02	rno-miR-449a-5p	MIMAT0001541	YP00204481
C03	rno-miR-181c-5p	MIMAT0000258	YP00204683	G03	rno-miR-495	MIMAT0002817	YP00206015
C04	rno-miR-181d-5p	MIMAT0002821	YP00204789	G04	rno-miR-497-5p	MIMAT0003453	YP00205164
C05	rno-miR-182	MIMAT0000211	YP00205089	G05	rno-miR-539-5p	MIMAT0003163	YP00205656
C06	rno-miR-183-5p	MIMAT0000261	YP00206030	G06	rno-miR-664-3p	MIMAT0012774	YP02106996
C07	rno-miR-186-5p	MIMAT0000456	YP00206053	G07	rno-miR-673-5p	MIMAT0005328	YP00205566
C08	rno-miR-195-5p	MIMAT0000461	YP00205869	G08	rno-miR-743b-3p	MIMAT0005280	YP00205570
C09	rno-miR-19a-3p	MIMAT0000073	YP00205862	G09	rno-miR-878	MIMAT0005286	YP00205576
C10	rno-miR-19b-3p	MIMAT0000074	YP00204450	G10	rno-miR-9a-5p	MIMAT0000441	YP00204513
C11	rno-miR-200a-3p	MIMAT0000682	YP00204707	G11	rno-miR-93-5p	MIMAT0000093	YP00204715
C12	rno-miR-200c-3p	MIMAT0000873	YP00205505	G12	rno-miR-98-5p	MIMAT0000096	YP00204640
D01	rno-miR-203a-3p	MIMAT0000264	YP00205914	H01	cel-miR-39-3p	MIMAT0000010	YP00203952
D02	rno-miR-205	MIMAT0000878	YP00205958	H02	cel-miR-39-3p	MIMAT0000010	YP00203952
D03	rno-miR-20a-5p	MIMAT0000075	YP00204292	H03	U6 snRNA (v2)	N/A	YP02119464
D04	rno-miR-20b-5p	MIMAT0001413	YP00204755	H04	5S rRNA	N/A	YP00203906
D05	rno-miR-21-5p	MIMAT0000076	YP00204230	H05	RNU5G	N/A	YP00203908
D06	rno-miR-221-3p	MIMAT0000278	YP00204532	H06	RNU1A1	N/A	YP00203909
D07	rno-miR-222-3p	MIMAT0000279	YP00204551	H07	UniSp2	N/A	YP00203950
D08	rno-miR-23a-3p	MIMAT0000078	YP00204772	H08	UniSp4	N/A	YP00203953
D09	rno-miR-23b-3p	MIMAT0000125	YP02119756	H09	UniSp5	N/A	YP00203955
D10	rno-miR-26a-5p	MIMAT0000082	YP00206023	H10	UniSp6	N/A	YP00203954
D11	rno-miR-26b-5p	MIMAT0000083	YP00204172	H11	UniSP3	N/A	YP02119288
D12	rno-miR-27a-3p	MIMAT0000084	YP00206038	H12	UniSP3	N/A	YP02119288

Functional Groupings

Inflammation: rno-let-7e-5p,rno-miR-106b-5p,rno-miR-140-5p,rno-miR-15b-5p,rno-miR-16-5p,rno-miR-17-5p,rno-miR-181a-5p,rno-miR-181b-5p,rno-miR-181c-5p,rno-miR-182,rno-miR-19a-3p,rno-miR-19b-3p,rno-miR-20a-5p,rno-miR-221-3p,rno-miR-222-3p,rno-miR-23b-3p,rno-miR-26a-5p,rno-miR-26b-5p,rno-miR-27b-3p,rno-miR-30e-5p,rno-miR-369-3p,rno-miR-98-5p.

Autoimmunity: rno-let-7d-5p,rno-let-7e-5p,rno-let-7f-5p,rno-let-7i-5p,rno-miR-128-3p,rno-miR-140-5p,rno-miR-16-5p,rno-miR-17-5p,rno-miR-181a-5p,rno-miR-181b-5p,rno-miR-182,rno-miR-186-5p,rno-miR-19a-3p,rno-miR-20a-5p,rno-miR-23a-3p,rno-miR-23b-3p,rno-miR-27b-3p,rno-miR-29a-3p,rno-miR-29c-3p,rno-miR-369-3p,rno-miR-93-5p.

Innate & Adaptive Immunity: rno-let-7i-5p,rno-miR-17-5p,rno-miR-181a-5p,rno-miR-182,rno-miR-19a-3p,rno-miR-19b-3p,rno-miR-20a-5p,rno-miR-322-5p.

Bioinformatics Predictions

Regulation of Acute-Phase Response: rno-let-7a-5p,rno-let-7b-5p,rno-let-7c-5p,rno-let-7d-5p,rno-let-7e-5p,rno-let-7f-5p,rno-let-7i-5p,rno-miR-106b-5p,rno-miR-125a-5p,rno-miR-125b-5p,rno-miR-141-3p,rno-miR-17-5p,rno-miR-200a-3p,rno-miR-200c-3p,rno-miR-20a-5p,rno-miR-20b-5p,rno-miR-21-5p,rno-miR-23a-3p,rno-miR-23b-3p,rno-miR-325-3p,rno-miR-34a-5p,rno-miR-34c-5p,rno-miR-351-5p,rno-miR-369-3p,rno-miR-374-5p,rno-miR-381-3p,rno-miR-410-3p,rno-miR-429,rno-miR-449a-5p,rno-miR-664-3p,rno-miR-93-5p,rno-miR-98-5p.

Regulation of Cytokine Metabolism: rno-let-7a-5p,rno-let-7b-5p,rno-let-7c-5p,rno-let-7d-5p,rno-let-7e-5p,rno-let-7f-5p,rno-let-7i-5p,rno-miR-125a-5p,rno-miR-125b-5p,rno-miR-128-3p,rno-miR-145-5p,rno-miR-148b-3p,rno-miR-152-3p,rno-miR-15b-5p,rno-miR-16-5p,rno-miR-181a-5p,rno-miR-181b-5p,rno-miR-181c-5p,rno-miR-181d-5p,rno-miR-186-5p,rno-miR-195-5p,rno-miR-19a-3p,rno-miR-19b-3p,rno-miR-205,rno-miR-221-3p,rno-miR-222-3p,rno-miR-26a-5p,rno-miR-26b-5p,rno-miR-27a-3p,rno-miR-27b-3p,rno-miR-291a-3p,rno-miR-29a-3p,rno-miR-29b-3p,rno-miR-29c-3p,rno-miR-30a-5p,rno-miR-30b-5p,rno-miR-30c-5p,rno-miR-30d-5p,rno-miR-30e-5p,rno-miR-322-5p,rno-miR-325-3p,rno-miR-34a-5p,rno-miR-34c-5p,rno-miR-351-5p,rno-miR-369-3p,rno-miR-374-5p,rno-miR-381-3p,rno-miR-384-5p,rno-miR-410-3p,rno-miR-448-3p,rno-miR-449a-5p,rno-miR-497-5p,rno-miR-539-5p,rno-miR-664-3p,rno-miR-673-5p,rno-miR-98-5p,rno-miR-9a-5p.

Regulation of Cytokine Receptors: rno-let-7a-5p,rno-let-7b-5p,rno-let-7c-5p,rno-let-7d-5p,rno-let-7e-5p,rno-let-7f-5p,rno-let-7i-5p,rno-miR-140-5p,rno-miR-142-3p,rno-miR-148b-3p,rno-miR-152-3p,rno-miR-15b-5p,rno-miR-16-5p,rno-miR-186-5p,rno-miR-195-5p,rno-miR-19a-3p,rno-miR-19b-3p,rno-miR-205,rno-miR-21-5p,rno-miR-23a-3p,rno-miR-23b-3p,rno-miR-29a-3p,rno-miR-29b-3p,rno-miR-29c-3p,rno-miR-322-5p,rno-miR-323-3p,rno-miR-325-3p,rno-miR-381-3p,rno-miR-497-5p,rno-miR-673-5p,rno-miR-743b-3p,rno-miR-98-5p,rno-miR-9a-5p.

Regulation of Cytokines: rno-let-7a-5p,rno-let-7b-5p,rno-let-7c-5p,rno-let-7d-5p,rno-let-7e-5p,rno-let-7f-5p,rno-let-7i-5p,rno-miR-101a-3p,rno-miR-101b-3p,rno-miR-106b-5p,rno-miR-125a-5p,rno-miR-125b-5p,rno-miR-128-3p,rno-miR-141-3p,rno-miR-144-3p,rno-miR-145-5p,rno-miR-148b-3p,rno-miR-152-3p,rno-miR-15b-5p,rno-miR-16-5p,rno-miR-17-5p,rno-miR-181a-5p,rno-miR-181b-5p,rno-miR-181c-5p,rno-miR-181d-5p,rno-miR-183-5p,rno-miR-186-5p,rno-miR-195-5p,rno-miR-19a-3p,rno-miR-19b-3p,rno-miR-200a-3p,rno-miR-203a-3p,rno-miR-205,rno-miR-20a-5p,rno-miR-20b-5p,rno-miR-21-5p,rno-miR-23a-3p,rno-miR-23b-3p,rno-miR-26a-5p,rno-miR-26b-5p,rno-miR-27a-3p,rno-miR-27b-3p,rno-miR-291a-3p,rno-miR-29a-3p,rno-miR-29b-3p,rno-miR-29c-3p,rno-miR-30a-5p,rno-miR-30b-5p,rno-miR-30c-5p,rno-miR-30d-5p,rno-miR-30e-5p,rno-miR-320-3p,rno-miR-322-5p,rno-miR-323-3p,rno-miR-325-3p,rno-miR-327,rno-miR-34a-5p,rno-miR-34c-5p,rno-miR-351-5p,rno-miR-369-3p,rno-miR-374-5p,rno-miR-384-5p,rno-miR-410-3p,rno-miR-448-3p,rno-miR-449a-5p,rno-miR-495,rno-miR-497-5p,rno-miR-539-5p,rno-miR-664-3p,rno-miR-673-5p,rno-miR-743b-3p,rno-miR-878,rno-miR-93-5p,rno-miR-98-5p,rno-miR-9a-5p.

Regulation of Humoral Immune Response: rno-let-7a-5p,rno-let-7b-5p,rno-let-7c-5p,rno-let-7d-5p,rno-let-7e-5p,rno-let-7f-5p,rno-let-7i-5p,rno-miR-145-5p,rno-miR-148b-3p,rno-miR-152-3p,rno-miR-15b-5p,rno-miR-16-5p,rno-miR-181a-5p,rno-miR-181b-5p,rno-miR-181c-5p,rno-miR-181d-5p,rno-miR-186-5p,rno-miR-195-5p,rno-miR-21-5p,rno-miR-23a-3p,rno-miR-23b-3p,rno-miR-27a-3p,rno-miR-27b-3p,rno-miR-291a-3p,rno-miR-29a-3p,rno-miR-29b-3p,rno-miR-29c-3p,rno-miR-322-5p,rno-miR-323-3p,rno-miR-325-3p,rno-miR-369-3p,rno-miR-374-5p,rno-miR-381-3p,rno-miR-410-3p,rno-miR-497-5p,rno-miR-539-5p,rno-miR-673-5p,rno-miR-98-5p.

Regulation of Other Immune Response Genes: rno-miR-101a-3p,rno-miR-101b-3p,rno-miR-106b-5p,rno-miR-136-5p,rno-miR-140-5p,rno-miR-142-3p,rno-miR-144-3p,rno-miR-145-5p,rno-miR-15b-5p,rno-miR-16-5p,rno-miR-17-5p,rno-miR-183-5p,rno-miR-186-5p,rno-miR-195-5p,rno-miR-19a-3p,rno-miR-19b-3p,rno-miR-200c-3p,rno-miR-20a-5p,rno-miR-20b-5p,rno-miR-221-3p,rno-miR-222-3p,rno-miR-23a-3p,rno-miR-23b-3p,rno-miR-291a-3p,rno-miR-322-5p,rno-miR-323-3p,rno-miR-34a-5p,rno-miR-34c-5p,rno-miR-369-3p,rno-miR-381-3p,rno-miR-410-3p,rno-miR-429,rno-miR-448-3p,rno-miR-449a-5p,rno-miR-497-5p,rno-miR-673-5p,rno-miR-93-5p.

Ordering Information

Product	Contents	Cat. no.
miRCURY LNA miRNA Focus PCR Panels	miRCURY LNA miRNA PCR Panels for application-based miRNome profiling, available in 96-well or 384-well format; for SYBR® Green-based detection	339325
miRCURY LNA miRNA miRNome PCR Panels	miRCURY LNA miRNA PCR Panels for PCR-based miRNome profiling, available in 384-well format; for SYBR® Green-based detection	339322
miRCURY LNA miRNA QC PCR Panel	miRCURY LNA miRNA PCR Panel of quality control assays, available in 96-well or 384-well format; for SYBR® Green-based detection	339331
miRCURY LNA miRNA Custom PCR Panels	8 identical, ready-to-use 96- or 384-well plates; each well contains primers sufficient for one 10 µl reaction; for SYBR® Green-based detection	339330
miRCURY LNA Custom PCR Panel Additional Plate	Additional miRCURY LNA Custom PCR Panel plates; set of 4 plates; only available in addition to the base plates ordered through the core product (cat. no. 339330)	339332
miRCURY LNA miRNA PCR Assays	Contains forward and reverse primers for 200 SYBR® Green-based, real-time qPCR reactions, 166 EvaGreen-based digital PCR reactions for Nanoplate 8.5k or 50 EvaGreen-based digital PCR reactions for Nanoplate 26k	339306
miRCURY LNA miRNA Probe PCR Assays	Complete premixed assays containing LNA-enhanced target-specific forward primer and probe. For 200 reactions.	339350
miRCURY LNA miRNA Custom Probe PCR Assays	Custom-designed, target-specific forward primer and probe for any user-defined miRNA target. Complete premixed assay for 200 reactions.	339351

Related Products

Product	Contents	Cat. no.
miRCURY LNA RT Kit	For 8–64 cDNA synthesis reactions: 5x RT SYBR Green Reaction Buffer, 5x RT Probe Reaction Buffer, 10x RT Enzyme Mix, UniSp6, RNA Spike-in template, RNase-Free Water	339340
RNA Spike-In Kit, For RT	Contains the UniSp2, UniSp4, and UniSp5 RNA Spike-in Template Mix and the cel-miR-39-3p RNA Spike-in Template	339390
miRCURY LNA SYBR® Green PCR Kits (200)	For 200 reactions: 2X miRCURY SYBR Green Master Mix, RNase-Free Water	339345
miRCURY LNA SYBR® Green PCR Kits (600)	For 600 reactions: 2X miRCURY SYBR Green Master Mix, RNase-Free Water	339346
miRCURY LNA SYBR® Green PCR Kits (4000)	For 4000 reactions: 2X miRCURY SYBR Green Master Mix, RNase-Free Water	339347
miRCURY LNA Probe PCR Kit (200)	For 200 reactions: 2X QuantiNova Probe Master Mix, 10X miRCURY Probe Univ. Primer, Rox Reference Dye, RNase-Free Water	339371
miRCURY LNA Probe PCR Kit (800)	For 800 reactions: 2X QuantiNova Probe Master Mix, 10X miRCURY Probe Univ. Primer, Rox Reference Dye, RNase-Free Water	339372
miRCURY LNA Probe PCR Kit (4000)	For 4000 reactions: 2X QuantiNova Probe Master Mix, 10X miRCURY Probe Univ. Primer, Rox Reference Dye, RNase-Free Water	339373
miRCURY LNA miRNA PCR Starter Kit	Two miRCURY LNA PCR Assays of your choice, spike-in control Assay (UniSp6), one candidate endogenous control assay (miR-103-3p) and all reagents for 20 reverse transcription reactions and 100 PCR amplifications; for SYBR® Green-based qPCR detection	339320

Important

The miRCURY LNA™ miRNA PCR Panels are Ready-to-Use and designed for optimal performance with the miRCURY LNA RT Kit and the miRCURY LNA SYBR® Green PCR Kit. The performance of the primer sets will be affected when used in combination with less than optimal reagents. Use the miRCURY LNA miRNA PCR Panel Handbook for experiment setup. RNA work requires specific handling and precautions should be taken to prevent RNase contamination and degradation of the RNA sample and reagents.

Shipping and storage

The Ready-to-Use miRCURY LNA miRNA PCR Panels are shipped at room temperature and can be stored at 4°C for at least 6 months. For long term storage, it is recommended to place the panels at -20°C. Under these conditions, the LNA PCR primers are stable for at least 6 months after receipt.

Intended use

miRCURY LNA miRNA PCR Assays are intended for molecular biology applications. This product is not intended for the diagnosis, prevention or treatment of a disease.

Safety information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). These are available online in convenient and compact PDF format at **www.qiagen.com/safety** where you can find, view and print the SDS for each QIAGEN kit and kit component.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, Sample to Insight®, GeneGlobe®, LNA™ (QIAGEN Group), SYBR® (Life Technologies Corporation). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law. 1107212 05/2017 HB-2376-001 © 2017 QIAGEN, all rights reserved.
